

Ti Application Note No. T- 49

Title:	Cationic surfactants in household cleaner by potentiometric two-phase titration
---------------	--

Summary:	Determination of cationic surfactants in household cleaner by potentiometric two-phase titration with sodium dodecylsulphate using the Surfactrode Resistant.
-----------------	---

Sample:	Household cleaner (containing ethanol and an excess of non-ionic surfactants)
----------------	---

Sample Preparation:	none
----------------------------	------

Instruments and Accessories:	702, 716, 736 or 751 Titrino or 726 Titroprocessor, 727 Titration Stand, 722 Propeller Rod Stirrer, 6.0507.130 Surfactrode Resistant
-------------------------------------	--

Analysis:	Weigh ca. 1 g sample (precision 0.1 mg) into a beaker, then add 10 mL each of ethanol and methyl isobutyl ketone (MIBK), ca. 70 mL dist. water and 10 mL buffer pH = 10.0. Make sure that the stirring speed is high enough to mix the solution thoroughly and form a stable emulsion, then titrate with sodium dodecylsulphate $c(\text{NaDDS}) = 0.004 \text{ mol/L}$.
------------------	---

Calculation:	$\text{RS1: } \text{meq cationics} / 100 \text{ g} = \text{EP1} * \text{C01} * \text{C02} / \text{C00}$ $\text{RS2: } \% \text{ cationics} = \text{RS1} * \text{C03}$ $\text{EP1} = \text{titrant consumption in mL}$ $\text{C00} = \text{ca. 1 (sample mass in g)}$ $\text{C01} = 0.004 \text{ (concentration of the titrant in mol/L)}$ $\text{C02} = 100 \text{ (conversion factor)}$ $\text{C03} = 0.350 \text{ (M(cationic surfactant) / 1000 in g/mol)}$
---------------------	--

Remarks:	Result: $\text{AVG}(4) = 1.320 \pm 0.01 \% \text{ cationics}$
-----------------	--