

Ti Application Note No. T- 66

Title:	Calcium in cheese
Summary:	Determination of calcium in cheese by potentiometric titration with EGTA using the Cu ISE.
Sample:	Any cheese
Sample Preparation:	If necessary grind the cheese before using it for analysis. Weigh approx. 1 g of it into a beaker, then add 50 mL dist. water and 10 mL $c(\text{H}_2\text{SO}_4) = 0.05 \text{ mol/L}$. Heat the mixture to 40°C and stir at this temperature for 10 min. After cooling down, pre-neutralise to about $\text{pH} = 7$ by adding $c(\text{NaOH}) = 0.1 \text{ mol/L}$.
Instruments and Accessories:	702, 716, 736, 751 or 785 Titrino or 726 Titroprocessor, 6.0502.140 Cu ISE, 6.0726.100 Ag/AgCl reference electrode
Analysis:	Add 4 mL CuEGTA complex solution (0.05 mol/L) and 10 mL $\text{NH}_3/\text{NH}_4\text{Cl}$ buffer $\text{pH} = 10$ to the prepared sample solution, then titrate with $c(\text{EGTA}) = 0.05 \text{ mol/L}$ using the MET mode.
Calculation:	$\% \text{Ca}^{2+} = \text{EP1} * \text{C01} * \text{C02} * \text{C03} / \text{C00}$ $\begin{aligned} \text{EP1} &= \text{titrant consumption in mL} \\ \text{C00} &= \text{approx. 1 (sample weight in g)} \\ \text{C01} &= 40.08 [\text{m}(\text{Ca}^{2+}) \text{ in g/mol}] \\ \text{C02} &= 0.05 (\text{concentration of the titrant in mol/L}) \\ \text{C03} &= 0.1 (\text{conversion factor for \%}) \end{aligned}$
Results:	$\text{AVG}(5) = 0.528 \pm 0.011 \% \text{Ca}^{2+}$ (for a Swiss cheese)
Remarks:	If the cheese sample is not heated and no sulphuric acid is added this leads to results that are too low.