

850 Professional IC



Anion – MCS – 2.850.2030

Manual

8.850.8046EN / 2019-11-28



Metrohm AG
CH-9100 Herisau
Switzerland
Phone +41 71 353 85 85
Fax +41 71 353 89 01
info@metrohm.com
www.metrohm.com

850 Professional IC

Anion – MCS – 2.850.2030

Manual

Technical Communication
Metrohm AG
CH-9100 Herisau
techcom@metrohm.com

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Table of contents

1	Introduction	1
1.1	Instrument description	1
1.2	Intended use	5
1.3	About the documentation	5
1.3.1	Symbols and conventions	5
1.4	Safety instructions	6
1.4.1	General notes on safety	6
1.4.2	Electrical safety	6
1.4.3	Tubing and capillary connections	7
1.4.4	Flammable solvents and chemicals	8
1.4.5	Recycling and disposal	8
2	Installation	9
2.1	About this chapter	9
2.2	Initial installation	9
2.3	Flow diagram	10
2.4	Setting up the instrument	12
2.4.1	Packaging	12
2.4.2	Checks	12
2.4.3	Location	12
2.5	Capillary connections in the IC system	13
2.6	Rear of the instrument	15
2.6.1	Rollers and handle	15
2.6.2	Positioning and connecting the detector	18
2.6.3	Transport locking screws	18
2.6.4	Leak sensor	18
2.6.5	Drainage tubings	19
2.7	Capillary and cable feed-throughs	21
2.8	Eluent	24
2.8.1	Connecting the eluent bottle	24
2.9	Eluent degasser	28
2.10	High pressure pump	30
2.10.1	Capillary connections high pressure pump/purge valve	30
2.10.2	Deaerating the high pressure pump	32
2.11	Inline filter	34
2.12	Pulsation absorber	35
2.13	Sample degasser	36

2.14	Injection valve	38
2.14.1	Connecting the injection valve	38
2.14.2	Mode of operation of the injection valve	39
2.14.3	Selecting the sample loop	40
2.15	Column thermostat	40
2.16	Metrohm Suppressor Module (MSM)	43
2.16.1	General information on the MSM	43
2.16.2	Connecting the MSM	43
2.17	Peristaltic pump	46
2.17.1	Principle of the peristaltic pump	46
2.17.2	Installing the peristaltic pump	47
2.18	Metrohm CO₂ Suppressor (MCS)	51
2.18.1	General information on the MCS	51
2.18.2	Connecting the MCS	51
2.18.3	Installing adsorption cartridges	52
2.19	Conductivity detector	55
2.20	Connecting the instrument to a computer	57
2.21	Connecting the instrument to the power grid	58
2.22	Guard column	59
2.23	Separation column	60
3	Start-up	63
3.1	Initial start-up	63
3.2	Conditioning	64
4	Operation and maintenance	66
4.1	General notes	66
4.1.1	Care	66
4.1.2	Maintenance by Metrohm Service	66
4.1.3	Operation	67
4.1.4	Shutting down	67
4.2	Capillary connections	67
4.2.1	Operation	67
4.3	Door	68
4.4	Eluent	68
4.4.1	Production	68
4.4.2	Operation	69
4.5	High pressure pump	69
4.5.1	Protection	69
4.5.2	Maintenance	70
4.6	Inline filter	80
4.6.1	Maintenance	80

4.7	Inline sample preparation	82
4.8	Rinsing the sample path	82
4.9	Sample degasser	84
4.9.1	Operation	84
4.10	Injection valve	84
4.10.1	Protection	84
4.11	Metrohm Suppressor Module (MSM)	84
4.11.1	Protection	84
4.11.2	Operation	84
4.11.3	Maintenance	85
4.12	Peristaltic pump	90
4.12.1	Operation	90
4.12.2	Maintenance	91
4.13	Metrohm CO₂ Suppressor (MCS)	93
4.13.1	Replacing the CO ₂ adsorption cartridge	93
4.13.2	Regenerating the H ₂ O adsorption cartridge	93
4.14	Conductivity detector	94
4.14.1	Maintenance	94
4.15	Separation column	95
4.15.1	Separating efficiency	95
4.15.2	Protection	95
4.15.3	Storage	95
4.15.4	Regeneration	96
5	Troubleshooting	97
5.1	Problems and their solutions	97
6	Technical specifications	102
6.1	Reference conditions	102
6.2	Device	102
6.3	Leak sensor	102
6.4	Ambient conditions	102
6.5	Housing	103
6.6	Eluent degasser	103
6.7	High pressure pump	103
6.8	Sample degasser	104
6.9	Injection valve	104
6.10	Column thermostat	105
6.11	Metrohm Suppressor Module (MSM)	105
6.12	Peristaltic pump	105

6.13	Metrohm CO ₂ Suppressor (MCS)	106
6.14	Conductivity measuring system	106
6.15	Power connection	107
6.16	Interfaces	107
6.17	Weight	108
7	Accessories	109
	Index	110

Table of figures

Figure 1	Front view of the instrument	2
Figure 2	Flow diagram with sequential suppression	11
Figure 3	Connection of capillaries with pressure screws	13
Figure 4	Rollers and handle	16
Figure 5	Handle as MPak holder	17
Figure 6	Connection of the leak sensor on the rear of the instrument	19
Figure 7	Drainage tubings	20
Figure 8	Capillary feed-throughs on the door	22
Figure 9	Capillary feed-throughs base tray/flask holder	23
Figure 10	Installing the eluent bottle cap	25
Figure 11	Mounting the aspiration filter	25
Figure 12	Installing the tubing weighting and aspiration filter	26
Figure 13	Eluent aspiration tubing fully equipped	26
Figure 14	Eluent bottle – connected	27
Figure 15	Eluent degasser	29
Figure 16	Capillary connections high pressure pump/purge valve	30
Figure 17	High pressure pump – Connect inlet	31
Figure 18	Deaerate the high pressure pump	33
Figure 19	Connecting the inline filter	35
Figure 20	Pulsation absorber – Connection	36
Figure 21	Sample degasser	37
Figure 22	Injection valve – connected	38
Figure 23	Injection valve – Positions	39
Figure 24	Column thermostat	41
Figure 25	MSM – connections	44
Figure 26	Peristaltic pump	46
Figure 27	Installing the pump tubing	47
Figure 28	Install pump tubing connection with filter	48
Figure 29	Install pump tubing connection without filter	49
Figure 30	MCS – Connection	51
Figure 31	Adsorption cartridge holder	53
Figure 32	Conductivity detector front	55
Figure 33	Conductivity detector rear	56
Figure 34	Connection detector – MCS	57
Figure 35	Pump head – removing the piston	71
Figure 36	Components of the piston cartridge	72
Figure 37	Tool for piston seal	73
Figure 38	Removing the piston seal	74
Figure 39	Inserting the piston seal into the tool	74
Figure 40	Inserting the piston seal into the pump head	75
Figure 41	Removing valves	76
Figure 42	Dismantling valve	77
Figure 43	Components of the inlet valve and outlet valve	78
Figure 44	Change filters (of the inline filter)	80
Figure 45	MSM – Components	87



Figure 46	Pump tubing connection – Changing the filter	92
-----------	--	----

1 Introduction

1.1 Instrument description

The instrument **850 Professional IC – Anion – MCS** (2.850.2030) is one of the model versions of the Professional IC line of instruments manufactured by the Metrohm Company. The Professional IC line of instruments is distinguished by

- the **intelligence** of its components, which are able to monitor and optimize all functions and to provide documentation according to FDA requirements.
- its **compactness**.
- its **flexibility**. A suitable model version exists for every application. Individual instruments can be converted, expanded or modified to create a different instrument as needed.
- its **transparency**. All components are easily accessible and arranged in a clear manner.
- its **safety**. Chemicals and electronics are separated and a leak sensor is integrated in the wet end.
- its **environmental compatibility**.
- **low noise emission**.

The instrument is operated with **MagIC Net** software. It is connected via a USB connection to a PC on which MagIC Net is installed. The software automatically recognizes the instrument and checks its functional readiness. MagIC Net controls and monitors the instrument, evaluates the measured data and administers it in a database. The operation of MagIC Net is described in the online help or in the tutorial for MagIC Net.




Figure 1 Front view of the instrument

1 Eluent bottle
See Chapter 2.8.1.

2 Conductivity detector
See Chapter 2.19.

3 MCS <i>See Chapter 2.18.</i>	4 Peristaltic pump <i>See Chapter 2.17.</i>
5 MSM <i>See Chapter 2.16.</i>	6 Injection valve <i>See Chapter 2.14.</i>
7 Eluent degasser <i>See Chapter 2.9.</i>	8 High pressure pump <i>See Chapter 2.10.</i>
9 Purge valve <i>See Chapter 2.10.1.</i>	10 Pulsation damper <i>See Chapter 2.12.</i>
11 Sample degasser <i>See Chapter 2.13. Use optional.</i>	12 Separation column <i>See Chapter 2.23.</i>
13 Column thermostat <i>See Chapter 2.15.</i>	

The instrument contains the following components:

Eluent degasser

The eluent degasser removes gas bubbles and dissolved gases from the eluent. For degassing, the eluent flows into a vacuum chamber through a special fluoropolymer capillary.

High pressure pump

The intelligent and low pulsation high pressure pump pumps the eluent through the system. It is equipped with a chip on which its technical specifications and "life history" (operating hours, service data, ...) are saved.

Inline filter

Inline filters protect the separation column securely against possible contamination from the eluent. Inline filters can however also just as well be used for the purpose of protecting other sensitive components against contaminations in the solutions used. The filter platelets with a pore size of 2 µm can be replaced quickly and easily. They remove particles like e. g. bacteria and algae from the solutions.

Pulsation absorber

The pulsation absorber protects the separation column from damage caused by pressure fluctuations when switching the injection valve, and reduces interfering pulsations during highly sensitive measurements.

Sample degasser

The sample degasser removes gas bubbles and dissolved gases from the sample. For degassing, the sample flows into a vacuum chamber through a special fluoropolymer capillary.



Injection valve

The injection valve connects the eluent and sample path through rapid and precise valve switchover. A precisely measured amount of sample solution is injected and rinsed with eluent onto the separation column.

Column thermostat

The column thermostat controls the temperature of the column and eluent channel and thus ensures stable measuring conditions. It provides space for 2 separation columns.

Metrohm Suppressor Module (MSM)

The MSM is used for chemical suppression during anion analysis. It is pressure-stable, robust and resistant to solvents.

Peristaltic pump

The Peristaltic pump is used for pumping sample and auxiliary solutions. It can rotate in both directions.

Metrohm CO₂ Suppressor (MCS)

The Metrohm CO₂ Suppressor (MCS) removes the CO₂ from the eluent stream. This lowers the background conductivity, improves detection sensitivity and minimizes the injection peak and carbonate peak.

Conductivity detector

The conductivity detector continuously measures the conductivity of the liquid passing through and indicates these signals in digital form (DSP – Digital Signal Processing). The conductivity detector exhibits outstanding thermal stability and thus guarantees reproducible measuring conditions.

Separation column

The intelligent separation column is the heart of the ion chromatographic analysis. It separates the different components corresponding to their interactions with the column. Metrohm separation columns are equipped with a chip on which their technical specifications and their history (first use / setting up, operating hours, injections, ...) are saved.

1.2 Intended use

The instrument **850 Professional IC – Anion – MCS** is used for ion chromatographic determination of anions or polar substances with **sequential suppression**:

- Chemical suppression with the Metrohm Suppressor Module (MSM) (*see chapter 2.16, page 43*) and subsequent
- CO₂ suppression with the Metrohm CO₂ Suppressor (MCS) (*see chapter 2.18, page 51*).

The use of sequential suppression reduces background conductivity to a minimum.

If required, the instrument can also be used for the determination of cations or anions without suppression.

The present instrument is suitable for processing chemicals and flammable samples. The usage of the 850 Professional IC – Anion – MCS therefore requires that the user has basic knowledge and experience in the handling of toxic and caustic substances. Knowledge with respect to the application of the fire prevention measures prescribed for laboratories is also mandatory.

1.3 About the documentation



CAUTION

Please read through this documentation carefully before putting the instrument into operation. The documentation contains information and warnings which have to be followed by the user in order to ensure safe operation of the instrument.

1.3.1 Symbols and conventions

The following symbols and styles are used in this documentation:

(5- 12)	Cross-reference to figure legend The first number refers to the figure number, the second to the instrument part in the figure.
1	Instruction step Carry out these steps in the sequence shown.

	Warning This symbol draws attention to a possible life hazard or risk of injury.
	Warning This symbol draws attention to a possible hazard due to electrical current.
	Warning This symbol draws attention to a possible hazard due to heat or hot instrument parts.
	Warning This symbol draws attention to a possible biological hazard.
	Caution This symbol draws attention to a possible damage of instruments or instrument parts.
	Note This symbol marks additional information and tips.

1.4 Safety instructions

1.4.1 General notes on safety



WARNING

Operate this instrument only according to the information contained in this documentation.

This instrument left the factory in a flawless state in terms of technical safety. To maintain this state and ensure non-hazardous operation of the instrument, the following instructions must be observed carefully.

1.4.2 Electrical safety

The electrical safety when working with the instrument is ensured as part of the international standard IEC 61010.



WARNING

Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.



WARNING

Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.

Supply voltage



WARNING

An incorrect supply voltage can damage the instrument.

Only operate this instrument with a supply voltage specified for it (see rear panel of the instrument).

Protection against electrostatic charges



WARNING

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Do not fail to pull the power cord out of the power socket before you set up or disconnect electrical plug connections at the rear of the instrument.

1.4.3 Tubing and capillary connections



CAUTION

Leaks in tubing and capillary connections are a safety risk. Tighten all connections well by hand. Avoid applying excessive force to tubing connections. Damaged tubing ends lead to leakage. Appropriate tools can be used to loosen connections.

Check the connections regularly for leakage. If the instrument is used mainly in unattended operation, then weekly inspections are mandatory.

1.4.4 Flammable solvents and chemicals




WARNING

All relevant safety measures are to be observed when working with flammable solvents and chemicals.

- Set up the instrument in a well-ventilated location (e.g. fume cupboard).
- Keep all sources of flame far from the workplace.
- Clean up spilled liquids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

1.4.5 Recycling and disposal



This product is covered by European Directive 2012/19/EU, WEEE – Waste Electrical and Electronic Equipment.

The correct disposal of your old instrument will help to prevent negative effects on the environment and public health.

More details about the disposal of your old instrument can be obtained from your local authorities, from waste disposal companies or from your local dealer.

2.1 About this chapter

- a list of all the steps to be carried out during the initial installation (*see chapter 2.2, page 9*).
- a schematic overview of the flow paths (*see chapter 2.3, page 10*).
- a description of the construction, connections and functionality of the instrument.
- step-by-step installation instructions. A part of these installation tasks has already been carried out before delivery. They are nevertheless described here in case they need to be carried out once again some time in the future (e.g. after maintenance of the high pressure pump).

 NOTICE


A large number of the capillary connections are already connected at the time the instrument is delivered.

- 1** Setting up the instrument (*see chapter 2.4, page 12*).
- 2** Remove handle and rollers (*see chapter 2.6.1, page 15*).
- 3** Place the detectors in the instrument and connect them (*see chapter 2.6.2, page 18*).
- 4** Remove transport locking screws (*see chapter 2.6.3, page 18*).
- 5** Connect the leak sensor (*see chapter 2.6.4, page 18*).
- 6** Connect drainage tubings (*see chapter 2.6.5, page 19*).
- 7** Connect the eluent bottle (*see chapter 2.8.1, page 24*).

- 8** Install sample path connections.
 - Connect the sample degasser (if required) (*see chapter 2.13, page 36*).
 - Connect sample path connections to the injection valve (*see chapter 2.14.2, page 39*).
- 9** Installing the MSM (*see chapter 2.16, page 43*) – with associated peristaltic pump (*see chapter 2.17, page 46*).
- 10** Connect the MCS (*see chapter 2.18.2, page 51*).
- 11** Connect detector capillaries (*see chapter 2.19, page 55*).
- 12** Mains connection .
- 13** Connect the instrument to the PC .
- 14** Initial start-up (*see chapter 3.1, page 63*).
- 15** Install the guard column (if used) (*see chapter 2.22, page 59*).
- 16** Install the separation column (*see chapter 2.23, page 60*).

2.3 Flow diagram

Figure 2 Flow diagram with sequential suppression shows the flow paths with the application of sequential suppression (MSM (2-**9**) and MCS (2-**10**)). The arrangement of the modules in the diagram corresponds to the front view of the instrument. Liquid containers (eluent bottle, sample vessel, waste container, auxiliary solutions container) and guard column (see chapter 2.22, page 59) are not shown in the diagram. Which pressure screws, connections and couplings are used, is described in the installation chapters of the individual modules.



2 Eluent degasser

See Chapter 2.9.

2.5 Capillary connections in the IC system

This chapter contains general information concerning the capillary connections in the IC instruments and systems.

Generally speaking, capillary connections between two components of an IC system are made up of one connection capillary and two pressure screws with which the capillary is connected to the respective components.

Pressure screws




Figure 3 Connection of capillaries with pressure screws

- | | |
|--|---|
| <p>1 PEEK pressure screw (6.2744.014)
Use on the injection valve.</p> | <p>2 Connection capillary</p> |
| <p>3 PEEK pressure screw, short (6.2744.070)
For use on the high pressure pump, the purge valve, the inline filter, the pulsation absorber, the guard column and the separation column.</p> | <p>4 PEEK pressure screw, long (6.2744.090)
Use on special components. Is not used on all instruments.</p> |



NOTICE

In order to keep the dead volume as low as possible, capillary connections should generally be as short as possible.



NOTICE

For an improved overview, capillary and tubing connections can be bundled with the 6.1815.010 spiral band.



Connection capillaries

PEEK capillaries and PTFE capillaries are used in the IC system.

PEEK capillaries (poly-
etheretherketone)

PEEK capillaries are temperature-resistant up to 100°C, stable under pressure up to 400 bar, flexible, chemically inert and exhibit an extremely smooth surface. They can be readily cut down to the desired length with the 6.2621.080 capillary cutter.

Usage:

- PEEK capillaries (6.1831.010) with an internal diameter of 0.25 mm for the entire high-pressure section.
- PEEK capillaries (6.1831.030) with an internal diameter of 0.75 mm for sample handling in the ultratrace range.



CAUTION

For the capillary connections between the injection valve and detector, PEEK capillaries with an internal diameter of 0.25 mm must be used. These are already connected to a newly delivered instrument.

PTFE capillaries (polytetrafluoroethylene)

PTFE capillaries are transparent and enable visual tracing of the liquids to be pumped. They are chemically inert, flexible and temperature-resistant up to 80°C.

Usage:

PTFE capillaries (6.1803.0x0) are used for the low-pressure section.

- PTFE capillaries with internal diameter of 0.5 mm for sample handling.
- PTFE capillaries with internal diameter of 0.97 mm for sample handling as well as for rinsing solutions (they do not have to be in the scope of delivery of the instrument).

Capillary connections

In order to achieve optimum analysis results, capillary connections in an IC system must be absolutely tight and free of dead volume. Dead volume occurs if two capillary ends connected to each other do not fit exactly, thus allowing liquid to escape. There are two possible reasons for this:

- The capillaries do not have exactly cut edges.
- The two capillary ends do not completely meet.

One prerequisite for dead volume free capillary connection is, that both capillary ends are cut exactly plane. Therefore we recommend only to cut PEEK capillaries with the capillary cutter (6.2621.080).




Figure 8 Capillary feed-throughs on the door

1 Luer connectors

For connecting a 6.2816.020 syringe. For manual sample feeding.

**3 PEEK pressure screws, short
6.2744.070**

2 Capillary feed-through

4 Door

Do not feed capillaries through the Luer connectors (8-**1**). The capillaries are fastened with PEEK pressure screws (8-**3**) from inside to the Luer connector. From outside, liquid can be aspirated or injected with a syringe.



1	Side panel (right) Right panel.	2	Rear of the instrument
3	Side panel (left) Left panel.	4	Front of the instrument
5	Capillary feed-through Upper. From front to right.	6	Capillary feed-through Upper. From front to back.
7	Capillary feed-through Upper. From front to back.	8	Capillary feed-through Upper. From front to left.
9	Capillary feed-through Lower. From front to right.	10	Capillary feed-through Lower. From front to back.
11	Capillary feed-through Lower. From front to back.	12	Capillary feed-through Lower. From front to left.

2.8 Eluent

2.8.1 Connecting the eluent bottle

The eluent is aspirated out of the eluent bottle via the eluent aspiration tubing (10-1).

The eluent aspiration tubing is connected to the eluent degasser (*see chapter 2.9, page 28*). The tubing must be threaded through a suitable capillary feed-through (*see chapter 2.7, page 21*) of the instrument before the other end can be equipped.


The following accessories are needed for equipping the eluent aspiration tubing:

- 6.1602.160 eluent bottle cap GL 45
- 6.2744.210 tubing adapter for aspiration filter
- 6.2821.090 aspiration filter

Proceed as follows for equipping the eluent aspiration tubing:

Assembling eluent aspiration tubing

- 1** Guide the free end of the eluent aspiration tubing (10-1) out of the instrument through a suitable capillary feed-through.
- 2 Installing the eluent bottle cap (6.1602.160)**
 - Slide tubing nipple (10-2) and O-ring (10-3) onto the eluent aspiration tubing (10-1).
 - Push eluent aspiration tubing (10-1) through the bottle cap (10-4) and screw tight.



- ### 3 Mounting the aspiration filter

-

Figure 11 Mounting the aspiration filter

- 25



4 Installing the tubing weighting and aspiration filter




Figure 12 Installing the tubing weighting and aspiration filter

1	Eluent aspiration tubing (6.1834.080)	2	Eluent bottle cap (6.1602.160)
3	Tubing weighting From accessory set (6.2744.210).	4	Clamping screw From accessory set (6.2744.210).
5	Aspiration filter (6.2821.090) With filter holder from accessory set (6.2744.210).		

- Slide the tubing weighting (12-**3**) onto the eluent aspiration tubing (12-**1**).
- Slide the clamping screw (12-**4**) onto the eluent aspiration tubing (12-**1**).
- Insert eluent aspiration tubing (12-**1**) into the aspiration filter (12-**5**). The end of the tubing should reach approximately to the center of the aspiration filter.
- Screw the clamping screw (12-**4**) to the filter holder (11-**1**).




Figure 13 Eluent aspiration tubing fully equipped

5 Mounting the eluent aspiration tubing to the eluent bottle

- Insert the eluent aspiration tubing into the eluent bottle (**14-10**).

- ## 6 Mounting the adsorber tube



- First, place a piece of cotton (14-3), then the CO₂ adsorber (14-4) in the large opening of the adsorber tube (14-2) and close with the plastic cover.
- Fasten the adsorber tube (14-2) on the bottle cap (14-11) using the SGJ clip (14-12).




Figure 14 Eluent bottle – connected

- 2 Adsorber tube (6.1609.000)**

3	Cotton	4	CO₂ adsorber Adsorbs CO ₂ from the air (e.g. Merck soda lime with indicator, no. 6839.1000).
5	Eluent	6	Aspiration filter (6.2821.090)
7	Filter holder From accessory set (6.2744.210).	8	Clamping screw From accessory set (6.2744.210).
9	Tubing weighting From accessory set (6.2744.210).	10	Eluent bottle (6.1608.070)
11	Bottle cap (6.1602.160)	12	SGJ clip (6.2023.020)
13	Tubing nipple	14	Threaded stopper

2.9 Eluent degasser

Gas bubbles in the eluent lead to an unstable baseline, as high pressure pumps can transport liquids, but not gases. This is why the eluent must be degassed before it reaches the high-pressure pump.

The eluent degasser removes gas bubbles and dissolved gases from the eluent. For degassing, the eluent flows into a vacuum chamber through a special fluoropolymer capillary.



NOTICE

The eluent degasser is already installed in the newly delivered instrument. The following installation instructions must only be followed, if the connections to the degasser had to be disconnected for maintenance.

The diagram illustrates the Eluent Degasser assembly. On the left, the degasser unit is shown with two ports: 'In' (top) and 'Out' (bottom). Port 1 is connected to the 'In' port, and Port 2 is connected to the 'Out' port. On the right, the HPLC system is shown with two ports: Port 3 (top) and Port 4 (bottom). Port 3 is connected to the 'In' port of the degasser, and Port 4 is connected to the 'Out' port of the degasser. The degasser unit is labeled 'Eluent Degasser'.

Figure 15 Eluent degasser

1	Eluent degasser input	2	Eluent degasser output
3	Tubing flare With tubing nipple.	4	Clamping screw
5	Eluent aspiration tubing (6.1834.080) For aspirating the eluent. The clamping screw (15- 4) is firmly mounted.	6	Connection tubing (6.1834.090) Connection from the eluent degasser to the high-pressure pump (<i>see chapter 2.10, page 30</i>). The clamping screw (15- 4) is firmly mounted.

CAUTION

The clamping screws (15-4) must be tightened carefully. Use the wrench (6.2621.050) to do this.

- Insert the eluent aspiration tubing (15-5) into the eluent degasser input (15-1).
- Carefully tighten the clamping screw (15-4).

- Insert connection tubing (15-6) (the end with the longer clamping screw (15-4)) into the eluent degasser output (15-2).
- Carefully tighten the clamping screw (15-4).
- Connect the other end of the connection capillary (15-6) (with the shorter clamping screw) to the high-pressure pump (16-9) (see "Connecting inlet to the high pressure pump", page 31).

2.10 High pressure pump

The intelligent and low pulsation high pressure pump pumps the eluent through the system. It is equipped with a chip on which its technical specifications and "life history" (operating hours, service data, ...) are saved.

The purge valve is used for deaerating (see chapter 2.10.2, page 32) the high pressure pump.

2.10.1 Capillary connections high pressure pump/purge valve



NOTICE

All of the capillary connections of the high pressure pump and the purge valve are already installed in the newly delivered instrument.




Figure 16 Capillary connections high pressure pump/purge valve

- | | |
|---|--|
| 1 Connection capillary
PEEK capillary, connects main piston and auxiliary piston. | 2 PEEK pressure screw, short (6.2744.070) |
| 3 Outlet valve holder | 4 Pump head (6.2824.110) |




Figure 19 Connecting the inline filter

1 Connection capillary

Connects the purge valve with the inline filter

3 Inline filter (6.2821.120)

Protects against particles.

2 PEEK pressure screws, short (6.2744.070)

4 Connection capillary

Connects the inline filter with the pulsation absorber.

- 1** Screw on the connection capillary running from the purge valve to the input side of the inline filter using a pressure screw (6.2744.070).
- 2** Screw on the connection capillary running to the pulsation absorber to the output side of the inline filter using a pressure screw (6.2744.070).

2.12 Pulsation absorber



NOTICE

The pulsation absorber is already installed in the newly delivered instrument.



CAUTION

The pulsation absorber is maintenance-free and may not be opened.

The pulsation absorber protects the separation column from damage caused by pressure fluctuations when switching the injection valve, and reduces interfering pulsations during highly sensitive measurements. In order to ensure these functionalities, it must be connected between the high pressure pump (see chapter 2.10, page 30) and injection valve (see chapter 2.14, page 38).

The pulsation absorber can be operated in both directions.

Figure 27 Installing the pump tubing

1	PEEK pressure screws, short (6.2744.070)	2	Tubing olive (6.2744.034)
3	Stopper The colors of the stopper indicate the inner diameter of the pump tubing.	4	Tubing cartridge (6.2755.000)
5	Contact pressure lever	6	Union nut
7	Adapter	8	Tubing olive Either with filter holder (6.2744.180) or without filter holder (6.2744.160).
9	Pump tubing (6.1826.xx0)	10	Snap-action lever

Mount the pump tubing as follows:

1 Removing the tubing cartridge

Release the tubing cartridge from the cartridge holder by pressing the snap-action lever and unhooking from the mounting pins (26-**1**).

2 Connecting the aspiration side

Place a 6.2744.034 tubing olive (27-**2**) on the aspiration side of the pump tubing.

3 Connecting the pressure side



NOTICE

Depending on the use of the peristaltic pump, on the pressure side you can either connect:

- **Case A:** a 6.2744.180 pump tubing connection **with filter**
(see figure 28, page 48) or
- **Case B:** a 6.2744.160 pump tubing connection **without filter**
(see figure 29, page 49).

For pumping the auxiliary solutions to the MSM or to the SPM, a 6.2744.180 pump tubing connection **with filter must** be used.

Case A: 6.2744.180 pump tubing connection with filter:




Figure 28 Install pump tubing connection with filter

1	Union nut	2	Adapter
3	Tubing olive with filter holder		

- Slide union nut (28-**1**) onto the pump tubing.
- Select a suitable adapter (28-**2**) and slide it onto the pump tubing. The type of adapter depends on the pump tubing (*see table 1, page 49*).
- Place the tubing olive with filter holder (28-**3**) onto the pump tubing.
- Screw the union nut (28-**1**) onto the tubing olive (28-**3**).

or

Case B: 6.2744.160 pump tubing connection without filter:



2 Adapter

- #### 4 Inserting the pump tubing

- ## 5 Inserting the tubing cartridge

- ## 6 Connecting the capillaries

- Table 1 Pump tubings and suitable adapters

850 Professional IC – Anion – MCS



- | | |
|---|---|
| <p>3 Aspiration capillary
For drawing in air with low CO₂ content (via the CO₂ adsorption cartridge (31-4)).</p> | <p>4 PEEK pressure screw, long (6.2744.090)</p> |
| <p>5 Connection capillary</p> | <p>6 Pressure screw, short (6.2744.070)
Installed on the air aspiration capillary.</p> |
| <p>7 Luer coupling (6.2744.120)
Mounted on the air aspiration capillary with a pressure screw (6.2744.070).</p> | |

Connecting the MCS

1 Connection from MSM

Use a long PEEK pressure screw (6.2744.090) (30-4) to connect the eluent outlet capillary (labeled **out**) to the input of the MCS (30-1).

2 Connection to the detector

Connect the detector input capillary (32-**3**) with a long PEEK pressure screw (6.2744.090) (30-**4**) to the MCS output (30-**2**).



CAUTION

If the MCS is not used, the inlet and outlet must be sealed with the stoppers (6.2744.220).

2.18.3 Installing adsorption cartridges

For an effective CO₂ removal, the air aspirated through the degassing cell should have as low a CO₂ content as possible. In order to achieve this, the air is aspirated through a CO₂ adsorption cartridge (6.2837.000) (31-4).

Moisture can block the CO₂ adsorption cartridge. To prevent this, a H₂O adsorption cartridge (6.2837.010) (31-7) is used first.



- 850 Professional IC – Anion – MCS

Installing the adsorption cartridges

1 Preparing the adsorption cartridge holder

Push the 4 clips (31-**3**) into the slots on the adsorption cartridge holder (31-**8**).

2 Removing the caps

- Remove the sealing caps at the top of both cartridges.
 - Replace the round sealing cap at the larger end of the H₂O adsorption cartridge with the star-shaped sealing cap.
- Important!** A small stopper is located in the center of the star-shaped sealing cap (on the air inlet (31-10)). It also has to be removed (*see leaflet for H₂O adsorption cartridge*).

3 Connecting the CO₂ adsorption cartridge

- Insert the CO₂ adsorber cartridge into the coupling (31-2) at the end of the MCS aspiration capillary .
- Click the CO₂ adsorption cartridge into both lower clips (31-3) of the adsorption cartridge holder (31-8).

4 Connecting the PVC tubing

- Insert the adapter (31-5) into the CO₂ adsorption cartridge.
- Fasten the PVC tubing (31-6) on the adapter (31-5).

5 Connecting the H₂O adsorption cartridge

- Insert the H₂O adsorption cartridge into the PVC tubing (31-**6**).
- Click the H₂O adsorption cartridge into both upper clips (31-**3**) of the adsorption cartridge holder (31-**8**).

6 Putting the adsorption cartridge holder into the instrument

- Place the adsorption cartridge holder with cartridges into the instrument's detector chamber.

2.21 Connecting the instrument to the power grid



WARNING

Electric shock from electrical potential

Risk of injury by touching live components or through moisture on live parts.

- Never open the housing of the instrument while the power cord is still connected.
- Protect live parts (e.g. power supply unit, power cord, connection sockets) against moisture.
- Unplug the power plug immediately if you suspect that moisture has gotten inside the instrument.
- Only personnel who have been issued Metrohm qualifications may perform service and repair work on electrical and electronic parts.

Connecting the power cord

Accessories

Power cord with the following specifications:

- Length: max. 2 m
- Number of cores: 3, with protective conductor
- Instrument plug: IEC 60320 type C13
- Conductor cross-section 3x min. 0.75 mm² / 18 AWG
- Power plug:
 - according to customer requirement (6.2122.XX0)
 - min. 10 A



NOTICE

Do not use a not permitted power cord!

1 Plugging in the power cord

- Plug the power cord into the instrument's power socket.
- Connect the power cord to the power grid.

2.22 Guard column

The use of guard columns serves for protecting the separation columns and increasing their service life considerably. The guard columns available from Metrohm represent either actual guard columns or are so-called guard column cartridges which are used together with a cartridge holder. The installation of a guard column cartridge in the associated holder is described in the leaflet of the guard columns.



Information regarding which guard column is suitable for your separation column can be found in the **Metrohm IC Column Program** (which is available from your Metrohm agent), the leaflet provided along with your separation column, the product information on the separation column at <http://www.metrohm.com> (product area Ion Chromatography), or obtained directly from your agent.



New guard columns are filled with solution and are sealed on both sides with stoppers or caps, respectively. Before using the guard column, you need to ensure that this solution is miscible with the eluents used (observe manufacturer's data).



The guard column may only be installed after the **initial start-up** (see *chapter 3.1, page 63*) of the instrument. Until then, use the coupling (6.2744.040) instead of the guard and separation column.



Metrohm recommends, always to work with guard columns. These protect the separation column and can be exchanged regularly.

Connecting and rinsing the guard column

1 Connecting the guard column



CAUTION

When inserting the guard column, always ensure that it is inserted correctly corresponding to the flow direction (if indicated).

- Remove sealing caps and/or stoppers from the guard column.
- Fasten the input of the guard column to the column inlet capillary using a short PEEK pressure screw (6.2744.070).
- In case the guard column is mounted on the separation column with a connecting capillary, connect the included connection capillary, which is included with the guard column, to the output of the guard column using the PEEK pressure screw, which is also included.

2 Rinsing the guard column

- Place beaker under the outlet capillary of the guard column.
- Set the flow rate of the high pressure pump according to the data given in the leaflet of the separation column.
- Start the high pressure pump and rinse the guard column approx. 5 minutes with eluent.
- Switch off the high pressure pump again.

2.23 Separation column

The intelligent separation column (iColumn) is the heart of the ion chromatographic analysis. It separates the different components corresponding to their interactions with the column. Metrohm separation columns are equipped with a chip on which their technical specifications and their history (first use / setting up, operating hours, injections, ...) are saved.



NOTICE

Information regarding which separation column is suitable for your application can be found in the **Metrohm IC Column Program**, the product information for your separation column at <http://www.metrohm.com> in the product area ion chromatography, or obtained directly from your agent.



2 Rinsing the separation column

- Place beaker under the outlet end of the separation column.
- Set the flow rate of the high pressure pump according to the data given in the leaflet of the separation column.
- Start the high pressure pump and rinse the separation column approx. 10 minutes with eluent.
- Switch off the high pressure pump again.

3 Mounting the separation column

- Fasten the column output capillary to the output of the separation column using a PEEK pressure screw (6.2744.070).
- Hang separation column with chip into the column holder.



NOTICE

The iColumns are equipped with a chip on which their operating data is saved. The chip has to be hooked into the chip holder provided for this so that the column recognition can function.

3 Start-up

The chapter *Start-up* is divided into 2 sections:

Initial start-up

The **initial start-up** is carried out during the **initial installation**.

Conditioning

Conditioning is carried out as a final installation step and each time after the system is started.

3.1 Initial start-up

The initial start-up is carried out during the initial installation. The entire system is rinsed before guard column and separation column are installed.



CAUTION

The separation column and guard column may not be installed for the initial start-up.

Make sure that the coupling (6.2744.040) is being used instead of the columns.

Perform the following steps during the initial start-up:

1 Preparing the software

- Start the PC program **MagIC Net™**.
- Open the **Equilibration** tab in MagIC Net™.
- Select (or create) a suitable method.

2 Preparing the instrument

- Ensure that the eluent aspiration tubing is immersed in the eluent and that there is enough eluent in the eluent bottle.
- Ensure that the aspiration tubings for the auxiliary solutions (regeneration solution and rinsing solution) are immersed into the respective solutions and that there is enough solution in the bottles.
- Switch on the instrument.

3 Starting equilibration

- In MagIC Net™, start the equilibration.



4 Deaerate the high pressure pump

- Deaerate the high pressure pump(s) via the purge valve (*see chapter 2.10.2, page 32*).

5 Set the contact pressure of the peristaltic pump



NOTICE

This work step needs to be performed only if a peristaltic pump is being used.

- If peristaltic pumps are used, set the contact pressure (see "Set flow rate", page 50).

6 Rinsing the instrument without columns

- Rinse the instrument (without columns) with eluent for 5 minutes.

The instrument is now ready for the installation of the columns (*see chapter 2.22, page 59*).

3.2 Conditioning

After the installation and after switching on the instrument, the system must be conditioned with eluent until a stable baseline is reached.



NOTICE

After a change of eluent (see chapter 4.4.2.3, page 69), the conditioning time can lengthen considerably.

Conditioning the system

1 Preparing the software



CAUTION

Ensure that the flow set is not higher than the flow permissible for the corresponding column (see column leaflet and chip data set).

- Start the PC program **MagIC Net™**.
- Open the **Equilibration** tab in MagIC Net™.

- Select (or create) a suitable method.

2 Preparing the instrument

- Ensure that the column is correctly mounted according to the flow direction indicated on the label (arrow must point in the direction of flow).
- Ensure that the eluent aspiration tubing is immersed in the eluent and that there is enough eluent in the eluent bottle.
- Ensure that the aspiration tubings for the auxiliary solutions (regeneration solution and rinsing solution) are immersed into the respective solutions and that there is enough solution in the bottles.

3 Checking leak-tightness

- In MagIC Net™, start the equilibration.
- Check all capillaries and their connections from the high pressure pump to the detector for signs of liquid escaping.
If eluent escapes anywhere, tighten the corresponding pressure screw or loosen the connection, check the end of the capillary, shorten it with a capillary cutter if necessary, and restore the connection.

4 Conditioning the system

Rinse the system with eluent until the required stability of the baseline is attained (normally 30 minutes).

During this time, step the MSM to the next position every 10 minutes.

The instrument is now ready for measuring samples.

4.4.2 Operation

4.4.2.1 Supply bottle

The supply bottle with the eluent must be connected as indicated in *chapter 2.8.1, page 24*. This is particularly important for eluents with volatile solvents (e.g. acetone).

Condensation in the eluent bottle must be avoided. Drop formation can change the concentration ratio in the eluent.

In the case of very sensitive measurements, we recommend that the eluent be stirred constantly with a magnetic stirrer (e. g. the 2.801.0010 with 6.2070.000).

4.4.2.2 Aspiration filter

To protect the IC system from foreign particles, we recommend aspirating the eluent via an aspiration filter (6.2821.090) (**11-2**). This aspiration filter must be replaced if it takes on a yellowish discoloration, and at least every 3 months.

4.4.2.3 Changing the eluent

Prevent any precipitations when changing the eluent. Immediately successive solutions must always be miscible. If the system has to be rinsed with organic solvents, use several solvents with rising or falling lipophilicity.

4.5 High pressure pump

4.5.1 Protection



CAUTION

The pump head is filled ex works with methanol/ultrapure water. It must be ensured that the eluent used is freely miscible with the solvent remaining in the pump head.

To protect the high pressure pump against **foreign particles**, we recommend that the eluent undergoes a **microfiltration** (filter 0.45 µm) before being aspirated via a 6.2821.090 aspiration filter (*see "Assembling eluent aspiration tubing", page 24*).

Salt crystals between the piston and seal cause abrasion particles which can find their way into the eluent. These lead to contaminated valves, a rise in pressure and in extreme cases scratched pistons. It is therefore essential to ensure that **no precipitates** can occur (*see chapter 4.4.2.3, page 69*).



- Remove the backup ring from the pump head and lay to the other parts.




Figure 36 Components of the piston cartridge

1	Piston cartridge screw	2	Retaining washer
3	Zirconium oxide piston with piston shaft Order number: 6.2824.070	4	Spring retainer
5	Spring Order number: 6.2824.060	6	Inner plastic sleeve Protects from metallic abrasion.
7	Piston cartridge	8	Backup ring

3 Cleaning the components of the piston

- Clean zirconium oxide pistons contaminated by abrasion or deposits with fine abrasive cleaning powder, rinse particle free with ultrapure water and dry.
Replace highly contaminated or scratched zirconium oxide pistons (spare part: 6.2824.070 zirconium oxide piston).
- Rinse the other parts of the piston and dry with a lint-free cloth.

4 Assembling the piston

- Insert the inner plastic sleeve, spring and spring retainer into the piston cartridge.
- Slide the zirconium oxide piston carefully into the piston cartridge until its tip emerges from the small opening of the piston cartridge.
- Attach screw and tighten by hand.

Replacing the piston seal

The special tool (6.2617.010) (see figure 37, page 73) is necessary to remove the piston seal from the pump head. It consists of two parts: a tip




Figure 43 Components of the inlet valve and outlet valve

1	Inlet valve (6.2824.170)	2	Outlet valve (6.2824.160)
3	Inlet valve housing	4	Outlet valve housing
5	Sealing ring (black)	6	Sleeve
7	Sapphire sleeve The shiny side must point to ruby ball.	8	Ruby ball
9	Ceramic holder for ruby ball	10	Seal The larger opening must point outwards.

4 Cleaning the components of the valve

Rinse the valve components with ultrapure water and/or acetone and dry with a lint-free cloth.

5 Reassembling the valve

Reassemble valve components *according to figure 43, page 78*.

- Insert the seal with the larger opening faced downwards into the recess of the tool.
- Lay the other valve components above another in the correct sequence (*see figure 43, page 78*).

- ## 6 Checking the flow direction

If this is not the case, the valve has to be dismantled again and be reassembled correctly (*see figure 43, page 78*).

If by mistake, the inlet valve is mounted instead of the outlet valve, an extreme pressure builds up within the working cylinder, which can destroy the piston seal!

- Insert the inlet valve into the inlet valve holder the way the seal is visible.
- Screw the inlet valve holder into the bottom of the pump head and tighten with a wrench (41-4).
- Insert the outlet valve into the outlet valve holder the way the seal is visible.
- Screw the outlet valve holder into the top of the pump head and tighten with a wrench (41-1).

Changing the filter

The flow must be stopped before changing the filter.

1 Removing the inline filter

- Unscrew the pressure screws (44-1) from the inline filter.

2 Unscrewing the filter screw

- Screw the filter screw (44-4) out of the filter housing (44-2) with the aid of two adjustable wrenches (6.2621.000).

3 Inserting the filter

- Remove the old filter (44-3) using tweezers.
- Place the new filter (44-3) flat in the filter housing (44-2) using tweezers.

4 Mounting filter screw

- Screw the filter screw (44-4) back into the filter housing (44-2) and tighten by hand. Then additionally tighten slightly with two adjustable wrenches (6.2621.000).

5 Remounting the inline filter

- Screw pressure screws (44-**1**) back onto the inline filter.

6 Rinsing the inline filter

- Dismantle the guard column (if present) and the separation column and replace with a coupling (6.2744.040).
- Rinse the instrument with eluent.

To ascertain the transfer time, proceed as follows:

Pump air through the sample path (pump tubing, tubing connections, capillary in the degasser, sample loop) for several minutes until all liquid is displaced by the air.

Aspirate a sample typical for the later application and use a stop watch to measure the time required by the sample to travel from the sample vessel to the end of the sample loop.

The time measured corresponds to the "transfer time". The rinsing time should be at least 3 times the transfer time.

It is possible to determine whether the rinsing time is adequate via a direct measurement of the sample carry-over. Proceed as follows:

- **Sample A:** A typical sample for the application.
- **Sample sample B:** Ultrapure water.

Let "Sample A" pass through the sample path for the duration of the rinsing time, then inject and measure.

Let "Sample B" pass through the sample path for the duration of the rinsing time, then inject and measure.

The degree of the sample carry-over corresponds to the ratio of the peak areas of the measurement for sample B to the measurement for sample A. The lower the ratio, the lower the sample carry-over. This ratio can be modified by varying the rinsing time – thus allowing the rinsing time required for the application to be ascertained.

In order to record every new chromatogram under comparable conditions, you should normally work with a freshly regenerated suppressor.



The MSM must never be switched over in a dry state, as there is a risk of jamming. If the MSM is in a dry state, the MSM must be rinsed for at least 5 minutes before it may be switched over.



In the case of reduced capacity or high backpressure, the MSM must be regenerated (see chapter 4.11.3.1, page 85), cleaned (see chapter 4.11.3.2, page 86) or replaced (see chapter 4.11.3.3, page 89).

4.11.3 Maintenance

4.11.3.1 Regenerating the MSM

If the suppressor units are loaded for a longer period with certain heavy metals (e.g. iron) or organic impurities, these can no longer be completely removed with the regeneration solution (50 mmol/L H₂SO₄) normally used. The capacity of the suppressor units is consequently affected, which can result in reduced phosphate sensitivity in less serious cases and a significant rise in the baseline in more serious cases. If such capacity problems occur at one or more positions, the suppressor units must be regenerated:

Regenerating the MSM

Regenerate the MSM as follows:

1 Disconnecting the MSM from the IC system

- Disconnect the MSM from the separation column and detector.

2 Regenerating the MSM



The pump tubing made of PVC must not be used for rinsing with solutions containing organic solvents. In this case, other pump tubing must be used for rinsing.



NOTICE

The high pressure pump can be used for regeneration. For this, remove the guard column and separation column and connect the capillary directly to the MSM (regenerate in opposite direction).

- Rinse the 3 suppressor units with the following solutions for about 15 minutes each:
 - **Contamination with heavy metals:**
1 mol/L H_2SO_4 + 0.1 mol/L oxalic acid
 - **Contamination with organic cationic complexing agents:**
0.1 mol/L H_2SO_4 / 0.1 mol/L oxalic acid / acetone 5%
 - **Severe contamination with organic substances:**
0.2 mol/L H_2SO_4 / acetone $\geq$ 20%

3 Connecting the MSM to the IC system

- Reconnect the MSM to the IC system. If the capacity problems remain, MSM rotor A must be replaced (*see chapter 4.11.3.3, page 89*).

4.11.3.2 Cleaning the MSM

It may be necessary to clean the MSM in the following cases:

- Increased backpressure onto the output capillaries of the MSM.
- Blockage of the MSM which cannot be eliminated (solutions can no longer be pumped through the MSM).
- Jamming of the MSM which cannot be eliminated (MSM can no longer be switched over).




Figure 45 MSM – Components

1	Union nut	2	MSM connecting piece 6.2832.010
3	MSM rotor A 6.2832.000	4	MSM housing
5	Slot in the MSM housing		

Cleaning the MSM

Clean the MSM as follows:

1 Disconnecting the MSM from the IC system

- Switch off the instrument.
- Disconnect the MSM from the separation column, peristaltic pump and detector.

2 Dismantling the MSM

- Unscrew union nut (45-**1**) from the MSM housing (45-**4**).
- Pull MSM connecting piece (45-**2**) and MSM rotor A (45-**3**) out of the MSM housing (45-**4**). The MSM connecting piece and MSM rotor A normally stick to one another - if this is not the case: Take a sharp object, insert into the slot (45-**5**) in the MSM housing, and pull out MSM rotor A (45-**3**) in this way.
- Detach the MSM connecting piece (45-**2**) from MSM rotor A (45-**3**).

3 Cleaning the supply and discharge lines

- Connect in turn each of the 6 capillary tubings fastened on the MSM connecting piece (45-**2**) on the high pressure pump (see chapter 2.10, page 30) and pump through ultra pure water.



- Check whether solution emerges at the MSM connecting piece (45-2). If one of the supply or discharge lines remained blocked, the SPM connecting piece (45-2) must be replaced (order number 6.2832.010).

4 Cleaning the MSM rotor A

- Clean sealing surface of the MSM rotor A (45-**3**) with ethanol using a lint-free cloth.

5 Inserting the MSM rotor A



CAUTION

An incorrectly inserted MSM rotor A (45-**3**) can be **destroyed** during start-up.

- Insert the MSM rotor A (45-3) into the MSM housing (45-4) in such a way that the tubing connections on the rear of MSM rotor A fit into the corresponding recesses inside the MSM housing and one of the three holes of the MSM rotor A is visible from below in the slot of the MSM housing (45-5).
- If the MSM rotor A (45-3) is correctly inserted, its sealing area will be approx. 4 mm within the MSM housing (45-4). If this is not the case, MSM rotor A must be moved into the right position from below using a sharp object (e.g. screwdriver).

6 Cleaning the MSM connecting piece

- Clean sealing surface of the MSM connecting piece (45-2) with ethanol using a lint-free cloth.

7 Inserting the MSM connecting piece

- Insert the MSM connecting piece (**45-2**) into the MSM housing (**45-4**) in such a way that the connector **1** is on top and the three pins of the MSM connecting piece fit into the corresponding recesses on the MSM housing (**45-4**).

8 Connecting and conditioning the MSM

- Reconnect the MSM to the IC system.
- Before switching the MSM over for the first time, rinse the three suppressor units with solution for 5 minutes.

4.11.3.3 Replacing parts of the MSM

It may be necessary to replace parts of the MSM in the following cases:

- Loss of suppression capacity which cannot be eliminated (reduced phosphate sensitivity and/or significant rise in the baseline).
- Blockage of the MSM which cannot be eliminated (solutions can no longer be pumped through the suppressor).

Both the MSM rotor A (45-3) and the MSM connecting piece (45-2) with its supply and discharge lines can be replaced.

Replacing parts of the MSM

Replace parts of the MSM as follows (see figure 45, page 87):

1 Disconnecting the MSM from the IC system

- Switch off the instrument.
- Disconnect the MSM from the separation column, peristaltic pump and detector.

2 Dismantling the MSM

- Unscrew union nut (45-1) from the MSM housing (45-4).
- Pull MSM connecting piece (45-2) and MSM rotor A (45-3) out of the MSM housing (45-4). The MSM connecting piece and MSM rotor A normally stick to one another - if this is not the case: Take a sharp object, insert into the slot (45-5) in the MSM housing, and pull out MSM rotor A (45-3) in this way.
- Detach the MSM connecting piece (45-2) from MSM rotor A (45-3).

3 Cleaning the new MSM rotor A

- Clean the sealing surface of the new MSM rotor A (45-3) with ethanol using a lint-free cloth.

4 Inserting the new MSM rotor A



CAUTION

An incorrectly inserted MSM rotor A (45-3) can be **destroyed** during start-up.



- Insert the new MSM rotor A (45-3) in the MSM housing (45-4) in such a way that the tubing connections on the rear of MSM rotor A fit into the corresponding recesses inside the MSM housing and one of the three holes of MSM rotor A is visible from below in the slot (45-5) of the MSM housing.
- If the MSM rotor A (45-3) is correctly inserted, its sealing area will be approx. 4 mm within the MSM housing (45-4). If this is not the case, MSM rotor A must be moved into the right position from below using a sharp object (e.g. screwdriver).

5 Cleaning the MSM connecting piece

- Clean the sealing surface of the MSM connecting piece (45-2) with ethanol using a lint-free cloth.

6 Inserting the new MSM connecting piece

- Insert the MSM connecting piece (45-2) into the MSM housing (45-4) in such a way that the connector 1 is on top and the three pins of the MSM connecting piece fit into the corresponding recesses on the MSM housing (45-4).

7 Connecting and conditioning the MSM

- Reconnect the MSM to the IC system.
- Before switching the MSM over for the first time, rinse the three suppressor units with solution for 5 minutes.

4.12 Peristaltic pump

4.12.1 Operation

The pumping capacity of the peristaltic pump depends on the drive speed (set via software), the contact pressure and, above all, the internal diameter of the pump tubing. Depending on the application, different pump tubings are used.



CAUTION

The service life of the pump tubings also depends on the contact pressure. Therefore fully lift the tubing cartridges by loosening the snap-action lever (27-**10**) on the right-hand side if the peristaltic pump is to be turned off for a longer period. Once set, the contact pressure remains unaffected.



 **CAUTION**

The 6.1826.xxx pump tubings consist of PVC or PP and therefore must not be used for rinsing with solutions containing acetone. In this case, use other pump tubings or use another pump for rinsing.

4.12.2 Maintenance

4.12.2.1 Pump tubings

Pump tubings used for the peristaltic pump are consumables with a limited lifetime.

LFL pump tubings with three stoppers are inserted into the tubing cartridge in such a way that the cartridge is located between two stoppers. This results in two possible positions for the tubing cartridge. If the pump tubing was to show pronounced signs of wear, you can insert it a second time in the other position.

Replace the pump tubings periodically, approx. every four weeks if used continuously.

Pump tubing selection

Pump tubing can differ in terms of material, diameter and thus flow rate. Different pump tubing is used depending on the application.

The following table shows the properties and areas of application of the pump tubings:

Table 2 Pump tubings

Order number	Name	Material	Inner diameter	Use
6.1826.020	Pump tubing (blue/blue), 2 stoppers	PVC (Tygon® ST)	1.65 mm	Pump tubing for online IC instruments and automation in voltammetry.
6.1826.310	Pump tubing LFL (orange/green), 3 stoppers	PVC (Tygon®)	0.38 mm	Pump tubing for bromate determination using the triiodide method.
6.1826.320	Pump tubing LFL (orange/yellow), 3 stoppers	PVC (Tygon®)	0.48 mm	For acceptor solutions for Inline Dialysis and for Inline Ultrafiltration.
6.1826.330	Pump tubing LFL (orange/white), 3 stoppers	PVC (Tygon®)	0.64 mm	No special applications.

Order number	Name	Material	Inner diameter	Use
6.1826.340	Pump tubing LFL (black/black), 3 stoppers	PVC (Tygon®)	0.76 mm	For sample solution in Inline Dialysis.
6.1826.360	Pump tubing LFL (white/white), 3 stoppers	PVC (Tygon®)	1.02 mm	For sample transfer.
6.1826.380	Pump tubing LFL (gray/gray), 3 stoppers	PVC (Tygon®)	1.25 mm	For Inline Dilution.
6.1826.390	Pump tubing LFL (yellow/yellow), 3 stoppers	PVC (Tygon®)	1.37 mm	For sample solution in Inline Ultrafiltration.

4.12.2.2 Pump tubing connection with filter

The 6.2821.130 filters (46-**2**) should be changed every 3 months, more frequently at higher backpressure.




Figure 46 Pump tubing connection – Changing the filter

- | | | | |
|----------|-----------------------|----------|--|
| 1 | Tubing olive | 2 | Filter 6.2821.130
Packaging contains 10 items. |
| 3 | Filter housing | | |

Replacing the filter

1 Unscrewing filter screw

- Screw the filter screw (46-3) out of the tubing olive (46-1) with the aid of two 6.2621.000 adjustable wrenches.

2 Replacing the filter

- Remove the old filter (46-2) with tweezers.
- Place the new filter (46-2) flat in the tubing olive (46-1) with tweezers.

3 Mounting filter screw

- Screw the filter screw (46-3) back into the tubing olive (46-1) and tighten by hand. Then additionally tighten with two 6.2621.000 adjustable wrenches.

4.13 Metrohm CO₂ Suppressor (MCS)

4.13.1 Replacing the CO₂ adsorption cartridge

The CO₂ adsorption cartridge (6.2837.000) (31-4) needs to be replaced at regular intervals, approximately every 6 months, to prevent blockage and capacity loss.

Blockage

Moisture blocks the CO₂ adsorption cartridge. This is indicated by a color change in the cartridge material (the orange part becomes colorless). As air flow is reduced, the vacuum sinks. As a protective measure, a H₂O adsorption cartridge (31-7) is installed prior to the CO₂ adsorption cartridge. A regular regeneration (*see chapter 4.13.2, page 93*) of the H₂O adsorption cartridge prolongs the service life of the CO₂ adsorption cartridge.

Capacity loss

The adsorption capacity of the CO₂ adsorption cartridge is limited. Depending on the operating time and the laboratory environment, the adsorption capacity decreases with time. This is indicated by an increasing baseline (since more CO₂ gets to the detector).

4.13.2 Regenerating the H₂O adsorption cartridge

The H₂O adsorption cartridge protects the CO₂ adsorption cartridge from moisture. The service life of the H₂O adsorption cartridge depends on the moisture content of the surrounding air. Moisture decreases the capacity of the H₂O adsorption cartridge which can be observed in the color change. Before the whole filling material changes its color (from orange to colorless, for Sigma-Aldrich art. no. 94098), the H₂O adsorption cartridge should be regenerated (*see leaflet*).

The filling material is replaced during regeneration.

Regenerating the H₂O adsorption cartridge

Proceed as follows to regenerate the H₂O adsorption cartridge:

- 1 Remove the material from the cartridge and let it dry at 140 °C over night. Then, refill the cartridge.

As an alternative, dispose of the old material and fill in new material.

- 2** Cover the packed material with cotton.

Two H₂O adsorption cartridges are provided so that you can continue to work without interruption while one of the H₂O adsorption cartridges is being regenerated.

4.14 Conductivity detector

4.14.1 Maintenance



CAUTION

The conductivity detector must not be opened!



WARNING

When rinsing the detector, the pressure must not exceed **5 MPa**. In order to ensure this, the maximum pressure of the high pressure pump must be set to **5 MPa** in MagIC Net.

If the conductivity detector is blocked, you should first of all check whether the blockage results from the capillary ends being pressed together too firmly. In this case, shorten the detector input capillary (32-**3**) or the detector output capillary (33-**2**) by several millimeters.

If this does not help, the conductivity detector can be rinsed against the normal flow direction. For this, connect the high pressure pump to the detector output capillary (33-2) and rinse - **the pressure may not exceed 5 MPa**.

4.15 Separation column

4.15.1 Separating efficiency

Which analysis quality can be attained, depends to a great extent on the separating efficiency of the separation column used. The separating efficiency of the selected separation column must be sufficient for the analysis problems present. If difficulties occur, you should always first check the quality of the separation column by recording a standard chromatogram.

You can find detailed information on the separation columns available from Metrohm in the leaflet provided along with your separation column, in the **Metrohm IC-Column Program** (available via your Metrohm agent) or in the Internet at <http://www.metrohm.com> in the product area Ion chromatography. You can request free information on special IC applications in the corresponding "**Application Bulletins**" or "**Application Notes**", which are available in the Internet at <http://www.metrohm.com> in the Applications area or via the Metrohm agent responsible.

4.15.2 Protection

To protect the separation column against foreign particles, which can affect the separating efficiency, we recommend that both the eluent and the samples undergo a microfiltration (filter 0.45 µm) before being aspirated via the aspiration filter (6.2821.090).

We recommend always to use a guard column (*see chapter 2.22, page 59*). This protects the separation column and considerably increases its service life. Information regarding which guard column is suitable for your separation column can be found in the **Metrohm IC Column Program** (which is available from your Metrohm agent), the leaflet provided along with your separation column, the product information on the separation column at <http://www.metrohm.com> (product area Ion Chromatography) or obtained directly from your agent.

The pulsation absorber (*see chapter 2.12, page 35*) must be installed in order to protect the column material from pressure concussion caused by injection.

4.15.3 Storage

Always store the separation columns sealed and filled according to the data of the column manufacturer when not using them.

4.15.4 Regeneration



NOTICE

The regeneration is considered as the last measure, and not to be carried out regularly.

If the separating properties of the column have deteriorated, the column can be regenerated according to the specifications of the column manufacturer. In the case of separation columns available from Metrohm, the specification for regeneration can be found on the leaflet provided along with each column.

5 Troubleshooting

5.1 Problems and their solutions

Problem	Cause	Remedy
Marked drop in pressure.	<i>Leak in the system.</i>	Check all capillary connections and seal leaks, if necessary (see chapter 2.5, page 13).
The baseline has a large amount of noise.	<i>High pressure pump – contaminated pump valves.</i>	Clean pump valves (see chapter 4.5.2, page 70).
	<i>Eluent – Leakage in eluent path.</i>	Check eluent path.
	<i>Eluent – Blockage in eluent path.</i>	Check eluent path.
	<i>High pressure pump – defective piston seals.</i>	Replace (see chapter 4.5.2, page 70) piston seals .
	<i>MCS – CO₂ adsorber cartridge depleted.</i>	Replace CO ₂ adsorber cartridge (see chapter 4.13.1, page 93).
	<i>Pulsation damper not connected.</i>	Connect the pulsation damper (see Chapter 2.12, Page 35).
	<i>Pulsation absorber not connected. or defective.</i>	Connect pulsation absorber (see chapter 2.12, page 35). or replace it.
The baseline is drifting.	<i>MCS – Vacuum pump defective.</i>	Contact Metrohm Service.
	<i>Thermal equilibrium not yet attained.</i>	Condition instrument with the column thermostat (see Chapter 2.15, Page 40) switched on .
	<i>Leak in the system.</i>	Check all capillary connections and seal leaks, if necessary (see chapter 2.5, page 13).
The pressure in the system markedly increases.	<i>Eluent – Evaporation of organic solvent in eluent.</i>	▪ Check the eluent bottle cap (see figure 12, page 26). ▪ Stir the eluent.
	<i>Inline filter (6.2821.120) blocked.</i>	Replace the filter (6.2821.130) (see chapter 4.6, page 80).

Problem	Cause	Remedy
	<i>Sample – gas bubbles in the sample.</i>	Use sample degasser (see chapter 2.13, page 36).
	<i>MCS – Not connected.</i>	Connect MCS.
The peristaltic pump is pumping too little.	<i>Peristaltic pump – contact pressure too weak.</i>	Correctly set contact pressure (see "Set flow rate", page 50).
	<i>Peristaltic pump – filter blocked.</i>	Replace the filter (see chapter 4.12.2.2, page 92).
	<i>Peristaltic pump – pump tubing defective.</i>	Replace pump tubing (see chapter 4.12.2.1, page 91).
Data of the separation column cannot be read.	<i>Column chip contaminated.</i>	Clean the contact surfaces of the column chip (with alcohol).
	<i>Column chip defective.</i>	1. Save column configuration in MagIC Net™. 2. Notify Metrohm Service.
Individual peaks are greater than expected.	<i>Sample – carry-over of the samples from previous measurement.</i>	Rinse system longer between two samples.
MSM – No (or insufficient) pumping of regeneration or rinsing solution	<i>Leak in the system.</i>	Check connections.
	<i>Peristaltic pump – contact pressure too weak.</i>	Correctly set contact pressure (see "Set flow rate", page 50).
	<i>Peristaltic pump – filter blocked (see figure 28, page 48).</i>	Replace the filter (see "Replacing the filter", page 92).
	<i>MSM – backpressure too high.</i>	Clean the MSM (see chapter 4.11.3.2, page 86) or replace parts (see chapter 4.11.3.3, page 89).
	<i>Peristaltic pump – pump tubing defective.</i>	Replace pump tubing (see figure 27, page 47).
The background conductivity is too high.	<i>MSM – not connected.</i>	Connect the MSM (see chapter 2.16, page 43).
	<i>MCS – Not connected.</i>	Connect MCS.
	<i>Wrong eluent.</i>	Change the eluent (see chapter 4.4.2.3, page 69).

Problem	Cause	Remedy
	<i>Injection valve – sample loop.</i>	Check installation of the sample loop (see chapter 2.14.1, page 38).
	<i>Sample – rinsing volume too low.</i>	Increase rinsing time (see chapter 4.8, page 82).
	<i>Injection valve – defective.</i>	Request Metrohm Service.
	<i>MCS – Vacuum too low.</i>	▪ Check connections. If they are OK: ▪ Contact Metrohm Service.
Vacuum is not being built	<i>Eluent degasser – Vacuum connector on the rear of the instrument not (tightly) sealed.</i>	▪ Seal the Vacuum connector tightly with a threaded stopper (6.1446.040).

6 Technical specifications

6.1 Reference conditions

The technical specifications listed in this chapter refers to the following reference conditions:

Ambient temperature	+25 °C (±3 °C)
Instrument status	> 40 minutes in operation (equilibrated)

6.2 Device

<i>IC system</i>	■ Metal-free IC system ■ Compact system with modular design ■ Up to two complete chromatographic systems in one housing
<i>Material</i>	Painted polyurethane hard foam without CFCs, fire class V0
<i>Operating pressure range</i>	■ 0...50 MPa (500 bar) high pressure pump ■ 0...35 MPa (350 bar) standard-PEEK system
<i>Intelligent components</i>	iPump, iDetector, iColumn, MagIC Net

6.3 Leak sensor

Type	Electronic, no calibration necessary
------	--------------------------------------

6.4 Ambient conditions

<i>Operation</i>	
<i>Ambient temperature</i>	+5...+45 °C
<i>Humidity</i>	20...80 % relative humidity
<i>Storage</i>	
<i>Ambient temperature</i>	-20...+70 °C
<i>Transport</i>	
<i>Ambient temperature</i>	-40...+70 °C

6.5 Housing

Dimensions

Width	365 mm
Height	642 mm
Depth	380 mm

Material of base tray, housing and flask holder	Polyurethane hard foam (PUR) with flame retardation for fire class V0, CFC-free, painted
---	--

Operating elements

Indicators	LED for power display
On/Off switch	On the rear of the instrument

6.6 Eluent degasser

Material	Fluoropolymer
Resistance to solvents	No restriction (except PFC)
Time to establish vacuum	< 60 s

6.7 High pressure pump

Type	▪ Serial dual-piston pump ▪ Intelligent pump head recognition ▪ Chemically inert ▪ Metal-free pump heads ▪ Materials in contact with the eluent: PEEK, ZrO₂, PTFE/PE ▪ Self-optimizing flow and pressure
------	--

Flow rate

Adjustable flow range	0.001...20.0 mL/min
Flow increment	1 µL/min
Reproducibility of the eluent flow	< 0.1 % deviation

6.10 Column thermostat

Type	Peltier technique thermostat for two intelligent separation columns
Adjustable temperature range	0...+ 80 °C, in increments of 0.1 °C
Heating	Ambient temperature +50 °C
Cooling	Ambient temperature –20 °C
Temperature reproducibility	± 0.2 °C
Stability	< 0.05 °C
Heating up time	< 30 minutes from 20 to 50 °C
Cooling time	< 40 minutes from 50 to 20 °C

6.11 Metrohm Suppressor Module (MSM)

<i>Resistance to sol-vents</i>	No restriction
<i>Switching duration</i>	typ.100 ms
<i>Operating pres-sure</i>	2.5 MPa (25 bar), valve function prevents damage at overpressure

6.12 Peristaltic pump

Type	2-channel peristaltic pump
Rotating direction	Counterclockwise/Clockwise rotation
Rotational speed	0...42 rpm in 7 stages at 6 rpm.
Pumping properties	0.3 mL/min at 18 rpm; with 6.1826.320 standard pump tubing.
Material of pump tubings	recommended: Tygon Long Flex Life

6.13 Metrohm CO₂ Suppressor (MCS)

<i>Material</i>	Fluoropolymer
<i>Resistance to solvents</i>	No restriction (except PFC)
<i>Vacuum</i>	
<i>Working range</i>	Microprocessor-controlled/stabilized
<i>Time to establish after starting</i>	< 30 s
<i>Capillary volume</i>	400 μ L
<i>Recommended flow range</i>	0.1 - 1.0 mL

6.14 Conductivity measuring system

Type	▪ Microprocessor-controlled Digital Signal Processing (DSP technology) ▪ Intelligent detector with 6 reference chromatograms
Measuring range	0...15000 $\mu\text{S}/\text{cm}$ single range
Noise	< 0.1 nS at 1 $\mu\text{S}/\text{cm}$
Deviations from the linearity	< 1 % for conductivity values greater than 1...16 $\mu\text{S}/\text{cm}$ (typical for analyses with sequential suppression)
Drift	< 0.2 nS/cm per hour
Measuring rate	10 measurements per second for optimum results without filtering
Resolution	0.0047 nS/cm
Baseline	Noise < 0.2 nS/cm typical for sequential suppression
Conductivity detector	
Cell volume	0.8 μL
Cell constant	▪ Individual calibration data saved in the detector ▪ adjustable range: 13.0...21.0 /cm
Electrodes	Ring-shaped electrodes made from stainless steel
Materials in contact with eluent	Chemically inert PCTFE

<i>Maximum operating pressure</i>	5.0 MPa (50 bar)
<i>Cell temperature</i>	20...50 °C in increments of 5 °C
<i>Temperature stability</i>	< 0.001 °C
<i>Temperature compensation</i>	0...5 %/K adjustable, default 2.3 %/K
<i>Heating up time</i>	< 30 minutes (40 °C)

6.15 Power connection

<i>Required supply voltage</i>	100 - 240 V ± 10% (autosensing)
<i>Required frequency</i>	50 - 60 Hz ± 3 Hz (autosensing)
<i>Power consumption</i>	65 W for typical analysis application 25 W standby (conductivity detector to 40 °C)
<i>Power supply unit</i>	- Up to 300 W maximum, electronically monitored - Internal fuse 3.15 A

6.16 Interfaces

USB

<i>Entry</i>	1 USB upstream, type B (for connection to the PC)
<i>Exit</i>	2 USB downstream, type A

<i>MSB</i>	2 MSB 8-pin Mini DIN (female) (for Dosino, stirrer, remote lines, etc.)
------------	---



CAUTION

When connecting an instrument to the MSB connector you **must** switch off the 850 Professional IC.

<i>Detector</i>	2 15-pin high-density DSUB (female)
<i>Column recognition</i>	3 (including 2 in the column thermostat (<i>see chapter 2.15, page 40</i>))
<i>Leak sensor</i>	1 jack plug
<i>Further connections</i>	1 15-pin DSUB (female)



6.17 Weight

1.850.2030	29.4 kg (without accessories)
1.850.9010 (<i>conductivity detector</i>)	2.3 kg (with accessories)
Transport carriage (<i>rollers and handle</i>)	1.8 kg

7 Accessories

Up-to-date information on the scope of delivery and optional accessories for your product can be found on the Internet. You can download this information using the article number as follows:

Downloading the accessories list

- 1 Enter <https://www.metrohm.com/> into your Internet browser.
- 2 Enter the article number (e.g. **2.850.2030**) into the search field.
The search result is displayed.
- 3 Click on the product.
Detailed information regarding the product is shown on various tabs.
- 4 On the **Included parts** tab, click on **Download the PDF**.
The PDF file with the accessories data is created.



NOTICE

Once you have received your new product, we recommend downloading the accessories list from the Internet, printing it out and keeping it together with the manual for reference purposes.

(6.2821.130) Filter	80
6.2821.090 aspiration filter	69

Adsorption cartridges	
Connection	52
Ambient conditions	102
Aspiration filter 6.2821.090	69
Aspiration tubing for eluent	24

Baseline	
Conditioning	65
Unstable	70
Blockage	
Conductivity detector	94
Blood	82

Cable feed-throughs	21
Capillaries	
Installation	13
Capillary feed-throughs	21
Carry-over	82
Cartridges	
Connection	52
Cleaning	
MSM	87
Valves of the high pressure pump	75
CO ₂ adsorption cartridge	53
Replace	93
Column	
see also "Separation column"	
.....	60
Column recognition	107
Column thermostat	
Column thermostat	105
Installation	40
Computer connection	57
Conditioning	65
Conductivity detector	
Capillary connector	55
Cell constant	106
Cell volume	106
Maintenance	94
Conductivity measuring system	
Technical data	106
Connect	
Power grid	58

To computer	57
Connection	
Mains	107
Connections	
Installation	13
Contamination	
High pressure pump	69
Valves in the high pressure pump	70
Contamination MSM	
Heavy metals	85
Crystallization	
High pressure pump	69

Deaerating	
High pressure pump	32
Purge valve	30
Degasser	
Eluent degasser	28
Sample degasser	36
Degassing	
Eluent	28
Detector	
Cable connector	18
Conductivity detector	55
Interface	107
Positioning	18
Diagram	10
Dilution	82
Dimensions	103
Door	68
Drainage tubings	
Installation	19

Electrostatic charge	7
Eluent	
Aspirate	24
Change	69
Production	68
Eluent bottle	
Figure	27
Installation	24
Operation	69
Eluent degasser	
Installation	28
Technical specifications	103
Equilibration	63, 65

Feed-throughs	
Capillaries	21
Fill	
Injection valve	40
Filter	
see also "inline filter"	34
Filter (6.2821.130)	80
Filter 6.2821.090	
Aspiration filter	69
Flow diagram	10
Flow fluctuations	70
Flow increment	103
Flow range	103
Flow rate	103
Frequency	107

Gas	28, 36
Guard column	
Installation	59
Rinsing	60

H ₂ O adsorption cartridge	53
Regenerate	93
Handle	15
Heating	
see also "Column thermostat"	
.....	40
Heavy metals	
Contamination of the MSM	85
High pressure pump	
Installation	30
Maintenance	69
Protection	18, 69
Technical specifications	103
Tubing connection	30
Valves	78
Housing	103
Humidity	102

IC column	
see also "Separation column"	
.....	60
Impurity MSM	
Organic	85
Initial installation	9

- Inject
 Injection valve 40
 Injection valve 4
 Fill 40
 Inject 40
 Installation 38, 104
 Maintenance 84
 Protection 84
 Inline filter 34
 Inline sample preparation 82
 Installation
 Column thermostat 40
 Conductivity detector 55
 Connections 13
 Drainage tubings 19
 Eluent bottle 24
 Eluent degasser 28
 Guard column 59
 High pressure pump 30
 Initial installation 9
 Injection valve 38, 104
 Leak sensor 18
 MCS 51
 MSM 43
 Peristaltic pump 47
 Pulsation absorber 35
 Pump tubings 47
 Sample degasser 36
 Separation column 60
 Interface
 MSB 107
 USB 107
 Interfaces 107
 Further connections 107
 Leak sensor 107
- L**
 Leak 70
 Leak sensor
 Installation 18
 Interface 107
 Technical specifications 102
 Leak-tightness 64, 65
 Leaking piston seals 70
 Loop
 see also "Sample loop" 40
- M**
 Maintenance
 Conductivity detector 94
 High pressure pump 69
 Injection valve 84
 MSM 84
 Peristaltic pump 90
 Pump head 70
 Material 103
 MCS
 Capillary connection 51
 Cartridge connection 52
 Installation 51
 Technical specifications 106
 Use 51
 Measuring range 106
 MPak
 Holder 17
 MSB 107
 MSM
 Cleaning 87
 Installation 43
 Maintenance 84
 Operation 84
 Protection 84
 Regeneration 85
 Replacing parts 89
 Switchover 85
 Technical data 105
- N**
 Noise 106
- O**
 Oil 82
 Operation
 MSM 84
 Peristaltic pump 90
 Sample degasser 84
 Organic impurities
 MSM 85
- P**
 Peristaltic pump
 Installation 47
 Maintenance 90
 Operation 90
 Principle 46
 Technical specifications 105
 Piston seal 70
 Pistons of the high pressure pump
 70
 Power connection 58, 107
 Power consumption 107
 Power supply unit 107
 Precipitates 69
 Pressure limit 104
 Pressure range 104
 Pressure screws
 Connection 13
 Protection
 Injection valve 84
 Inline filter 34
 MSM 84
 Pulsation 70
 Pulsation absorber
 Installation 35
 Pump head
 Maintenance 70
 Pump tubings
 Installing 47
 Overview 91
 Service life 90
 Purge valve 30
- R**
 Reference conditions 102
 Regeneration 66
 MSM 85
 Rinsing
 Conductivity detector 94
 Guard column 60
 Pump tubings 91
 Sample path 82
 Separation column 62
 Rinsing time 83
 Rise in pressure 69
 Rollers 15
- S**
 Safety instructions 6
 Safety shutdown 104
 Sample
 Carry-over 82
 Sample loop 40
 Transfer time 83
 Sample degasser
 Installation 36
 Operation 84
 Technical specifications 104
 Sample loop 40
 Sample path
 Rinsing 82
 Sample preparation 82
 Screws
 Connection 13
 Separation column
 Installation 60
 Protection 3, 35, 95
 Regeneration 96
 Rinsing 62
 Separating efficiency 95
 Storage 95
 Service 6, 66
 Shutting down 67
 Start-up 63

Index

Storage 102
Supply voltage 7, 107
Suppressor
 Maintenance 84
 Operation 84
 see also "MSM" 43

T

Technical data
 Column thermostat 105
 Conductivity measuring system
 106
 MSM 105
Technical specifications
 Detector 107
 Eluent degasser 103

High pressure pump 103
Interfaces 107
Leak sensor 102
MCS 106
Peristaltic pump 105
Reference conditions 102
Sample degasser 104
Temperature 102
Thermostat
 see also "Column thermostat"
 40
Transfer time 83
Transport 102
 Rollers 15
Transport locking screws 18

Tubing 10
Tubing pump
 See also "Peristaltic pump" . 46
Tubings
 Installation 13

U

USB 107

V

Vacuum pump
 Protection 18
Valve
 see also "Injection valve" 38
Valves of the high pressure pump
 78