

850 Professional IC



AnCat – 2.850.3010

Manual

8.850.8051EN / 2019-10-10



Metrohm AG
CH-9100 Herisau
Switzerland
Phone +41 71 353 85 85
Fax +41 71 353 89 01
info@metrohm.com
www.metrohm.com

850 Professional IC

AnCat – 2.850.3010

Manual

Technical Communication
Metrohm AG
CH-9100 Herisau
techcom@metrohm.com

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Table of contents

1	Introduction	1
1.1	Instrument description	1
1.2	Intended use	4
1.3	About the documentation	5
1.3.1	Symbols and conventions	5
1.4	Safety instructions	6
1.4.1	General notes on safety	6
1.4.2	Electrical safety	6
1.4.3	Tubing and capillary connections	7
1.4.4	Flammable solvents and chemicals	7
1.4.5	Recycling and disposal	7
2	Installation	9
2.1	About this chapter	9
2.2	Initial installation	9
2.3	Flow diagram	10
2.4	Setting up the instrument	14
2.4.1	Packaging	14
2.4.2	Checks	14
2.4.3	Location	14
2.5	Capillary connections in the IC system	14
2.6	Rear of the instrument	17
2.6.1	Rollers and handle	17
2.6.2	Positioning and connecting the detector	20
2.6.3	Transport locking screws	20
2.6.4	Leak sensor	20
2.6.5	Drainage tubings	21
2.7	Capillary and cable feed-throughs	23
2.8	Eluent	26
2.8.1	Connecting the eluent bottle	26
2.9	Eluent degasser	30
2.10	High pressure pump	32
2.10.1	Capillary connections high pressure pump/purge valve	32
2.10.2	Deaerating the high pressure pump	34
2.11	Inline filter	36
2.12	Pulsation absorber	37
2.13	Sample degasser	38

4.9	Sample degasser	82
4.9.1	Operation	82
4.10	Injection valve	82
4.10.1	Protection	82
4.11	Metrohm Suppressor Module (MSM)	82
4.11.1	Protection	82
4.11.2	Operation	82
4.11.3	Maintenance	83
4.12	Peristaltic pump	88
4.12.1	Operation	88
4.12.2	Maintenance	89
4.13	Servicing the detector	91
4.14	Separation column	91
4.14.1	Separating efficiency	91
4.14.2	Protection	91
4.14.3	Storage	92
4.14.4	Regeneration	92
5	Troubleshooting	93
5.1	Problems and their solutions	93
6	Technical specifications	98
6.1	Reference conditions	98
6.2	Device	98
6.3	Leak sensor	98
6.4	Ambient conditions	98
6.5	Housing	99
6.6	Eluent degasser	99
6.7	High pressure pump	99
6.8	Sample degasser	100
6.9	Injection valve	100
6.10	Column thermostat	101
6.11	Metrohm Suppressor Module (MSM)	101
6.12	Peristaltic pump	101
6.13	Detector	102
6.14	Power connection	102
6.15	Interfaces	102
6.16	Weight	103



7 Accessories	104
Index	105

Table of figures

Figure 1	Front view of the instrument	2
Figure 2	Flow diagram for a 2-channel system	12
Figure 3	Connection of capillaries with pressure screws	15
Figure 4	Rollers and handle	18
Figure 5	Handle as MPak holder	19
Figure 6	Connection of the leak sensor on the rear of the instrument	21
Figure 7	Drainage tubings	22
Figure 8	Capillary feed-throughs on the door	24
Figure 9	Capillary feed-throughs base tray/flask holder	25
Figure 10	Installing the eluent bottle cap	27
Figure 11	Mounting the aspiration filter	27
Figure 12	Installing the tubing weighting and aspiration filter	28
Figure 13	Eluent aspiration tubing fully equipped	28
Figure 14	Eluent bottle – connected	29
Figure 15	Eluent degasser	31
Figure 16	Capillary connections high pressure pump/purge valve	32
Figure 17	High pressure pump – Connect inlet	33
Figure 18	Deaerate the high pressure pump	35
Figure 19	Connecting the inline filter	37
Figure 20	Pulsation absorber – Connection	38
Figure 21	Sample degasser	39
Figure 22	Injection valve – connected	40
Figure 23	Injection valve – Positions	41
Figure 24	Column thermostat	43
Figure 25	MSM – connections	46
Figure 26	Peristaltic pump	48
Figure 27	Installing the pump tubing	49
Figure 28	Install pump tubing connection with filter	50
Figure 29	Install pump tubing connection without filter	51
Figure 30	Front of conductivity detector	53
Figure 31	Rear of conductivity detector	54
Figure 32	Connection detector–MSM	55
Figure 33	Pump head – removing the piston	69
Figure 34	Components of the piston cartridge	70
Figure 35	Tool for piston seal	71
Figure 36	Removing the piston seal	72
Figure 37	Inserting the piston seal into the tool	72
Figure 38	Inserting the piston seal into the pump head	73
Figure 39	Removing valves	74
Figure 40	Dismantling valve	75
Figure 41	Components of the inlet valve and outlet valve	76
Figure 42	Change filters (of the inline filter)	78
Figure 43	MSM – Components	85
Figure 44	Pump tubing connection – Changing the filter	90




Figure 1 Front view of the instrument

1	Conductivity detectors <i>See Chapter 2.18.</i>	2	Peristaltic pump <i>See Chapter 2.17.</i>
3	MSM <i>See Chapter 2.16.</i>	4	Injection valve <i>See Chapter 2.14.</i>
5	Eluent degasser <i>See Chapter 2.9.</i>	6	High pressure pump <i>See Chapter 2.10.</i>
7	Purge valve <i>See Chapter 2.10.1.</i>	8	Injection valve <i>See Chapter 2.14.</i>
9	Eluent degasser <i>See Chapter 2.9.</i>	10	High pressure pump <i>See Chapter 2.10.</i>
11	Purge valve <i>See Chapter 2.10.1.</i>	12	Pulsation damper <i>See Chapter 2.12.</i>

13 Sample degasser*See Chapter 2.13. Use optional.***14 Separation columns***See Chapter 2.22.***15 Column thermostat***See Chapter 2.15.*

The instrument contains the following components:

Eluent degasser

The eluent degasser removes gas bubbles and dissolved gases from the eluent. For degassing, the eluent flows into a vacuum chamber through a special fluoropolymer capillary.

High pressure pump

The intelligent and low pulsation high pressure pump pumps the eluent through the system. It is equipped with a chip on which its technical specifications and "life history" (operating hours, service data, ...) are saved.

Inline filter

Inline filters protect the separation column securely against possible contamination from the eluent. Inline filters can however also just as well be used for the purpose of protecting other sensitive components against contaminations in the solutions used. The filter platelets with a pore size of 2 µm can be replaced quickly and easily. They remove particles like e. g. bacteria and algae from the solutions.

Pulsation absorber

The pulsation absorber protects the separation column from damage caused by pressure fluctuations when switching the injection valve, and reduces interfering pulsations during highly sensitive measurements.

Sample degasser

The sample degasser removes gas bubbles and dissolved gases from the sample. For degassing, the sample flows into a vacuum chamber through a special fluoropolymer capillary.

Injection valve

The injection valve connects the eluent and sample path through rapid and precise valve switchover. A precisely measured amount of sample solution is injected and rinsed with eluent onto the separation column.

Column thermostat

The column thermostat controls the temperature of the column and eluent channel and thus ensures stable measuring conditions. It provides space for 2 separation columns.



More details about the disposal of your old instrument can be obtained from your local authorities, from waste disposal companies or from your local dealer.

2 Installation

2.1 About this chapter

The *Installation* chapter contains

- a list of all the steps to be carried out during the initial installation (*see chapter 2.2, page 9*).
- a schematic overview of the flow paths (*see chapter 2.3, page 10*).
- a description of the construction, connections and functionality of the instrument.
- step-by-step installation instructions. A part of these installation tasks has already been carried out before delivery. They are nevertheless described here in case they need to be carried out once again some time in the future (e.g. after maintenance of the high pressure pump).

2.2 Initial installation



NOTICE

A large number of the capillary connections are already connected at the time the instrument is delivered.

The following steps still have to be carried out after delivery:



NOTICE

Each of the steps numbered 3, 7, 8, 10, 14 and 15 are to be performed for the anion and cation channels, respectively.

- 1** Setting up the instrument (*see chapter 2.4, page 14*).
- 2** Remove handle and rollers (*see chapter 2.6.1, page 17*).
- 3** Place the detectors in the instrument and connect them (*see chapter 2.6.2, page 20*).
- 4** Remove transport locking screws (*see chapter 2.6.3, page 20*).




Figure 2 Flow diagram for a 2-channel system

1 Eluent input

Connection to the eluent bottle (see chapter 2.8.1, page 26).

2 Eluent degasser

See Chapter 2.9.

A diagram of a three-part assembly. It consists of two identical components, each with a conical tip and a cylindrical body, connected by a central rod. The components are labeled 1 and 3. The central rod is labeled 2. A fourth label, 4, points to a cylindrical component with a flange, which is shown separately above the main assembly.

Figure 3 Connection of capillaries with pressure screws

- | | |
|--|---|
| <p>1 PEEK pressure screw (6.2744.014)
Use on the injection valve.</p> | <p>2 Connection capillary</p> |
| <p>3 PEEK pressure screw, short (6.2744.070)
For use on the high pressure pump, the purge valve, the inline filter, the pulsation absorber, the guard column and the separation column.</p> | <p>4 PEEK pressure screw, long (6.2744.090)
Use on special components. Is not used on all instruments.</p> |



In order to keep the dead volume as low as possible, capillary connections should generally be as short as possible.



For an improved overview, capillary and tubing connections can be bundled with the 6.1815.010 spiral band.

Connection capillaries

PEEK capillaries and PTFE capillaries are used in the IC system.

PEEK capillaries (poly-
etheretherketone)

PEEK capillaries are temperature-resistant up to 100°C, stable under pressure up to 400 bar, flexible, chemically inert and exhibit an extremely smooth surface. They can be readily cut down to the desired length with the 6.2621.080 capillary cutter.

Usage:

- PEEK capillaries (6.1831.010) with an internal diameter of 0.25 mm for the entire high-pressure section.



- PEEK capillaries (6.1831.030) with an internal diameter of 0.75 mm for sample handling in the ultratrace range.



CAUTION

For the capillary connections between the injection valve and detector, PEEK capillaries with an internal diameter of 0.25 mm must be used. These are already connected to a newly delivered instrument.

PTFE capillaries (polytetrafluoroethylene)

PTFE capillaries are transparent and enable visual tracing of the liquids to be pumped. They are chemically inert, flexible and temperature-resistant up to 80°C.

Usage:

PTFE capillaries (6.1803.0x0) are used for the low-pressure section.

- PTFE capillaries with internal diameter of 0.5 mm for sample handling.
- PTFE capillaries with internal diameter of 0.97 mm for sample handling as well as for rinsing solutions (they do not have to be in the scope of delivery of the instrument).

Capillary connections

In order to achieve optimum analysis results, capillary connections in an IC system must be absolutely tight and free of dead volume. Dead volume occurs if two capillary ends connected to each other do not fit exactly, thus allowing liquid to escape. There are two possible reasons for this:

- The capillaries do not have exactly cut edges.
- The two capillary ends do not completely meet.

One prerequisite for dead volume free capillary connection is, that both capillary ends are cut exactly plane. Therefore we recommend only to cut PEEK capillaries with the capillary cutter (6.2621.080).

Creating dead volume free capillary connections

To create dead volume free capillary connections, proceed as follows:

- 1** Slide the pressure screw over the capillary. Ensure that the capillary protrudes 1–2 mm from the tip of the pressure screw.
- 2** Plug the capillary all the way into the connection or coupling until the stop.
- 3** Only then start turning the pressure screw, while keeping the capillary pressed in space.

The enclosed set of varicolored sleeves for PEEK capillaries (6.2251.000) serves to easily differentiate the various flows of liquid in the system through color coding. Each capillary leading a given liquid (e. g. eluent) can be highlighted with sleeves of the same color.

- 1 Slide a sleeve of a selected color over a capillary and move it to an easily visible position.

2.6 Rear of the instrument

In order to make transport easier, the instrument is equipped with rollers and a handle.

Section of the 6.1816.020 silicon tubing.
Guides escaped liquid into a waste container.

For draining the escaped liquid from the base tray through the connected drainage tubing.

For supplying escaped liquid through the connected drainage tubing to the leak sensor.

Installing drainage tubings

- 1 Plug drainage tubing (7-2) into the drainage tubing connector (7-1) of the flask holder and shorten to required length.
- 2 Plug drainage tubing (7-4) into the drainage tubing connector (7-3) of the detector chamber and shorten to required length.
- 3 Connect drainage tubing (7-2) from the flask holder and drainage tubing (7-4) from the detector chamber to the Y connector (7-5).
- 4 Connect drainage tubing (7-6) to the Y connector (7-5), shorten to required length and plug the other end of the drainage tubing into the drainage tubing connector (7-9) of the base tray.
- 5 Plug drainage tubing (7-7) into the drainage tubing connector (7-8) of the base tray and guide the other end into a waste container.

Several openings have been integrated for feeding through capillaries and cables. They are located on the doors (*see figure 8, page 24*), on the rear panel or below the flask holder or above the base tray (*see Figure 9, Page 25*).




Figure 8 Capillary feed-throughs on the door

1 Luer connectors

For connecting a 6.2816.020 syringe. For manual sample feeding.

**3 PEEK pressure screws, short
6.2744.070**

2 Capillary feed-through

4 Door

Do not feed capillaries through the Luer connectors (8-**1**). The capillaries are fastened with PEEK pressure screws (8-**3**) from inside to the Luer connector. From outside, liquid can be aspirated or injected with a syringe.



1	Side panel (right) Right panel.	2	Rear of the instrument
3	Side panel (left) Left panel.	4	Front of the instrument
5	Capillary feed-through Upper. From front to right.	6	Capillary feed-through Upper. From front to back.
7	Capillary feed-through Upper. From front to back.	8	Capillary feed-through Upper. From front to left.
9	Capillary feed-through Lower. From front to right.	10	Capillary feed-through Lower. From front to back.
11	Capillary feed-through Lower. From front to back.	12	Capillary feed-through Lower. From front to left.

2.8 Eluent

2.8.1 Connecting the eluent bottle

The eluent is aspirated out of the eluent bottle via the eluent aspiration tubing (10-1).

The eluent aspiration tubing is connected to the eluent degasser (*see chapter 2.9, page 30*). The tubing must be threaded through a suitable capillary feed-through (*see chapter 2.7, page 23*) of the instrument before the other end can be equipped.

The following accessories are needed for equipping the eluent aspiration tubing:

- 6.1602.160 eluent bottle cap GL 45
- 6.2744.210 tubing adapter for aspiration filter
- 6.2821.090 aspiration filter

Proceed as follows for equipping the eluent aspiration tubing:

Assembling eluent aspiration tubing

- 1** Guide the free end of the eluent aspiration tubing (10-1) out of the instrument through a suitable capillary feed-through.
- 2 Installing the eluent bottle cap (6.1602.160)**
 - Slide tubing nipple (10-2) and O-ring (10-3) onto the eluent aspiration tubing (10-1).
 - Push eluent aspiration tubing (10-1) through the bottle cap (10-4) and screw tight.



- ### 3 Mounting the aspiration filter

-

Figure 11 Mounting the aspiration filter

- 27



4 Installing the tubing weighting and aspiration filter




Figure 12 Installing the tubing weighting and aspiration filter

1	Eluent aspiration tubing (6.1834.080)	2	Eluent bottle cap (6.1602.160)
3	Tubing weighting From accessory set (6.2744.210).	4	Clamping screw From accessory set (6.2744.210).
5	Aspiration filter (6.2821.090) With filter holder from accessory set (6.2744.210).		

- Slide the tubing weighting (12-**3**) onto the eluent aspiration tubing (12-**1**).
- Slide the clamping screw (12-**4**) onto the eluent aspiration tubing (12-**1**).
- Insert eluent aspiration tubing (12-**1**) into the aspiration filter (12-**5**). The end of the tubing should reach approximately to the center of the aspiration filter.
- Screw the clamping screw (12-**4**) to the filter holder (11-**1**).




Figure 13 Eluent aspiration tubing fully equipped


5 Mounting the eluent aspiration tubing to the eluent bottle

- Insert the eluent aspiration tubing into the eluent bottle (**14-10**).

- ## 6 Mounting the adsorber tube



- First, place a piece of cotton (14-3), then the CO₂ adsorber (14-4) in the large opening of the adsorber tube (14-2) and close with the plastic cover.
- Fasten the adsorber tube (14-2) on the bottle cap (14-11) using the SGJ clip (14-12).



- 2 Adsorber tube (6.1609.000)**




Figure 19 Connecting the inline filter

1 Connection capillary

Connects the purge valve with the inline filter

3 Inline filter (6.2821.120)

Protects against particles.

2 PEEK pressure screws, short (6.2744.070)

4 Connection capillary

Connects the inline filter with the pulsation absorber.

- 1** Screw on the connection capillary running from the purge valve to the input side of the inline filter using a pressure screw (6.2744.070).
- 2** Screw on the connection capillary running to the pulsation absorber to the output side of the inline filter using a pressure screw (6.2744.070).

2.12 Pulsation absorber



NOTICE

The pulsation absorber is already installed in the newly delivered instrument.



CAUTION

The pulsation absorber is maintenance-free and may not be opened.

The pulsation absorber protects the separation column from damage caused by pressure fluctuations when switching the injection valve, and reduces interfering pulsations during highly sensitive measurements. In order to ensure these functionalities, it must be connected between the high pressure pump (see chapter 2.10, page 32) and injection valve (see chapter 2.14, page 40).

The pulsation absorber can be operated in both directions.



Connecting the sample degasser

- 

If the sample degasser is not used, the input and output **must** be sealed with the threaded stoppers (6.2744.220).

The sample loop can be replaced, depending on requirements. For additional information concerning selection of the appropriate sample loop, *see Chapter 2.14.3, page 42*.



1 Removing existing sample loop

- ## 2 Mounting new sample loop

- Fasten one end of the sample loop (22-**2**) with a 6.2744.010 PEEK pressure screw (22-**7**) to connector 3.
- Fasten the other end of the sample loop (22-**2**) with a second 6.2744.010 PEEK pressure screw (22-**7**) to connector 6.

The injection valve (see figure 23, page 41) can adopt two valve positions - **FILL** and **INJECT**. Switching back and forth between the two valve positions determines whether the sample path or the eluent path is guided through the sample loop. The following figure provides a schematic display of the flow paths of the two valve positions.



B Position INJECT



NOTICE

The column input capillary is already threaded into the capillary recesses of the column thermostat at the time the instrument is first delivered. The following installation instructions need **not** be carried out at the time of initial installation.

Threading capillaries

- 1 Introduce column input capillary into the column thermostat via a suitable capillary feed-through (24-**1**).
- 2 Slide the column input capillary from below into the outer of the two capillary recesses (24-**2**). Slide it through under the holder plate until it emerges again at the top.
- 3 Carefully bend the column input capillary downward and slide from from above to below through the inner capillary recess until it emerges at the lower edge of the holder plate.



NOTICE

The columns (guard column and separation column) may only be installed after the initial start-up (see chapter 3.1, page 61).

- **Before initial start-up:**
Fasten the 6.2744.040 coupling with a 6.2744.010 pressure screw to the end of the column input capillary.
- **After initial start-up:**
Fasten the guard column (if used) or the separation column with a 6.2744.010 pressure screw to the end of the column input capillary.




Figure 25 MSM – connections

1	Union nut	2	MSM connecting piece 6.2832.010
3	Eluent input capillary Labeled with <i>Eluent</i> .	4	Eluent output capillary Labeled with <i>Detector</i> .
5	Rinsing solution input capillary Labeled with <i>H2O</i> .	6	Rinsing solution output capillary Labeled with <i>Waste</i> .
7	Regeneration solution output capillary Labeled with <i>Waste</i> .	8	Regeneration solution input capillary Labeled with <i>H2SO4</i> .

The PTFE capillaries firmly mounted on the MSM are connected to the other components of the IC system as follows:



CAUTION

As the PTFE capillaries are very soft, the pressure screws should not be overtightened.

Flattened capillaries must be shortened with the help of a 6.2621.080 capillary cutter.




Figure 31 Rear of conductivity detector

1 Detector cable

With installed plug.

2 Detector output capillary

Permanently installed.

3 Type plate

With serial number.



NOTICE

In order to prevent unnecessary peak widening after separation, the connection between the outlet of the separation column and the inlet to the detector should be kept as short as possible.

Connecting the detector input capillary to the suppressor

- 1 ■ Connect the detector input capillary (32-**1**) and the capillary of the MSM labeled *out* (32-**2**) to each other using a coupling (6.2744.040) (32-**3**) and two short pressure screws (6.2744.070) (32-**4**).

Connecting and rinsing the guard column

1 Connecting the guard column



CAUTION

When inserting the guard column, always ensure that it is inserted correctly corresponding to the flow direction (if indicated).

- Remove sealing caps and/or stoppers from the guard column.
- Fasten the input of the guard column to the column inlet capillary using a short PEEK pressure screw (6.2744.070).
- In case the guard column is mounted on the separation column with a connecting capillary, connect the included connection capillary, which is included with the guard column, to the output of the guard column using the PEEK pressure screw, which is also included.

2 Rinsing the guard column

- Place beaker under the outlet capillary of the guard column.
- Set the flow rate of the high pressure pump according to the data given in the leaflet of the separation column.
- Start the high pressure pump and rinse the guard column approx. 5 minutes with eluent.
- Switch off the high pressure pump again.

2.22 Separation column

The intelligent separation column (iColumn) is the heart of the ion chromatographic analysis. It separates the different components corresponding to their interactions with the column. Metrohm separation columns are equipped with a chip on which their technical specifications and their history (first use / setting up, operating hours, injections, ...) are saved.



NOTICE

Information regarding which separation column is suitable for your application can be found in the **Metrohm IC Column Program**, the product information for your separation column at <http://www.metrohm.com> in the product area ion chromatography, or obtained directly from your agent.



You can find the separation columns and guard columns currently available from Metrohm in the Metrohm IC Column Program, or in the Internet at <http://www.metrohm.com> in the product area Ion Chromatography. A test chromatogram and an leaflet are provided along with each column. You can request detailed information on special IC applications in the corresponding "**Application Bulletins**" or "**Application Notes**", available in the Internet at <http://www.metrohm.com> in the Applications area or via the Metrohm agent responsible free of charge.



Connecting and rinsing the separation column

1 Connect the separation column



- Remove stoppers from the separation column.
- Attach the guard column to the input of the separation column.
OR
Connect the input of the separation column to the outlet capillary of the guard column, using the PEEK pressure screw (6.2744.070) included.
OR
If no guard column is used (not recommended), connect the column input capillary to the input of the separation column, using a PEEK pressure screw (6.2744.070).

3 Start-up

The chapter *Start-up* is divided into 2 sections:

Initial start-up	The initial start-up is carried out during the initial installation .
Conditioning	Conditioning is carried out as a final installation step and each time after the system is started.

3.1 Initial start-up

The initial start-up is carried out during the initial installation. The entire system is rinsed before guard column and separation column are installed.



CAUTION

The separation column and guard column may not be installed for the initial start-up.

Make sure that the coupling (6.2744.040) is being used instead of the columns.

Perform the following steps during the initial start-up:

1 Preparing the software

- Start the PC program **MagIC Net™**.
- Open the **Equilibration** tab in MagIC Net™.
- Select (or create) a suitable method.

2 Preparing the instrument

- Ensure that the eluent aspiration tubing is immersed in the eluent and that there is enough eluent in the eluent bottle.
- Ensure that the aspiration tubings for the auxiliary solutions (regeneration solution and rinsing solution) are immersed into the respective solutions and that there is enough solution in the bottles.
- Switch on the instrument.

3 Starting equilibration

- In MagIC Net™, start the equilibration.

- Select (or create) a suitable method.

2 Preparing the instrument

- Ensure that the column is correctly mounted according to the flow direction indicated on the label (arrow must point in the direction of flow).
- Ensure that the eluent aspiration tubing is immersed in the eluent and that there is enough eluent in the eluent bottle.
- Ensure that the aspiration tubings for the auxiliary solutions (regeneration solution and rinsing solution) are immersed into the respective solutions and that there is enough solution in the bottles.

3 Checking leak-tightness

- In MagIC Net™, start the equilibration.
- Check all capillaries and their connections from the high pressure pump to the detector for signs of liquid escaping.
If eluent escapes anywhere, tighten the corresponding pressure screw or loosen the connection, check the end of the capillary, shorten it with a capillary cutter if necessary, and restore the connection.

4 Conditioning the system

Rinse the system with eluent until the required stability of the baseline is attained (normally 30 minutes).

During this time, step the MSM to the next position every 10 minutes.

The instrument is now ready for measuring samples.



- Remove the backup ring from the pump head and lay to the other parts.




Figure 34 Components of the piston cartridge

1	Piston cartridge screw	2	Retaining washer
3	Zirconium oxide piston with piston shaft Order number: 6.2824.070	4	Spring retainer
5	Spring Order number: 6.2824.060	6	Inner plastic sleeve Protects from metallic abrasion.
7	Piston cartridge	8	Backup ring

3 **Cleaning the components of the piston**

- Clean zirconium oxide pistons contaminated by abrasion or deposits with fine abrasive cleaning powder, rinse particle free with ultrapure water and dry.
Replace highly contaminated or scratched zirconium oxide pistons (spare part: 6.2824.070 zirconium oxide piston).
- Rinse the other parts of the piston and dry with a lint-free cloth.

4 **Assembling the piston**

- Insert the inner plastic sleeve, spring and spring retainer into the piston cartridge.
- Slide the zirconium oxide piston carefully into the piston cartridge until its tip emerges from the small opening of the piston cartridge.
- Attach screw and tighten by hand.

Replacing the piston seal

The special tool (6.2617.010) (see figure 35, page 71) is necessary to remove the piston seal from the pump head. It consists of two parts: a tip

The diagram shows two tubes. Tube 2 is a standard cylindrical tube. Tube 1 is a tapered tube, wider at one end and narrower at the other, with a label '1' pointing to its wider section.

2 Sleeve
Sleeve for inserting the new piston seal.



CAUTION

Screwing the special tool for the piston seal (6.2617.010) into the piston seal destroys this completely!

1 Removing the piston seal



CAUTION

Avoid touching the sealing surface in the pump head (16-4) with the tool.

Screw the special tool for the piston seal (35-**1**) with the narrow side just as far into the piston seal as the same can be removed.




Figure 38 Inserting the piston seal into the pump head

4 Replacing the piston cartridge

Screw the assembled piston cartridge back into the pump head and tighten, first by hand, then additionally by approx. 15° with a wrench.

Cleaning the inlet valve and outlet valve

1 Removing valves

- Unscrew the connection capillary for the auxiliary piston (16-1) from the outlet valve holder.
- Unscrew the holders for the inlet and outlet valves and remove the valves (39-3) and (39-2).




Figure 39 Removing valves

1	Outlet valve holder	2	Outlet valve Order number: 6.2824.160
3	Inlet valve Order number: 6.2824.170	4	Inlet valve holder

2 Cleaning undissected valve

Clean contaminated or blocked valves initially **without** dismantling them completely.

- Rinse the valve in eluent flow and counterflow direction using a spray bottle filled with ultrapure water, RBS solution or acetone.
- The rinsing effect is further increased through a short treatment (lasting for a maximum of 20 s) in an ultrasonic bath.



NOTICE

Longer lasting ultrasonic baths can damage the ruby ball of the valve.

Only if this cleaning is useless, dismantle the valves separately and clean the components.

3 Dismantling valve

Dismantle every valve separately.



NOTICE

For dismantling the valve the tool for valve cartridges (6.2617.020) is required.




Figure 41 Components of the inlet valve and outlet valve

1	Inlet valve (6.2824.170)	2	Outlet valve (6.2824.160)
3	Inlet valve housing	4	Outlet valve housing
5	Sealing ring (black)	6	Sleeve
7	Sapphire sleeve The shiny side must point to ruby ball.	8	Ruby ball
9	Ceramic holder for ruby ball	10	Seal The larger opening must point outwards.

4 Cleaning the components of the valve

Rinse the valve components with ultrapure water and/or acetone and dry with a lint-free cloth.

5 Reassembling the valve

Reassemble valve components *according to figure 41, page 76*.

- Insert the seal with the larger opening faced downwards into the recess of the tool.
- Lay the other valve components above another in the correct sequence (*see figure 41, page 76*).

- ## 6 Checking the flow direction

If this is not the case, the valve has to be dismantled again and be reassembled correctly (*see figure 41, page 76*).

If by mistake, the inlet valve is mounted instead of the outlet valve, an extreme pressure builds up within the working cylinder, which can destroy the piston seal!

- Insert the inlet valve into the inlet valve holder the way the seal is visible.
- Screw the inlet valve holder into the bottom of the pump head and tighten with a wrench (39-4).
- Insert the outlet valve into the outlet valve holder the way the seal is visible.
- Screw the outlet valve holder into the top of the pump head and tighten with a wrench (39-1).

Changing the filter

The flow must be stopped before changing the filter.

1 Removing the inline filter

- Unscrew the pressure screws (42-**1**) from the inline filter.

2 Unscrewing the filter screw

- Screw the filter screw (42-4) out of the filter housing (42-2) with the aid of two adjustable wrenches (6.2621.000).

3 Inserting the filter

- Remove the old filter (42-**3**) using tweezers.
- Place the new filter (42-**3**) flat in the filter housing (42-**2**) using tweezers.

4 Mounting filter screw

- Screw the filter screw (42-4) back into the filter housing (42-2) and tighten by hand. Then additionally tighten slightly with two adjustable wrenches (6.2621.000).

5 Remounting the inline filter

- Screw pressure screws (42-**1**) back onto the inline filter.

6 Rinsing the inline filter

- Dismantle the guard column (if present) and the separation column and replace with a coupling (6.2744.040).
- Rinse the instrument with eluent.



4.7 Inline sample preparation

To protect the separation column (see chapter 2.22, page 58) against foreign particles which can affect the separating efficiency, we recommend that all samples undergo a microfiltration (filter 0.45 µm). The ultrafiltration cell can be used for **filtration** (see manual of the *IC Equipment for Ultrafiltration*).

Samples containing significant amounts of **gas** should be degassed. The sample degasser (*see chapter 2.13, page 38*) is used for degassing.

Matrix-loaded samples (e.g. blood, oil) should be prepared for the measurement by means of dialysis (see manual on the *IC Equipment for Dialysis*).

If the concentration of the sample is too high, the sample should be **diluted** before feeding (see documentation on the *IC Equipment for Sample Dilution*).

For the sample preparation methods **Neutralization** (replacement of e.g. Na^+ with H^+) and **cation exchange** (replacement of e.g. heavy metals with H^+), a sample preparation module (SPM) is used.

For an overview of all Metrohm inline sample preparation methods go to the following website: <http://misp.metrohm.com>

4.8 Rinsing the sample path

Before a new sample can be measured, the sample path must be rinsed with it so that the measuring result is not falsified by the previous sample (**Sample carry-over**).

In the case of automated sample feeding, the rinsing time should be at least 3 times the **transfer time**.

The transfer time is the time required by the sample to flow from the sample vessel to the end of the sample loop. The transfer time depends on the pump capacity of the peristaltic pump or the Dosino used, the total capillary volume and the volume of the gas removed by the sample degasser (if used) - in other words the amount of gas in the sample.

To ascertain the transfer time, proceed as follows:

Pump air through the sample path (pump tubing, tubing connections, capillary in the degasser, sample loop) for several minutes until all liquid is displaced by the air.

Aspirate a sample typical for the later application and use a stop watch to measure the time required by the sample to travel from the sample vessel to the end of the sample loop.

The time measured corresponds to the "transfer time". The rinsing time should be at least 3 times the transfer time.

It is possible to determine whether the rinsing time is adequate via a direct measurement of the sample carry-over. Proceed as follows:

- **Sample A:** A typical sample for the application.
- **Sample sample B:** Ultrapure water.

Let "Sample A" pass through the sample path for the duration of the rinsing time, then inject and measure.

Let "Sample B" pass through the sample path for the duration of the rinsing time, then inject and measure.

The degree of the sample carry-over corresponds to the ratio of the peak areas of the measurement for sample B to the measurement for sample A. The lower the ratio, the lower the sample carry-over. This ratio can be modified by varying the rinsing time – thus allowing the rinsing time required for the application to be ascertained.

4.9 Sample degasser

4.9.1 Operation

If you are working with sample degassing, the longer "transfer time" (see *Ascertaining the transfer time, page 81*) means that it is also necessary to rinse for longer (with the subsequent sample). The rinsing time should be at least 3 times the "transfer time" so as to minimize the carry-over effects. The "transfer time" depends on the pump capacity, the total capillary volume and the volume of the removed gas (i.e. the gas amount in the sample).



NOTICE

The rinsing time increases by at least 2 minutes when the sample degasser is used.

4.10 Injection valve

4.10.1 Protection

To prevent contamination of the injection valve, a 6.2821.120 inline filter (see chapter 2.11, page 36) should be mounted between the high pressure pump and the pulsation damper.

4.11 Metrohm Suppressor Module (MSM)

4.11.1 Protection

To protect the MSM against foreign particles or bacterial growth, a pump tubing connection with filter must be mounted between the peristaltic pump (*see chapter 2.17, page 48*) and the input capillaries of the MSM (*see figure 28, page 50*).

4.11.2 Operation



NOTICE

The suppressor units must never be regenerated in the same flow direction in which the eluent is pumped. Therefore always mount the input and output capillaries according to diagram outlined in *Figure 25, Page 46*.

The MSM consists of 3 suppressor units, which are used for suppression, regenerated with sulfuric acid or rinsed with ultra pure water in rotation.

In order to record every new chromatogram under comparable conditions, you should normally work with a freshly regenerated suppressor.



The MSM must never be switched over in a dry state, as there is a risk of jamming. If the MSM is in a dry state, the MSM must be rinsed for at least 5 minutes before it may be switched over.



In the case of reduced capacity or high backpressure, the MSM must be regenerated (see chapter 4.11.3.1, page 83), cleaned (see chapter 4.11.3.2, page 84) or replaced (see chapter 4.11.3.3, page 87).

4.11.3 Maintenance

4.11.3.1 Regenerating the MSM

If the suppressor units are loaded for a longer period with certain heavy metals (e.g. iron) or organic impurities, these can no longer be completely removed with the regeneration solution (50 mmol/L H_2SO_4) normally used. The capacity of the suppressor units is consequently affected, which can result in reduced phosphate sensitivity in less serious cases and a significant rise in the baseline in more serious cases. If such capacity problems occur at one or more positions, the suppressor units must be regenerated:

Regenerating the MSM

Regenerate the MSM as follows:

1 Disconnecting the MSM from the IC system

- Disconnect the MSM from the separation column and detector.

2 Regenerating the MSM



The pump tubing made of PVC must not be used for rinsing with solutions containing organic solvents. In this case, other pump tubing must be used for rinsing.



NOTICE

The high pressure pump can be used for regeneration. For this, remove the guard column and separation column and connect the capillary directly to the MSM (regenerate in opposite direction).

- Rinse the 3 suppressor units with the following solutions for about 15 minutes each:
 - **Contamination with heavy metals:**
1 mol/L H_2SO_4 + 0.1 mol/L oxalic acid
 - **Contamination with organic cationic complexing agents:**
0.1 mol/L H_2SO_4 / 0.1 mol/L oxalic acid / acetone 5%
 - **Severe contamination with organic substances:**
0.2 mol/L H_2SO_4 / acetone $\geq$ 20%

3 Connecting the MSM to the IC system

- Reconnect the MSM to the IC system. If the capacity problems remain, MSM rotor A must be replaced (*see chapter 4.11.3.3, page 87*).

4.11.3.2 Cleaning the MSM

It may be necessary to clean the MSM in the following cases:

- Increased backpressure onto the output capillaries of the MSM.
- Blockage of the MSM which cannot be eliminated (solutions can no longer be pumped through the MSM).
- Jamming of the MSM which cannot be eliminated (MSM can no longer be switched over).




Figure 43 MSM – Components

1	Union nut	2	MSM connecting piece 6.2832.010
3	MSM rotor A 6.2832.000	4	MSM housing
5	Slot in the MSM housing		

Cleaning the MSM

Clean the MSM as follows:

1 Disconnecting the MSM from the IC system

- Switch off the instrument.
- Disconnect the MSM from the separation column, peristaltic pump and detector.

2 Dismantling the MSM

- Unscrew union nut (43-**1**) from the MSM housing (43-**4**).
- Pull MSM connecting piece (43-**2**) and MSM rotor A (43-**3**) out of the MSM housing (43-**4**). The MSM connecting piece and MSM rotor A normally stick to one another - if this is not the case: Take a sharp object, insert into the slot (43-**5**) in the MSM housing, and pull out MSM rotor A (43-**3**) in this way.
- Detach the MSM connecting piece (43-**2**) from MSM rotor A (43-**3**).

3 Cleaning the supply and discharge lines

- Connect in turn each of the 6 capillary tubings fastened on the MSM connecting piece (43-**2**) on the high pressure pump (see chapter 2.10, page 32) and pump through ultra pure water.



- Check whether solution emerges at the MSM connecting piece (43-2). If one of the supply or discharge lines remained blocked, the SPM connecting piece (43-2) must be replaced (order number 6.2832.010).

4 Cleaning the MSM rotor A

- Clean sealing surface of the MSM rotor A (43-**3**) with ethanol using a lint-free cloth.

5 Inserting the MSM rotor A



CAUTION

An incorrectly inserted MSM rotor A (43-**3**) can be **destroyed** during start-up.

- Insert the MSM rotor A (43-3) into the MSM housing (43-4) in such a way that the tubing connections on the rear of MSM rotor A fit into the corresponding recesses inside the MSM housing and one of the three holes of the MSM rotor A is visible from below in the slot of the MSM housing (43-5).
- If the MSM rotor A (43-3) is correctly inserted, its sealing area will be approx. 4 mm within the MSM housing (43-4). If this is not the case, MSM rotor A must be moved into the right position from below using a sharp object (e.g. screwdriver).

6 Cleaning the MSM connecting piece

- Clean sealing surface of the MSM connecting piece (43-2) with ethanol using a lint-free cloth.

7 Inserting the MSM connecting piece

- Insert the MSM connecting piece (43-2) into the MSM housing (43-4) in such a way that the connector 1 is on top and the three pins of the MSM connecting piece fit into the corresponding recesses on the MSM housing (43-4).

8 Connecting and conditioning the MSM

- Reconnect the MSM to the IC system.
- Before switching the MSM over for the first time, rinse the three suppressor units with solution for 5 minutes.

4.11.3.3 Replacing parts of the MSM

It may be necessary to replace parts of the MSM in the following cases:

- Loss of suppression capacity which cannot be eliminated (reduced phosphate sensitivity and/or significant rise in the baseline).
- Blockage of the MSM which cannot be eliminated (solutions can no longer be pumped through the suppressor).

Both the MSM rotor A (43-**3**) and the MSM connecting piece (43-**2**) with its supply and discharge lines can be replaced.

Replacing parts of the MSM

Replace parts of the MSM as follows (see figure 43, page 85):

1 Disconnecting the MSM from the IC system

- Switch off the instrument.
- Disconnect the MSM from the separation column, peristaltic pump and detector.

2 Dismantling the MSM

- Unscrew union nut (43-1) from the MSM housing (43-4).
- Pull MSM connecting piece (43-2) and MSM rotor A (43-3) out of the MSM housing (43-4). The MSM connecting piece and MSM rotor A normally stick to one another - if this is not the case: Take a sharp object, insert into the slot (43-5) in the MSM housing, and pull out MSM rotor A (43-3) in this way.
- Detach the MSM connecting piece (43-2) from MSM rotor A (43-3).

3 Cleaning the new MSM rotor A

- Clean the sealing surface of the new MSM rotor A (43-**3**) with ethanol using a lint-free cloth.

4 Inserting the new MSM rotor A



CAUTION

An incorrectly inserted MSM rotor A (43-**3**) can be **destroyed** during start-up.



- Insert the new MSM rotor A (43-3) in the MSM housing (43-4) in such a way that the tubing connections on the rear of MSM rotor A fit into the corresponding recesses inside the MSM housing and one of the three holes of MSM rotor A is visible from below in the slot (43-5) of the MSM housing.
- If the MSM rotor A (43-3) is correctly inserted, its sealing area will be approx. 4 mm within the MSM housing (43-4). If this is not the case, MSM rotor A must be moved into the right position from below using a sharp object (e.g. screwdriver).

5 Cleaning the MSM connecting piece

- Clean the sealing surface of the MSM connecting piece (43-2) with ethanol using a lint-free cloth.

6 Inserting the new MSM connecting piece

- Insert the MSM connecting piece (43-2) into the MSM housing (43-4) in such a way that the connector 1 is on top and the three pins of the MSM connecting piece fit into the corresponding recesses on the MSM housing (43-4).

7 Connecting and conditioning the MSM

- Reconnect the MSM to the IC system.
- Before switching the MSM over for the first time, rinse the three suppressor units with solution for 5 minutes.

4.12 Peristaltic pump

4.12.1 Operation

The pumping capacity of the peristaltic pump depends on the drive speed (set via software), the contact pressure and, above all, the internal diameter of the pump tubing. Depending on the application, different pump tubings are used.



CAUTION

The service life of the pump tubings also depends on the contact pressure. Therefore fully lift the tubing cartridges by loosening the snap-action lever (27-**10**) on the right-hand side if the peristaltic pump is to be turned off for a longer period. Once set, the contact pressure remains unaffected.



4.12.2.1 Pump tubings

LFL pump tubings with three stoppers are inserted into the tubing cartridge in such a way that the cartridge is located between two stoppers. This results in two possible positions for the tubing cartridge. If the pump tubing was to show pronounced signs of wear, you can insert it a second time in the other position.

Pump tubing selection

The following table shows the properties and areas of application of the pump tubings:

Order number	Name	Material	Inner diameter	Use
6.1826.020	Pump tubing (blue/blue), 2 stoppers	PVC (Tygon® ST)	1.65 mm	Pump tubing for online IC instruments and automation in voltammetry.
6.1826.310	Pump tubing LFL (orange/green), 3 stoppers	PVC (Tygon®)	0.38 mm	Pump tubing for bromate determination using the triiodide method.
6.1826.320	Pump tubing LFL (orange/yellow), 3 stoppers	PVC (Tygon®)	0.48 mm	For acceptor solutions for Inline Dialysis and for Inline Ultrafiltration.
6.1826.330	Pump tubing LFL (orange/white), 3 stoppers	PVC (Tygon®)	0.64 mm	No special applications.

3 Mounting filter screw

- Screw the filter screw (44-**3**) back into the tubing olive (44-**1**) and tighten by hand. Then additionally tighten with two 6.2621.000 adjustable wrenches.

4.13 Servicing the detector

Follow the maintenance instructions in the the detector's manual.

4.14 Separation column

4.14.1 Separating efficiency

Which analysis quality can be attained, depends to a great extent on the separating efficiency of the separation column used. The separating efficiency of the selected separation column must be sufficient for the analysis problems present. If difficulties occur, you should always first check the quality of the separation column by recording a standard chromatogram.

You can find detailed information on the separation columns available from Metrohm in the leaflet provided along with your separation column, in the **Metrohm IC-Column Program** (available via your Metrohm agent) or in the Internet at <http://www.metrohm.com> in the product area Ion chromatography. You can request free information on special IC applications in the corresponding "**Application Bulletins**" or "**Application Notes**", which are available in the Internet at <http://www.metrohm.com> in the Applications area or via the Metrohm agent responsible.

4.14.2 Protection

To protect the separation column against foreign particles, which can affect the separating efficiency, we recommend that both the eluent and the samples undergo a microfiltration (filter 0.45 µm) before being aspirated via the aspiration filter (6.2821.090).

We recommend always to use a guard column (*see chapter 2.21, page 57*). This protects the separation column and considerably increases its service life. Information regarding which guard column is suitable for your separation column can be found in the **Metrohm IC Column Program** (which is available from your Metrohm agent), the leaflet provided along with your separation column, the product information on the separation column at <http://www.metrohm.com> (product area Ion Chromatography) or obtained directly from your agent.

The pulsation absorber (*see chapter 2.12, page 37*) must be installed in order to protect the column material from pressure concussion caused by injection.



4.14.3 Storage

Always store the separation columns sealed and filled according to the data of the column manufacturer when not using them.

4.14.4 Regeneration



NOTICE

The regeneration is considered as the last measure, and not to be carried out regularly.

If the separating properties of the column have deteriorated, the column can be regenerated according to the specifications of the column manufacturer. In the case of separation columns available from Metrohm, the specification for regeneration can be found on the leaflet provided along with each column.

5 Troubleshooting

5.1 Problems and their solutions

Problem	Cause	Remedy
Marked drop in pressure.	<i>Leak in the system.</i>	Check all capillary connections and seal leaks, if necessary (see chapter 2.5, page 14).
The baseline has a large amount of noise.	<i>High pressure pump – contaminated pump valves.</i>	Clean pump valves (see chapter 4.5.2, page 68).
	<i>Eluent – Leakage in eluent path.</i>	Check eluent path.
	<i>Eluent – Blockage in eluent path.</i>	Check eluent path.
	<i>High pressure pump – defective piston seals.</i>	Replace (see chapter 4.5.2, page 68) piston seals .
	<i>Pulsation damper not connected.</i>	Connect the pulsation damper (see Chapter 2.12, Page 37).
	<i>Pulsation absorber not connected. or defective.</i>	Connect pulsation absorber (see chapter 2.12, page 37). or replace it.
The baseline is drifting.	<i>Thermal equilibrium not yet attained.</i>	Condition instrument with the column thermostat (see Chapter 2.15, Page 42) switched on .
	<i>Leak in the system.</i>	Check all capillary connections and seal leaks, if necessary (see chapter 2.5, page 14).
	<i>Eluent – Evaporation of organic solvent in eluent.</i>	- Check the eluent bottle cap (see figure 12, page 28). - Stir the eluent.
The pressure in the system markedly increases.	<i>Inline filter (6.2821.120) blocked.</i>	Replace the filter (6.2821.130) (see chapter 4.6, page 78).
	<i>MSM – blocked.</i>	Regenerate the MSM (see chapter 4.11.3.1, page 83). Note: 6.2821.180 pump tubing connection with filter must be used (28-3).

Problem	Cause	Remedy
	<i>Peristaltic pump – filter blocked.</i>	Replace the filter (see chapter 4.12.2.2, page 90).
	<i>Peristaltic pump – pump tubing defective.</i>	Replace pump tubing (see chapter 4.12.2.1, page 89).
Data of the separation column cannot be read.	<i>Column chip contaminated.</i>	Clean the contact surfaces of the column chip (with alcohol).
	<i>Column chip defective.</i>	1. Save column configuration in MagIC Net™. 2. Notify Metrohm Service.
Individual peaks are greater than expected.	<i>Sample – carry-over of the samples from previous measurement.</i>	Rinse system longer between two samples.
MSM – No (or insufficient) pumping of regeneration or rinsing solution	<i>Leak in the system.</i>	Check connections.
	<i>Peristaltic pump – contact pressure too weak.</i>	Correctly set contact pressure (see "Set flow rate", page 52).
	<i>Peristaltic pump – filter blocked (see figure 28, page 50).</i>	Replace the filter (see "Replacing the filter", page 90).
	<i>MSM – backpressure too high.</i>	Clean the MSM (see chapter 4.11.3.2, page 84) or replace parts (see chapter 4.11.3.3, page 87).
	<i>Peristaltic pump – pump tubing defective.</i>	Replace pump tubing (see figure 27, page 49).
The background conductivity is too high.	<i>MSM – not connected.</i>	Connect the MSM (see chapter 2.16, page 45).
	<i>Wrong eluent.</i>	Change the eluent (see chapter 4.4.2.3, page 67).
	<i>MSM – regeneration or rinsing solution flow problems.</i>	Check the flow of the regeneration and rinsing solution (see chapter 2.16.2, page 45).
The retention times are poorly reproducible.	<i>Eluent – Leakage in eluent path.</i>	Check eluent path.
	<i>Eluent – Blockage in eluent path.</i>	Check eluent path.

6 Technical specifications

6.1 Reference conditions

The technical specifications listed in this chapter refers to the following reference conditions:

Ambient temperature	+25 °C (±3 °C)
Instrument status	> 40 minutes in operation (equilibrated)

6.2 Device

IC system	■ Metal-free IC system ■ Compact system with modular design ■ Up to two complete chromatographic systems in one housing
Material	Painted polyurethane hard foam without CFCs, fire class V0
Operating pressure range	■ 0...50 MPa (500 bar) high pressure pump ■ 0...35 MPa (350 bar) standard-PEEK system
Intelligent components	iPump, iDetector, iColumn, MagIC Net

6.3 Leak sensor

Type	Electronic, no calibration necessary
------	--------------------------------------

6.4 Ambient conditions

<i>Operation</i>	
<i>Ambient temperature</i>	+5...+45 °C
<i>Humidity</i>	20...80 % relative humidity
<i>Storage</i>	
<i>Ambient temperature</i>	-20...+70 °C
<i>Transport</i>	
<i>Ambient temperature</i>	-40...+70 °C

6.5 Housing

Dimensions

Width	365 mm
Height	642 mm
Depth	380 mm

Material of base tray, housing and flask holder	Polyurethane hard foam (PUR) with flame retardation for fire class V0, CFC-free, painted
---	--

Operating elements

Indicators	LED for power display
On/Off switch	On the rear of the instrument

6.6 Eluent degasser

Material	Fluoropolymer
Resistance to solvents	No restriction (except PFC)
Time to establish vacuum	< 60 s

6.7 High pressure pump

Type	▪ Serial dual-piston pump ▪ Intelligent pump head recognition ▪ Chemically inert ▪ Metal-free pump heads ▪ Materials in contact with the eluent: PEEK, ZrO₂, PTFE/PE ▪ Self-optimizing flow and pressure
------	--

Flow rate

Adjustable flow range	0.001...20.0 mL/min
Flow increment	1 µL/min
Reproducibility of the eluent flow	< 0.1 % deviation

Pressure range

<i>Pump</i>	0...50.0 MPa (0...500 bar)
<i>Pump head</i>	0...35.0 MPa (0...350 bar) (applies for the standard PEEK pump head)
<i>Residual pulsa- tion</i>	< 1 %

Safety shutdown

<i>Function</i>	Automatic shutdown upon reaching the pressure limit values
<i>Maximum pressure limit</i>	▪ Adjustable from 0.1...50 MPa (1...500 bar) ▪ The pump is automatically shut down at the first piston stroke above the maximum limit value
<i>Minimum pressure limit</i>	▪ Adjustable from 0...49 MPa (0...490 bar) ▪ The shutdown mechanism is inactive at 0 MPa ▪ The shutdown mechanism only becomes active 2 minutes after system start ▪ The pump is automatically shut down after 3 piston strokes below the minimum pressure limit

<i>Gradient capacity</i>	Isocratic or gradient (extendable to quaternary)
<i>Profile</i>	Step, linear, convex and concave
<i>Resolution</i>	< 1 nL/min flow increments

6.8 Sample degasser

<i>Material</i>	Fluoropolymer
<i>Resistance to solvents</i>	No restriction (except PFC)
<i>Time to establish vacuum</i>	< 60 s

6.9 Injection valve

<i>Actuator time</i>	typ.100 ms
<i>Max. operating pressure</i>	35 MPa (350 bar)
<i>Material</i>	PEEK

6.10 Column thermostat

Type	Peltier technique thermostat for two intelligent separation columns
Adjustable temperature range	0...+ 80 °C, in increments of 0.1 °C
Heating	Ambient temperature +50 °C
Cooling	Ambient temperature –20 °C
Temperature reproducibility	± 0.2 °C
Stability	< 0.05 °C
Heating up time	< 30 minutes from 20 to 50 °C
Cooling time	< 40 minutes from 50 to 20 °C

6.11 Metrohm Suppressor Module (MSM)

<i>Resistance to sol-vents</i>	No restriction
<i>Switching duration</i>	typ.100 ms
<i>Operating pres-sure</i>	2.5 MPa (25 bar), valve function prevents damage at overpressure

6.12 Peristaltic pump

Type	2-channel peristaltic pump
Rotating direction	Counterclockwise/Clockwise rotation
Rotational speed	0...42 rpm in 7 stages at 6 rpm.
Pumping properties	0.3 mL/min at 18 rpm; with 6.1826.320 standard pump tubing.
Material of pump tubings	recommended: Tygon Long Flex Life

6.13 Detector

You can find the technical specifications for the detector in the manual for the detector.

6.14 Power connection

<i>Required supply voltage</i>	100 - 240 V $\pm$ 10% (autosensing)
<i>Required frequency</i>	50 - 60 Hz $\pm$ 3 Hz (autosensing)
<i>Power consumption</i>	■ 65 W for typical analysis application ■ 25 W standby (conductivity detector to 40 °C)
<i>Power supply unit</i>	■ Up to 300 W maximum, electronically monitored ■ Internal fuse 3.15 A

6.15 Interfaces

USB

Entry	1 USB upstream, type B (for connection to the PC)
Exit	2 USB downstream, type A
USB	2 MSB 8-pin Mini DIN (female) (for Dosino, stirrer, remote lines, etc.)

 MSB 

CAUTION

When connecting an instrument to the MSB connector you **must** switch off the 850 Professional IC.

<i>Detector</i>	2 15-pin high-density DSUB (female)
<i>Column recognition</i>	3 (including 2 in the column thermostat (<i>see chapter 2.15, page 42</i>))
<i>Leak sensor</i>	1 jack plug
<i>Further connections</i>	■ 1 15-pin DSUB (female)

Index

(6.2821.130) Filter 78
6.2821.090 aspiration filter 67

A

Ambient conditions 98
Anion channel 10
Aspiration filter 6.2821.090 67
Aspiration tubing for eluent 26

B

Baseline
 Conditioning 63
 Unstable 68
Blood 80

C

Cable feed-throughs 23
Capillaries
 Installation 14
Capillary feed-throughs 23
Carry-over 80
Cation channel 10
Cleaning
 MSM 85
 Valves of the high pressure pump 73
Column
 see also "Separation column"
 58
Column recognition 102
Column thermostat
 Column thermostat 101
 Installation 42
Computer connection 55
Conditioning 63
Conductivity detector
 Capillary connection 53
Connect
 Power grid 56
 To computer 55
Connection
 Mains 102
Connections
 Installation 14
Contamination
 High pressure pump 67
 Valves in the high pressure pump 68
Contamination MSM
 Heavy metals 83

Crystallization
 High pressure pump 67

D

Deaerating
 High pressure pump 34
 Purge valve 32
Degasser
 Eluent degasser 30
 Sample degasser 38
Degassing
 Eluent 30
Detector
 Cable connector 20
 Conductivity detector 53
 Interface 102
 Positioning 20
Diagram 10
Dilution 80
Dimensions 99
Door 66
Drainage tubings
 Installation 21

E

Electrostatic charge 7
Eluent
 Aspirate 26
 Change 67
 Production 66
Eluent bottle
 Figure 29
 Installation 26
 Operation 67
Eluent degasser
 Installation 30
 Technical specifications 99
Equilibration 61, 63

F

Feed-throughs
 Capillaries 23
Fill
 Injection valve 42
Filter
 see also "inline filter" 36
Filter (6.2821.130) 78
Filter 6.2821.090
 Aspiration filter 67
Flow diagram 10

Flow fluctuations 68
Flow increment 99
Flow range 99
Flow rate 99
Frequency 102

G

Gas 30, 38
Guard column
 Installation 57
 Rinsing 58

H

Handle 17
Heating
 see also "Column thermostat"
 42
Heavy metals
 Contamination of the MSM 83
High pressure pump
 Installation 32
 Maintenance 67
 Protection 20, 67
 Technical specifications 99
 Tubing connection 32
 Valves 76
Housing 99
Humidity 98

I

IC column
 see also "Separation column"
 58
Impurity MSM
 Organic 83
Initial installation 9
Inject
 Injection valve 42
Injection valve 3
 Fill 42
 Inject 42
 Installation 40, 100
 Maintenance 82
 Protection 82
Inline filter 36
Inline sample preparation 80
Installation
 Column thermostat 42
 Conductivity detector 53
 Connections 14

Drainage tubings	21
Eluent bottle	26
Eluent degasser	30
Guard column	57
High pressure pump	32
Initial installation	9
Injection valve	40, 100
Leak sensor	20
MSM	45
Peristaltic pump	49
Pulsation absorber	37
Pump tubings	49
Sample degasser	38
Separation column	58
Interface	
MSB	102
USB	102
Interfaces	102
Further connections	102
Leak sensor	102

L

Leak	68
Leak sensor	
Installation	20
Interface	102
Technical specifications	98
Leak-tightness	62, 63
Leaking piston seals	68
Loop	
see also "Sample loop"	42

M

Maintenance	
High pressure pump	67
Injection valve	82
MSM	82
Peristaltic pump	88
Pump head	68
Material	99
MPak	
Holder	19
MSB	102
MSM	
Cleaning	85
Installation	45
Maintenance	82
Operation	82
Protection	82
Regeneration	83
Replacing parts	87
Switchover	83
Technical data	101

O

Oil	80
Operation	
MSM	82
Peristaltic pump	88
Sample degasser	82
Organic impurities	
MSM	83

P

Peristaltic pump	
Installation	49
Maintenance	88
Operation	88
Principle	48
Technical specifications	101
Piston seal	68
Pistons of the high pressure pump	68
Power connection	56, 102
Power consumption	102
Power supply unit	102
Precipitates	67
Pressure limit	100
Pressure range	100
Pressure screws	
Connection	15
Protection	
Injection valve	82
Inline filter	36
MSM	82
Pulsation	68
Pulsation absorber	
Installation	37
Pump head	
Maintenance	68
Pump tubings	
Installing	49
Overview	89
Service life	88
Purge valve	32

R

Reference conditions	98
Regeneration	64
MSM	83
Rinsing	
Guard column	58
Pump tubings	89
Sample path	80
Separation column	60
Rinsing time	81
Rise in pressure	67
Rollers	17

S

Safety instructions	6
Safety shutdown	100
Sample	
Carry-over	80
Sample loop	42
Transfer time	81
Sample degasser	
Installation	38
Operation	82
Technical specifications	100
Sample loop	42
Sample path	
Rinsing	80
Sample preparation	80
Screws	
Connection	15
Separation column	
Installation	58
Protection	3, 37, 91
Regeneration	92
Rinsing	60
Separating efficiency	91
Storage	92
Service	6, 64
Shutting down	65
Start-up	61
Storage	98
Supply voltage	6, 102
Suppressor	
Maintenance	82
Operation	82
see also "MSM"	45

T

Technical data	
Column thermostat	101
MSM	101
Technical specifications	
Detector	102
Eluent degasser	99
High pressure pump	99
Interfaces	102
Leak sensor	98
Peristaltic pump	101
Reference conditions	98
Sample degasser	100
Temperature	98
Thermostat	
see also "Column thermostat"	
.....	42
Transfer time	81
Transport	98
Rollers	17



Transport locking screws 20

Tubing 10

Tubing pump

 See also "Peristaltic pump" . 48

Tubings

 Installation 14

U

USB 102

V

Vacuum pump

 Protection 20

Valve

 see also "Injection valve" 40

Valves of the high pressure pump

..... 76