

916 Ti-Touch



Manual

8.916.8005EN / v11 / 2024-10-30



Metrohm AG
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

916 Ti-Touch

Program version 5.916.0046

Manual

This documentation is protected by copyright. All rights reserved.

This documentation is an original document.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g. conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

Melody for the BEEP command: excerpt from "En Altfrentsche", with kind permission of the Laseyer Quartett, Appenzell.

Table of contents

1	Introduction	1
1.1	Device description	1
1.1.1	Titration and measuring modes	1
1.1.2	Connectors	3
1.2	Displaying accessories	3
1.3	About the documentation	4
1.3.1	Symbols and conventions	4
2	Safety	6
2.1	Intended use	6
2.2	Responsibility of the operator	6
2.3	Requirements for operating personnel	7
2.4	Safety instructions	7
2.4.1	General notes on safety	7
2.4.2	Electrical safety	7
2.4.3	Tubing and capillary connections	8
2.4.4	Flammable solvents and chemicals	9
2.4.5	Recycling and disposal	9
3	Overview of the instrument	10
3.1	Front of the instrument	10
3.2	Rear of the instrument	11
4	Installation	12
4.1	Setting up the instrument	12
4.1.1	Packaging	12
4.1.2	Checks	12
4.1.3	Location	12
4.2	Connecting the power supply unit	12
4.3	Mounting the electrode holder	14
4.4	Connecting the propeller stirrer (depending on the product version)	15
4.5	Connecting MSB devices	16
4.5.1	Connecting a dosing device	17
4.5.2	Connecting an additional stirrer or titration stand	18
4.5.3	Connecting a Remote Box	18
4.6	Connecting USB devices	19
4.6.1	General	19
4.6.2	Connecting a USB hub	19

7.4	Measured value display	60
7.5	Acoustic signals	60
8	Titrants	61
8.1	Adding a new titrant	62
8.2	Editing titrant data	63
8.3	Monitoring the working life	65
8.4	Dosing unit	66
8.4.1	Parameters for preparing (PREP) and emptying (EMPTY)	67
8.4.2	Tubing parameters	69
8.4.3	Shift direction of the valve disk	72
8.5	Exchange unit	72
8.5.1	Parameters for the preparation (PREP)	74
8.5.2	Tubing parameters	75
8.6	GLP test for exchange unit and dosing unit	76
8.7	Titer determination options and data	78
8.7.1	Titer validity	78
8.7.2	Properties of the previous titer determinations	79
9	Reagents	81
9.1	Editing reagent data	81
9.2	Reagent monitoring	82
10	Sensors	85
10.1	Adding a new sensor	86
10.2	Editing the sensor data	87
10.3	Monitoring the working life	88
10.4	Calibration data (for pH electrodes only)	89
10.4.1	Properties of the previous calibrations	91
10.5	Limit values for the calibration data	93
10.6	Monitoring the calibration interval (for pH electrodes only)	94
11	Device manager	96
11.1	Adding a new device	97
11.2	Configuring the device	97
11.3	Ti-Touch	98
11.3.1	E-mail	99
11.3.2	PC/LIMS report	100
11.3.3	Shared memory	102
11.3.4	TCP/IP settings	103

12 File manager 137

14.2	Properties of common variables	159
14.3	Monitoring validity	160
14.4	Assigning a result automatically to a common variable	161
15	Templates	163
15.1	Sample data	163
15.1.1	Sample identification list	164
15.1.2	Sample assignment table	165
15.2	Custom result templates	167
15.2.1	Editing result templates	168
15.3	Input lines	171
15.3.1	Editing the input signal	172
15.4	Output lines	173
15.4.1	Editing the output signal	175
15.5	Custom calibration buffers	176
15.5.1	Defining calibration buffers	176
15.6	Report header	178
15.7	Custom electrode type	179
15.7.1	Limit values for the electrode rating	180
16	Methods	183
16.1	Creating a new method	183
16.2	Saving a method	184
16.3	Loading a method	185
16.4	Editing a method	186
16.4.1	Inserting a command	187
16.5	Method options	188
16.5.1	Start options	189
16.5.2	Stop options	191
16.5.3	Sample data	192
16.5.4	Method properties	196
16.5.5	Note	197
16.5.6	Saving a determination automatically	198
17	Electronic signatures	200
17.1	Signing methods/determinations electronically	200
17.2	Deleting electronic signatures	201
18	Control	202

25 Statistics	247
25.1 Displaying details for a result	249
25.2 Deleting statistical data	250
25.3 Adding a determination to a determination series	250
26 Result table	252
26.1 Properties	253
26.2 Saving the result table	257
26.3 Loading the result table	257
27 Printing	258
27.1 General report options	260
27.2 Settings of the individual reports	261
27.3 List of all printable reports	262
28 Manual control	266
28.1 Opening and closing the manual control	267
28.2 Measuring	268
28.2.1 Parameter description	269
28.3 Dosing	270
28.3.1 Continuous dosing	273
28.3.2 Dosing fixed volumes	274
28.3.3 Preparing	277
28.3.4 Emptying	278
28.3.5 Filling	278
28.3.6 Replacing reagent	278
28.4 Stirring	279
28.5 Remote	281
28.6 Sample Processor	282
28.6.1 Moving the lift	284
28.6.2 Moving to a rack position	286
28.6.3 External positions	287
29 Parameters	292
29.1 Dynamic equivalence point titrations (DET)	292
29.1.1 Start conditions	292
29.1.2 Titration parameters	295
29.1.3 Stop conditions	299
29.1.4 Potentiometric evaluation	300
29.1.5 Control instrument	305
29.1.6 Sensor	305
29.1.7 Dosing device	307

X ■■■■■■

29.8 Electrode test for pH electrodes (ELT)	374
29.8.1 Electrode test parameters	375
29.8.2 Control instrument	376
29.8.3 Sensor	376
29.8.4 Stirrer	377
29.8.5 Additional information	378
29.9 Evaluations (EVAL)	380
29.9.1 Fixed endpoint evaluation (EVAL FIX-EP)	381
29.9.2 pK value and half neutralization potential evaluation (EVAL pK/HNP)	382
29.9.3 Minimum and maximum evaluation (EVAL MIN/MAX)	383
29.9.4 Break point evaluation (EVAL BREAK)	385
29.9.5 Rate evaluation (EVAL RATE)	387
29.10 Calculations	388
29.10.1 Calculations (CALC)	388
29.10.2 Calculations (CALC LIVE)	396
29.10.3 Formula editor	398
29.11 Reports (REPORT)	404
29.11.1 General report options	404
29.11.2 Settings of the individual reports	405
29.11.3 List of reports	406
29.12 Dosing and Liquid Handling	409
29.12.1 Preparing an exchange or dosing unit (PREP)	409
29.12.2 Emptying a dosing unit (EMPTY)	410
29.12.3 Dosing a specified volume (ADD)	411
29.12.4 Liquid Handling (LQH)	416
29.13 Communication	419
29.13.1 Scanning remote lines (SCAN)	419
29.13.2 Setting remote lines (CTRL)	421
29.13.3 Scanning the RS-232 interface (SCAN RS)	422
29.13.4 Defining RS-232 commands (CONTROL RS)	423
29.14 Automation	424
29.14.1 Rotating sample rack (MOVE)	424
29.14.2 Moving the lift (LIFT)	426
29.14.3 Controlling pumps (PUMP)	427
29.14.4 Resetting the rack (RACK)	428
29.14.5 Defining the sample variable (SAMPLE)	428
29.14.6 Creating a subsequence (SUBSEQ)	429
29.15 Miscellaneous commands	432
29.15.1 Controlling a stirrer (STIR)	432
29.15.2 Pausing the method run (WAIT)	433
29.15.3 Scan data (REQUEST)	434
29.15.4 Defining an acoustic signal (BEEP)	435
29.15.5 Signing a determination (SIGN)	435
29.15.6 Canceling the method run (END)	435

XII ■■■■■■

32.9.11	Precisa	465
32.9.12	Merck CertiPUR 20 / Titrisol	466
32.9.13	Merck CertiPUR 25	467
32.10	Using AuditTrailViewer	469
32.10.1	Installing AuditTrailViewer	469
32.10.2	Opening the Audit Trail	470
32.10.3	Contents of the Audit Trail	471
32.10.4	Filtering the Audit Trail	472
32.10.5	Exporting the Audit Trail	473
32.11	Diagnosis	473
32.11.1	LCD test	474
32.11.2	Formatting an external storage medium	474
32.11.3	Removing an external storage medium	475
32.11.4	Adjusting the touch screen	475
32.11.5	Testing the touch screen	476
32.11.6	Software update (loading program versions and language files)	477
32.11.7	Service	480
32.12	Arithmetic algorithms in the Ti-Touch	481
33	Technical specifications	483
33.1	Touch screen	483
33.2	Measuring inputs	484
33.2.1	Potentiometry	484
33.2.2	Polarizer	484
33.2.3	Temperature	485
33.3	Built-in magnetic stirrer (depending on the product version)	486
33.4	Interfaces	486
33.5	Energy supply	487
33.6	Declarations of Conformity	487
33.7	Ambient conditions	488
33.8	Reference conditions	488
33.9	Dimensions	488
33.10	Storage capacity	489
Index		490

1 Introduction

1.1 Device description

The 916 Ti-Touch is a compact titration system for volumetric titration. This titrator combines in a single instrument the touch-sensitive color monitor for convenient and efficient operation, the titration unit and as well as the magnetic stirrer or an integrated stirrer interface for a propeller stirrer, depending on the product version. The upper side of the housing offers space for the titrant and the titration vessel. The titrator is standard-equipped for operation with an external dosing drive of the *800 Dosino* type with a dosing unit. You can, however, also use a *805 Dosimat* with an exchange unit. Thanks to its compact construction, you can use the 916 Ti-Touch in a small space as a stand-alone titrator.

You manage titrants, sensors, methods, etc. conveniently in the internal memory of the 916 Ti-Touch. You can also save your files externally through the USB connector, e.g. on a USB flash drive. On this storage medium you can not only store your methods and determinations, but also create a backup together with all of the data and settings of your system.

The integrated Ethernet connector is available to you should you wish to connect your 916 Ti-Touch to a network. The network connection offers you the following advantages:

- Saving data to a computer within the network
- Printing reports on a network printer
- Sending displayed messages as e-mails

1.1.1 Titration and measuring modes

The 916 Ti-Touch supports the following titration and measuring modes.

- **DET**

Dynamic equivalence point titration. The reagent is added in variable volume steps.

Measuring modes:

- **pH** (pH measurement)
- **U** (potentiometric voltage measurement)
- **Ip_{ol}** (voltametric measurement with selectable polarization current)
- **Up_{ol}** (amperometric measurement with selectable polarization voltage)



- **MET**

Monotonic equivalence point titration. The reagent is added in constant volume steps.

Measuring modes:

- **pH** (pH measurement)
- **U** (potentiometric voltage measurement)
- **I_{pol}** (voltametric measurement with selectable polarization current)
- **U_{pol}** (amperometric measurement with selectable polarization voltage)

- SET

Endpoint titration at one or two specified endpoints.

Measuring modes:

- **pH** (pH measurement)
- **U** (potentiometric voltage measurement)
- **I_{pol}** (voltametric measurement with selectable polarization current)
- **U_{pol}** (amperometric measurement with selectable polarization voltage)

- **STAT**

Titration with constant maintenance of the measured value.

Measuring modes:

- **pH** (pH measurement)
- **U** (potentiometric voltage measurement)

- MAT

Manual titration with manual dosing and manual ending.

Measuring modes (optional):

- **pH** (pH measurement)
- **U** (potentiometric voltage measurement)

- **MEAS**

Measuring modes:

- **pH** (pH measurement)
- **U** (potentiometric voltage measurement)
- **I_{pol}** (voltametric measurement with selectable polarization current)
- **U_{pol}** (amperometric measurement with selectable polarization voltage)
- **T** (temperature measurement)

- CAL

pH electrode calibration.

Measuring mode:

- **pH** (calibration of pH electrodes)
- **ELT** (Electrode test for pH electrodes)

1.1.2 Connectors

The 916 Ti-Touch is equipped with the following connectors:

- **Power connection**
For connecting to the power grid using the power supply unit provided.
- **Two MSB connectors (Metrohm Serial Bus)**
For connecting dosing devices, one additional stirrer or a Remote Box.
- **USB connector**
For connecting peripheral devices (printer, PC keyboard, etc.), a USB flash drive, a Sample Processor or a USB hub.
- **Sensor connectors**
One connector each for:
 - Potentiometric electrodes (pH, ISE, metal)
 - Reference electrodes
 - Polarizable electrodes
 - Intelligent electrodes (iTrodes)
 - Temperature sensor (Pt1000 or NTC)
- **Stirrer connector (depending on the product version)**
For connecting a propeller stirrer (only available for the version without magnetic stirrer).
- **Ethernet connector**
For connecting the Ti-Touch to a network.

1.2 Displaying accessories

Up-to-date information on the scope of delivery and on optional accessories can be found on the Metrohm website.

1 Searching for a product on the website

- Go to <https://www.metrohm.com>.
- Click on .
- Enter the article number of the product (e.g. **2.1001.0010**) into the search field and press **[Enter]**.

The search result is displayed.

2 Displaying product information

- To display the products matching the search term, click on **Product models**.
- Click on the desired product.

Detailed information regarding the product is displayed.

3 Displaying accessories and downloading the accessories list

- To display the accessories, scroll down to **Accessories and more**.
 - The **scope of delivery** is displayed.
 - Click on **[Optional parts]** for the optional accessories.
- To download the accessories list, click on **[Download accessories PDF]** under **Accessories and more**.



NOTE

Metrohm recommends keeping the accessories list for reference purposes.

1.3 About the documentation



CAUTION

Please read through this documentation carefully before putting the device into operation. The documentation contains information and warnings which the user must follow in order to ensure safe operation of the device.

1.3.1 Symbols and conventions

The following symbols and formatting may appear in this documentation:

(5-12)	Cross-reference to figure legend
	The first number refers to the figure number, the second to the instrument part in the figure.
1	Instruction step Perform the steps one after the other.
Method	Dialog text, parameter in the software
File ► New	Menu or menu item
[Continue]	Button or key
	WARNING This symbol draws attention to a possible life-threatening hazard or risk of injury.



WARNING

This symbol draws attention to a possible hazard due to electrical current.



WARNING

This symbol draws attention to a possible hazard due to heat or hot instrument parts.



WARNING

This symbol draws attention to a possible biological hazard.



WARNING

Warning of optical radiation



CAUTION

This symbol draws attention to possible damage to instruments or instrument parts.



NOTICE

This symbol highlights additional information and tips.

2.3 Requirements for operating personnel

Only qualified personnel may operate the product. Qualified personnel are persons who meet the following requirements:

- Basic regulations on occupational safety and accident prevention for chemical laboratories are known and complied with.
- Knowledge of handling hazardous chemicals is present. Personnel have the ability to recognize and avoid potential dangers.
- Knowledge regarding the application of fire prevention measures for laboratories is available.
- Safety-relevant information is communicated and understood. The personnel can operate the product safely.
- The user documentation has been read and understood. The personnel operate the product according to the instructions in the user documentation.

2.4 Safety instructions

2.4.1 General notes on safety



WARNING

Operate this instrument only according to the information contained in this documentation.

This instrument left the factory in a flawless state in terms of technical safety. To maintain this state and ensure non-hazardous operation of the instrument, the following instructions must be observed carefully.

2.4.2 Electrical safety

The electrical safety when working with the instrument is ensured as part of the international standard IEC 61010.



WARNING

Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.



WARNING

Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.

Supply voltage



WARNING

An incorrect supply voltage can damage the instrument.

Only operate this instrument with a supply voltage specified for it (see rear panel of the instrument).

Protection against electrostatic charges



WARNING

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Do not fail to pull the power cord out of the power socket before you set up or disconnect electrical plug connections at the rear of the instrument.

2.4.3 Tubing and capillary connections



CAUTION

Leaks in tubing and capillary connections are a safety risk. Tighten all connections well by hand. Avoid applying excessive force to tubing connections. Damaged tubing ends lead to leakage. Appropriate tools can be used to loosen connections.

Check the connections regularly for leakage. If the instrument is used mainly in unattended operation, then weekly inspections are mandatory.

2.4.4 Flammable solvents and chemicals




WARNING

All relevant safety measures are to be observed when working with flammable solvents and chemicals.

- Set up the instrument in a well-ventilated location (e.g. fume cupboard).
- Keep all sources of flame far from the workplace.
- Clean up spilled liquids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

2.4.5 Recycling and disposal



Properly dispose of chemicals and of the product to reduce negative effects on the environment and public health. Local authorities, waste disposal companies or dealers provide more detailed information on disposal. Observe the WEEE EU directive (WEEE = Waste Electrical and Electronic Equipment) for the proper disposal of waste electronic equipment within the European Union.

3 Overview of the instrument

3.1 Front of the instrument




Figure 1 Front 916 Ti-Touch

1 Display Touch screen.	2 Fixed key [Home] Opens the main dialog.
3 Fixed key [Back] Saves the entry and opens the next-higher dialog page.	4 Fixed key [Help] Opens the online help for the dialog displayed.
5 Fixed key [Print] Opens the print dialog.	6 Fixed key [Manual] Opens the manual control.
7 Fixed key [STOP] Cancels the running determination.	8 Fixed key [START] Starts a determination.
9 Bottle holder With holding clips, for reagent bottle.	10 Support rod (lower part) For mounting the support rod (upper part).
11 Titration stand For placing the titration vessel (with built-in magnetic stirrer, depending on the product version).	

3.2 Rear of the instrument





Figure 2 Rear 916 Ti-Touch

- | | |
|---|---|
| 1 Electrode connector (Ind.)
For connecting pH, metal or ion-selective electrodes with integrated or separate reference electrode. Socket F. | 2 Electrode connector (Ref.)
For connecting reference electrodes, e.g. Ag/AgCl reference electrode. Socket B, 4 mm. |
| 3 USB connector (type A)
For connecting a printer, USB flash drive, USB hub, Sample Processor etc. | 4 Ethernet connector (RJ-45)
For connecting to a network. |
| 5 Temperature sensor connector (Temp.)
For connecting temperature sensors (Pt1000 or NTC). Two B sockets, 2 mm. | 6 Electrode connector (Pol.)
For connecting polarizable electrodes, e.g. Pt wire electrodes. Socket F. |
| 7 MSB connector (MSB 1 and MSB 2)
Metrohm Serial Bus. For connecting external dosing devices, one additional stirrer or a Remote Box. Mini DIN, 8-pin. | 8 Power socket (Power)
For connecting the external power supply unit. |
| 9 Power switch
Switches the instrument on/off. | 10 Electrode connector (iConnect)
For connecting electrodes with integrated data chip (iTrodes). |
| 11 Stirrer connector (present, depending on the product version)
For connecting the propeller stirrer (802 Stirrer) to instruments without built-in magnetic stirrer. | 12 Type plate
Contains the serial number. |



- 1 Connect the plug of the external power supply unit with the power socket of the Ti-Touch (*see figure 3, page 13*).



The plug of the power supply unit is protected against accidental disconnection of the cable by means of a pull-out protection feature. If you wish to pull out the plug, you first need to pull back the outer plug sleeve marked with arrows.

- 2 Connect the power cord with the external power supply unit of the Ti-Touch and with the power grid.



Switch off the Ti-Touch correctly by pressing the power switch before you disconnect the electricity supply. If this is not done, then there is a danger of data loss.

4.4 Connecting the propeller stirrer (depending on the product version)

The 916 Ti-Touch is equipped with a built-in magnetic stirrer or a built-in stirrer connector, depending on the product version. You can connect the 802 Stirrer propeller stirrer to the stirrer connector.



NOTE

Before the rod stirrer is connected, the electrode holder with the clamping ring must be mounted on the support rod, see (see chapter 4.3, page 14).

Proceed as follows:

- 1 Plug in the connection cable of the 802 Stirrer propeller stirrer to the stirrer connector of the Ti-Touch (2-**11**).




Figure 4 Connecting the propeller stirrer

- 2 Insert the propeller stirrer 802 Stirrer without the stirring propeller from above into the center opening of the electrode holder.
- 3 Plug the stirring propeller (6.1909.010) from below to the propeller stirrer.



- No additional stirrer may be connected to the MSB 1!
- Only one device of the same type may be used at a single MSB connector at one time.
- When making the connection, take care to ensure that the flat part of the MSB plug marked with arrows is pointing in the direction of the marking on the MSB connector (*see figure 6, page 17*).



MSB connections can be extended with the 6.2151.010 cable. The maximum connection length permitted is 6 m.

4.5.1 Connecting a dosing device

You can connect two dosing devices to the Ti-Touch.

The types of dosing devices that are supported are:

- 800 Dosino
- 805 Dosimat

Proceed as follows:

- 1** Switch off the Ti-Touch.

- 2 Connect the dosing device connection cable to an MSB connector (2-7) on the rear side of the Ti-Touch.
- 3 Switch on the Ti-Touch.

4.5.2 Connecting an additional stirrer or titration stand

Depending on the product version, in addition to the built-in magnetic stirrer or the built-in stirrer connector for the propeller stirrer 802, you can also use the magnetic stirrers 801 Stirrer, 803 Ti Stand or the 804 Ti Stand with the propeller stirrer 802 Stirrer .

Proceed as follows:

- 1 Switch off the Ti-Touch.
- 2 Connect the connection cable of the magnetic stirrer or of the titration stand to MSB 2 (2-**7**) on the rear of the Ti-Touch.
- 3 Switch on the Ti-Touch.

4.5.3 Connecting a Remote Box

Instruments that are controlled via remote lines and/or that send control signals via remote lines can be connected via the 6.2148.010 Remote Box. In addition to Metrohm, other instrument manufacturers also use similar connectors that make it possible to connect different instruments together. These interfaces are also frequently given the designations "TTL Logic", "I/O Control" or "Relay Control" and they generally have a signal level of 5 volts.

Control signals are understood to be electrical line statuses or brief (> 200 ms) electrical pulses which display the operating status of a device or which trigger or report an event. Sequences on a variety of instruments can thus be coordinated in a single complex automation system. However, no exchange of data is possible.

Proceed as follows:

- 1 Switch off the Ti-Touch.
- 2 Connect the Remote Box connection cable to an MSB connector (2-7) on the rear side of the Ti-Touch.
- 3 Switch on the Ti-Touch.

You can connect an 869 Compact Sample Changer. The Remote Box also has an MSB socket at which a further MSB device, e.g. a dosing device, can be connected.

You will find precise information concerning the pin assignment of the interface on the Remote Box in Appendix (*see chapter 32.6, page 447*).

4.6 Connecting USB devices

4.6.1 General

The 916 Ti-Touch has a USB connector (type A socket) for peripheral devices with USB interface and for Sample Processors. If you wish to connect more than one device to the USB, you can use a commercially available USB hub.



NOTE

We recommend that the Ti-Touch be switched off while you set up or disconnect connections between the devices.

4.6.2 Connecting a USB hub

Use a USB hub with its own power supply.

Connect the USB hub as follows:

- 1 With the help of the 6.2151.030 cable (length 0.6 m) or the 6.2151.020 cable (length 1.8 m), connect the USB connector of the Ti-Touch (Type A) with the USB connector of the hub (Type B, see manual for the USB hub).

The USB hub is recognized automatically.

4.6.3 Connecting a printer

Printers connected to the 916 Ti-Touch must meet the following requirements:

- Printer languages: HP-PCL, HP-PCL-GUI, Canon BJL Commands or Epson ESC P/2
- Paper size: A4 or Letter.

Connect the printer as follows:

- 1 With the aid of the 6.2151.020 cable, connect the USB connector of the Ti-Touch (type A) with the USB connector of the printer (type B, see manual for the printer).

- 2 Configure the printer in the device manager of the Ti-Touch (see *chapter 11.7, page 125*).

4.6.4 Connecting a balance

If you want to connect a balance to the Ti-Touch you need a USB/RS-232 adapter (6.2148.050).

The following table offers an overview of the balances that you can use together with the Ti-Touch and of which cable you will need for connection to the RS-232 interface:

Balance	Cables
AND ER, FR, FX with RS-232 interface (OP-03)	6.2125.020 + 6.2125.010
Mettler AB, AG, PR (LC-RS9)	In the scope of delivery for the balance
Mettler AM, PM, PE with interface option 016 or Mettler AJ, PJ with interface option 018	6.2146.020 + 6.2125.010 Also from Mettler: ME 47473 adapter and either ME 42500 hand switch or ME 46278 foot switch
Mettler AT	6.2146.020 + 6.2125.010 Also from Mettler: ME 42500 hand switch or ME 46278 foot switch
Mettler AX, MX, UMX, PG, AB-S, PB-S, XP, XS	6.2134.120
Mettler AE with interface option 011 or 012	6.2125.020 + 6.2125.010 Also from Mettler: ME 42500 hand switch or ME 46278 foot switch
Ohaus Voyager, Explorer, Analytical Plus	Cable AS017-09 from Ohaus
Precisa balances with RS-232-C interface	6.2125.080 + 6.2125.010
Sartorius MP8, MC, LA, Genius, Cubis	6.2134.060
Shimadzu BX, BW	6.2125.080 + 6.2125.010

Connect the balance as follows:

- 1** Connect the USB plug of the USB/RS-232 adapter with the USB connector of the Ti-Touch (Type A).
The USB/RS-232 adapter will be recognized automatically and entered in the device manager of the Ti-Touch.
- 2** Connect the RS-232 interface of the USB/RS-232 adapter with the RS-232 interface of the balance (see table for cable).
- 3** Switch on the balance.
- 4** If necessary, switch on the RS-232 interface of the balance.
- 5** Configure the RS-232 interface of the USB/RS-232 adapter in the device manager of the Ti-Touch (*see chapter 11.9, page 131*).
- 6** Enter and configure the balance in the device manager of the Ti-Touch (*see chapter 11.8, page 129*).
- 7** Make sure that the parameters of the USB/RS-232 adapter configured in the device manager match those of the balance.

4.6.5 Connecting a PC keyboard

The PC keyboard is used as an aid for text and numerical input.

Connect the PC keyboard as follows:

- 1** Connect the USB plug of the keyboard with the USB connector of the Ti-Touch (Type A).
- 2** Enter and configure the keyboard in the device manager of the Ti-Touch (*see chapter 11.10, page 133*).

4.6.6 Connecting a barcode reader

The barcode reader is used as an aid for text and numerical input. You can connect a barcode reader with USB interface.

Connect the barcode reader as follows:

- 1** Connect the USB plug of the barcode reader with the USB connector of the Ti-Touch (Type A).

- 2 Enter and configure the barcode reader in the device manager of the (see chapter 11.11, page 134).

Settings on the barcode reader:

Program the barcode reader as follows (also see manual for the barcode reader):

- 1 Switch the barcode reader to programming mode.
- 2 Specify the desired layout for the keyboard (USA, Germany, France, Spain, German-speaking Switzerland).

This setting must match the setting in the device manager.
- 3 Make sure that the barcode reader is set in such a way that Ctrl characters (ASCII 00 to 31) can be sent.
- 4 Program the barcode reader in such a way that the ASCII character 02 (STX or Ctrl B) is sent as the first character. This first character is normally referred to as the "Preamble" or "Prefix Code".
- 5 Program the barcode reader in such a way that the ASCII character 04 (EOT or Ctrl D) is sent as the last character. This last character is normally referred to as the "Postamble", "Record Suffix" or "Postfix Code".
- 6 Exit the programming mode.

4.6.7 Connecting a Sample Processor

If you wish to integrate your Ti-Touch in an automation system, then you can connect the following Sample Processors to the USB connector:

- 814 USB Sample Processor
- 815 Robotic USB Sample Processor XL
- 810 Sample Processor (from firmware version 5.916.0040)

To connect a Sample Processor, you need the 6.2151.000 controller cable.

Connect the Sample Processor as follows:

- 1 Connect the Sample Processor to the power grid.
- 2 Connect the Sample Processor to the Ti-Touch with the controller cable.

 **NOTE**

The plug on the controller cable is protected against accidental disconnection by means of a pull-out protection feature. If you wish to pull out the plug, you first need to pull back the outer plug sleeve marked with arrows.

- 3 Configure the Sample Processor in the device manager (see *chapter 11.5, page 108*).

4.7 Connecting sensors

4.7.1 General

The measuring interface includes one high-ohm measuring input (**Ind.**) for a pH, metal or ion-selective electrode, one input (**Ref.**) for a separate reference electrode, one measuring input (**Temp.**) for a temperature sensor (Pt1000 or NTC), one measuring input (**Pol.**) for a polarizable electrode and one connector (**iConnect**) for an iConnect. This contains a measuring interface for electrodes with integrated data chip, referred to as iTrodes.

 **NOTE**

Before the sensor is connected, the electrode holder must be mounted on the support rod with the clamping ring, see (*see chapter 4.3, page 14*).

4.7.2 Connecting a pH, metal or ion-selective electrode

Connect the pH, metal or ion-selective electrode as follows:

- 1** Plug the electrode plug into the **Ind.** socket of the Ti-Touch.




Figure 7 Connecting a pH, metal or ion-selective electrode




NOTE

The electrode cable is protected against accidental disconnection of the cable by means of a pull-out protection. If you wish to pull out the plug again, you first need to pull back the outer plug sleeve.

4.7.3 Connecting a reference electrode


Connect the reference electrode as follows:

- 1** Plug the electrode plug into the **Ref.** socket of the Ti-Touch.



Connect the polarizable electrode as follows:

- 1





NOTE

The electrode cable is protected against accidental disconnection of the cable by means of a pull-out protection. If you wish to pull out the plug again, you first need to pull back the outer plug sleeve.

4.7.5 Connecting the temperature sensor or electrode with integrated temperature sensor

A temperature sensor of the Pt1000 or NTC type can be connected to the **Temp.** connector.

Connect the temperature sensor or the electrode with integrated temperature sensor as follows:

- 1 Plug the temperature sensor plugs into the **Temp.** sockets of the Ti-Touch.




Figure 10 Connecting the temperature sensor or electrode with integrated temperature sensor



NOTE

The red plug must always be plugged into the red socket for the purpose of shielding against disruptions.

If you use an electrode with an integrated NTC sensor, then you must plug the red plug into the red socket.

4.7.6 Connecting the iConnect

The external measuring interface 854 iConnect can be connected to the **iConnect** connector.

Connecting the 6.2168.000 mini USB adapter cable to an instrument

If the mini USB adapter cable is not connected to the instrument yet, proceed as follows:

- 1** Plug the mini USB adapter cable (2) into the **iConnect** connector of the instrument (1). Observe the correct orientation (markings).




Figure 11 Plugging in the mini USB adapter cable

Depending on the instrument's model version, there may be several **iConnect** connectors.

- 2** Leave the adapter cable plugged in to protect the connector inside the instrument (1) against mechanical impact.

Connecting the 854 iConnect to the adapter cable

Make sure that the mini USB adapter cable is connected to the instrument.

The 854 iConnect can also be connected while the instrument is switched on.

- 1 Plug the plug of the 854 iConnect (3) into the socket of the mini USB adapter cable (2). Observe the correct orientation (markings).




Figure 12 Plugging in the 854 iConnect

As soon as the instrument is switched on, the 854 iConnect is detected automatically and entered as measuring input into the device properties.

Connecting the electrode

The 854 iConnect is used as measuring input for iTrodes (electrodes with integrated memory chip).

Connect the electrode as follows:

- 1** Remove the protective cap of the 854 iConnect.




Figure 13 Removing the protective cap

- Align the guide pin (5) of the electrode with the recess in the 854 iConnect (4).




Figure 14 Aligning the guide pin

- 3** Attach the electrode to the 854 iConnect.




Figure 15 Attaching the electrode

The guide pin guarantees correct connection in such a way that the contact pins cannot be damaged.

- 4** Tighten the screw cap by hand.

If there is an electrode in the sensor list of the firmware or software, the electrode is detected automatically when connecting it.



NOTE

Mount the protective cap again as soon as the 854 iConnect is not in use anymore and no electrode is connected.

Removing the 854 iConnect

The 854 iConnect can also be removed while the instrument is switched on.

- 1 Unplug the 854 iConnect (3) from the socket of the mini USB adapter cable (2).




Figure 16 Unplugging the 854 iConnect

- 2** Leave the mini USB adapter cable (2) plugged into the socket of the instrument (1).



NOTE

Position the mini USB adapter cable in such a way that it cannot be removed by accident.



NOTE

For further information on the 854 iConnect, see *Manual iConnect 854*, article number 8.854.8002.

4.7.7 Differential potentiometry

Potentiometric measurements with high-ohm measuring chains can be disrupted by electrostatic and electromagnetic fields in media with low conductivity. Use our 6.0229.100 Solvotrode or other special electrodes for pH measurements in organic solvents. If no reliable measurements are possible with these, then a 6.5104.030 (230 V) or 6.5104.040 (115 V) differential amplifier can be used. The differential amplifier is connected to the high-ohm measuring input (**Ind.**).

4.8 Connecting the Ti-Touch to a network

The 916 Ti-Touch has a network connection (Ethernet). This can be used to integrate your Ti-Touch in your network. You can, for example, store data on a PC within the network or print reports on a network printer. In *Chapter 11, page 96*, you will find information as to which settings are necessary for the network connection.

5 Titrations

5.1 Dynamic equivalence point titration (DET)

Dynamic equivalence point titration is a titration mode for all standard titrations with an s-shaped curve progression. The reagent is added in variable volume steps. The volume steps vary as a function of the slope of the curve. An attempt is made to reach constant measured value changes with each dosing. The optimal volume for dosing is determined from the measured value changes of the previous dosings. Measured value acceptance is measured value drift-controlled (equilibrium titration) or after a waiting time. Equivalence points are evaluated automatically.




Figure 18 Reagent dosing for DET

5.2 Monotonic equivalence point titration (MET)

Monotonic equivalence point titration is a robust titration mode for titrations with any curve shape and for slow titrations or slow-response electrodes. The reagent is added in constant volume steps. Measured value acceptance is measured value drift-controlled (equilibrium titration) or after a waiting time. Equivalence points are evaluated automatically.




Figure 19 Reagent dosing for MET

5.3 Endpoint titration (SET)

Endpoint titration is a titration mode for rapid routine determinations to a preset endpoint (e.g. titrations in accordance with special standards) and titrations for which reagent overflow must be avoided. The titration termination at the endpoint takes place either volume drift-controlled or after a waiting time. The volume dosed until the endpoint can be used for further calculations (e.g. the content of the sample).




Figure 20 Reagent dosing for SET

5.4 Titration with constant maintenance of the measured value (STAT)

In STAT titrations, reagent is added until a preset measured value (control point) has been reached, which is then maintained at a constant level. Because the substance that is released as a result of the reaction is immediately titrated off with the titrant, the STAT controller has to be able to maintain a preset control point at a constant level until the preset stop condition has been reached, i.e. without oscillating around the control point.

STAT titrations are used, for example, in enzyme analysis where the enzyme activity can be determined from the resulting dosing rate (in mL/min).

Reagent dosing takes place in two phases during the titration:

- **Dosing outside the control range**

The dosing rate increases continuously during this phase. It starts with the Min. rate and increases to the Max. rate. The dosing continues to be carried out at the Max. rate until the Control range is reached. In this range, STAT mode and SET mode control equally.

- **Dosing inside the control range**

The STAT controller compensates minimal measured value deviations from the control point quickly by dosing in the control range faster than for the SET controller.




Figure 21 Reagent dosing for STAT

5.5 Manual titration (MAT)

Manual titration is a titration mode that is suitable in particular for titrations that are carried out conventionally manually with a Schellbach buret or a Dosimat.

The reagent is added by manual dosing. The endpoint is detected by a color change or with an electrode connected over a measuring signal. The command always has to be canceled manually with the **Cancel command** button or with the **[FILL]** key on the Manual Dosing Controller.

An advantage of the MAT titration mode is that the manual reagent addition as well as the measuring signal can be recorded optionally. These data points (time, volume, measured value) can be evaluated after the titration, e.g. to assess the titration quality, especially the accuracy of reaching the transition point.

The system is FDA compliant and enables a manual titration according to the FDA guidelines with this mode.

6 Operation

6.1 Switching the instrument on and off

Switching on the instrument



CAUTION

Peripheral devices (e.g., printers) must be connected and switched on before you switch on the 916 Ti-Touch.



NOTE

English is set as the default dialog language when the instrument is switched on for the first time.


Proceed as follows:

- 1 ■ Press the power switch on the left-hand side of the back panel of the 916 Ti-Touch.
The 916 Ti-Touch is initialized. A system test is performed. This process takes some time.



NOTE

If a buret unit is connected, then a request appears to carry out the **Prepare** function:



All tubings and the cylinder are rinsed with the **Prepare** function. The preparing of the buret unit is described in *chapter 28.3.3, page 277*.

- Confirm the message with **[OK]**.

The main dialog is displayed:

Switching off the instrument



CAUTION

The 916 Ti-Touch must be switched off by pressing the power switch on the rear of the instrument before the electricity supply is interrupted. If this is not done, then there is a danger of data loss.

Proceed as follows:

- 1 Press the power switch on the left-hand side of the back panel of the 916 Ti-Touch.

The current data is saved and the system is shut down. This process takes just a short time. At the same time, all other instruments connected to the 916 Ti-Touch via a USB cable are also being switched off.

6.2 Fundamentals of operation

6.2.1 Touch-sensitive screen

The entire 916 Ti-Touch user interface is touch-sensitive. Simply touch a few of the buttons on the interface to learn how a touch-sensitive screen reacts. You can always return to the main dialog by touching .

In order to enable an element on the 916 Ti-Touch user interface, just touch the screen with your fingertip, the eraser of a pencil or a stylus (special tool for operating instruments with touch-sensitive screens).




CAUTION

Never touch the touch screen with a pointed or sharp object such as a ballpoint pen.

In the default setting, the software is configured in such a way that an acoustic signal will be generated every time an active control is touched. This setting can be deactivated in the system settings (*see chapter 7.5, page 60*).

6.2.2 Display elements and controls



The following display elements and controls are available:

Table 1 Fixed keys which are always available



	[Home] always opens the main dialog.
	[Back] saves the entry and opens the next-higher dialog page.
	[Help] opens the online help for the dialog displayed.
	[Print] opens the printing dialog.
	[Manual] opens the manual control.
	[Stop] cancels the running determination.
	[Start] starts a determination.

The file name of the loaded method, the time and the system status are displayed in the main dialog in the **Title bar**.

In the other dialogs, the title bar shows the headings of the next upper level and of the displayed dialog. This is an aid for orientation during navigation through the user dialog.






Table 2 Screen elements

Buttons open a new dialog when they are tapped.



The current status of the system is displayed in the upper right-hand corner of the title bar.

Table 3 Status displays



Number editor

Edit command / Stop conditions

Stop time

off^s

Input:
1 ... 999999

Default value:
off

7

8

9

off

4

5

6

R1 ▾

1

2

3

0

+/-

.

Cancel

Delete entry

OK

Table 5 Editing functions

Editing function	Description
[OK]	The modification is applied and the editing dialog is exited.
[Cancel]	The editing dialog is exited without applying the modification.
[Delete entry]	The content of the input field is deleted completely.
[off]	If not only numbers but also special values (e.g., off) can be entered, then the corresponding buttons will be shown to the right of the numerical keypad.
[R1]	For many parameters, a result previously defined in the method can also be entered in place of a number (<i>see chapter 32.7, page 450</i>). You can select the result variable by touching [R1] .



NOTE

A commercially available USB keyboard can be connected to make text and numerical input easier. The key assignment is described in *chapter 11.10, page 133*.

- Tap on the **Dialog language** list box and select the desired language.

Tap on the fixed keys [] or [].

The main dialog is displayed in the respective dialog language.

7.1.2 Setting the date, time and local time

Proceed as follows to set the date and time:

- In the main dialog, tap on **[System]**.
- Tap on **[System settings]**.

The **System / System settings** dialog is displayed.

- Tap on the input field for the date.
The editor opens.
- Enter the current date in the format **YYYY-MM-DD** and confirm with **[OK]**.
The arrow keys **[←]** and **[→]** are used to move the cursor to the left or to the right by one character.

The input is saved and the editor is closed.

- Tap on the input field for the time.
The editor opens.
- Enter the current time in the format **hh:mm:ss** (24-hour format) and confirm with **[OK]**.
The arrow keys **[←]** and **[→]** are used to move the cursor to the left or to the right by one character.

The input is saved and the editor is closed.

4 Enter the local time

- Tap on the list box **Local time - UTC** and select the difference from the UTC (Coordinated Universal Time).
The selection **off** means that the time is saved with no difference from the UTC.


5 Saving the settings

Tap on the fixed keys [

The time settings are saved.

7.2 System-specific dialog options

Main dialog: **System** ► **System settings** ► **Dialog options**



If you work without a login function, in this dialog you can define whether the system should generally be operated in expert mode or in routine mode. If you work with the login function activated, you must define this setting separately for each user (*see chapter 7.3.1, page 49*).

Dialog

Dialog mode in which the user may operate the system.

Selection	Expert dialog Routine dialog
Expert dialog	All functions that are supported by the system are available.
Routine dialog	The user dialog can be limited for routine operations. Only selected functions are available (<i>see Configuring the routine dialog, page 46</i>).

[Fixed keys]

Block unneeded fixed keys (see *"Blocking unneeded fixed keys"*, page 45).

[Routine dialog]

Configure functions for the routine dialog (see *"Configuring the routine dialog"*, page 46).

Selecting the dialog mode

Proceed as follows to change the dialog mode:

1 Select the dialog mode

Open the selection list **Dialog** and select either **Expert dialog** or **Routine dialog**.

2 Saving the settings

Tap on the fixed keys [↔] or [🏠].

The setting will apply to all dialogs.



NOTE

If you have selected **Routine dialog** and if the routine dialog was configured in such a way that the dialog **System settings / Dialog options** is blocked, then you can switch back over to the expert dialog as follows:

- Operation without login function:
In the main dialog, enter **User = Metrohm**.
- Operation with login function:
A user who works with expert dialog must log in.

Blocking unneeded fixed keys

This following configurations apply for **both** dialog modes: routine dialog and expert dialog.

Blocking fixed keys

Proceed as follows to block unneeded fixed keys:

1 Display fixed keys which can be blocked

Tap on the **[Fixed keys]** button.

All fixed keys which can be blocked are displayed.

2 Deactivate buttons

Deactivate those buttons which are not permitted to be used.

All deactivated buttons will appear grayed-out, i.e. they are inactive.

3 Deactivate other functions

Many additional buttons and parameters can be disabled in the dialogs of **[System]**, **[Load method]**, **[Control]**, **[Edit parameters]** and **[Results]**. The corresponding option must be activated in order to enable these buttons.

4 Saving the settings

Tap on the fixed keys [] or [].

All deactivated functions will appear grayed-out, i.e. they are inactive.

7.3 User administration

Main dialog: **System** ► **System settings** ► **User admin.**

All of the functions of the user administration are described in this chapter. No user administration is mandatory for the simple operation of a titration system. If, however, you wish to be in compliance with *FDA Guidance 21 CFR Part 11*, then you must use the functions of the user administration. Additional information regarding the requirements of the FDA Guidance can be found in the document *Compliance Guide Touch Control*.



NOTE

If you work with the login function activated, then the user administration is accessible only for users with administrator rights. This means that you must ensure that at least two users have administrator rights so at least one of them will be available. Keep the access rights for a user with administrator rights in a safe place so that they are accessible in an emergency.



NOTE

Once users have worked with the login function and password protection, they can no longer be deleted, even if the password protection is disabled again. The status of these users must be set to **inactive** (requirement of FDA Guidance 21 CFR Part 11).

The last user with administrator rights cannot be deleted.

[Edit]

Editing the data of the selected user (*see chapter 7.3.1, page 49*).

7.3.1 Editing the user configuration

User list: **User ► New / Edit**

User

The designation of the user is used for unambiguous identification, e.g. the company internal shorthand symbol or the personal number. The user name is printed out in all reports containing determination data and stored in the determination file. Each file always contains the name of the user who created it and the name of the last user to edit it.

Entry	24 characters maximum
Default value	empty

Full name

Complete name of the user.

Entry	24 characters maximum
Default value	empty

Dialog

Dialog mode in which the user may operate the system.

Selection **Expert dialog** | **Routine dialog**

Expert dialog

All functions that are supported by the system are available.

Routine dialog

The user dialog can be limited for routine operations. Only selected functions are available (*see Configuring the routine dialog, page 46*).

Status

Status of the user. Users can be deactivated. This function is useful, for instance, if the user is no longer authorized to operate the system or no longer works for the company.

Selection	active inactive
Default value	active

active

The user is authorized to operate the system.

inactive

The user is not authorized to operate the system and can no longer log in.

Admin. rights

on | off (Default value: **off**)

If this option is activated, then the user has administrator rights.



NOTE

The last user with administrator rights cannot be deleted anymore.

[Favorites]

Saving methods and sample tables as user-specific favorites (see *chapter 19, page 206*).

[Signature method]

Assign the permission to use and sign methods (see *"Edit user / Signature method" dialog*, page 50).

[Signature determ.]

Assign the permission to use and sign determinations (see *"Edit user / Signature determination" dialog*, page 51).

"Edit user / Signature method" dialog

User: **Edit** ► **Signature method**

Use only released methods

If this option is activated, then the user may only start methods that have been released.

on | off (Default value: **off**)

Release methods (signature level 2)

on | off (Default value: **off**)

If this option is activated, then the user may only sign methods which have the status **reviewed**. The method is assigned the status **released**.



If the two options **Review methods** and **Release methods** are selected for a user, then this user can sign different methods at either level 1 or level 2, but not the same method at both level 1 and level 2.

on | off (Default value: **off**)

If this option is activated, then the user may delete the signatures of a method which has been released. The method is assigned the status **saved**. The signatures of a method can only be deleted if it has the status **released**.

"Edit user / Signature determination" dialog

User: **Edit ► Signature determ.**

Authorizations for signing determinations can be defined in this dialog. These settings are only effective if you are working with the login function enabled and password protection.

on | off (Default value: **off**)

If this option is activated, then the user may only sign determinations at the first level. The determination is assigned the status **reviewed**.

Release determinations (signature level 2)

on | off (Default value: **off**)

If this option is activated, then the user may only sign determinations at the second level. The determination is assigned the status **released**.

Delete signatures

on | off (Default value: **off**)

If this option is activated, then the user may delete the signatures of a determination which has been released. The signatures of a determination can only be deleted if it has the status **released**.

7.3.2 Creating an identification profile

User list: **User** ► **Create ID profile**

If you plan to carry out the login with an identification profile (*see chapter 7.3.3, page 53*), then you must first create an identification profile on a storage medium for each user. A check can then be made at the time of login as to whether or not the user does exist and whether or not he or she is working in the expert dialog or in the routine dialog. After a successful login the routine dialog settings stored on the card are loaded.



NOTE

In addition to the user name, the current routine dialog settings are also stored in this identification profile. This means you can define user-specific routine dialog settings for each user. However, you must configure them (see *"Configuring the routine dialog"*, page 46), **before** you create the identification profile.

Before you create the identification profile, check whether the desired routine dialog settings are active.

1 Plug in a storage medium

Plug in the USB storage medium.

2 Select a user

In the user list, select the user for whom the profile is to be created.

3 Create the identification profile

Tap on **[Create ID profile]**.

The user configuration and the current routine dialog settings are saved.

7.3.3 Defining login options

Main dialog: **System ► System settings ► User admin. ► Login options**

There are a number of different ways to log onto the system:

- Without login
- Login via user name
- Login via user name and password
- Login via identification profile
- Login via identification profile and password



NOTE

If you work with the login function activated, then the user administration is accessible only for users with administrator rights. This means that you must ensure that at least two users have administrator rights so at least one of them will be available. Keep the access rights for a user with administrator rights in a safe place so that they are accessible in an emergency.



NOTE

If you exit this dialog with [↩] or [🏠], and if you have selected one of the login variants **Login via user name** or **Login via identification profile**, then the login dialog will open automatically and you must also log in to the system.

This means that you must make sure that you have first defined all the users and created the identification profiles before you activate the login function.

User administration / Login options

☒ Login via user name

☐ Login via identification profile

☒ Password required

☐ Logout automatically

after min

☐ Login only for the same user

Audit Trail Reasons Modific. options Password options

Login via user name

on | off (Default value: **off**)

If this option is activated, then the user must log in with his or her unambiguous identification.

Login via identification profile

on | off (Default value: **off**)

If this option is activated, then the login will take place via USB storage medium with the identification profile stored on it.

Password required

on | off (Default value: **off**)

If this option is activated, then the user must enter a password in addition to his or her user name or identification profile.

Logout automatically

on | off (Default value: **off**)

If this option is activated, then the user will be logged out automatically after the specified time.

Input range	1 to 60 min
-------------	--------------------

Login only for the same user

on | off (Default value: **off**)

If this option is activated, then only the same user may log in again after he or she has logged out. Users with administrator rights can, however, log in at any time.

[Audit Trail]

Define the settings for recording an Audit Trail (see chapter 7.3.7, page 58).

[Reasons]

Predefining a list of reasons which can be selected for the modification/
signature of a method or determination (*see chapter 7.3.6, page 57*).

[Modific. options]

Defining the modifications for which a reason is required (*see chapter 7.3.5, page 57*).

[Password options]

Define the settings for the password, see following chapter.

7.3.4 Password options

Main dialog: **System ► System settings ► User admin. ► Login options ► Password options**

You can make various settings for password entry in the password options.

Login options / Password options

Minimum password length

1

No. of entry attempts

off

☐

Special characters required

☐

Password expires

every

365

days

Minimum password length

Minimum number of characters of the passwords.

Input range	1 to 10
Default value	1

No. of entry attempts

If the user has logged in incorrectly this many times, then it will automatically be deactivated. It can only be reactivated by a user with administrator rights.

Input range	2 to 5
Selection	off
Default value	off

Special characters required

on | off (Default value: **off**)

If this option is activated, then the password must contain one of the following special characters: ° § + ! @ * # \$ % & ¬ () = ' ^ ` ~ [] { } - _ . : ; , > < £ !

Password expires

on | off (Default value: **off**)

If this option is activated, then the user must define a new password after the time specified. A password that has already been used cannot be used again.

Input range	1 to 999 days
Default value	365 days

Forgotten password



CAUTION

If a user has forgotten his password, a new user name must be defined. The same user name can only be used again after a re-installation and re-creating the user list.

Proceed as follows:

1 Create a backup

Create a backup (see chapter 12.3, page 143).

2 Executing a system initialization

Switch off the Ti-Touch and execute a system initialization (see *chapter 30.1, page 436*).

3 Restore backed-up data

Use the function **Restore** to reload the data from the backup into your system (see chapter 12.3.1, page 143).

Deactivate the options **User list** and **System settings / User admin..**

4 Switch on the Ti-Touch

Switch the Ti-Touch back on again after a few seconds.

5 Restore the user list

Create the user list again and redefine the login options.

7.3.5 Modification options

Main dialog: **System ► System settings ► User admin. ► Login options ► Modific. options**

In the dialog **Login options / Modification options**, you can define for which actions a reason must be entered. These reasons are documented in the Audit Trail (*see chapter 7.3.7, page 58*) together with the modification. The reason for the last modification is shown in the properties of the method or determination.



NOTE

The reasons are only requested when working with activated login function and password.

Saving modified method

on | off (Default value: **off**)

If this option is activated, then a reason must be given when saving a method modification.

Recalculating determination

on | off (Default value: **off**)

If this option is activated, then a reason must be given when determinations are being recalculated.

7.3.6 Reasons

Main dialog: **System ► System settings ► User admin. ► Login options ► Reasons**

In the dialog **Login options / Reasons**, you can create a selection list containing reasons from which a selection can be made when signing and modifying methods and determinations. Some reasons have already been provided.



[New]

Add a new reason to the list.

[Delete]

Delete the selected reason from the list.

[Edit]

Change the designation of the selected reason.

7.3.7 Audit Trail

Main dialog: **System** ► **System settings** ► **User admin.** ► **Login options** ► **Audit Trail**

An Audit Trail is an automatically generated log of all user activities. An Audit Trail contains precise logs of user actions (date, time, user, action, etc.). Recording an Audit Trail is important for compliance with *FDA Guidance 21 CFR Part 11* when using PC programs. A step-by-step description of how you must proceed to remain in compliance with the FDA Guidance is included in the document *Compliance Guide Touch Control*. The Audit Trail is saved to the internal memory.

You can also use the Audit Trail function specifically to record the data which are of interest to you.

You can view, filter and export the Audit Trail on a PC with the supplied software program *AuditTrailViewer*. Details concerning utilization of the *AuditTrailViewer* can be found in the Appendix (*chapter 32.10, page 469*).

If necessary, you can delete all of the entries in the Audit Trail with **[Delete Audit Trail]**. Do not fail to create a backup before doing so, however.

Security log

on | off (Default value: **off**)

If this option is activated, then the following events will be recorded: user login/logout, password changes, automatic deactivation of users and messages displayed during the login procedure.

User administration log

on | off (Default value: **off**)

If this option is activated, then all changes with respect to user administration will be recorded (modifying user list/user data, changing login options, etc.).

Method log

on | off (Default value: **off**)

If this option is activated, then all changes with respect to methods and determinations will be recorded (saving, deleting, renaming, copying and loading methods; editing method parameters; starting, stopping, pausing determinations).

Data log

on | off (Default value: **off**)

If this option is activated, then the following events will be recorded: changing sample data; settings in the **Control** dialog, recalculating determinations, changing statistics data, changing the titer/concentration of a titrant, changing the calibration data of a sensor, changing the value of a common variable.

System log

on | off (Default value: **off**)

If this option is activated, then the program start and the displayed messages of the following types will be recorded:

- : General warning messages
- : Error messages

7.4 Measured value display

Main dialog: **System ► System settings ► Meas. value display**

The number of decimal places can be defined for pH values and voltages.

This setting refers exclusively to the display of the measured values in the live display and in the manual control. The values will, however, always be stored with their full accuracy.

7.5 Acoustic signals

Main dialog: **System** ► **System settings** ► **Acoustic signals**

You can define acoustic signals in order to direct attention to particular events. You can define signals for the following events:

- **Wrong manipulation**
An acoustic signal will sound each time an invalid action is attempted (e.g. pressing [?] again when Help is open).
- **Display a message**
A short beep will sound each time a message appears on the display. This informs the user that the message must be confirmed.
- **Button contact**
Each time a button is touched on the touch screen, this will be confirmed by an acoustic signal.
- **External data input**
An acoustic signal will provide confirmation each time data is received from external devices (e.g. balance, barcode reader).

8 Titrants

Main dialog: **System ► Titrants**

This chapter describes how you can create a list of titrants used in the system. Titrants can be used in intelligent buret units or in nonintelligent buret units. Intelligent buret units have a built-in data chip on which the data for the titrant is stored. This data is automatically read out during attachment and entered in the titrant list.

System / Titrants				
Titrant	Cyl.	Type	Dos.device	
c(AgNO ₃) = 0.1 mol/L	10 mL	IDU	D1/Ti-Touch	
c(HCl) = 0.1 mol/L	10 mL	IDU		
c(NaOH) = 0.1 mol/L	5 mL	IDU		

The titrant list can contain a maximum of 30 titrants. The following data is specified for each titrant:

- Designation
- Cylinder volume
- Type
 - **IDU**: dosing unit with integrated data chip
 - **IEU**: exchange unit with integrated data chip
 - **EU**: exchange unit without data chip
- MSB connector of the dosing device/control instrument (only when exchange unit/dosing unit is attached)

Titrants in exchange units/dosing units with integrated data chips are depicted in green lettering.

The following titrant data is stored in the titrant list:

- Name
Each titrant in the system is identified by its unambiguous name.
- Concentration
- Current titer
- Working life



- Data on exchange/dosing unit:
 - Parameters for the function **PREP**
 - Length and diameter of the tubings
 - Port assignment of the dosing unit
 - Cylinder volume
 - Serial number
 - etc.
- etc.



NOTE

If data is read out from the data chip, then a check is made whether the titrant list already contains a titrant of the same type with the identical serial number. If this is the case, then the older data set will **always** be overwritten by the new data set, no matter whether the data set in the titrant list or the data set on the data chip is the most recent one.

[New]

Adding a new titrant to the list(see chapter 8.1, page 62).

[Delete]

Delete the selected titrant from the list.

[Edit]

Editing the data of the selected titrant (*see chapter 8.2, page 63*).

8.1 Adding a new titrant

Before you can use a titrant, you must add it to the titrant list. To do this, use the button **[New]**.

- Exchange unit/dosing unit with data chip:
The exchange unit or dosing unit must be attached. All dosing devices on which non-configured exchange or dosing units have been detected are included in a selection list. Tapping on the button **[Edit]** opens the properties dialog, see following chapter.
- Exchange unit without data chip:
The properties dialog is opened after the dosing unit type has been selected, see following chapter.

8.2 Editing titrant data

Titrant list: **Titrants** ► **New / Edit**

All of the data for the selected titrant is displayed in the dialog **Titrants / Edit**.

Titrant

The designation of the titrant is used for unambiguous identification.

Entry	24 characters maximum
Selection	Selection of frequently used titrants

Concentration

Concentration of the titrant.

Input range	–999999999 to 9999999999
Default value	1.000

Unit of the concentration.

Entry	10 characters maximum
Selection	μmol/mL mmol/L mol/L g/L mg/L mg/mL μg/L ppm % mEq/L
Default value	mol/L

Comment

Entry	24 characters maximum
-------	------------------------------

Titer

Titer of the titrant.

Input range	–999999999 to 9999999999
Default value	1.000

8.3 Monitoring the working life

Titrant: **Edit ► Working life**

In the dialog **Edit titrant / Working life**, you can define the time interval after which the titrant must be replaced. This is particularly important if your titrant has a limited working life. If you do not wish to monitor the working life, then you can enter only the date of manufacture for documentation purposes.

Preparation date

Date on which the reagent was manufactured or the bottle was opened. For new titrants, the time the preparation was made will be specified.

Format: YYYY:MM:DD

Monitoring

on | off (Default value: **off**)

If this parameter is activated, then the working life will be monitored.

Working life

If you define a time interval for the working life, then the **Expiry date** will be tracked automatically.

Input range	1 to 999 days
Default value	999 days

Expiry date

If you define an expiry date, then the **Working life** will be tracked automatically.

Format: YYYY:MM:DD

Action

Selection of the action which is carried out when the time interval has expired.

Selection	Display message Document message Cancel determination
Default value	Display message For all three options it is documented in the determination data (see dialog More determination data / Messages), that the time interval has been expired.

Display message

A message is displayed. You can select whether you want to continue with the determination or cancel the run.

Document message

In the determination data it will be documented, that the time interval has been expired.

Cancel determination

The determination is stopped.

8.4 Dosing unit

Titrant: **Edit** ► **Dosing unit**

Edit titrant / Dosing unit 

Titrant: c(12) = 0.05 mol/L

Name

Order number 6.3032.250

Serial number 10170552

Cylinder volume 50 mL

Serial no. cyl.

You can edit data for the dosing unit of the titrant in this dialog.

Name

Designation of the exchange or dosing unit.

Entry	24 characters maximum
Default value	empty

Order number

Order number of the exchange or dosing unit. It will be read out automatically on units with integrated data chips.

Entry	24 characters maximum
Default value	empty

Serial number

Serial number of the exchange or dosing unit. It will be read out automatically on units with integrated data chips.

Entry	8 digits maximum
-------	------------------

Cylinder volume

Cylinder volume of the dosing unit. It will be read out automatically on dosing units with integrated data chips.

Selection	2 5 10 20 50
Default value	20

Serial no. cyl.

Serial number of the dosing cylinder. It will be read out automatically on exchange or dosing units with integrated data chips. The number can be changed manually at any time, e.g. when a cylinder is replaced.

Entry	8 digits maximum
-------	-------------------------

[Valve disk]

Specifying the shift direction of the valve disk (*see chapter 8.4.3, page 72*).

[GLP test]

Defining the time interval for the GLP test (*see chapter 8.6, page 76*).

[PREP param.]

Entering the parameters for the preparation (*see chapter 8.4.1, page 67*).

[Tubing param.]

Entering the parameters for the connected tubing (*see chapter 8.4.2, page 69*).

8.4.1 Parameters for preparing (PREP) and emptying (EMPTY)

Titrant: **Edit ► Dosing unit ► PREP param.**

In the dialog **Dosing unit / PREP parameters**, you can adjust the parameters for the execution of the **Prepare** (command PREP) and **Empty** (command EMPTY) functions. The **Prepare** function is used to rinse the cylinder and tubing of the dosing unit and fill it air bubble-free. You should carry out this function before the first determination or once a day. The **EMPTY** function empties the cylinder and the tubings of the dosing unit.

Dosing port PREP/EMPTY

Dosing port through which the cylinder contents are ejected.

Selection	Dosing port 1 Dosing port 2 Fill port Special port
Default value	Dosing port 1

Rate used for the aspiration and ejection of the reagent via dosing port 1.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

The maximum dosing rate depends on the cylinder volume (see chapter 32.2, page 444).

Rate used for the aspiration and ejection of the reagent via dosing port 2.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

The maximum dosing rate depends on the cylinder volume (*see chapter 32.2, page 444*).

Rate used for the aspiration and ejection of the reagent via fill port.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

The maximum dosing rate depends on the cylinder volume (see chapter 32.2, page 444).

Rate used for the aspiration and ejection of the reagent via the special port.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

The maximum dosing rate depends on the cylinder volume (see *chapter 32.2, page 444*).

8.4.2 Tubing parameters

Titrant: **Edit ► Dosing unit ► Tubing param.**

You can enter the length and diameter of the connected tubings in the dialog **Dosing unit / Tubing parameters**. The values which have already been entered correspond to the dimensions of the supplied standard tubings. In addition, the port assignment can be modified.



NOTE

These parameters are important for the correct execution of the **Preparing** (PREP command) and **Emptying** (EMPTY command) functions, because the volumes of the tubing connections are taken into account.

Dosing port 1

Port

Port to be used as dosing port 1 for the **PREP** and **EMPTY** (see figure 22, page 71) functions.

Selection	Port 1 Port 2 Port 3 Port 4
Default value	Port 1

Length

Length of the tubing.

Input range	0.0 to 999.9 cm
Default value	40.0 cm
The setting 0.0 means that this tubing will neither be rinsed nor emptied.	

Diameter

Diameter of the tubing.

Input range	0.0 to 9.9 mm
Default value	2.0 mm

Dosing port 2

Port

Port to be used as dosing port 2 for the **PREP** and **EMPTY** functions (see figure 22, page 71).

Selection	Port 1 Port 2 Port 3 Port 4
Default value	Port 3

Input range	0.0 to 999.9 cm
Default value	0.0 cm
The setting 0.0 means that this tubing will neither be rinsed nor emptied.	

Diameter

Diameter of the tubing.

Input range	0.0 to 9.9 mm
Default value	2.0 mm




Figure 22 Dosing unit – port assignment

1 VENT This port is set up for the deaeration of the reagent bottle. An adsorber tube (filled with desiccant) is usually mounted here.	2 Port 1 The default definition of this port is dosing port 1.
3 Port 2 The default definition of this port is the filling port. A riser tube is usually mounted to it.	4 Port 4 Air is suctioned through this port during the Emptying function.
5 Port 3 The default definition of this port is dosing port 2.	

Monitoring

on | off (Default value: **off**)

If this parameter is activated, then the time interval after which a GLP test has to be carried out again will be monitored.

GLP test interval

If you define a time interval for the GLP test, then the date in **Next GLP test** will be tracked automatically.

Input range	1 to 999 days
Default value	999 days

Next GLP test

If you define a date for the next GLP test, then the **GLP test interval** will be tracked automatically.

Format: YYYY:MM:DD

Action

Selection of the action which is carried out when the time interval has expired.

Selection	Display message Document message Cancel determination
Default value	Display message For all three options it is documented in the determination data (see dialog More determination data / Messages), that the time interval has been expired.

Display message

A message is displayed. You can select whether you want to continue with the determination or cancel the run.

Document message

In the determination data it will be documented, that the time interval has been expired.

Cancel determination

The determination is stopped.

Titer validity

If you define a time interval for the validity of the titer, then the date in **Next titer determ.** will be tracked automatically.

Input range	1 to 999 days
Default value	999 days

Next titer determ.

If you define a date for the next titer determination, then the time interval for the **Titer validity** will be tracked automatically.

Format: YYYY:MM:DD

Action

Selection of the action which is carried out when the time interval has expired.

Selection	Display message Document message Cancel determination
Default value	Display message For all three options it is documented in the determination data (see dialog More determination data / Messages), that the time interval has been expired.

Display message

A message is displayed. You can select whether you want to continue with the determination or cancel the run.

Document message

In the determination data it will be documented, that the time interval has been expired.

Cancel determination

The determination is stopped.

8.7.2 Properties of the previous titer determinations**Dialog "Titer options / History"**

Titration: **Edit ► Titer options ► History**

The date, time and titer of the last ten titer determinations are displayed in tabular form in the dialog **Titer options / History**. Titrations that were determined automatically will be displayed in green; manually entered titer values will be displayed in black with the designation **(m)**. You can delete these entries, e.g. if you have opened a new bottle.



NOTE

If multiple determinations have been carried out in order to determine the titer, then only one entry will be made in the history.

[Delete History]


Delete the entire history.

[Graph]

Open the diagram of the titer values, see following chapter.

Dialog "History / Graph"

Titrant: **Edit** ► **Titer options** ► **History** ► **Graph**



In this diagram, titer values are plotted against the date of the titer determination. You can define warning limits (blue dashed lines) and intervention limits (red dashed lines). These limits will not, however, be monitored.

[Limits]

Define warning and intervention limits.

9 Reagents

Main dialog: **System ► Reagents**

This chapter describes how you can create a list of the reagents used in the system.

System / Reagents	
Reagent	Type
Reagent A	Volumetric
Reagent B	Volumetric
Reagent C	Volumetric

New

Delete

Edit

The list of reagents gives the designation and type of each reagent.

[New]

Add a new reagent to the list, see following chapter.

[Delete]

Delete the selected reagent from the list.

[Edit]

Edit the data of the selected reagent, see following chapter.

9.1 Editing reagent data

Main dialog: **System ▶ Reagents ▶ New / Edit**



Reagents / Edit

Reagent

Reagent A

Comment

Reagent monitoring

Reagent

The designation of the reagent is used for unambiguous identification.

Entry	24 characters maximum
-------	-----------------------

Comment

Entry	24 characters maximum
-------	-----------------------

[Reagent monitoring]

Set the parameters for the reagent monitoring, see following chapter.

9.2 Reagent monitoring

The conditions for the monitoring of the reagent are defined in the dialog **Edit reagent / Reagent monitoring**.

Edit reagent / Reagent monitoring

Number of determ.

off

days

Working life

off

days

Volume

off

mL

Status

Reagent exchange

If one of the following values is reached, then the reagent must be replaced. The values are checked in the following cases:

- at the start of the determination.

- at the end of the determination.

Number of determ.

The number of determinations to be carried out with a certain amount of reagent depends on the type of sample and its amount.

Input range	1 to 999
Selection	off
Default value	off

Working life

Working life of the reagent.

Input range	1 to 999 days
Selection	off
Default value	off

Volume

Volume of titrant dosed.

Input range	1.0 to 999.9 mL
Selection	off
Default value	off

[Status]

Display the status overview of the current values of the reagent monitoring.

[Reagent replacement]

Edit the parameters for the reagent exchange.

"Reagent monitoring / Status" dialog

The current reagent monitoring values are displayed in this dialog.

[Reset]

Reset the values to zero.

Dialog "Reagent monitoring / Reagent replacement"

The parameters for the reagent replacement are defined in this dialog.

Reagent replacement

The reagent can either be exchanged manually or automatically.

Selection	manual auto
Default value	manual

manual

If a monitored parameter has reached the limit set, a message is being displayed. Then the reagent has to be exchanged manually.

10 Sensors

Main dialog: **System ► Sensors**

This chapter describes how you can create a list of the sensors used in the system.

System / Sensors		
Sensor	Sensor type	Meas. input
Ecotrode Plus	pH	
iUnitrode with Pt1000	pH IS	I1/Ti-Touch
pH electrode	pH	
Ag Titrode	Metal	
Metal electrode	Metal	
Temperature sensor	Temp.	
New Delete Edit		

Three standard sensors are defined in the sensor list: **pH electrode**, **Metal electrode** and **Temperature sensor**. These sensors cannot be deleted or renamed. A maximum of 25 additional sensors can be added to these sensors.

For each sensor, the following data is displayed in the sensor list:

- Designation
- Type
 - **pH**: pH electrode
 - **Metal**: Metal electrode
 - **Temp.**: Temperature sensor
 - **Other**: Other sensor
- Measuring input/control instrument (only for intelligent sensors if they are connected)

Intelligent sensors are also indicated by **IS** and are depicted in green lettering.

The following sensor data is stored in the list of sensors:

- Name
Each sensor in the system is identified by its unambiguous name.
- Calibration data (for pH sensors only)
- Calibration interval (for pH sensors only)
- Working life
- etc.



NOTE

If data is read out from the data chip of an intelligent sensor, then a check is made whether the sensor list already contains a sensor of the same serial number. If this is the case, then the older data set will **always** be overwritten by the new data set, no matter whether the data set in the sensor list or the data set on the data chip is the most recent one.

[New]

Adding a new sensor to the list (see chapter 10.1, page 86).

[Delete]

Delete the selected sensor from the list.

[Edit]

Editing the data of the selected sensor (see chapter 10.2, page 87).

10.1 Adding a new sensor

Before you can use a sensor, you must add it to the sensor list. To do this, use the button **[New]**.

- Conventional sensors:
The properties dialog is opened after the sensor type has been selected, see following chapter. The following sensor types can be selected:
 - pH electrode
 - Metal electrode (Pt electrode, Ag Titrode, Ag/AgCl reference electrode, etc.)
 - Other sensor, e.g. Spectrosense
 - Temperature sensor
- Intelligent sensors (also known as iTrodes):
If the 854 iConnect with iTrode is connected to the Ti-Touch, then the sensor is automatically entered in the sensor list and can be configured, see following chapter.

10.2 Editing the sensor data

Sensor list: **Sensor ► New / Edit**

All of the data for the selected sensor is displayed in the dialog **Sensors / Edit**.

Sensor

The designation of the sensor is used for unambiguous identification.

Entry	24 characters maximum
-------	------------------------------

Order number

Order number of the sensor. With intelligent sensors it is read out automatically.

Entry	24 characters maximum
Default value	empty

Serial number

Serial number of the sensor. With intelligent sensors it is read out automatically.

Entry	8 digits maximum
-------	-------------------------

Comment

Entry	24 characters maximum
-------	------------------------------

[Working life]

Defining the working life of the sensor (*see chapter 10.3, page 88*).

[Limit values]

This button is displayed only for pH electrodes.

Action

Selection of the action which is carried out when the time interval has expired.

Selection	Display message Document message Cancel determination
Default value	Display message For all three options it is documented in the determination data (see dialog More determination data / Messages), that the time interval has been expired.

Display message

A message is displayed. You can select whether you want to continue with the determination or cancel the run.

Document message

In the determination data it will be documented, that the time interval has been expired.

Cancel determination

The determination is stopped.

10.4 Calibration data (for pH electrodes only)

Sensor: **Edit ► Calibration data**

Dialog: Edit sensor / Calibration data

Sensor: iUnitrode with Pt1000

Slope: 99.8 % pH(0): 7.033

Electrode test: Excellent electrode

Calibration temp: 24.7 °C (Pt1000)

Calibration date: 2011-08-31 15:09:21

Cal. method: CAL pH

User: SL

Measuring input: 1.854.0010 9318

Buttons: Initial data, GLP test, History

Detailed information concerning calibration is displayed in this dialog **Edit sensor / Calibration data**.

Slope

Slope of the electrode.

10.4.1 Properties of the previous calibrations

Dialog "Calibration data / History"

Sensor: **Edit ► Calibration data ► History**

The date, time and calibration data of the last ten calibrations are displayed in tabular form in the dialog **Calibration data / History**. Calibrations that were carried out automatically will be displayed in green; manually entered calibration data will be displayed in black with the designation **(m)**.

[Delete History]

Delete the entire history.

[Graph slope]


Open the diagram of the electrode slopes, see following chapter.

[Graph pH(0)]

Open the diagram of the electrode zero points, see following chapter.

Dialog "History / Graph"

Sensor: **Edit ► Calibration data ► History ► Graph Slope / pH(0)**.



In this diagram, either slope or pH(0) is plotted against the date of the calibration. You can define warning limits (blue dashed lines) and intervention limits (red dashed lines). These limits will not, however, be monitored.

[Limits]

Define warning and intervention limits.

Dialog "Calibration data / GLP test"

Sensor: **Edit ► Calibration data ► GLP test**

In the dialog **Calibration data / GLP test**, you can define the time interval after which a GLP test must be carried out again for the sensor.

10.5 Limit values for the calibration data

Sensor: **Edit ► Limit values**

Dialog: Edit sensor / Limit values

Sensor: iUnitrode with Pt1000

☒ Monitoring slope

Lower limit: 96.0 %

Upper limit: 101.0 %

☒ Monitoring pH(0)

Lower limit: 6.750

Upper limit: 7.250

You can define the following limit values in the dialog **Edit sensor / Limit values**:

- Slope (pH electrodes)
- Electrode zero point (pH electrodes)

These values are monitored during the calibration. If these limits are infringed, then a message will be displayed and you can decide whether to accept the calibration data or not.

Monitoring slope

on | off (Default value: **off**)

If this parameter is activated, then the slope will be monitored.

Lower limit

pH electrodes:

Input range	–999.9 to 999.9 %
Default value	96.0 %

Upper limit

pH electrodes:

Input range	–999.9 to 999.9 %
Default value	101.0 %

Control Remote Box

Shows to which control device and MSB connector the Remote Box is connected.

The "Control Remote Box" is the interface via which the system can be started and stopped externally. If multiple Remote Boxes are connected, then the one that is recognized first when the program starts will be used as the "Control Remote Box."

Selection	Name of the control device / Number of the MSB connector Not available
-----------	--

11.3.1 E-mail

916 Ti-Touch: **Edit ► E-mail**

The system allows you to send displayed messages as e-mails. The Ti-Touch must be connected to a network for this to function. The following types of messages can be sent:

- : general warning messages
- : error messages

Configuring e-mail dispatch

Proceed as follows so that messages can be sent as e-mails:

1 Activating the option

- In the instrument properties of the 916 Ti-Touch, tap on the **[E-mail]** button.
- In the **Edit device / E-mail** dialog, activate the **Send the following messages as e-mail:** option.

2 Configuring e-mail addresses

- Tap on the **[E-mail settings]** button.
The **E-mail / Settings** dialog is displayed.
- Enter the addresses of the mail server, the sender and the desired recipient.

Parameter description

Send the following messages as e-mail:

on | off (Default value: **off**)

If this parameter is activated, then messages with the following symbols will be sent as e-mails:

- : General warning messages

Memory

Memory location where the PC/LIMS report is stored as a file. The report will be saved in the directory *pc_lims_report*. This directory will be created the first time a PC/LIMS report is generated. All three memory locations are listed as possible selections, even if they cannot be accessed at the moment.

Selection	off External memory 1 External memory 2 Shared memory
Default value	off

off

The report will not be saved as a file.

Shared memory

The report will be saved in a shared directory on the network. The shared directory is selected in the **Edit device / Shared memory** dialog (see chapter 11.3.3, page 102).

RS-232

The RS-232 interface via which the PC/LIMS report is sent. The interface parameters are adjusted in the **Edit device / Port parameters** dialog (see chapter 11.9, page 131).

Selection	off COM 1 COM 2
Default value	off

off

The report will not be sent via an RS-232 interface.

COM 2

This interface is inactive.

Coding

Format in which the PC/LIMS report is coded and stored.

Selection	ISO 8859-1 UTF-8
Default value	ISO 8859-1

ISO 8859-1

This format is recommended for all languages that use the extended ASCII code (e.g. German, English, Spanish, etc.).

UTF-8

This format is required for all languages that do not use the extended ASCII code (e.g. Russian, Chinese, Korean, etc.).

User

Password

[Connect]

11.3.4 TCP/IP settings


If you have connected your Ti-Touch to your network, then you will have to define the network-relevant settings in this dialog. The Ti-Touch requires an IP address so that it is identifiable on the network. The Ti-Touch can acquire the IP address either dynamically from a DHCP server or you can enter the address directly.

on | off (Default value: **on**)

If this parameter is activated, then the Ti-Touch will obtain its IP address directly from a DHCP server. In this case, the remaining parameters can no longer be edited.

IP address for the Ti-Touch. IP addresses are 32-bit numbers and are written as sequences of four decimals, each separated by a period, e.g. "10.157.212.8".

916 Ti-Touch



The dialog shows the control instrument with its measuring inputs, the MSB connectors and connected peripheral devices (dosing device, stirrer, Remote Box).

[Properties]


Open the properties dialog of the highlighted entry.

The following properties can be displayed and edited to a certain extent:

- Properties of the measuring inputs (*see chapter 11.4.1, page 105*)
- Properties of the MSB connectors (*see chapter 11.4.2, page 107*)
- Properties of the peripheral devices on the MSB connector (*see chapter 11.4.3, page 108*)

11.4.1 Properties – Measuring input

Instrument list: **Control instrument ► Edit ► Measuring input 1 ► Properties**



In this dialog, you define the type of temperature sensor you have connected to the selected measuring input.

11.4.2 Properties – MSB connector

Instrument list: **Control instrument** ► **Edit** ► **MSB connector**
1/2 ► **Properties**

Device manager / Edit

Ti-Touch / MSB 1

When would you like to be asked to prepare the dosing device on MSB 1?

☒ Switch on

☒ Attach an exchange or dosing unit

☐ Time interval h

In this dialog, you can define when the request to carry out the **Prepare** function for connected dosing devices is to be displayed. This setting applies for all dosing devices of the selected MSB connector. The manual control contains a description of how to prepare the exchange unit and/or the dosing unit (*see chapter 28.3.3, page 277*).

Switch on

on | off (Default value: **on**)

If this parameter is activated, then you will be requested to prepare the dosing device when the Ti-Touch is switched on.

Attach an exchange or dosing unit

on | off (Default value: **on**)

If this parameter is activated, then you will be requested to prepare the dosing device when the exchange/dosing unit is attached.

Time interval

on | off (Default value: **off**)

Activate this parameter if you wish to receive a regular request to prepare the dosing unit.

Input range	0.1 to 999.9 h
Default value	12.0 h

[Properties]

11.5.1 Properties – Sample Processor

Device name

Entry	24 characters maximum
-------	-----------------------

Entry	24 characters maximum
-------	-----------------------

Program version of the instrument software.

Shows the serial number of the device.

Name of the currently attached sample rack.

Display the internal adjustment data (EEPROM data) of the Sample Processor.



■■■■■■■ 109

off

No check takes place.

Rotation increment

Amount by which the sample rack can be rotated forward or backward relative to its current position. This parameter is required in the command **MOVE** for the settings **Destination = rotate +** and **rotate -**.

Input range	0 to 353.99 °
Default value	5.00 °

Axial distance

Distance between the axis of rotation of the sample rack and the swing axis of the robotic arm.

for 814 USB Sample Processor:

Input range	100.0 to 300.0 mm
Default value	166.0 mm

for 815 Robotic USB Sample Processor XL:

Input range	100.0 to 300.0 mm
Default value	196.0 mm

for 810 Sample Processor:

Input range	100.0 to 300.0 mm
Default value	166.0 mm

11.5.3 Properties – Swing Head

Instrument list: **Sample Processor ► Edit ► Swing Head ► Properties**

Edit device / Properties		
Sample Processor / Tower 1 / Swing Head 786		
Serial number 01108		
Ext. pos.	Angle [°]	Work [mm]
1	60.00	100
2	60.00	0
3	60.00	0
4	60.00	0
Ext. pos. 1-4	Robotic arm	Adjustment data
Edit		

The dialog **Edit device / Properties** shows a list of all external positions with the assigned swing angle and the specific work position for each.

[Ext. pos. 1-4]

Edit the settings that apply for all external positions (see "Properties – External positions 1-4", page 112).

[Robotic arm]

Edit the robotic arm settings (see *"Properties – Robotic arm"*, page 113).

[Adjustment data]

Display the overview of the EEPROM data of the Swing Head. This dialog cannot be edited.

[Edit]

Edit specific settings of the selected external position (see "Properties – External position", page 114).

Properties – External positions 1-4

Swing Head: **Properties** ► **Ext. pos. 1-4**

Properties / Edit external positions 1-4

Sample Processor / Tower 1 / Swing Head 786

External positions 1-4

Shift position	0	mm
Rinse position	0	mm
Swing increment	10.00	°

Properties which apply to all four external positions can be defined in the dialog **Properties / Edit external positions 1-4**.

Shift position

Lift position at which the robotic arm rotates to the external positions.

Input range	0 to 'Maximum stroke path' mm
Default value	0 mm A lift position of 0 mm corresponds to the "home position", i.e. the lift is located at the upper stop position. The maximum stroke path is defined in the properties of the tower (see " <i>Maximum stroke path</i> ", page 110).

Rinse position

Lift position used for rinsing.

Input range	0 to 'Maximum stroke path' mm
Default value	0 mm A lift position of 0 mm corresponds to the "home position", i.e. the lift is located at the upper stop position. The maximum stroke path is defined in the properties of the tower (see " <i>Maximum stroke path</i> ", page 110).

Swing increment

Amount by which the robotic arm can be swung relative to its current position. This parameter is used in the command **MOVE** for the settings **Destination = swing +** and **swing -**.

Input range	0.00 to 180.00 °
Default value	10.00 °

Properties – Robotic arm

Swing Head: **Properties** ► **Robotic arm**

Properties / Robotic arm

Robotic arm Tower 1

Swing offset

0.00

°

Maximum swing angle

84.00

°

Swing radius

110.00

mm

Rotation offset

0.00

°

Swing direction

-



Properties of the robotic arm can be defined in the dialog **Properties / Robotic arm**.


Swing offset

The swing offset is a physical angle offset of a specific robotic arm model. The required values can be found in the leaflet for the robotic arm.

Input range	-270.00 to 270.00 °
Default value	0.00 °

Maximum swing angle

Useable swing range of the robotic arm. Each robotic arm model displays a different value on the basis of its construction. The range can also be reduced if necessary. You will find the required values in the leaflet for the robotic arm.



Properties which apply only to the selected external position can be defined in the dialog **Properties / Edit external position**.

Angle

Swing angle for the selected external position.

Input range	(Offset) to (Offset + max. swing range) °
Default value	60.00 °
	The offset is made up of a design-dependent angle (approx. 8...9°) together with the robotic arm offset from the robotic arm properties. The maximum swing range is also defined under the robotic arm properties (see "<i>Properties – Robotic arm</i>", page 113).

Work position

Work position for the selected external position.


Input range	0 to 'maximum stroke path' mm A lift position of 0 mm corresponds to the "home position", i.e. the lift is located at the upper stop position. The maximum stroke path is defined in the properties of the tower (see "<i>Maximum stroke path</i>", page 110).
-------------	--

Loading the sample rack

When you receive the file for a new sample rack from your Metrohm representative, you can easily import this file into your existing system. Proceed as follows:

1 Copying the file to an external storage medium

- Copy the file to a subdirectory of "Files".



If this structure is not maintained, the new rack will not be found because the software directly accesses subdirectories of "Files".

- Plug in the external storage medium at the Ti-Touch.

2 Display the list of saved sample racks

- In the **Edit device / Sample racks** dialog, tap on the **[Load]** button.

The selection of file groups on the external storage medium is displayed. If only one group is available, then the list of the saved sample rack files will be displayed directly.

- Select the group with the desired sample rack.
- Tap on **[Show files]**.

The list with the saved sample rack files is opened.

3 Load the sample rack file

- Select the desired file.
- Tap on **[Load]**.

The new sample rack is now loaded and appears in the list of available racks.

Input range	1.0 to 100.0 mm
Selection	off
Default value	off

No check takes place.

Beaker sensor

Selection	Tower Robotic arm off
Default value	off

A Swing Head with beaker sensor must be mounted. In addition, a suitable work position must be defined for the lift, so that the robotic arm touches the sample vessel. The work position is moved to for the purpose of beaker detection.

No check takes place.

1. Calibration pos.

Input range	1 to "highest rack position" It is imperative that the remaining buffers/standards be placed on the rack positions directly following.
Selection	Special beaker 1...n off
Default value	off

If a special beaker is selected as the first calibration position, then the number of buffers used for calibration will determine the number of special beakers to be defined. Any rack position can be defined as a special beaker. It is, however, preferable to set them at high rack positions in order to be able to begin sample series at rack position 1. The special beakers are moved to in ascending order.

Rack offset

■■■■■■■■ 119

Input range	-10.00 to 10.00 °
Default value	0.00 °

[Adjust rack]

Adjust the rack (see chapter 11.6.2, page 124).

[Lift pos. tower 1]/[Lift pos. tower 2]

Define rack specific lift positions (see "Lift positions Tower 1/2", page 120).

[Special beakers]

Define rack positions as special beakers (see *"Editing special beakers"*, page 121).

Lift positions Tower 1/2

Sample rack list: **Rack ▶ Edit ▶ Lift pos. tower 1 / Lift pos. tower 2**

Edit rack data / Lift positions tower 1

Rack 6.2041.410, Code 001010, 142 Positions

Work position

120

mm

Rinse position

0

mm

Shift position

0

mm

Special position

0

mm

You can define rack-specific lift positions in the dialog **Edit rack data / Lift positions tower 1/2**. These then apply for all rack positions except those that have been defined as special beakers.

These positions can be moved to directly under manual control and with the LIFT command. Only lift positions within the maximum stroke path can be entered. This is defined in the device properties of the tower.



NOTE

These lift positions can also be assigned directly in manual control after moving to the desired lift height (*see chapter 28.6.1, page 284*).

Work position

At this lift position the electrodes, stirrer and buret tips are optimally positioned for work.

Rinse position

Input range	0 to 'Maximum stroke path' mm
Default value	0 mm A lift position of 0 mm corresponds to the "home position", i.e. the lift is located at the upper stop position. The maximum stroke path is defined in the properties of the tower (<i>see "Maximum stroke path", page 110</i>).

The lift is raised to this position with each rotational movement of the rack whenever it is located at a lower lift position.

Special position

Input range	0 to 'maximum stroke path' mm A lift position of 0 mm corresponds to the "home position", i.e. the lift is located at the upper stop position. The maximum stroke path is defined in the properties of the tower (see "<i>Maximum stroke path</i>", page 110).
-------------	--

Sample rack list: **Rack** ► **Edit** ► **Special beakers**



Edit rack data / Special beakers	
Rack 6.2041.410, Code 001010, 142 Positions	
Special beaker	Rack position
Special beaker 1	142
Special beaker 2	0
Special beaker 3	0
Special beaker 4	0
Special beaker 5	0
Special beaker 6	0
Special beaker 7	0

A maximum of 16 rack positions can be defined as special beakers for each sample rack. The dialog **Edit rack data / Special beakers** shows a list of all special beakers with their assigned rack positions.

[Edit]

Edit the data of the selected special beaker, see the following.

Special beaker / Edit

Special beaker 1

Rack position

142

Work pos. tower 1

100

mm

Work pos. tower 2

0

mm

Beaker radius

off

mm

Beaker sensor

Tower

▼

You can edit the data of the selected special beaker in the dialog **Special beaker / Edit**.

Rack position

Number of rack position for selected special beaker. Any rack position can be defined as a special beaker. It is, however, preferable to set them at high rack positions in order to be able to begin sample series at rack position 1. Rack positions defined as special beakers can no longer be used as sample positions.

Input range	0 to 'maximum number of rack positions'
-------------	---

Work pos. tower 1/2

Work position for the selected special beaker. One specific work position can be defined for tower 1 and tower 2.

Input range	0 to 'maximum stroke path' mm A lift position of 0 mm corresponds to the "home position", i.e. the lift is located at the upper stop position. The maximum stroke path is defined in the properties of the tower (see " <i>Maximum stroke path</i> ", page 110).
-------------	---

Beaker radius

Actual radius of the selected special beaker.

This beaker radius may not be less than the minimum beaker radius defined in the tower properties (see "*Min. beaker radius*," page 110). If the lift is to be moved to the work position, then these two values will be compared with one another.

Input range	1.0 to 100.0 mm
Selection	off
Default value	off

off

No check takes place.

Beaker sensor

Each time this special beaker is to be moved to with the **MOVE** command, the beaker sensor checks whether a vessel is present. In the **MOVE** command, you define the action that takes place if the beaker sensor does not detect a vessel at the position being moved to.

Selection	Tower Robotic arm off
Default value	off

Robotic arm

A Swing Head with beaker sensor must be mounted. In addition, a suitable work position must be defined for the lift, so that the robotic arm touches the sample vessel. The work position is moved to for the purpose of beaker detection.

off

No check takes place.

The lift is moved to the uppermost stop.

11.7 Printer

Instrument list: **Printer ▶ Edit**

Device manager / Edit

Device type: Printer

Device name

Printer

Comment

Printer

PCL Inkjet ▼

Connector

USB ▼

PDF settings

Network printer

More parameters

The list of devices always includes a printer, even if no corresponding device is connected. In this dialog, you also define when you would like to have a report generated as a PDF file.

Device name

Here you can enter a device name of your choice.

Entry	24 characters maximum
-------	------------------------------

Comment

Entry	24 characters maximum
-------	------------------------------

Printer

Selection of the printer type.

Selection	PCL Inkjet PCL Laser HP Epson Epson new Canon CUSTOM NEOS off
-----------	--

PCL Inkjet
for HP DeskJet printers.

PCL Laser
for HP LaserJet printers.

CUSTOM NEOS
Due to the narrow paper size, only reports marked with # in the "List of all printable reports" (*see chapter 27.3, page 262*) can be printed.

off
The reports will not be printed out on paper.

Connector

Selection of the connection type for the printer.

Selection	USB Ethernet
Default value	USB

USB

Local printer at a USB interface.

Ethernet

Network printer.



NOTE

The **Ti-Touch** prints the reports with a fixed resolution of 300 dpi. If you are using a printer with a resolution of 360 dpi (or a multiple thereof, e.g. an Epson), then the text will be printed out somewhat smaller than with printers with a resolution of 300 dpi (or a multiple thereof, e.g. a Canon or HP).

11.7.1 PDF settings

Printer: **Edit** ► **PDF settings**

The settings for saving a report as a PDF file are defined in this dialog.

Memory

Memory location where the PDF file will be saved. The report will be saved in the directory *PDF_Report*. This directory is created when a PDF file is generated. The file name is appended with the date and time when the determination was started and the date and time when the PDF was generated (e.g.: **Identification 1**-20130320-141044-20130320-141220). The first date and time set indicates the start of the determination, the second date and time set indicates the generation of the PDF file. This prevents previously created PDF files from being overwritten each time a new file is created, e.g. when a result is recalculated and the report is generated again.

All three memory locations are listed as possible selections, even if they cannot be accessed at the moment.

Selection	off External memory 1 External memory 2 Shared memory
Default value	off

off

The report will not be saved as a PDF file.

Shared memory

The report will be saved in a shared directory on the network. The shared directory is selected in the **Edit device / Shared memory** dialog (see chapter 11.3.3, page 102).

Copy or extract content allowed

on | off (Default value: **on**)

If this option is activated, then content can be copied or deleted from the PDF file.

Printing allowed

on | off (Default value: **on**)

If this option is activated, then the PDF file can be printed.

Change the document allowed

on | off (Default value: **off**)

If this option is activated, then the PDF file can be edited.

Add or change comments allowed

on | off (Default value: **on**)

If this option is activated, then comments can be added to the PDF file.

11.7.2 Network printer

Printer: **Edit** ► **Network printer**

If you have connected your Ti-Touch to your network, then you can specify a network printer for your reports in this dialog.



NOTE

If you have shared a memory location on a PC within your network and if you have configured the Ti-Touch accordingly (*see chapter 11.3.3, page 102*), then the settings for the parameters **Domain**, **User** and **Password** will be applied and can then no longer be modified in this dialog.



CAUTION

If the computer on which the network printer is configured does not have the same subnet as the Ti-Touch, then a WINS server must be present. The computer must be entered in this server. In order for the communication to work, the SMB protocol version 1 (SMBv1) must be activated. Since this function is not standard anymore in Windows 10, it has to be activated manually. From firmware version 5.916.0044, the SMB protocol version 2 and 3 is also supported.

Print server

Host name of the print server or the computer on which the network printer is configured. If you are working with Windows, then you will find the host name as follows: In the Windows **command prompt** window, enter the command **ipconfig -all**. The host name is listed together with other parameters of the print server or the computer.

Do not under any circumstances enter an IP address in this input field.

Entry	max. 32 characters
Default value	empty

Share name

Share name of the network printer on the above-specified computer. Please note that the Share name of a network printer often does not match the name of the printer. The share name can be found on the Release tab in the Properties dialog of the released network printer.

Entry	max. 32 characters
Default value	empty

Domain

Network domain in which the above-specified computer is located. If you are working with Windows, then you will find the domain name as follows: In the Windows **command window**, enter the command **ipconfig -all**. The domain is listed together with other parameters of the computer. Leave the field empty if the computer is not located within a domain.

Entry	max. 32 characters
Default value	empty

User

User name of the user authorized to access the shared network printer.

Entry	max. 32 characters
Default value	empty

Device manager / Edit

Device type: Balance

Device name

Balance

Comment

Balance type

Sartorius

▼

RS-232

COM1

▼

You define the type of balance and its connector in the dialog **Device manager / Edit**.

Device name

Here you can enter a device name of your choice.

Entry	24 characters maximum
-------	-----------------------

Comment

Entry	24 characters maximum
-------	-----------------------

Balance type

If you have connected a balance, then you have to specify the balance type here.

Selection	AND Mettler Mettler AT Mettler AX Ohaus Precisa Sartorius Shimadzu
Default value	Sartorius

RS-232

RS-232 interface to which the balance is connected. The interface parameters are adjusted in the **Edit device / Port parameters** dialog (see chapter 11.9, page 131).

Selection	COM 1 COM 2
Default value	COM 1

COM 2

The interface is disabled.

The following table indicates the balance type that needs to be selected for the balance model:

Balance	Balance type
AND	AND
Mettler AB, AE, AG, AM, AJ, PE, PM, PJ, PR, XP, XS	Mettler
Mettler AT	Mettler AT
Mettler AX, MX, UMX, PG, AB-S, PB-S	Mettler AX
Ohaus Voyager, Explorer, Analytical Plus	Ohaus
Precisa	Precisa
Sartorius	Sartorius
Shimadzu BX, BW	Shimadzu

11.9 USB/RS-232 adapter

Instrument list: **USB/RS-232 adapter** ► **Edit** ► **COM interface** ► **Edit**

As a rule, balances are equipped with a serial RS-232 interface. To connect a balance, you will require the 6.2148.050 cable. You can adjust the interface parameters in the dialog **Edit device / Port parameters**. These settings must match the settings on the connected device.



NOTE

If you make changes to parameter settings in this dialog, then you must switch the Ti-Touch off and back on in order for the changes to take effect.

Edit device / Port parameters

Baud rate

9600

▼

Data bits

8

▼

Parity

none

▼

Stop bits

1

▼

Handshake

Hardware (DTR/CTS)

▼

Baud rate

Transfer rate in characters per second.

Selection	1200 2400 4800 9600 19200 38400 57600 115200
Default value	9600

Data bits

Number of data bits.

Selection	7 8
Default value	8

Parity

Type of parity testing.

Selection	even odd none
Default value	none

Stop bits

Number of stop bits.

Selection	1 2
Default value	1

Handshake

Type of data transfer protocol.



NOTE

In case of communication problems, try the software handshake (**Software (XON/XOFF)**).

Selection	none Software (XON/XOFF) Hardware (DTR/CTS)
Default value	Hardware (DTR/CTS)

Software (XON/XOFF)

Use the software handshake when you send a PC/LIMS report via an RS-232 interface.

11.10 PC keyboard

Instrument list: **PC keyboard** ► **Edit**

A commercially available USB keyboard can be connected to make text and numerical input easier. It will be recognized automatically and entered in the list of devices with default settings.

Device manager / Edit

Device type: PC keyboard

Device name

PC keyboard

Comment

Keyboard layout

English US



Device name

Here you can enter a device name of your choice.

Entry	24 characters maximum
-------	-----------------------

Comment

Entry	24 characters maximum
-------	-----------------------

Keyboard layout

Define the country-specific keyboard layout here.

Selection	English US German DE French FR Spanish ES German CH
Default value	English US

To enter texts and numbers with the PC keyboard, the appropriate text and number input dialog must be opened on the Ti-Touch. Only the following keys on the PC keyboard have a function:

Table 6 Key assignment

Function in the editing dialog	PC keyboard key
[Cancel]	[Esc]

A barcode reader can be connected to read in sample data or other texts. A connected barcode reader will be recognized automatically and entered in the list of devices with default settings.

You will hear an acoustic signal as confirmation that a character string has been transmitted by the barcode reader and accepted.

NOTE

The data will be read only if the system is in normal status, i.e. when no determination is running.

Here you can enter a device name of your choice.

Selection of the input field for the character string read in by the barcode reader.

The character string is entered in the input field of the opened text- or number-input dialog.

The character string is entered in the input field **Method**.

The character string is entered in the input field **Identification 1**.



Identification 2

The character string is entered in the input field **Identification 2**.

Sample size

The character string is entered in the input field **Sample size**. Character strings containing characters other than numbers and decimal separators will be ignored.

Keyboard layout

Specify the country-specific keyboard layout for the emulation of the PC keyboard. This setting must match the setting on the barcode reader (see documentation for the barcode reader).

Selection	English US German DE French FR Spanish ES German CH
Default value	English US



12 File manager

Main dialog: **System ► File manager**

The saved methods, determinations, sample tables, etc. are managed in the file manager. You can also create a backup of your system (all data and settings). Also, an existing backup can be reloaded.

The file memory is organized as follows:

- **Internal memory**

The following files can be stored in the internal memory:

- Methods

- **External memory**

You can use a USB flash drive as an auxiliary storage medium, for instance. The following files can be stored on an external memory:

- Backup
- Methods
- Determinations
- Sample tables
- Result tables
- Reports as *.pdf
- PC/LIMS report

12.1 Managing files

File manager: **Internal memory / External memory 1 / External memory 2**

The saved files can be organized into groups. These groups are comparable to file directories on your PC, although unlike your PC, only one level is possible.




NOTE

The file names must be unique for each memory location, i.e. you cannot save two files with the same name, even if they are in different groups.



NOTE

If you use an external storage medium with the FAT or FAT32 file system, then you can save a maximum of 999 files per group. If you find it necessary to store more than 999 files in a single group, then you must reformat the storage medium with the file system **ExFAT** (see *chapter 32.11.2, page 474*).




[Show all]

Display the list of all the files in the selected memory location.

[Show files]

Display the list of all the files in the selected group.



[Rename]

Rename the selected file (see chapter 12.1.2, page 139).

[Properties]

Display the file properties (see chapter 12.1.3, page 140).

[Copy]

Copy the selected file to a different memory location (*see chapter 12.1.1, page 139*).

[Delete]

Delete the selected file.

[Load]

Load the selected file.

12.1.1 Copying a file

Proceed as follows to copy a file:

1 Selecting the file

- Select the desired file.
- Tap on **[Copy]**.

The selection of memory locations is displayed.

2 Copy the file

- Tap on the desired memory location.
Only memory locations currently being accessible are active.

The file is copied. The file group is retained, i.e. the group will be created again in the new memory location if it does not already exist there.

12.1.2 Renaming a file

Proceed as follows to rename a file:

1 Selecting the file

- Select the desired file.
- Tap on **[Rename]**.


The **Show files / Rename** dialog is displayed.

2 Change the file name

- Tap on the input field **File name**.
The text editor opens.
- Enter a new file name (max. 32 characters) and confirm with **[OK]**.
- Tap on **[OK]**.

The file is saved under the new name.

12.2 External storage medium



[Backup]

Creating a backup of all data and settings on this storage medium (*see chapter 12.3, page 143*).

[Restore]

Load the backup. This function is active only if a backup is available (see *chapter 12.3, page 143*).

[Memory info]

Display information on the storage medium, see following chapter.

Storage medium – Information

Detailed information about the storage medium is displayed in the dialog **External memory 1/2 / Memory info**:

- **Name**
Designation of the storage medium.
- **Write protection**
If this option is activated, then the file cannot be saved, deleted, moved or renamed. This is only an internal write protection feature and is independent of the write protection of the operating system on your computer.
- **Storage capacity**
Total capacity of the storage medium.
- **Used memory**
Storage capacity used on the storage medium.
- **Free memory**
Free storage capacity on the storage medium.

Directory structure

The directory structure appears as follows on the external storage medium:




Figure 24 Directory structure on the external storage medium

Backup	All of the files of the backup are stored in this directory. The directory will be created the first time a backup is created.
Files	This directory contains all of the groups as well as the files contained in them (methods, determinations, etc.).
pc_lims_report	PC/LIMS reports are stored in this directory as TXT files. The directory will be created the first time a PC/LIMS report is printed out.
PDF_Report	PDF reports are stored in this directory. The directory will be created the first time a report is saved as a PDF file.
Profiles	The identification profile is stored in this directory. The directory will be created the first time an identification profile is generated.
SwUpdates	Files for software updates must be saved in this directory.

12.3 Creating backups / Restoring data

File manager: **External memory 1 / External memory 2**

You can use the **[Backup]** function to easily create a backup containing all the data and settings of your system. You should make a backup at regular intervals in order to avoid data loss.



NOTE

Only **one** backup can be created on a storage medium.

If a backup has already been stored on the medium, then it will be overwritten when this function is carried out again.

12.3.1 Restoring data

You can use the function **[Restore]** to restore either a complete backup or only certain data.



CAUTION

All of the methods in the internal memory will be deleted.

Backups are not backward compatible, i.e. backups of the latest version cannot be restored using earlier versions.

Proceed as follows to restore backed-up data:

1 Select the storage medium

- Connect the storage medium on which the backup is located.
- Select the storage medium under **System ► File manager**.

2 Select data

- Tap on **[Backup]**.
The selection of the data that can be restored separately is displayed (see the following).
- Deselect the data sets which are not to be restored.
- Tap on **[Load data]**.

The selected data blocks are restored.

3 Complete the restoration

- Switch the Ti-Touch off and on again.

13 GLP manager

Main dialog: **System ► GLP manager**

In the GLP manager you can document data for various GLP tests. The results of the automatic system tests carried out after switching on are also documented.

The following tests can be documented:

- GLP test "Measurement" (*see chapter 13.3, page 147*)
- GLP test "Titration" (*see chapter 13.3, page 147*)
- System validation (*see chapter 13.4, page 151*)

In addition, you can also:

- Create a list with your test tools (*see chapter 13.2, page 146*).
- Define a service interval for having routine maintenance carried out by a Metrohm service technician (*see chapter 13.5, page 155*).
- Define an interval for the regular performance of backups (*see chapter 13.5, page 155*).

Further information about quality management and validation is also given in the documentation series **Quality Management with Metrohm** and **Application Bulletin AB 252** (Validation of Metrohm titrators (potentiometric) in accordance with GLP/ISO 9001).

System / GLP manager		
Test	Last test	Next test
Automatic system test	2010-08-30	Power on
System validation	2010-08-02	2011-01-04
GLP test "Measurement"	2010-07-01	2011-01-03
GLP test "Titration"	2010-07-02	2011-01-04
Test tools GLP tests HW/SW System validation Monitoring View test data		

The table shows the last time that each test was carried out and when the next test is to be carried out. A test is entered in the list the first time that it is documented.

[Test tools]

Configuring the test tools for GLP tests (*see chapter 13.2, page 146*).

[GLP tests HW/SW]

Document the GLP tests "Measurement" and "Titration" (*see chapter 13.3, page 147*).

[System validation]

Document the system validation (see chapter 13.4, page 151).

[Monitoring]

Defining the time intervals for system maintenance and backups (*see chapter 13.5, page 155*).

[\[View test data\]](#)

Display the data of the selected test.

13.1 Automatic system test

The system test is carried out automatically when the Ti-Touch is switched on. The result of each individual test is shown in green if no error has occurred. If a result is shown in red, then an error occurred when the respective test was carried out. Switch the Ti-Touch off and back on again. If the error still occurs please notify Metrohm Service.

Print report at system start

on | off (Default value: **off**)

If this parameter is activated, then the result of the system test will be printed out automatically when the Ti-Touch is switched on.

13.2 Test tools

Main dialog: **System ► GLP manager ► Test tools**

You can create a list of test tools to be used in the tests in the GLP manager.

The following test tools have already been defined:

- **767 Calibrated Reference:** Device for checking measuring inputs and electrode cable.
- **822 Titr.Curve Simulator:** Curve simulator for checking the hardware and software.
- **773 pH/mV Simulator:** Device for checking measuring inputs and electrode cable.
- **868 UR Generator:** Device for checking measuring inputs and electrode cable.

[New]

Add a new test tool to the list.

[Delete]

Delete the selected test tool from the list.

[Edit]

Change the designation of the selected test tool.

13.3 GLP tests for measurement and titration

Main dialog: **System ► GLP manager ► GLP tests HW/SW ► GLP test "Measurement" / GLP test "Titration"**

You can document the GLP tests for measurements and titrations in the dialog **GLP manager / GLP tests Hardware/Software** and its subdialogs. The following procedure describes the GLP test "Measurement," but also applies for the GLP test "Titration".

Proceed as follows:

1 Open the properties dialog

- In the **System / GLP manager** dialog, tap on the **[GLP tests HW/SW]** button and then on the **[GLP test "Measurement"]** button.

2 Edit data

- Define the test method, the test results, etc. (see "Dialogs "GLP tests HW/SW / GLP test "Measurement"" and "GLP tests HW/SW / GLP test "Titration"", page 148).

3 Define the hardware used

- Tap on **[Hardware]**.

- Select the hardware used for the test (see "Dialogs "GLP test "Measurement" / Hardware" and "GLP test "Titration" / Hardware", page 150).
- Tap on the fixed key .

The dialog **GLP tests HW/SW / GLP test "Measurement"** appears again.

4 Define the test interval

- Tap on **[GLP test interval]**.
- Enter the time interval or the date for the next GLP test (see "Diagnosis "GLP test "Measurement" / Test interval" and "GLP test "Titration" / Test interval", page 149).
- Tap on the fixed key **[↩]**.

13.3.1 Parameter description

Dialogs "GLP tests HW/SW / GLP test "Measurement"" and "GLP tests HW/SW / GLP test "Titration""

Method

Method the GLP test has been carried out with.

Entry	32 characters maximum
Selection	Selection of methods stored in the internal memory

User

User who carried out the GLP test.

Entry	24 characters maximum
Selection	Selection of configured users

Test date

Date on which the GLP was carried out.

Format: YYYY:MM:DD

Test result

Result of the test.

Selection	Test OK Test not OK
Default value	Test OK

Comment

Entry	24 characters maximum
-------	-----------------------

[GLP test interval]

Defining the time interval for the GLP test (see "Dialogs "GLP test "Measurement" / Test interval" and "GLP test "Titration" / Test interval", page 149).

[Hardware]

Documenting the hardware with which the GLP test was carried out (see "Dialogs "GLP test "Measurement" / Hardware" and "GLP test "Titration" / Hardware", page 150).

Dialogs "GLP test "Measurement" / Test interval" and "GLP test "Titration" / Test interval"

Test date

Date on which the GLP was carried out.

Format: YYYY:MM:DD

Monitoring

on | off (Default value: **off**)

If this parameter is activated, then the time interval after which a GLP test has to be carried out again will be monitored.

GLP test interval

If you define a time interval for the GLP test, then the date in **Next GLP test** will be tracked automatically.

Input range	1 to 999 days
Default value	999 days

Next GLP test

If you define a date for the next GLP test, then the **GLP test interval** will be tracked automatically.

Format: YYYY:MM:DD

Action

Selection of the action which is carried out when the time interval has expired.

Selection	Display message Document message Cancel determination
Default value	Display message For all three options it is documented in the determination data (see dialog More determination data / Messages), that the time interval has been expired.

13.4 System validation

Main dialog: **System ► GLP manager ► System validation**

In the dialog **GLP manager / System validation** and its subdialogs, you can document the system validation results and define the time interval after which validation must be carried out again.

Proceed as follows:

1 Open the properties dialog

- In the **System / GLP manager** dialog, tap on the **[System validation]** button.

2 Edit data

- Define the method, the result, etc. (see *"GLP manager / System validation" dialog*, page 152).

3 Define the validation interval

- Tap on **[Validation interval]**.
- Enter the time interval or the date for the next system validation (see *"Dialog "System validation / Validation interval"*, page 153).
- Tap on the fixed key [**↩**].

The dialog **GLP manager / System validation** appears again.

4 Enter statistical data

- Tap on **[Test statistics]**.
- Enter the statistical data for the last system validation (see *"System validation / Test statistics" dialog*, page 154).
- Tap on the fixed key [**↩**].

The dialog **GLP manager / System validation** appears again.

5 Insert a note

- Tap on **[Note (SOP)]**.
- Enter a brief description, e.g. a summary of the SOP (standard operating procedure) according to which the system validation was carried out (*see "System validation / Note (SOP)" dialog", page 153*).
- Tap on the fixed key **[↩]**.

13.4.1 Parameter description

"GLP manager / System validation" dialog

Method

Method with which the system validation has been carried out.

Entry	32 characters maximum
Selection	Selection of methods stored in the internal memory

User

User who carried out the system validation.

Entry	24 characters maximum
Selection	Selection of configured users

Test date

Date on which the system validation was carried out.

Format: YYYY:MM:DD

Test result

Result of the test.

Selection	Test OK Test not OK
Default value	Test OK

Comment

Entry	24 characters maximum
-------	-----------------------

[Note (SOP)]

Enter a brief description, e.g. a summary of the SOP (standard operating procedure) according to which the system validation has been carried out.

[Validation interval]

Defining the time interval for the system validation (see "Dialog "System validation / Validation interval"", page 153).

14 Common variables

Main dialog: **System ► Common variables**

You can save 25 **method-independent variables**, or common variables. These variables can be used in future calculations (as variables **CV01... CV25**). Common variables are useful, e.g. for the following applications:

- Determination of a blank value which will be taken into account during the content determination of the sample.
- Determination of the content of a standard solution, which will be taken into account during the content determination of the sample.

System / Common variables		
CV	Name	Value
01	Blank value	0.0143 mL
02	Factor	1.059
03		
04		
05		
06	Density	0.986 g/mL
07		
08		
		Delete Edit

The designation and the value (including the unit) are displayed for every common variable in the list.

[Delete]

Delete the selected common variable from the list.

[Edit]

Edit the data of the selected common variable, see following chapter.

14.1 Editing common variables

List of common variables: **Common Variable** ► **Edit**

The common variables can be modified as follows:

- Edit manually, see the following.
- Automatic assignment from the determination run. A calculation result must be configured accordingly for this purpose (*see chapter 14.4, page 161*).

Common variables / Edit

CV01: Blank value

Name

Blank value

Value

0.0143

mL

▼

Validity

Properties

Name

Designation of the common variable.

Entry	24 characters maximum
Default value	empty

Value

Value of the common variable.

Input range	–999999999 to 999999999
Default value	empty

Unit of the common variable.

Entry	10 characters maximum
Default value	empty
Selection	% mol/L mmol/L g/L mg/L mg/mL mg/100 g ppm g mg µg mL µL mg/ piece °C µg/min mL/min µL/min
Default value	%

Document message

In the determination data it will be documented, that the time interval has been expired.

Cancel determination

The determination is stopped.

14.4 Assigning a result automatically to a common variable

**NOTE**

This instruction is based on the assumption that the method contains a calculation command with a calculation.

Proceed as follows to assign a result to a common variable:

1 Opening the editing dialog of the result

- In the command list, select the command **CALC**.
- Tap on the **[Edit command]** button.
- Select the result whose value is to be assigned to a common variable and tap on **[Edit]**.

The editing dialog of the result is displayed.

2 Define result options

- Tap on the **[Result options]** button.

- Activate the parameter **Save as common variable**.
- Define **Variable** = CV01...CV25.



3 Saving the settings

Tap on the fixed keys [.

In the future, the calculated result will be assigned to the selected common variable (result name, value and unit).

15 Templates

Main dialog: **System ► Templates**

You have the option of defining **system-specific templates**. You can use these templates when editing the respective data.

You can create the following templates:


- **Sample data** (see chapter 15.1, page 163)
Create sample identifications and sample assignments.
- **Custom result templates** (see chapter 15.2, page 167)
Define formulas for result calculations.
- **Input lines** (see chapter 15.3, page 171)
Define input signals on the remote interface.
- **Output lines** (see chapter 15.4, page 173)
Defining output signals on the remote interface.
- **Custom calib. buffers** (see chapter 15.5, page 176)
Define a buffer series for the calibration of pH electrodes.
- **Report header** (see chapter 15.6, page 178)
Create a system-specific report header that contains information concerning the laboratory, for example.
- **Custom electrode type** (see chapter 15.7, page 179)
Define customized limit values for the electrode test of pH electrodes.

15.1 Sample data

Main dialog: **System ► Templates ► Sample data**

The screenshot shows the 'Templates / Sample data' dialog box. It has a title bar with the text 'Templates / Sample data' and a green icon. The main area is light blue and contains several sections. The first section is 'Sample identification list' with a checkbox 'Use sample identification list' that is checked. The second section is 'Sample assignment table' with a checkbox 'Use sample assignment table' that is checked. Below this is an 'Assignment ident.' section with a dropdown menu showing 'Identification 1'. The last section is 'Request assignment identification' with an unchecked checkbox.

You can create the following templates in the dialog **Templates / Sample data**:



You can create a system-specific list with sample identifications in the dialog **Sample data / Sample identification list**. This list makes it easier for you to enter frequently used sample identifications, i.e. this list is available to you for entering them (main dialog, command **REQUEST**, etc.). It may be a good idea to define the unchanging part of the identification as a template and to add the variable part during sample data input.

[New]

Add a new sample identification to the list.

[Delete]

Delete the selected sample identification from the list.

[Edit]

Change the designation of the selected sample identification.

15.1.2 Sample assignment table

Main dialog: **System ► Templates ► Sample data ► Sample assignment table**

Sample data / Sample assignment table	
Sample ident.	Method
Apple juice	Me2115
Mineral water	Me3901
Orange juice	Me4612
Tap water *	Me4155

New
Delete
Edit

In the sample assignment table, a particular method is assigned to a sample identification. In this way, you ensure that your samples will be processed with the correct method; mix-ups are no longer possible. When you start a determination, you need only enter the sample identification; the method is loaded automatically.

[New]

Add a new sample assignment to the list, see following chapter.

[Delete]

Delete the selected sample assignment.

[Edit]

Edit the selected sample assignment, see following chapter.

Editing the sample assignment

Identification

Identification of the sample.



NOTE

You can place an `*` as a wildcard at the beginning or end of the character string. Doing so allows you to prefix or suffix a sequential number, for example, to the identification, which will be ignored when methods are being assigned.

Entry	24 characters maximum
Selection	Selection of defined identifications in the sample identification list

Memory

Memory location the method is loaded from. All memory locations are selectable, even if they are currently not accessible.

Selection	Internal memory External memory 1 External memory 2 Shared memory
Default value	Internal memory

Shared memory

Shared directory in the network.


Method

Method that is loaded when a sample with the above-specified identification is processed. You can also enter a method that does not yet exist in the selected memory. When a determination is started there is a check whether the method is available.

Entry	32 characters maximum
Selection	Selection of methods stored in the selected memory

15.2 Custom result templates

Main dialog: **System ► Templates ► Custom result templates**



In the dialog **Templates / Custom result templates**, you can define the result calculations which can be loaded in the calculation command **CALC**. With the exception of the following points, creating a result template is identical to editing a calculation in the command CALC (*see chapter 29.10.1.2, page 392*):

- No result variable can be assigned to the template.
- The definition of limit values is not possible.



NOTE

Up to nine **wildcards** can be inserted in the calculation formula, the **Variables F1...F9**. Use these variables, e.g. for the molar mass of your analyte. If you have loaded a result template with these wildcards in the calculation command, then you will automatically be prompted to enter the numerical values.

[New]

Add a new result template to the list, see following chapter.

[Delete]

Delete the selected result template from the list.

[Edit]

Edit the data of the selected result template, see following chapter.

Variable for mean value

If the statistic calculations have been activated (see method options), the mean value of the single results will be saved as variable SMN1 to SMN9. As default value, always the first free variable is being displayed.

Selection	SMN1 ... SMN9	off
-----------	----------------------	------------

off

For the result not statistic calculations are carried out.

Save as titer

The result can be saved as titer for the selected titrant.



NOTE

The titer is assigned to the titrant defined in the last titration command before the calculation. Please note that the CALC command containing the titer assignment is inserted after the titration command with which the titer is determined.

Selection	Single value Mean value off
Default value	off

Single value

The result of the current determination is saved as titer.

Mean value

If statistic calculations are carried out for the result, the current mean value of the determination series is saved.

off

The result is not saved as titer.

Save as common variable

on | off (Default value: **off**)

The calculated result can be saved as a method-independent variable, a so-called common variable. The result is then also available in other methods for calculations. In all cases, only the single value will be saved, even if the statistics function is enabled.

Variable

Selection of the common variable to which the result is being assigned.

Selection **CV01 ... CV25**

15.3 Input lines

Main dialog: **System** ► **Templates** ► **Input lines**

Templates / Input lines	
Signal name	Input signal
Cond OK	*****1*
End1	*****1***
End2	*1*****
EndMeter	***11***
Ready*	**1*****1
Ready1	*****1
Ready2	**1*****
Sample ready	***1****

In the dialog **Templates / Input lines**, you can define the input signals at the remote interface as a template. You can select these templates in the command **SCAN**. The list can contain a maximum of 20 templates.

Table 7 List of predefined input signals

Signal name	Input signal	Function
Cond OK	*****1*	queries "Cond OK" condition of device.
End1	****1***	waits for the EOD pulse from Device 1 (Titrimo, Titrimo).
End2	*1*****	waits for the EOD pulse from Device 2.
EndMeter	***11***	waits for the EOD pulse from 780/781 pH/Ion meter (stirrer 1 will be switched on during the waiting period).
Ready*	**1****1	queries the "Ready" condition of Devices 1 and 2 (Titrimo, Titrimo). The status of devices working in parallel can be scanned with this signal. In this case the "Ready" line of both devices must be statically set (i.e. constantly) at the end of a determination. Devices which only transmit a brief impulse cannot be controlled in parallel.
Ready1	*****1	queries the "Ready" condition of Device 1.
Ready2	**1*****	queries the "Ready" condition of Device 2.

15.4 Output lines

Main dialog: **System** ► **Templates** ► **Output lines**

Templates / Output lines	
Signal name	Output signal
Continue	*****p****
Init	000000000000000
Meter Cal C	*****1001*
Meter Cal pH	*****0101*
Meter enter	*****1111*
Meter Mode C	*****1000*
Meter Mode I	*****0100*
Meter Mode pH	*****0001*

In the dialog **Templates / Output lines**, you can define the output signals at the remote interface as a template. You can select these templates in the command **CTRL**. The list can contain a maximum of 20 templates.

Table 8 List of predefined output signals

Signal name	Output signal	Function
Continue	*****p****	sends a stepping pulse to the connected Sample Processor.
Init	00000000000000	initializes the remote interface.
Meter Cal C	*****1001*	switches the 781 pH/Ion meter to concentration calibration.
Meter Cal pH	*****0101*	switches the 780/781 pH/Ion meter to pH calibration and starts the calibration.
Meter enter	*****1111*	simulates the [Enter] key of the 780/781 pH/Ion meter; mandatory for pH calibration in order to start the measurement of the second buffer.
Meter Mode C	*****1000*	switches the 781 pH/Ion meter to concentration measurement.
Meter Mode I	*****0100*	switches the 780/781 pH/Ion meter to voltametric measurement with polarization current and starts the measurement.
Meter Mode pH	*****0001*	switches the 780/781 pH/Ion meter to pH measurement and starts the measurement.

15.4.1 Editing the output signal

List of input signals: **Output signal** ► **New / Edit**

Signal name

Designation of the template.

Entry	24 characters maximum
-------	-----------------------

Output signal

Entering the desired bit pattern:

Entering the bit pattern:

- 0 = line inactive
- 1 = line active
- * = retain line status
- p = set pulse

The output lines are always numbered from right to left, i.e. with the signal *******1** line 0 is set.



NOTE

We recommend masking lines that are of no interest or for which no defined condition can be predicted with an asterisk (*).

Entry	Bit pattern of exactly 14 characters
Default value	*****

Pulse length

Duration of the sent pulses.

Input range	100 to 1000 ms
Default value	200 ms

16 Methods


16.1 Creating a new method

Proceed as follows to create a new method:

1 Opening the method table

- In the main dialog, tap on **[Load method]**.
- Tap on **[New method]**.

The method table with the stored templates opens:



2 Loading the method

- Tap on **[Empty method]**.

or

- Select the desired template and tap on **[Load template]**.

The method is now loaded. **New method** is displayed in the main dialog in the title line.

If a new method has been created, then the individual parameters can be modified with **[Edit parameters]**.

- The method will be saved and the command sequence is displayed.

- Tap on the input field **File name**.
The text editor opens.
- Enter a new file name (max. 32 characters) and confirm with **[OK]**.
- Tap on **[Save]**.

The method will be saved and the command sequence is displayed.


16.3 Loading a method

Proceed as follows to load a method:

1 Opening the method table

- In the main dialog, tap on **[Load method]**.
The selection of memory locations is displayed. This selection will be skipped if only the internal memory is accessible.
- Select the memory location where the desired method is saved.
The selection of file groups is displayed (*see chapter 12.1, page 137*). If only one group is available, then the method table will be displayed directly.
- Select the group with the desired method.
- Tap on **[Show files]**.

The method table with the stored methods opens:



2 Selecting the method

- Select the desired method.

[Insert command]

Insert a new method command. It is inserted before the selected command.

**NOTE**

Not all commands are available for subsequences. Commands which cannot be inserted into subsequences are disabled.

[Delete command]

Delete the selected method command.

[Edit command]

Edit the selected method command.

16.4.1 Inserting a command

A method is comprised of individual commands. When you start a determination, the commands are executed one after the other.

The following constraints apply:


- Maximum of 99 commands
- Maximum of nine commands for titrations, measurements and calibrations

Inserting a command

To insert a command in a method, proceed as follows:

1 Open the command overview

- In the main dialog, tap on **[Edit parameters]**.
- Tap on the **[Insert command]** button.



This dialog contains all of the commands which can be inserted into a method run, organized into thematic groups.

2 Select a command group

- Tap on the desired command group.

In the case of **[Report]**, the command is inserted directly into the command list. For all other command groups, the available commands will be displayed (titration / measuring commands, etc.).


3 Inserting a command

Tap on the desired command.

The dialog is closed and the command is displayed in the command list.

16.5 Method options

Main dialog: **Edit parameters** ► **Method options**



Method options are settings which apply to the method as a whole and not to an individual command.

The following settings are available:

- Activate/deactivate statistics calculations
- Save determination automatically
- Sign a method electronically
- Check the method when the determination starts
- Define a fixed sample size
- etc.

Method check at start

on | off (Default value: **on**)

If this option is activated, then the following tests are carried out at the method start.

- Are the required control and peripheral devices connected?
- Are the titrant and the sensor on hand?
- Have the monitoring intervals expired?
- etc.



NOTE

We recommend that the method check only be deactivated if the titration must begin immediately after the start of the determination (e.g. with very rapid reactions such as enzyme-catalyzed reactions). If the method check is deactivated, it could happen that the determination run will be interrupted, e.g. because devices, titrants or sensors are not on hand.

The following settings can be implemented for a connected Sample Processor:

Rack reset

on | off (Default value: **off**)

If this option is activated, then the sample rack is initialized at the time of the method start. The following actions are carried out:

- The lift is moved upwards.
- The sample rack is rotated to the starting position.
- The rack code of the rack attached is being read out and the corresponding rack data is transferred to the Sample Processor.
- The sample variable is reset to the value 1.

Check rack

Definition of the sample rack which must be attached at the time of the method start. This ensures that the method will be carried out only with this sample rack.



NOTE

The sample rack can, however, only be checked if the **Rack reset** option is activated.

Selection	no	Selection of the configured sample racks
Default value	no	

Increase sample variable automatically

on | off (Default value: **on**)

The sample variable describes the current position of the sample on the rack of the Sample Processor. If this option is activated, then the sample variable is automatically increased by 1 at the end of each determination.

This option must be deactivated if the sample variable is to be altered specifically with a **SAMPLE** command.

Request rack position at start

on | off (Default value: **off**)

If this option is activated, then the rack position of the first sample to be processed is requested.

16.5.2 Stop options

Main dialog: **Edit parameters ► Method options ► Start/Stop options ► Stop options**

In the dialog **Method options / Stop options**, you can define the actions which are carried out when a method is canceled.

The method can be canceled as follows:

- Manual stop with the fixed key []
- Stop because of an error
- Stop by remote signal via the Control Remote Box

Switch off pumps

on | off (Default value: **on**)

If this option is activated, then all of the connected pumps will be switched off.

Switch off stirrers

on | off (Default value: **on**)

If this option is activated, then all of the connected stirrers will be switched off.

Set remote lines

Selection of the signal out of the templates or entering the required bit pattern. Templates are defined under **System ► Templates ► Output lines**.



NOTE

A line set active is not being reset automatically, not even at the end of the determination.

Entering the bit pattern:

- 0 = line inactive
- 1 = line active
- * = retain line status
- p = set pulse

The output lines are always numbered from right to left, i.e. with the signal *******1** line 0 is set. With a pulse, the length is set to 200 ms. If you wish to set pulses with other lengths, you have to define a corresponding template.

Entry	Bit pattern containing exactly 14 characters or a max. of 24 characters for the name of the template
Default value	*****
Selection	Selection of the templates defined

16.5.3 Sample data

Main dialog: **Edit parameters** ► **Method options** ► **Sample data**

You can define the following method-specific settings in the dialog

Method options / Sample data:

- Modify the designation for the sample identifications
- Define a fixed sample size
- Define limits for the sample size

If you define limits, these will be monitored at the start and at the end of the determination.

Below you will find an example demonstrating the effects on the main dialog of the settings in this dialog (see *"Modifying sample data", page 194*).

Identification 1, Identification 2

on | off (Default value: **on**)

If this option is activated, then the input field for the sample identification in the main dialog will be displayed.

Designation

A designation of your own choosing can be defined for the input field for each method.

Entry	16 characters maximum
-------	------------------------------

Fixed sample size

on | off (Default value: **off**)

Activate this parameter if the same sample amount is always to be used for all determinations. If this is case, then you can define the sample size here. It will be displayed in the main dialog, but it can no longer be changed there.

Sample size

This parameter can only be edited when **Fixed sample size** is activated.

Input range	–999999999 to 999999999
Default value	1.0
Unit of sample size.	
Selection	g mg µg mL µL pieces
Default value	g



NOTE

You can also use the parameter **Fixed sample size** if you would like to establish a default for the sample size.

Example: For a particular determination, you require 10 mL of sample. Now and again, however, you have too little sample available and you must perform the determination with 5 mL, for instance. In this case, proceed as follows:

1. Activate the parameter **Fixed sample size**.
2. Define the value and unit for the sample size, in the above example, for instance, **[10 mL]**.
3. Deactivate the parameter **Fixed sample size** again.
4. Save method.

⇒ Each time you load this method, 10 mL is entered as the sample size, but in the example mentioned you can enter the lesser quantity of 5 mL.

[Sample size limits]

This button is accessible only if **Fixed sample size** is deactivated.

Define the limit values for the sample size (see "Defining limit values for the sample size", page 195).

Modifying sample data

This instruction is for the purpose of clarifying the settings which are possible in this dialog. The following settings should be made:

- Change the designation for the input field of the first sample identification.
- Hide the input field for the second sample identification.
- Define a fixed sample size.

Proceed as follows:

1 Change the designation for the input field

- Tap on the input field **Identification 1**.
The text editor opens.
- Change the designation in **Batch** and confirm with **[OK]**.
- Deactivate the parameter **Identification 2**.

2 Define a fixed sample size

- Activate the parameter **Fixed sample size**.
- Enter the value **10 mL** as **Sample size**.

3 Displaying the main dialog

Tap on the fixed key [].

The data is saved and the adjusted main dialog appears.

New method

09:08:13

User

jb

Batch

Sample size

10

mL

System

Load method

Control/ Logout

Edit parameters

Results

Only Identification 1 is shown with the title "Batch." The fixed sample size "10 mL" is shown, but cannot be edited.

Defining limit values for the sample size

In the dialog **Sample data / Sample size limits**, you can define the limit values for the sample size. These values are not monitored during sample data input, but rather:

- when the determination is started.
- when the determination is finished.
- upon automatic sample data query via the command **REQUEST**.
- when the determination is being recalculated.

An entry is made in the message list automatically when these limits are infringed (see chapter 24.2, page 234).

Monitoring sample size limits

on | off (Default value: **off**)

If this parameter is activated, then the sample size will be monitored.

Lower limit

Input range	–999999999 to 999999999
Default value	0

Upper limit

Input range	–999999999 to 999999999
Default value	999999999

Display message

on | off (Default value: **off**)

If this parameter is activated, a message is displayed in case the limits are infringed. You can select whether you want to continue with the determination or cancel the run.

If this parameter is deactivated, then the message will be entered only in the message list of the determination.

Timeout

The message is displayed during this time. Afterwards, the determination is continued automatically.

Input range	0 to 999999 s
Default value	30 s
Selection	off

off

The determination will not be continued until after the message has been confirmed.

- **Reviewed on**
Date and time at which the method has been signed at level 1.
- **Reason**
Reason for the signing of the method.
- **Comment**
Comment for the signing of the method.

The following information is only shown if the method has been signed at level 2:

- **Released by**
User who signed the method at level 2.
- **Released on**
Date and time at which the method has been signed at level 2.
- **Reason**
Reason for the signing of the method.
- **Comment**
Comment for the signing of the method.

[Delete signature]

Delete all of the signatures for the method or the determination. This button is inactive if level 2 has not yet been signed (*see chapter 17.2, page 201*).

[Sign]

Sign the method or determination. This button is active only when working with activated login function and password (*see chapter 17.1, page 200*).

16.5.5 Note

Main dialog: **Edit parameters ► Method options ► Note**

In this dialog, you can enter a short text, e.g. important information about carrying out the determination.

[Display options]

Definition when the note is being displayed.

Automatically after loading the method

on | off (Default value: **off**)

If this parameter is activated, the note is displayed when loading the method. It can otherwise only be read in this dialog.

16.5.6 Saving a determination automatically

Main dialog: **Edit parameters** ► **Method options** ► **Save automat.**

In the dialog **Method options / Save automatically**, you can specify whether the determination is to be saved automatically and/or whether a PC/LIMS report is to be created.

Save determination automatically

on | off (Default value: **off**)

If this parameter is activated, then the determination will be saved automatically.

Memory

Memory location where the determination will be saved.

Selection	External memory 1 External memory 2
	Shared memory
Default value	External memory 1

Shared memory

The determination will be saved in a shared directory on the network. The shared directory is selected in the **Edit device / Shared memory** dialog (see chapter 11.3.3, page 102).

Group

Directory where the determination will be saved.

Entry	32 characters maximum
Selection	Selection of available directories empty
Default value	empty

File name

File name of the determination. The file name is always extended to include the date and the time of day (YYYYMMDD-hhmmss) in order to ensure that the file names are unique for all determinations.

Entry	max. 16 characters
Selection	Identification 1 Identification 2 Method
Default value	Identification 1

Identification 1

The first 16 characters of the text which was entered in the main dialog for Identification 1 + YYYYMMDD-hhmmss.

Identification 2

The first 16 characters of the text which was entered in the main dialog for Identification 2 + YYYYMMDD-hhmmss.

Method

The first 16 characters of the text, of the method name + YYYYMMDD-hhmmss.

Write protection

on | off (Default value: **on**)

If this parameter is activated, then the file cannot be saved, deleted, or renamed. This is only an internal write protection feature and is independent of the write protection of the operating system on your computer. This write protection feature protects the saved determination data against accidental modifications or modifications by unauthorized persons.

Create PC/LIMS report

on | off (Default value: **off**)

If this parameter is activated, then a machine-readable report with all of the important data for the determination is created, a so-called PC/LIMS report.

The settings are implemented in the device manager (*see chapter 11.3.2, page 100*).

When you save the PC/LIMS report as a file, the file name will be generated automatically as follows: **PC_LIMS_Report**-first 16 characters of the value of the parameter **File name-Date-Time.txt**. If the value of the parameter **File name** contains one or more of the following characters, then these will be automatically converted in the file name into the character "_": / \ : * ? " < > |.

Selection	Selection of configured reasons
-----------	---------------------------------

Entry	24 characters maximum
-------	-----------------------


Sign the method or determination with the data entered above.

Selection	Selection of configured reasons
-----------	---------------------------------

Entry	24 characters maximum
-------	-----------------------

18 Control

Main dialog: **Control**



In the dialog **Control**, the settings for the execution of a single determination or of one sample series are defined.

If you work with the login function activated, you can log off the system in this dialog. The login dialog will then be shown immediately.

Statistics

on | off (Default value: **off**)

With this option, you can enable or disable the statistics calculations for individual determinations (see *"Statistics"*, page 204).



NOTE

The parameter **Statistics** in the method options will not be changed as a result. This setting here will be adjusted automatically to match the one in the method options.

Sample table

on | off (Default value: **off**)

If this option is activated, the sample data for a sample series can be entered in a table (see *chapter 21*, page 213).

Autostart

on | off (Default value: **off**)

If this option is activated, a new determination is started automatically at the end of a determination. This continues until the number specified has been reached (siehe **Number of autostarts**).



NOTE

If you deactivate this option while a determination is running, then that determination will be ended, although the next one will however also not be started. This can be used to cancel a sample series.

Number of autostarts

Number of automatic starts.

Input range	1 to 9999
Selection	Table
Default value	Table

Table

The number of automatic starts corresponds to the number of samples in the sample table.

Sample number

The sample number is primarily used for documentation. Each time a determination is started, it is increased by one. The sample number can, however, also be entered manually. Each time the Ti-Touch is switched on, the value is reset to zero.

Input range	0 to 9999
Default value	0

[Logout]

User logout (visible only when work is done with activated login function).

[Favorites]

Save methods/sample tables as a favorite (*see chapter 19, page 206*). Depending on whether you are working with or without login function, user-specific or common favorites can be created. On the main dialog, a button is created for each favorite. This way, a method or a sample table can be started with one single mouse click.

[Delete table]

Delete the sample table completely.

[Reset autostart]

Reset the value of the autostart counter to zero.

Continuing the sample series

19.1 Creating favorites

Control / Favorites	
Name	Position
PREP	1
EMPTY	6
Me4155	5
Me3901	10
New Delete Edit	

The list of favorites gives the designation and position of the button on the main dialog of each favorite configured.

[New]

Add a new favorite to the list, see following chapter.

[Delete]

Delete the selected favorite from the list.

[Edit]

Edit the properties of the selected favorite, see following chapter.

19.1.1 Editing favorites

Favorites / Edit	
Position	5 ▼
Name	Me4155
Type	Method ▼
Memory	Internal memory ▼
Method / Sample table	Me4155 ▼
More options	

on | off (Default value: **off**)

Number of autostarts

Number of automatic starts.

Input range	1 to 9999
Default value	1
Selection	Table

The number of automatic starts corresponds to the number of samples in the sample table.

Sample identification. The sample identification can be used in calculations as the variable **C12**. If you have defined a sample identification list (see *chapter 15.1.1, page 164*) then the entries can be selected here.

Sample size. The value of the sample size can be used in calculations as the variable **C00**.



Me2115



01 REQUEST Data request

Identification 1

Sample size

10

mL

Continue

To accomplish this, you must insert and configure the command **REQUEST** in the method (see chapter 29.15.3, page 434).

If the **Hold sequence** parameter is activated, then the run will be paused and must be continued with **[Continue]** after the sample data has been entered. If the **Hold sequence** parameter is deactivated, then the determination will continue in the background until the measured data is required for further processing. This dialog will be displayed until the entering of the sample data is confirmed with **[Continue]**, even if the titration/measurement is already completed. This ensures that the sample data is available for calculations.

If a fixed sample size has been defined in the method (see *chapter 16.5.3, page 192*), then it will be displayed at the time of the entry but it will not be editable.

21 Sample table

21.1 General

Main dialog: **Sample table**

The sample table is a table in which the sample data for up to 999 samples can be entered. If you are processing larger-sized sample series, it makes sense to enter the sample data in this table.


Activating the sample table

Proceed as follows to enable the sample table:

1 Opening the control dialog

- In the main dialog, tap on **[Control]**.


The control dialog is displayed:



2 Activating the sample table

- Activate the **Sample table** option.
- Tap on the fixed key **[]**.

The main dialog is displayed:



The display will show the number of determinations already carried out and the total number of sample lines containing data. The sample table is still empty in this example.

The sample table contains numbered lines; each line represents one determination. In the standard settings, the first sample identification and the sample size of each sample are displayed. In the properties, you define the columns to be displayed (*see chapter 21.3, page 218*).

Sample table		
No.	Identification 1	Sample size
1	#2370015	1.0 g
2	#2370015	1.0 g
3	#2370016	1.0 g
4	#2370016	1.0 g
5	...	

Load/Save

Properties

Insert line

Delete

Edit

[Load/ Save]

Load a saved sample table (see "Loading a sample table", page 216) or save the current sample table (see "Saving a sample table", page 215).

[Properties]

Edit the properties of the sample table (see chapter 21.3, page 218).

[Insert line]

Insert a new line above the line selected.

[Delete]

Delete the selected line from the sample table.

[Edit]

Edit the data of the selected line (see chapter 21.2, page 216).

Saving a sample table



NOTE

Sample tables can only be saved on external storage media or on a shared file system.

Proceed as follows to save a sample table:

1 Opening the dialog

- In the **Sample table** dialog, tap on the **[Load/ Save]** button.

The **Sample table / Load/Save** dialog is displayed.

If no external storage medium is plugged in or if no file system is shared, then **[Save]** is disabled.

2 Modify/apply a name

- Tap on **[Save]**.
A name is suggested for new sample tables. If the sample table has already been saved, the previous name will be displayed.

Applying the name:

- Tap on **[Save]**.

The sample table is saved.

Entering a new name:

- Tap on the input field **File name**.
The text editor opens.
- Enter a new file name (max. 32 characters) and confirm with **[OK]**.
- Tap on **[Save]**.

The sample table is saved.

At the very top you will see the line number of the selected line. In this example, the sample data of the first line is displayed. You can scroll between individual data sets with the keys **[-]** and **[+]**.

Method

Method used for processing the sample.

If you use the sample assignment table (*see chapter 15.1.2, page 165*), then the method can no longer be defined here.

Selection	Selection of stored methods empty
Default value	empty

empty

The currently loaded method is used.

Identification 1

Sample identification. The sample identification can be used in calculations as the variable **CI1**. If you have defined a sample identification list (*see chapter 15.1.1, page 164*) then the entries can be selected here.

Entry	24 characters maximum
Default value	empty

Identification 2

Sample identification. The sample identification can be used in calculations as the variable **CI2**. If you have defined a sample identification list (*see chapter 15.1.1, page 164*) then the entries can be selected here.

Entry	24 characters maximum
Default value	empty

Sample size

Sample size. The value of the sample size can be used in calculations as the variable **C00**.

Input range	-99999999 to 999999999
Default value	1.0

Unit of sample size.

Entry	6 characters maximum
Selection	g mg µg mL µL pieces
Default value	g

Unit defined in method:

The unit defined in the method is ignored. The unit entered in the sample table is used.

Method

Only the assigned method is displayed for each sample.

Identification 1

The first sample identification and the sample size are displayed for each sample.

Identification 2

The second sample identification and the sample size are displayed for each sample.

Number of lines

Maximum number of lines which can be present in the sample table.

Input range	2 to 999
Default value	99

Lock editing

on | off (Default value: **off**)

If this option is activated, then the sample data can no longer be altered.

Display options

Configure the sample data to be displayed in the dialog **Sample table / Edit** (see *"Properties / Display options" dialog*, page 219).

Editing options

Specify the settings for the processing of the sample data (see *"Properties / Editing options" dialog*, page 219).

Identifications

Define the designations for the input fields of Identification 1 and Identification 2 (see *"Dialog "Properties / Identifications"*, page 220).

"Properties / Display options" dialog

Sample table: **Properties ► Display options**

In this dialog, you can configure which of the following sample data can be edited in the sample table:

- **Method**
- **Identification 1**
- **Identification 2**
- **Sample size**

"Properties / Editing options" dialog

Sample table: **Properties ► Editing options**

In this dialog, you can define the settings so that the input of the sample data becomes simpler and faster.



Properties / Editing options

End line with

Sample size

▼

Copy automatically into next line:

☒ Method
 ☒ Identification 2

☒ Identification 1
 ☒ Sample size

End line with

Selection as to with which input the editing dialog for the next sample is automatically displayed.

Selection	Sample size Identification 1 Identification 2 manual
Default value	Sample size

manual

The editing dialog for the next sample can be displayed with the **[New sample]** button in the editing dialog.

If, for example, you have to enter similar data for each sample, then you can automatically copy the data which is identical for each sample into the next line:

- Method
- Identification 1
- Identification 2
- Sample size

Dialog "Properties / Identifications"

Sample table: **Properties** ► **Identifications**

In the dialog **Properties / Identifications**, you can alter the designations for the sample identifications in accordance with specific methods.

22 Determination run

This chapter provides a description of how to carry out a determination, i.e. how to process a single sample or a sample series. The current sample data is used for calculations.

22.1 Carrying out a single determination

When you carry out a determination, you can enter the sample data in various ways (see chapter 20, page 210). The following instructions describe how to enter the sample data in the main dialog.

To carry out a single determination, proceed as follows:

1 Loading the method

- See Chapter "Loading a method", Page 185.

2 Preparing the sample



NOTE

Calculate the amount of the sample in such a way that a titrant consumption of 10...90% of the cylinder volume will result.

- Weigh in or measure the sample in a sample vessel.
- Add solvent.
- Place the sample vessel on the stirrer.
- Immerse the electrode and buret tip in the solution.

3 Entering the sample size

- Tap on the input field **Sample size** in the dialog. The number editor opens.
- Enter the sample size and confirm with **[OK]**.

4 Starting a determination

- Tap on the fixed key [▶].

Sequence of the sample series

- After a determination is started, a check is made as to
 - whether all of the devices required are connected and ready for operation.
 - whether all titrants and sensors are present.
 -
 - whether the monitoring intervals for titrants, sensors and common variables are still valid.
 -
 - if a sample assignment table is activated (*see chapter 15.1.2, page 165*):
 - whether the entered sample identification appears in the sample assignment table. If this is the case, then the assigned method will be loaded automatically and the actual determination will be started. If a different method has been entered in the sample table, then it will be ignored.
- At the end of the determination,
 - the line in the sample table is deleted.
 - the next determination is started (if you are working with the autostart function enabled (*see chapter 18, page 202*)).

22.3 Canceling determinations manually

You can cancel a determination at any time with the fixed key []. The command being carried out at that moment is canceled and only calculations and reports are executed or created.

If you are working with the autostart function activated (*see chapter 18, page 202*), the entire series will be canceled.



NOTE

If you do not wish to cancel the entire determination, but only the current command, use the function **[Cancel command]** in the live display. The subsequent method commands are carried out normally.

23.2 Editing the sample table while a determination is running

You can insert new lines or delete existing ones or edit sample data while a determination is running.



NOTE

To ensure that no problems arise during the run and that the current data is always available for the calculation, we recommend that the input dialogs always be closed.

Editing the sample table

Proceed as follows to edit the sample table:

1 Display the main dialog

- Tap on the fixed key [].

The main dialog is displayed. The determination continues to run in the background.

2 Select the sample data

- Tap on **Sample table**.
The sample table is displayed. The first line always contains the sample data of the running determination.
- Select the desired line.
- Tap on the button **[Edit]**.

3 Edit the sample data

- Edit the sample data and apply with **[OK]**.



NOTE

In addition to the sample data, the method can also be modified, except in cases where the determination is running.


4 Display the live dialog

- Tap on the fixed key [].
The main dialog is displayed.
- Tap on **[Live display]**.

The live dialog is displayed once again.

23.3 Live display

Main dialog: **Live display**



The current method and the system status are displayed in the title line. The currently running command is displayed directly underneath the title line.

The following functions are available to you while a determination is running:

- **[Hold]**
Pause the determination. The label changes to **[Continue]**. This can be used to continue the run.
- **[Cancel command]**
Cancel the current method command. The next command is started.
- **[View]** (only for titrations/measurements)
To modify the settings for the curve and measured value display, see below.
- **[Sample data]** (visible only if a titration or a measurement is not currently running)
Switch to the main dialog in order to edit the sample data.
- **[Edit parameters]**
Edit the method options and live parameters (*see chapter 23.4, page 228*).
- **[Results]**
Display the results dialog (*see chapter 24, page 230*).

Dialog "View"

Main dialog: **Live display ► View**



Selection	Curve Curve & meas. values Measured values Curve & sample data
Default value	Curve & meas. values

Only the curve is displayed.

The measured values are displayed to the right next to the curve.

Only the measured values are displayed.

The sample data is displayed to the right next to the curve.

Define settings for the display of the measured values.

A maximum of three different values can be displayed at the same time. The selection depends on the titration/measuring mode. The settings are saved separately for each mode.

Define the settings for the curve display (see "Curve options", page 244). The settings are saved separately for each titration/measuring mode and apply not only for the live display but also for the curve display in the results dialog.

- Tap on the fixed key [].
The main dialog is displayed.
- Tap on **[Live display]**.

- Tap on **[Live display]**.

24 Results and more determination data

Main dialog: **Results**

Detailed specifications concerning the current determination are displayed in the **Results** dialog:

The screenshot shows a software window titled "Results". It contains the following information:

- Current method: Me7625
- Determination: Sample 43C-20100903-155708
- Determ. time: 2011-09-03 15:57:08
- Total acid: 0.53 g/L
- 01 MET pH Monotonic pH titration
 - EP1 pH 5.383 2.2667 mL
- 02 EVAL FIX-EP Evaluation fixed endpoints
 - FP1 pH 5.500 2.2799 mL

At the bottom, there are five buttons: "Result table", "Statistics", "Recalculate", "Curve", and "More data".

- **Current method**
Method with which the determination was carried out.
- **Determination**
Name of the determination. Displayed only if the determination has been saved.
- **Determination time**
Date and time at which the determination was carried out.
- **Results**
Calculated results. If monitoring of the result limits is activated, then the result will be highlighted in color:
 - Green lettering: The result is within the limits
 - Red lettering: The result is outside the limits
- **Commands**
List of the data-generating commands (titrations, measurements, etc.). The most important data is displayed for each command.

[Result table]

Open result table (*see chapter 26, page 252*). The result table is a table in which the results for several determinations can be saved.

[Statistics]

Display the statistical overview of a determination series (*see chapter 25, page 247*).

[Recalculate]

Recalculate the determination. The procedure will be executed immediately. Recalculation cannot be undone (*see chapter 24.8, page 245*).

[Sign]

Sign the current determination.

If the method contains a command **SIGN**, then this button will be displayed instead of **[Recalculate]**. The determination can be signed directly on level 1 (with a second command also on level 2).

[Curve]


Display the curve of the current determination (*see chapter 24.7, page 243*).

[More data]

Display detailed specifications for the determination (determination properties, measuring point list, endpoint details, etc.), save determination (*see chapter 24.1, page 231*).

24.1 More determination data

Results dialog: **More data**



All of the commands for which you can display other data (titrations, measurements, evaluations, etc.) are listed in the **Results / More determination data** dialog.

[Load/ Save]

Load saved determinations (*see chapter 24.5, page 238*) and save current determination (*see chapter 24.6, page 243*).

[Properties]

Display the properties of the current determination (*see chapter 24.4, page 235*).

[Local com. var.]

(available only with loaded determinations)

Display the list of the common variables defined at the time of the determination (see *chapter 24.3, page 235*).

[Messages]

Display list with all of the messages which occurred during the determination (see chapter 24.2, page 234).

[\[View data\]](#)

Display all of the details for the selected command.

24.1.1 Details

Results dialog: **More data** ► **View data**

More determination data / View data		
01 MET pH Monotonic pH titration		1M
Stop criterion	Stop volume reached	MTS
Titration duration	48.9 s	MCD
Start volume	1.0000 mL	MSV
End volume	3.0000 mL	MCV
Temp. measurement	Pt1000	MTM
Initial measured value	3.507	MIM
Initial temperature	27.8 °C	MIT
Start duration	6.5 s	MSD
Start measured value	3.682	MSM
Start temperature	27.8 °C	MST
End measured value	6.812	MCM
Endpoint details		Measuring point list

All of the variables generated by the determination are listed in the **More determination data / View data** dialog. These can be used in calculations.

[Calibration data]

(available only with loaded determinations)

Display calibration data of the sensor used which was valid at the time of the determination (see "Calibration and titrant data", page 234).

[Titrant data]

(available only with loaded determinations)

Display titrant data which was valid at the time of the determination (see "Calibration and titrant data", page 234).

[Endpoint details]

Display details for each endpoint found (see "Endpoint details", page 233).

[Measuring point list]

Display measuring point list of titrations and measurements (see "Measuring point list", page 233).

Measuring point list

Results dialog: **More data ► View data ► Measuring point list**

Results dialog: **Curve ► Display curve ► Measuring point list**

View data / Measuring point list			
01 MET pH Monotonic pH titration			1M
Volume [mL]	Meas.value [pH]	Delta MV [pH]	
1.00000	3.682	0.000	
1.10000	3.713	0.032	
1.20000	3.749	0.036	
1.30000	3.790	0.040	
1.40000	3.836	0.046	
1.50000	3.889	0.053	
1.60000	3.950	0.061	
			View

The individual measured values are listed in the **View data / Measuring point list** dialog. You can configure which columns are displayed (**[View]** button). These settings are saved separately for each mode (DET, MET, etc.).

Column 1 / Column 2 / Column 3

Quantity displayed in the columns 1, 2 and 3. Both the selection of the quantities and the default setting depend on the mode.

Selection	Volume Measured value ERC Delta meas.value Time Temperature Volume drift Signal drift Index
-----------	--

Endpoint details

Results dialog: **More data ► View data ► Endpoint details**

Results dialog: **Curve ► Display curve ► Endpoint details**

Detailed information concerning each equivalence point or endpoint determined is displayed in the **View data / Endpoint details** dialog.

The volume, the measured value, the equivalence point criterion ERC (only DET, MET), the time and the temperature are specified for each equivalence point or endpoint determined.

The equivalence point criterion determined (ERC = Equivalence point Recognition Criterion) is a measure for the size of the jump in the titration

curve. The ERC determined is compared with the EP criterion defined in the method for the evaluation of the titration curve *see "Evaluation and equivalence point criterion with DET", page 304 and "Evaluation and equivalence point criterion with MET", page 320*). The ERC is thus an important quantity when you need to adjust the parameters for the evaluation of the equivalence points.

Calibration and titrant data

Results dialog: **More data** ► **View data** ► **Calibration data**

Results dialog: **More data** ► **View data** ► **Titrant data**

The calibration data or titrant data which was current at the time the determination was performed is displayed in the dialogs **View data / Calibration data** and **View data / Titrant data**.



NOTE

Recalculation:

The concentration and titer of the titrant can be modified retroactively. The values displayed here are used for the recalculation of a loaded determination. The values saved in the dialog **System / Titrants** are not changed as a result.

24.2 Messages

Results dialog: **More data** ► **Messages**

If messages occurred during a determination, they will be displayed in the **More determination data / Messages** dialog. The following data for the message is indicated under **[Details]**:

- **Message number**
Unique identification number.
- **Time**
Moment that the message occurred.
- **Message**
Message text.

24.3 Local common variables

Results dialog: **More data ► Local com. var.**

The list of common variables as it existed at the time of the determination is displayed in the **More determination data / Local common variables** dialog.



NOTE

Recalculation:

The common variables can be modified retroactively. The values displayed here are used for the recalculation of a loaded determination. The list of common variables in the **System / Common variables** dialog is not changed as a result.

24.4 Determination properties

Results dialog: **More data ► Properties**

Detailed information concerning the determination is displayed in the **More determination data / Properties** dialog:

- **User**
User who carried out the determination.
- **Sample number**
Number which was entered in the dialog **Control** at the start of the determination.
- **Start counter**
This value, just like the sample number, is increased by one each time a determination starts. In contrast to the sample number, this value cannot be modified by the user.
- **Duration determ.**
Duration of the entire determination.

- **Method status**

Current status of the method.

- **new**
The method has been newly created and not yet saved.
- **saved**
The current method has been saved.
- **modified**
The current method has been modified but not saved yet (e.g. live modification).
- **reviewed**
The current method has been signed at level 1.
- **released**
The current method has been signed at level 2.

- **Determination**

Name of the determination. Displayed only with saved determinations.

- End of determin.

Way in which the determination was ended.

- **Regular without errors**
The determination was automatically ended at the end of the method.
- **Manual stop**
The determination has been canceled with the fixed key [☐].
- **Error**
The determination has been canceled due to an error.

- **Status of deter.**

- **original**
The determination has not been recalculated.
- **recalculated**
The determination has been recalculated.
- **loaded**
The determination has been loaded.
- **saved**
The determination has been saved. The version number indicates how many times the file has been saved.
- **reviewed**
The determination has been signed at level 1.
- **released**
The determination has been signed at level 2.

The following data is displayed only if the option **Recalculating determination** under **Login options / Modification options** is activated:

- **Recalculated by**
User who recalculated the determination.
- **Recalculated on**
Date and time at which the determination was recalculated.

- **Reason**
Reason for the recalculation of the determination.
- **Comment**
Comment for the recalculation of the determination.

The following information is displayed only if the determination has been signed at level 1:

- **Reviewed by**
User who signed the determination at level 1.
- **Reviewed on**
Date and time at which the determination was signed at level 1.
- **Reason**
Reason for the signing of the determination.
- **Comment**
Comment for the signing of the determination.

The following information is displayed only if the determination has been signed at level 2:

- **Released by**
User who signed the determination at level 2.
- **Released on**
Date and time at which the determination was signed at level 2.
- **Reason**
Reason for the signing of the determination.
- **Comment**
Comment for the signing of the determination.
- **Used devices**
Hardware that was used for the determination.
 - Ti-Touch with program version of the software and the serial number.
 - Measuring inputs with ADC type und serial number.
 - MSB connector with peripheral devices connected to it (dosing device, stirrer).

[Delete signature]

Delete all of the signatures for the method or the determination. This button is inactive if level 2 has not yet been signed (*see chapter 17.2, page 201*).

[Sign]

Sign the method or determination. This button is active only when working with activated login function and password (*see chapter 17.1, page 200*).

**NOTE**

When you load a determination, the method with which the determination was carried out will be loaded automatically.

For this reason, save any modifications which may have been made to the method currently loaded before you load a determination.

24.5.1 Determination list

Results dialog: **More data ► Load/ Save ► Load ► Show files ► View**

Show files / View			
Sorted by Determ. time		10	
No.	Identification 1	Result 1	
1	Blank methanol	0.01 %	
2	Blank methanol	0.02 %	
3	Calcium hardness	2.015 mmol/L	
4	Calcium hardness	1.981 mmol/L	
5	Titer of NaOH	0.99615	
6	Titer of NaOH	0.99417	
7	Titer of NaOH	0.98729	
Previous 100			
Next 100			
Properties			
Delete			
Load			

In the **Show files / View** dialog you can depict the individual determinations in greater detail. The sorting criterion and the number of saved determinations are displayed at the very top. The list contains numbered lines; each line represents one determination. In the default settings, the first sample identification and the first calculated result of each determination are displayed. In the properties, you define which data is displayed in the columns and the criteria according to which the determinations are sorted (*see "Properties", page 240*). In addition, you can define search filters to aid you in finding the desired determinations (*see "Filter criteria", page 241*).

[Previous 100]

Display the previous 100 determinations.

[Next 100]

Display the next 100 determinations.

[Properties]

Configure the display of the columns and define filter criteria.

[Delete]

Delete the selected file.

[Load]

Load the selected file.

Properties

Determination list: **Properties**

View / Properties

Column 1

Identification 1

▼

Column 2

Result 1

▼

Sort by

Determ. time

▼

Sort sequence

descending

▼

Filter

In the **View / Properties** dialog, you can define which data is displayed in the columns and the criteria according to which the determinations are to be sorted.

Column 1

Selection of what is to be displayed in the first column.

Selection	Identification 1 Identification 2 Determ. time Result name Method Sample size User
Default value	Identification 1

Column 2

Selection of what is to be displayed in the second column.

Selection	Result 1...9 The numbering of the results corresponds to the order in which the results were calculated in the determination run.
-----------	--

Sort by

Selection according to which criterion the list is to be sorted.

Selection	Identification 1 Identification 2 Determ. time Result Method Sample size User
Default value	Determ. time

Identification 1

Sort alphabetically according to the first sample identification.

Identification 2

Sort alphabetically according to the second sample identification.

Determ. time

Sort chronologically according to the determination time.

Result

Sort according to the numerical value of the result.

Method

Sort alphabetically according to the method name.

Sample size

Sort according to the numerical value of the sample size.

User

Sort alphabetically according to the user who carried out the determination.

Sort sequence

Selection	ascending descending
Default value	descending

[Filter]

Specify the filter criteria for the list.

Filter criteria

Determination list: **Properties** ► **Filter**

Properties / Filter

Filter

Date

is equal to

Date from

2011-07-01

to

2011-07-08

Value from

to

You can set filters in the **Properties / Filter** dialog. A filter is a rule defining what is to be displayed in the determination list. If a filter has been set,

this is indicated with a corresponding symbol in the **Show files / View** dialog.

Filter

Define the criterion according to which the list is to be filtered. The only lines which will be displayed are those which fulfil the filter criterion.

Selection	no filter Identification 1 Identification 2 Date Result Method Sample size User
Default value	no filter

no filter

The list itself will not be filtered.

Identification 1

In the input field **is equal to**, select or enter the first sample identification according to which filtering is to be carried out.

Identification 2

In the input field **is equal to**, select or enter the second sample identification according to which filtering is to be carried out.

Date _____

In the input fields **Date from ... to**, specify the period of time according to which filtering is to be carried out.

Result

In the input fields **Value from ... to**, specify the value range according to which filtering is to be carried out.

Method

In the field **is equal to**, select or enter the method name according to which filtering is to be carried out.

Sample size

In the input fields **Value from ... to**, specify the value range according to which filtering is to be carried out.

User

In the field **is equal to**, select or enter the user according to whom filtering is to be carried out.



NOTE

The filter can be applied to empty fields as well, for example, if no sample identification has been entered. In this case, select only the desired filter criterion and leave the fields below empty.

24.6 Saving a determination

Results dialog: **More data ► Load/ Save ► Save**

Determinations can be saved as follows:


- Manually in this dialog.
The suggested default setting for file names is the first 16 characters of Identification 1, followed by the determination time in the format YYYYMMDD-hhmmss. Extending the file name with the determination time has the advantage that the file names of all determinations are unique.
- Automatically at the end of the determination (*see chapter 16.5.6, page 198*).

The following data is stored for each determination:

- Calculated results
- Measuring point lists, equivalence points, curves and variables which were generated during the determination
- Determination properties
- Statistics (if carried out)
- Method with which the determination was carried out
- Common variables (values and properties of all common variables at the time of the determination)
- Titrant data of the titrant used
- Calibration data of the sensor used

24.7 Curves

Results dialog: **Curve**



The curve of the current determination is displayed in the **Results / Curve** dialog. You can configure the curve display.

[Endpoint details]

Display details for each endpoint found (see "Endpoint details", page 233).

[Measuring point list]

Display measuring point list of titrations and measurements (see "*Measuring point list*", page 233).

[Curve options]

Define the settings for the curve display (see "Curve options", page 244). The settings are saved separately for each titration/measuring mode and apply not only for the live display but also for the curve display in the results dialog.

Curve options

Results dialog: **Curve ► Curve options**

Curve / Options

Display of the curve for mode MET

x axis

Volume

▼

y1 axis

Measured value

▼

y2 axis

none

▼

Color

■

▼

Color

■

▼

☒ Grid

☐ Display measuring points

You can make the settings for the curve display in the **Curve / Options** dialog. These settings are saved separately for each mode (DET, MET, etc.) and also apply for the live display. Both the selection of the quantities and the default setting depend on the mode.

x axis

Quantity which is shown on the x axis of the graph.

Selection	Volume	Measured value	ERC	Delta meas.value	Time	Temperature	Volume drift	Signal drift
-----------	--------	----------------	-----	------------------	------	-------------	--------------	--------------

y1 axis

Quantity which is shown on the y1 axis of the graph.

y2 axis

Selection	none Volume Measured value ERC Delta meas.value Time Temperature Volume drift Signal drift
Default value	none

Selection of a color for the curves.

on | off (Default value: **on**)

Display measuring points

on | off (Default value: **off**)

If this parameter is activated, the curve is shown with the single measured values.

Results dialog: **Recalculate**

You can use the **[Recalculate]** function to recalculate and reevaluate determinations. This is necessary in the following cases:

- If you have modified the evaluation parameters in a way that changes the recognized endpoints or equivalence points.
- If you have modified the calculations.
- If you have modified the variables used in calculations, e.g. sample size, titer or common variables.



Neither the recalculation nor the reevaluation can be undone.

We therefore recommend that the original determination be saved beforehand (*see chapter 24.6, page 243*). The best way of doing this is to use the **Save determination automatically** function (*see chapter 16.5.6, page 198*).

The following method commands are executed once again with the **[Recalculate]** function:

- DET and MET titrations:
The potentiometric evaluation will be carried out again.
- EVAL commands:
All evaluations are carried out again. You can also insert new evaluations retroactively.
- CALC command:
All calculations are carried out again. You can also define new calculations retroactively.
If you subsequently modify the titrant in a method, then the variables TITER and CONC of the "new" titrant will be used. If a result is assigned to the variable TITER in a calculation, then the titer of the "new" titrant will also be overwritten.



NOTE

If you modify the calibration data of the sensor used in the **System / Sensors** dialog, these modifications will not be taken into account. Changing the sensor in the method also has no influence on the measured data.

If a determination has been recalculated, this will be documented in the determination properties (see chapter 24.4, page 235).

Recalculation/reevaluation of loaded determinations

With loaded determinations, it is assumed that they were carried out a long time ago and that the **titrant data** and **common variables** of the system no longer correspond to those of the determination. This is why the most important titrant data of the titrants used and the list of **common variables** are stored with each determination. This data will be used if you recalculate a loaded determination. If you wish to recalculate the result with a corrected **common variable** or with a corrected titer, you must therefore modify the respective values in the **More determination data / Local common variables** or **More determination data / View data** dialog.

If a result is assigned to the variable TITER or to a **common variable**, then the titer or **common variable** which is currently valid in the system will not be overwritten until after a query has been confirmed.

If you modify the titrant retroactively in the method, this will not be taken into account.

The recalculated results of loaded determinations are not saved in the result table.

25 Statistics

Main dialog: **Results ► Statistics**


In the dialog **Results**, you can use **[Statistics]** to display the statistical overview of a determination overview. You can statistically evaluate a maximum of nine results calculated in a determination. A statistical series can contain a maximum of 20 determinations.

Activating the statistics function

Proceed as follows to activate the statistics function:

1 Editing the method options

- In the main dialog, tap on **[Edit parameters]**.
- Tap on **[Method options]**.



- Activate the **Statistics** option.
- Under **Number of determinations**, enter the desired number of determinations.

The results for which statistical calculations are performed are defined in the calculation command **CALC** (see "Variable for mean value", page 169).

25.1 Displaying details for a result

Results dialog: **Statistics ► Details**

Statistics / Details		
Result name: Content (mmol/L)		
Mean value	100.31 mmol/L	SMN1 n=03
s abs	0.236 mmol/L	
s rel	0.24 %	
No.	Sample size	Result
1	5.0 mL	100.56 mmol/L
2	4.5 mL	100.09 mmol/L
3	5.5 mL	100.28 mmol/L
Sample data Determ. on/off Result on/off		

In the dialog **Statistics / Details**, the following details are displayed:

- Result name
- Mean value
On the right-hand side of the dialog the number of single results from which the mean value was calculated is displayed. In this example, it is 3.
- Absolute standard deviation (**s abs**)
- Relative standard deviation (**s rel**)
- Result and sample size of each single determination
Results that could not be calculated are indicated with **invalid** and ignored.

[Sample data]

Display the sample data of the selected determination.

[Determ. on/off]

Remove all of the results of the selected determination from the statistics. All of the entries of this line are marked with an asterisk (*), the statistics will be recalculated automatically.

[Result on/off]

Remove the selected result from the statistics. The result is marked with an asterisk (*), the statistics will be recalculated automatically.



NOTE

If the mean value has been assigned to a common variable or to the variable TITER (see command **CALC**), then the determination must be recalculated manually (*see chapter 24.8, page 245*).

25.2 Deleting statistical data

In the dialog **Control**, you can delete all of the statistics data manually with the function **[Delete statistics]**. The statistics data is deleted automatically in the following cases:

- when all of the determinations of the determination series have been carried out and a new determination has been started afterwards.
 - when a new method is loaded (even if it is the same method)
- If the **sample table** or the **sample assignment table** is used, the statistics data is only deleted when another method is loaded.

25.3 Adding a determination to a determination series

You can use the function **[Add determ.]** to add an additional sample to a determination series, e.g. because a determination was faulty and had to be removed from the statistics.

Results / Statistics	
Current method: Me3901 3 of 3 (+2)	
Result name	Mean value
Content (mmol/L)	100.31 mmol/L
Volume of HCl	5.0 mL

The number of determinations by which the set statistics counter has been increased is shown above the list in brackets. This modification is also displayed in the main dialog and in the statistics report.



The number of samples defined in the method options does not, however, change as a result.

[Delete]

Delete the selected line from the result table.

[Details]

Display details concerning the selected determination.


"Result table / Details" dialog

The following determination data is displayed in this dialog:

- **User**
User who carried out the determination.
- **Method**
Method with which the determination was carried out.
- **Determ. time**
Date and time at which the determination was carried out.
- Sample data (Identification 1 and 2 and sample size).
- All of the results of this determination saved in the result table.

26.1 Properties

Result table: **Properties**



In the **Result table / Properties** dialog, you can select the data to be displayed with each determination as well as the sorting criterion for the list.

Column 1

Selection of what is to be displayed in the first column.

Selection	Identification 1 Identification 2 Determ. time Result name Method Sample size User
Default value	Identification 1

Column 2

Selection of what is to be displayed in the second column.

Selection	Result 1...9 The numbering of the results corresponds to the order in which the results were calculated in the determination run.
-----------	--

Sort by

Selection according to which criterion the list is to be sorted.

Selection	Identification 1 Identification 2 Determ.
	time Result Method Sample size User
Default value	Determ. time

Identification 1

Sort alphabetically according to the first sample identification.

Identification 2

Sort alphabetically according to the second sample identification.

Determ. time

Sort chronologically according to the determination time.

Result

Sort according to the numerical value of the result.

Method

Sort alphabetically according to the method name.

Sample size

Sort according to the numerical value of the sample size.

User

Sort alphabetically according to the user who carried out the determination.

Sort sequence

Selection	ascending descending
Default value	descending

Max. number of lines

If the result table contains the quantity of lines defined here, then no further results can be saved. If this is the case, save the result table and create a new one.

Input range	10 to 200
Default value	100

[Filter]

Specify the filter criteria for the list.

Save result table

on | off (Default value: **on**)

If this option is activated, the result table will be deleted as soon as you have saved it.

26.2 Saving the result table



NOTE

When you save the result table, the complete determination will not be saved with each entry. If you wish to save the determination data for each determination, use the **Save determination automatically** function (see chapter 16.5.6, page 198).

26.3 Loading the result table

Proceed as follows to load a result table:

1 Opening the current result table

- In the results dialog, tap on **[Result table]**.
- Tap on **[Load/ Save]**.
- Tap on **[Load]**.
The selection of memory locations is displayed. This selection will be skipped if only one external storage medium is available.
- Select the memory location where the desired result table is saved.
The selection of file groups is displayed (see chapter 12.1, page 137). If only one group is available, then the list of the saved result tables will be displayed directly.
- Select the group with the desired result table.
- Tap on **[Show files]**.

The list with the saved result tables opens.

2 Selecting the result table

- Select the desired result table.

3 Loading the result table


- Tap on **[Load]**.

The result table is now loaded and will be displayed immediately.

More reports

Fixed key : **More reports**

The dialog **Print / More reports** displays a list of all of the available reports.



NOTE

If you are using a **CUSTOM Neos** printer, you can print only the reports marked with a #. The designation on the reports is only visible once the **CUSTOM NEOS** printer is selected (see screenshot above).

[System reports]

Opening the list of reports for the system settings.



NOTE

If you are using a **CUSTOM NEOS** printer, you cannot print system reports.

[Method reports]

Opening the list of reports for the current method.



NOTE

If you are using a **CUSTOM NEOS** printer, you cannot print method reports.

[Send]

Create a PC/LIMS report.

[Edit]

Edit the settings of the selected report.

[Print]

Print the selected report.

27.1 General report options

Fixed key []: **Report options**

The settings for the report printout can be adjusted in the dialog **Print / Report options** to indicate, for example, whether a report header, a signature line or a frame is to be printed.

Report header

Defining the output of the report header. The report header contains general information such as device type, serial number, program version and the printing date.



NOTE

You can also define your own report header, which is printed in addition to this report header (**System ► Templates ► Report header**).

Selection	off once on each page
Default value	on each page

off

The report header will not be printed.

once

The report header will only be printed on the first page.

on each page

The report header will be printed on every page.

Signature line

Output of a special line for date and signature. This line will be printed at the very bottom of every page.

Selection	off once on each page
Default value	off

off

No signature line will be printed.

once

The signature line will only be printed on the last page.

on each page

The signature line will be printed on every page.

Frame

on | off (Default value: **on**)

If this parameter is activated, a frame is printed as lateral marking.

27.2 Settings of the individual reports

Fixed key []: **More reports ► Report ► Edit**

Command REPORT: **Report ► Edit**


Settings can be edited for the following reports (for details, see online help):

- **Result report**
- **Curve**
Definition of the curve size, display of the individual measuring points and grid lines, etc.
- **Measuring point list**
Definition of the method command to which the measuring point list is to be printed.
- **Calculations**
Definition of the accuracy with which the variables used as well as the results calculated are printed.
- **Statistics short**
Definition of the time at which the report is to be printed (with each determination or only at the end of a sample series).
- **Statistics overview**
Definition of the time at which the report is to be printed (with each determination or only at the end of a sample series).
- **Sample table**
- **Result table**
Definition of the time at which the report is to be printed (with each determination, at the end of a sample series or at the end of the sample table).
- **Used devices**
- **Variables**
- **Monitoring**

Report	Contents
PC/LIMS report	Machine-readable report with all of the important data for a determination. This report can be saved as a TXT or UTF-8 file or sent to a terminal program or a LIMS via an RS-232 interface. The definition is carried out in the device manager (<i>see chapter 11.3.2, page 100</i>).
Sample data	
Current sample data	Sample data of the current determination, <i>context-sensitive only</i> from the main dialog.
Sample table	List of all determinations in the sample table with the respective sample data, as entered in the sample table.
Result table	
Result table	List of all determinations in the result table with results and with the determination data, as saved in the result table.
More reports ► Method reports	
Method sequence	Method properties and the list of all method commands.
Parameters full	Method properties and options, all method commands with all parameters. All parameters which no longer have their default settings will be printed in bold . All parameters which have been modified in comparison to the stored version of the method are indicated by an *.
Titration & measuring param.	Method properties; titration, measurement and calibration commands with all parameters. All parameters which no longer have their default settings will be printed in bold . All parameters which have been modified in comparison to the stored version of the method are indicated by an *.
Modified parameters	Method properties, all method parameters which have been modified in comparison to the stored version of the method.
Non default parameters	Method properties, all method parameters which no longer have any default settings.
More reports ► System reports	
System settings	
System settings	Settings for acoustic signals and accuracy of the measured value display.

Report	Contents
Common variable list	List of all the common variables defined in the system, together with their most important data (name, value, status).
All common variable properties	Properties of all common variables (name, value, validity, status).
Common variable properties	Properties of the selected common variable (name, value, validity, status), <i>context-sensitive only</i> from the editing dialog.
Local Common variable list	List with the common variables available at the time of the determination, together with their most important data (name, value, status), <i>context-sensitive only</i> from the dialog "More determination data / Local common variables".
Local common var. properties	Properties of the selected common variable (name, value, validity, status), <i>context-sensitive only</i> from the dialog "More determination data / Local common variables".
Templates	
Templates sample data	Sample identification list and sample assignment table.
Result template list	List with all of the user-generated result templates.
All result templates details	Details of the all of the user-generated result templates (calculation formula, result options, note).
Result template details	Details of the selected result template, <i>context-sensitive only</i> from the editing dialog.
Input/Output lines	List with all of the defined input and output lines at the remote interface (name, bit pattern).
Custom calibration buffers	Temperature tables for all defined custom calibration buffers.
Rack tables	
Sample rack list	List with all of the sample racks configured in the system, together with designation, number of positions and rack code.

28 Manual control



The following functions are available in the manual control:

[Measure]

Carrying out manual measurements (*see chapter 28.2, page 268*).

Measuring modes:

- **pH** (potentiometric pH measurement)
- **U** (potentiometric voltage measurement)
- **Ipol** (voltametric measurement with selectable polarization current)
- **Upol** (amperometric measurement with selectable polarization voltage)
- **T** (temperature measurement)

[Dosing]



NOTE

Only active if a Dosino or dosing unit is connected.

Manual dosing (*see chapter 28.3, page 270*).

The following dosing functions are available:

- Preparing the exchange unit or dosing unit
- Emptying the dosing unit
- Filling the dosing cylinder of the exchange unit/dosing unit
- Dosing a specified volume
- Dosing continuously

[Stir]

Controlling the stirrer manually (*see chapter 28.4, page 279*).

[Remote]

Scanning the input lines on the remote interface and setting the output lines (see chapter 28.5, page 281).

[Sample Processor]

Controlling the Sample Processor manually (see chapter 28.6, page 282).


28.1 Opening and closing the manual control

Opening the manual control

1 Opening the dialog window

- Tap on the fixed key [].

Manual control opens:



2 Select the function

- Tap on the button for the desired function.

If the function can be carried out by several devices, a dialog to select the device will appear first; otherwise the dialog for the selected function will appear immediately.

Closing the manual control

- 1 ▪ Tap on the fixed key [].



NOTE

Manual control can also be exited when a manual sequence has been started but has not yet ended. The fact that a sequence has been started in manual control can be recognized by the corresponding symbol in the title line (*see chapter 6.2.3, page 39*).

28.2 Measuring

Manual control ► Measure

Manual control / Measure

Device: Ti-Touch

Sensor: pH electrode

Measuring input: 1

Measuring mode: pH

Temperature: 25 °C

Info sensor Start

Manual measurements can be carried out with the function **[Measure]**.

Proceed as follows:

1 Selecting a sensor

- Select the desired sensor from the sensor list.
All of the sensors available in the sensor list are displayed. Sensors are defined under **System ► Sensors**.
- Select the measuring input at which the sensor is connected.

2 Select a measuring mode


- Select the desired measuring mode.
Only those measuring modes are displayed which are advisable for the selected sensor.

- Enter the measuring temperature if no temperature sensor is connected. If a temperature sensor is connected, then the temperature will be measured automatically. This temperature is used for automatic temperature compensation with pH measurements.

4 Starting the measurement

- Tap on **[Start]**.

The measured value display is opened:



The ongoing measurement is canceled with **[Stop]**. The label switches once again to **[Start]**.

28.2.1 Parameter description

"Manual control / Measure" dialog

Sensor


Selection of the sensor from the sensor list. Sensors are defined under **System ► Sensors**.

Selection	Selection of configured sensors
-----------	---------------------------------

Measuring input

Selection of the measuring input the sensor is connected to. The selection depends on whether the control device has one or two measuring interfaces.

Selection	1 2
Default value	1



[Dosing device X]

Selection of the dosing device. The MSB connector and the type of the connected dosing device are displayed.

Various information is shown next to the button:

- Titrant
For dosing devices of the type 8XX with an intelligent exchange unit or dosing unit the titrant name is read off the data chip.
- Cylinder volume
For dosing devices of the type 8XX with an exchange unit without data chip only the cylinder volume is displayed.
- ---
For exchange or dosing units which are not or not correctly attached.
- Manual busy
When a dosing function is being carried out.

[All dosing devices]

Carrying out the functions **[Prepare]**, **[Empty]** and **[Fill]** at the same time with several dosing devices of the control device.

Selecting the dosing function

If only one dosing device is connected, then the properties of this dosing device will be shown directly. This information varies according to the type of dosing device:

- Type 8XX with integrated data chip:
If an exchange or dosing unit with integrated data chip is attached, the data stored on the data chip is displayed.

Manual control / Dosing device 1

Device: Ti-Touch

Dosing dev. type 800

Cylinder volume 10 mL

Titrant $c(\text{NaOH}) = 0.1 \text{ mol/L}$

Titer 1.077

Empty Prepare Fill Dosing fixed vol. Dosing

- Type 8XX without integrated data chip:
If an exchange unit without integrated data chip is attached, then the titrant can be selected from the titrant list. We recommend always selecting the titrant. This ensures that the parameters defined for the titrant for the preparing are used. Titrants are defined under **System ► Titrants**.

Manual control / Dosing device 2

Device: Ti-Touch

Dosing dev. type 805

Cylinder volume 20 mL

Titrant 

Prepare Fill Dosing fixed vol. Dosing

[Empty]

Empty the cylinder and the tubings of the dosing unit (see chapter 28.3.4, page 278).

[Prepare]

Rinse the cylinder and the tubings of the exchange / dosing unit (see chapter 28.3.3, page 277).

[Fill]

Fill the cylinder of the exchange / dosing unit (see chapter 28.3.5, page 278).

[Dosing fixed vol.]


Dose a specified volume (*see chapter 28.3.2, page 274*).

[Dosing]

Dose continuously (*see chapter 28.3.1, page 273*).

28.3.1 Continuous dosing

Manual control: **Dosing ► Dosing**



With the function **[Dosing]**, dosing is continuous for as long as you hold the button **[Dose]** pressed down.

Proceed as follows:

1 Configure the dosing function**NOTE**

- The dosing and filling rates should be decreased for viscous liquids.
- The maximum dosing and filling rates depend on the cylinder volume (*see chapter 32, page 443*).

- Enter the dosing rate.
- Enter the filling rate.

2 Start dosing

- Tap on **[Dose]**.

The volume dosed is displayed. When the volume of one cylinder has been dosed, the dosing cylinder will be refilled automatically.

3 Fill the cylinder

- Tap on **[Fill]**.

The dosing cylinder is filled. The displayed volume value will be reset to 0.0000 mL.



NOTE

If you exit the dialog with the fixed key [, the dosing cylinder is filled automatically.

Parameter description

Dosing rate

Rate at which it is dosed.

Input range	0.01 to 166.00 mL/min
Selection	maximum dynamic
Default value	maximum

maximum

The maximum dosing rate depends on the cylinder volume (see chapter 32, page 443).

dynamic

This setting can only be selected when the dosing cylinder is filled. The dosing is being carried out faster and faster until the maximum dosing rate is reached (starting with 1 mL/min and then doubling the dosing rate every 1.5 s).

Filling rate

Rate at which the dosing cylinder is filled.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum filling rate depends on the cylinder volume (see chapter 32, page 443).

28.3.2 Dosing fixed volumes

Manual control: **Dosing** ► **Dosing fixed vol.**

Dosing device 1 / Dosing fixed volume

Device: Ti-Touch
Titrant: c(NaOH) = 0.1 mol/L

Volume mL

Dosing rate mL/min

Filling rate mL/min

Fill automatic. ☒

You can dose a particular volume with the **[Dosing fixed vol.]** function.
Proceed as follows:

1 Configure the dosing function



NOTE


- The dosing and filling rates should be decreased for viscous liquids.
- The maximum dosing and filling rates depend on the cylinder volume (*see chapter 32, page 443*).

- Enter the desired volume.
- Enter the dosing rate.
- Enter the filling rate.

2 Start dosing

- Tap on **[Start]**.

The volume display is shown. When the volume of one cylinder has been dosed, the dosing cylinder will be refilled automatically.



Continuous dosing is paused with **[Hold]**. The label changes to **[Continue]**. It can be used to resume the run.

3 Fill the cylinder

With the default settings (see "Fill automatic.", page 277) the dosing cylinder is filled automatically.

Otherwise:

- Tap on **[Fill]**.

The dosing cylinder is filled. The displayed volume value will be reset to 0.0000 mL.

Parameter description

Volume

Volume which is dosed.

Input range	0.00000 to 99999.9 mL
Default value	0.10000 mL

Dosing rate

The rate at which dosing takes place.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum dosing rate depends on the cylinder volume (see chapter 32, page 443).

Filling rate

Rate at which the dosing cylinder is filled.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum filling rate depends on the cylinder volume (*see chapter 32, page 443*).

Fill automatic.

on | off (Default value: **on**)

If this parameter is activated, then the dosing cylinder is being automatically filled at the end of the dosing.

28.3.3 Preparing

Manual control: **Dosing ► Prepare**

The **Prepare** function is used to rinse and fill the cylinder and tubings of the exchange or dosing unit air bubble-free. You should carry out this function before the first determination or once a day.

If the titrant is selected (*see "Selecting the dosing function", page 271*), the parameters defined for the titrant for preparing/emptying and the tubing parameters will be used (dosing unit only). If the titrant is not selected, default parameters will be used (*see Chapter 32.2.2, page 444 and Chapter 32.3.2, page 445*).

The following two figures show the live display of an exchange unit and a dosing unit:




Figure 25 Live display "Preparing the dosing unit"




Figure 26 Live display "Preparing the exchange unit"

28.3.4 Emptying

Manual control: **Dosing ► Empty**



NOTE

The function **Empty** is possible only with dosing units.

The cylinders and the tubings of the dosing unit are emptied with this function.

If the titrant is selected (see *"Selecting the dosing function"*, page 271), the parameters defined for the titrant for preparing/emptying and the tubing parameters will be used. If the titrant is not selected, default parameters will be used (see *Chapter 32.3.2*, page 445 and *Chapter 32.2.2*, page 444).

28.3.5 Filling

Manual control: **Dosing** ► **Fill**

You can use the function **[Fill]** to fill the dosing cylinder of the exchange unit/dosing unit manually. The filling rate cannot be configured; the maximum filling rate will be used.

28.3.6 Replacing reagent



NOTE


You can easily change the reagent in a dosing unit without no contact with the chemicals using the functions **Prepare** and **Empty**.

Proceed as follows:

- 1 Execute the function **[Empty]**.
- 2 Attach the dosing unit to the bottle with a suitable solvent.
- 3 Execute the function **[Prepare]**.
The last traces of the reagent in the cylinder and tubings are removed with the solvent.
- 4 Execute the function **[Empty]** again.
- 5 Attach the dosing unit to the bottle with the new reagent.
- 6 Execute the function **[Prepare]**.
The dosing unit can now be used.

28.4 Stirring

Manual control ► Stir



The MSB connector and type of stirrer is indicated for each connected stirrer.

You can control a connected stirrer manually with the function **[Stir]**.

**NOTE**

Note that the rack has to be set to a valid position for all lift movements! Otherwise a corresponding error message will be displayed.

Move to a valid rack position with **[Previous]** or **[Next]**. Alternatively, you can carry out a MOVE command.

[Home position]

Move the lift to the home position (corresponds to lift position 0 mm).

[Work position]

Move the lift to the work position of the current rack position or the external position.

The work position can directly be defined here (see *"Assigning lift positions"*, page 285) or in the device manager (see *"Lift positions Tower 1/2"*, page 120).

[↑]

Continuously move the lift upwards as long as the button is pressed down. The lift rate used is defined in the dialog **Sample Processor / Lift**.

[↓]

Continuously move the lift downwards as long as the button is pressed down. The lift rate used is defined in the dialog **Sample Processor / Lift**.

[Previous]

Move to the previous rack position.

[Next]

Move to the next rack position.

[Reset rack]

Initialize the rack. This function is identical with the command **RACK**. The following actions are carried out:

- The lift is moved upwards.
- The sample rack is rotated to the starting position.
- The rack code of the rack attached is being read out and the corresponding rack data is transferred to the Sample Processor.
- The sample variable is being reset to the value 1.

[Pump]

Switch on and off the pumps mounted or connected to the tower.

[Lift]

Move the lift to any or a specific lift position and assigning lift positions (see chapter 28.6.1, page 284).

Input range	0 to 'maximum stroke path' mm A lift position of 0 mm corresponds to the "home position", i.e. the lift is located at the upper stop position. The maximum stroke path is defined in the properties of the tower (<i>see "Maximum stroke path", page 110</i>). If a higher value is entered, an error message will be displayed.
Selection	Work position Shift position Rinse position Special position
Default value	Work position

Lift rate

Rate at which the lift is moved in the manual control.

Input range	5 to 25 mm/s
Default value	25 mm/s

[Assign lift pos.]

Assign the current lift position to a preset lift position.

Assigning lift positions

Specific lift position are stored separately for every sample rack used as well as for tower 1 and tower 2. The following specific lift positions are available:

- **General rack positions**
For general rack positions you can define one work position, one shift position and one rinse position as well as one special position.
- **Special beaker positions**
A specific work position can be defined for each special beaker. The definitions of the general rack positions are used for the shift position, rinse position and special position of the affected tower.
- **External positions** (only with robotic arm and mounted Swing Head, *see Chapter 28.6.3, page 287*)
For each of the four possible external positions, a specific work position can be defined. Shift position and rinse position can only be defined for all four external positions commonly. A special position is not possible.

Proceed as follows to assign the current lift position to a specific lift position:

1 Move to lift position

- Enter the desired lift position in mm and tap on **[Start]**.

The lift moves to the desired lift position.

2 Assign lift position

- Tap on **[Assign lift pos.]**.

The **Lift / Assign lift position** dialog is displayed:

Lift / Assign lift position	
Current lift pos. 100 mm	
Position name	Position height
Work position	130 mm
Shift position	50 mm
Rinse position	100 mm
Special position	0 mm
Assign	

- Select the desired specific lift position and tap on **[Assign]**.

The current lift position is assigned to the specific lift position.

28.6.2 Moving to a rack position

Manual control: **Sample Processor** ► **Rack position**

Sample Processor / Rack position

Tower 1

Current rack pos. 1

Rack position

Shift rate %/s

Shift direction

Swing rate %/s

Start

In the dialog **Sample Processor / Rack position**, you can move to any position on the attached rack.

Current rack pos.

Current rack position. If the rack is not set to a valid position, e.g. after a rack reset, -- is displayed to indicate this.

Rack position

Desired rack position.

Special beaker 1...16

The selection depends on how many rack positions are defined as special beakers (see *"Editing special beakers"*, page 121).

Rate at which the sample rack is moved.

ection in which the rack is moved.

Direction in which the rack is moved.

auto
A shift direction, with which the shorter way has to be passed, is automatically selected.

+
Counterclockwise rotation.

-
Clockwise rotation.

Rate at which the robotic arm is swung in the manual control.

28.6.3 External positions

Manual control: **Sample Processor** ► **External position**



Sample Processor / External position

Tower 1

Lift position 100 mm Angle 105.0 °

°/s

Move to ext. pos. External position 1

NOTE

This dialog is accessible only if you have a Swing Head with robotic arm attached to the Sample Processor.

In the dialog **Sample Processor / External position**, you can swivel the robotic arm to an external position outside of the sample rack. As with the rack positions, here you can also define up to four frequently used external positions as specific positions. You can assign specific lift positions to these positions, analogous to the rack positions.

External positions are required if you have, for example, a measuring vessel or a titration cell mounted in addition to the rack.

Lift position

Current lift position.

Angle

Current angle position of the robotic arm.

[↑]

Only enabled, when the robotic arm has been moved to a specific external position with **[Start]**.

Continuously move the lift upwards as long as the button is pressed down. The lift rate used is defined in the dialog **Sample Processor / Lift**.

[↓]

Only enabled, when the robotic arm has been moved to a specific external position with **[Start]**.

Continuously move the lift downwards as long as the button is pressed down. The lift rate used is defined in the dialog **Sample Processor / Lift**.

[⇐]

Swing the robotic arm in clockwise direction as long as the button is pressed down.

[⇒]

Swing the robotic arm in counterclockwise direction as long as the button is pressed down.

Swing rate

Rate at which the robotic arm is swung in the manual control.

Input range	10 to 55 °/s
Default value	55 °/s

Move to ext. pos.

Selecting a predefined position or entering any swing angle.

Input range	(Offset) to (Offset + max. swing range) ° The offset is made up of a design-dependent angle (approx. 8...9°) together with the robotic arm offset from the robotic arm properties. The maximum swing range is also defined under the robotic arm properties (see "Properties – Robotic arm", page 113).
Selection	External position 1...4

[Assign Ext. pos.]

Assign the current angle position of the robotic arm to an external position as swing angle.

[Assign lift pos.]

Assign the current lift position to a preset lift position.

Assigning swing angles and lift positions

For each of the four possible external positions, a specific work position can be defined. Shift position and rinse position can only be defined for all four external positions commonly. A special position is not possible.

Proceed as follows to assign a swing angle to an external position and to assign the current lift position to a specific lift position:

1 Move to the external position

- Swing the robotic arm to the desired position with the arrow keys [⇐] or [⇒].

2 Assign the external position

- Tap on **[Assign Ext. pos.]**.

The **External position / Assign swing angle** dialog is displayed:

External position / Assign swing angle	
Current robotic arm angle 98.0 °	
External position	Angle
1	105.0 °
2	60.0 °
3	60.0 °
4	60.0 °

- Select the desired external position and tap on **[Assign]**.

The current angle position of the robotic arm is assigned to the external position.

3 Move to the external position

- Tap on [←].

The next higher-level dialog is displayed.

- Select **Move to ext. pos. = External position X** (X = 1...4) and tap on **[Start]**.

The robotic arm swings to the desired position.

4 Move to lift position

- Move the lift to the desired height with the arrow keys [↑] or [↓].

5

- Tap on **[Assign lift pos.]**.

The **Lift / Assign lift position** dialog is displayed:

Lift / Assign lift position	
Current lift pos. 100 mm	
Position name	Position height
Work position	130 mm
Shift position	50 mm
Rinse position	100 mm
Special position	0 mm

Assign

- Select the desired specific lift position.
- Tap on **Assign**.

The current lift position is assigned to the specific lift position.

29 Parameters



NOTE

For most of the numerical parameters, the result of a calculation can also be utilized in place of a number. A description of this can be found in the appendix (*see chapter 32.7, page 450*).

29.1 Dynamic equivalence point titrations (DET)

29.1.1 Start conditions

The parameters that are carried out before the start of titration are defined under **[Start conditions]**. The start conditions are processed in the following sequence:

1. Start volume
2. Start measured value
3. Start slope
4. Pause

Start volume

Volume that is dosed prior to the start of the titration.

Input range	0.00000 to 9999.99 mL
Default value	0.00000 mL

Dosing rate

Rate at which the start volume is dosed.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum dosing rate depends on the cylinder volume (see chapter 32, page 443).

Pause

Waiting time, e.g. for the electrode to settle down after the start or a reaction time after the dosing of a start volume. The time entered here is only running when all start conditions have been fulfilled.

Input range	0 to 999999 s
Default value	0 s

More start conditions

Start meas. value

Before starting the titration, it will be dosed until this measured value is reached. If the start measured value is reached by the dosing of a start volume, then the titration starts directly.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measuring mode pH:

Input range	-20.000 to 20.000
Selection	off
Default value	off

Measuring mode U, Ipol:

Input range	-2000.0 to 2000.0 mV
Selection	off
Default value	off

Measuring mode Upol:

Input range	-200.0 to 200.0 µA
Selection	off
Default value	off

Start slope

Before starting the titration, it will be dosed until this slope (measured value per volume) is reached. If the start slope is reached by the dosing of a start volume, then the titration starts directly.

Measuring mode pH:

Input range	00.000 to 9,999 pH/mL
Selection	off
Default value	off

Measuring mode U, Ipol:

Input range	0 to 999 mV/mL
Selection	off
Default value	off

Measuring mode Upol:

Input range	0 to 99 µA/mL
Selection	off
Default value	off

Dosing rate

Rate to be dosed with until the start measured value or the start slope is reached.

Input range	0.01 to 166.00 mL/min
Default value	5.00 mL/min
Selection	maximum

maximum

The maximum dosing rate depends on the cylinder volume (see chapter 32, page 443).

Initial measured value

The **[Initial meas. value]** is measured before carrying out the **start conditions**. This value is used for the **titration direction determination**, which is necessary for the **start measured value** and **stop measured value** parameters.

Signal drift

Maximum permissible drift for the measured value acceptance, i.e. maximum change of the measured value per minute. For this value, a suitable maximum waiting time is automatically calculated except you already have modified this waiting time.

Measuring mode pH, U and I_{pol}:

Input range	0.1 to 999.0 mV/min
Selection	off
Default value	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed.

Measuring mode Upol:

Input range	0.01 to 99.90 μA/min
Selection	off
Default value	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed.

Min. waiting time

The measured value is not accepted until the minimum waiting time has elapsed, even if the signal drift has already been reached. The minimum waiting time is only important for drift-controlled measurements.

Input range	0 to 999999 s
Default value	0 s

Max. waiting time

If the signal drift has been switched off or has not yet been reached, then the measured value will be accepted when the maximum waiting time has elapsed. As default value, a suitable waiting time for the signal drift is automatically calculated.

Input range	0 to 999999 s
Default value	1 s

29.1.2 Titration parameters

Under **[Titration parameters]**, the parameters influencing the run of the entire titration are defined.

Titration rate

Three predefined sets of parameters can be selected for the titration rate.

Selection	slow optimal fast user
Default value	optimal

slow

For titrations in which the finest details are also to be visible. This could however also lead to an increase in noise, which could result in unwanted equivalence points.

optimal

For all standard titrations. The parameters have been optimized for the most frequent applications.

fast

For fast and less critical titrations.

user

The individual titration parameters can be modified.



NOTE

Select **optimal** as titration rate when you are developing a new titration method. This parameter is suitable for virtually all titrations and only needs adaptation in special cases.

The settings of the individual titration rates are listed in *table 9, page 298*.

Temperature

Temperature entered manually. The temperature is being continuously measured when a temperature sensor is connected and when **Temp. measurement = automatic** or **continuous** is defined (see sensor dialog

of the command). This value is used for temperature correction in pH measurements.

Input range	–20.0 to 150.0 °C
Default value	25.0 °C

User-defined parameters

These parameters are only accessible when **Titration rate = user**.

Meas. point density

A small value means small volume increments, i.e. a high measuring point density. The curve then shows all the finest details which also include noise; this could cause unwanted equivalence points to be found. A larger value, i.e. a smaller measuring point density, permits quicker titrations. If you are using a dosing device with a small cylinder volume then a smaller measuring point density value may be beneficial. However, you should also set a smaller signal drift and a higher EP criterion at the same time.

Input range	0 to 9
Default value	4

Min. increment

This smallest permitted volume increment is added at the start of the titration and with steep curves in the region of the equivalence point. Very small values should only be used if a low titrant consumption is expected; otherwise unwanted equivalence points could be evaluated.

Input range	0.05 to 999.90 µL
Default value	10.00 µL

Max. increment

A maximum volume increment should be selected in the following cases:

- when titration consumption is very low up until the equivalence point is reached.
- when a start volume is dosed up until shortly before the equivalence point is reached.
- when the change of direction in the jumping range is very abrupt, because otherwise it could easily happen that an excessively large volume will be dosed in the region of the equivalence point.

The value should not be less than 1/100 cylinder volume.

Input range	0.1 to 9999.9 μL
Selection	off
Default value	off

**NOTE**

It is not advisable to select similar volumes for the minimum and the maximum increment. Monotonic equivalence point titration (MET) is appropriate for these applications.

Dosing rate

Rate at which the volume increments are dosed.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum dosing rate depends on the cylinder volume (*see chapter 32, page 443*).

Signal drift

Maximum permissible drift for the measured value acceptance, i.e. maximum change of the measured value per minute. This type of titration is often referred to as equilibrium titration. For this value, a suitable maximum waiting time is automatically calculated except you already have modified this waiting time.

**NOTE**

A constant measured value is often only reached after a certain time, as mixing and the reaction itself require a certain time. The response time of an electrode can also increase with time, i.e. reaching a constant measured value takes longer and longer. Drift-controlled measured value acceptance is particularly advisable in such cases, as the measured values are only accepted when equilibrium has almost been reached.

Measuring mode pH, U and Ipol:

Input range	0.1 to 999.0 mV/min
Default value	50.0 mV/min
Selection	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed. This can be useful when the titration reaction proceeds slowly or the electrode is slow to respond.

29.1.3 Stop conditions

The conditions for canceling the titration are defined under **[Stop conditions]**.

Stop volume

The titration is stopped when the specified volume has been dosed since the start of the titration. This volume should be adjusted to the size of the titration vessel in order to prevent the contents from running over.

Input range	0.00000 to 9999.99 mL
Default value	100,000 mL
Selection	off

Stop meas. value

The titration is canceled when the specified measured value has been reached since the start of the titration.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measuring mode pH:

Input range	-20.000 to 20.000
Selection	off
Default value	off

Measuring mode U, Ipol:

Input range	-2000.0 to 2000.0 mV
Selection	off
Default value	off

Measuring mode Upol:

Input range	-200.0 to 200.0 µA
Selection	off
Default value	off

Stop EP

The titration is stopped when the specified number of equivalence points has been found.

Input range	1 to 9
Default value	9
Selection	off

Volume after EP

This volume will be added when the number of equivalence points defined under **Stop EP** has been found. The curve shape after the equivalence point can also be seen this way.

Input range	0.01000 to 9999.99 mL
Selection	off
Default value	off

Stop time

The titration is stopped when the specified time has elapsed since the start of the titration.

Input range	0 to 999999 s
Selection	off
Default value	off

Filling rate

Rate at which the dosing cylinder is filled.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum filling rate depends on the cylinder volume (*see chapter 32, page 443*).

29.1.4 Potentiometric evaluation

The parameters for the evaluation of the titration curve are defined under **[Potentiometr. evaluation]**.

EP criterion

The Equivalence point Recognition Criterion (ERC) found is compared with this value. Equivalence points whose ERC is less than the value defined here will be ignored.

Input range	0 to 200
Default value	5

EP recognition

This parameter allows you to filter out only the equivalence points that are being sought.

Selection	all greatest last ascending descending off
Default value	all

all

All equivalence points will be recognized.

greatest

Only the equivalence point with the greatest ERC value, i.e. the steepest jump, will be recognized.

last

Only the last equivalence point will be recognized.

ascending

Only equivalence points with a positive slope of the titration curve will be recognized.

descending

Only equivalence points with a negative slope of the titration curve will be recognized.

off

No evaluation takes place.

Set windows

This parameter allows to recognize equivalence points only in a particular range (so-called window) of the curve. Equivalence points outside this window will not be recognized. A maximum of 9 windows can be defined. Setting windows is advisable when interference and unnecessary equivalence points are to be suppressed.

Selection	Measured value Volume off
Default value	off

Measured value

The windows are defined on the measured value axis.

Volume

The windows are defined on the volume axis.

off

The entire titration curve is being evaluated.

Setting windows

The list of defined windows is opened with **[Set window]**.

When the list is opened for the first time, a window over the entire measured value range or volume range is already defined. The windows must not overlap, they only may lie next to one another. For each window, own criteria can be defined for the equivalence point recognition.



NOTE

When you define a result variable instead of a numerical value as limit value (see *chapter 32.7, page 450*), then just one window can be defined.

[New]

Define a new window. This is only possible when not the entire range is covered yet.

[Delete]

Delete the selected window.

[Edit]

Edit the settings of the selected window.

Lower limit

Measured value or volume for the lower limit.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measured value (measuring mode pH):

Input range	-20.000 to 20.000
Default value	-20.000

Measured value (measuring mode U, Ipol):

Input range	−2000.0 to 2000.0 mV
Default value	−2000.0 mV

Measured value (measuring mode Upol):

Input range	−200.00 to 200.00 μA
Default value	−200.00 μA

Volume:

Input range	0.00000 to 9999.99 mL
Default value	0.00000 mL

Upper limit

Measured value or volume for the upper limit.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measured value (measuring mode pH):

Input range	–20.000 to 20.000
Default value	20.000

Measured value (measuring mode U, I_{pol}):

Input range	–2000.0 to 2000.0 mV
Default value	2000.0 mV

Measured value (measuring mode U_{pol}):

Input range	–200.00 to 200.00 µA
Default value	200.00 µA

Volume:

Input range	0.00000 to 9999.99 mL
Default value	9999.99 mL

EP criterion

The Equivalence point Recognition Criterion (ERC) found is compared with this value. Equivalence points whose ERC is less than the value defined here will be ignored.

Input range	0 to 200
Default value	5

EP recognition

This parameter allows you to filter out only the equivalence points that are being sought.

Selection	first greatest last ascending descending
Default value	first

first

Only the first equivalence point will be recognized.

greatest

Only the equivalence point with the greatest ERC value, i.e. the steepest jump, will be recognized.

last

Only the last equivalence point will be recognized.

ascending

Only equivalence points with a positive slope of the titration curve will be recognized.

descending

Only equivalence points with a negative slope of the titration curve will be recognized.

Only one equivalence point will be recognized in a window. The numbering of the equivalence points (EP) is defined by the numbering of the win-

For the recognition of the EPs found, the set EP criterion is compared to the ERC (Equivalence point Recognition Criterion) found. The ERC is the first derivative of the titration curve combined with a mathematical function that is more sensitive for flat jumps than for steeper ones. EPs whose ERC is smaller than the defined EP criterion will not be recognized. The ERC is displayed in the results dialog for each discovered and recognized EP. If you adjust the EP criterion later in order to recognize more or fewer EPs, then you can reevaluate the determination (**[Recalculate]** function in the results dialog).

The control instrument the determination is carried out with is selected under **[Control device]**. Control instruments are defined under **System ► Device manager**.

NOTE

Selection of the control instrument from the list of devices. Only those devices are displayed which are able to carry out the command.

Selection of the measuring input the sensor is connected to. The selection is not dependent on whether the control device has one or two measuring interfaces.

Selection of the sensor from the sensor list. The selection depends on the measuring mode. Sensors are defined under **System ► Sensors**. You can also enter a sensor name which is not contained in the sensor list. When a determination is started there is a check whether the sensor is contained in the sensor list.

Selection	Selection of configured sensors
-----------	---------------------------------

I(pol)

The polarization current is the current that is applied to a polarizable electrode during voltametric measurement. This parameter is available only with I(pol) determinations.

The values that are actually adjustable may vary from the values specified here. The binding values can be found in the manual, in the *Technical specifications* chapter.

Input range	-125.0 to 125.0 μA (Increment: 2.5)
Default value	5.0 μA
Selection	-1.0 μA 1.0 μA

U(pol)

The polarization voltage is the voltage applied to the polarizable electrode during an amperometric measurement. This parameter is available only with U(pol) determinations.

The values that are actually adjustable may vary from the values specified here. The binding values can be found in the manual, in the *Technical specifications* chapter.

Input range	-1,250 to 1,250 mV (Increment: 25)
Default value	400 mV

Electrode check

on | off (Default value: **off**)

For the following electrodes, an electrode check can be carried out:

- pH electrodes
- Metal electrodes
- Ion-selective electrodes

A check is made that the electrode is properly connected and that no short-circuit is present. The electrode check is carried out when this command is started. Note that this electrode check has nothing to do with the electrode test (command **ELT**).

Temp. measurement

Type of temperature measurement.

Selection	continuous automatic off
Default value	automatic

continuous

A temperature sensor must be connected. The temperature is measured continuously.

automatic

If a temperature sensor is connected then the temperature will be measured continuously. Otherwise, the temperature entered manually will be used (see dialog of the titration and measuring parameters).

off

The temperature will not be measured. The temperature entered manually is used (see dialog of the titration and measuring parameters).

29.1.7 Dosing device

The parameters for the dosing device are edited under **[Dosing device]**.

Dosing device

Selection of the MSB connector to which the dosing device is connected. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4
Default value	1

Titrant

Selection of the titrant from the titrant list. We recommend always selecting the titrant. This ensures that the correct data (titer, concentration, etc.) is always used for the calculation. Titrants are defined under **System ► Titrants**. You can also enter a name which is not contained in the titrant list. When a determination is started there is a check whether the solution is contained in the list.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct titrant has been attached and whether the type of dosing drive matches. For buret units without an integrated data chip, the cylinder volume and the type of dosing drive are checked. For the selected titrant, the validity of the titer, the working life of the titrant and the time interval for the GLP test of the buret unit are checked at the start of the determination.

Selection	Selection of configured titrants not defined
Default value	not defined

not defined

No check takes place.

29.1.8 Stirrer

The parameters for the stirrer are edited under **[Stirrer]**.

Selection of the MSB connector the stirrer is connected to. The selection does not depend on the number of MSB connectors of the control instrument.

off

Setting the stirring rate. It can be set in steps of -15 to $+15$. The default setting **8** corresponds to 1000 rpm. The formula for calculating the rotational speed is specified in *chapter 32.4, page 445*. The optimum stirring rate can be tested in the manual control.

- "+": counterclockwise rotation
- "-": clockwise rotation

on | off (Default value: **on**)

The parameters that are carried out before the start of titration are defined under **[Start conditions]**. The start conditions are processed in the following sequence:

Volume that is dosed prior to the start of the titration.

Dosing rate

Rate at which the start volume is dosed.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum dosing rate depends on the cylinder volume (*see chapter 32, page 443*).

Pause

Waiting time, e.g. for the electrode to settle down after the start or a reaction time after the dosing of a start volume. The time entered here is only running when all start conditions have been fulfilled.

Input range	0 to 999999 s
Default value	0 s

More start conditions

Start meas. value

Before starting the titration, it will be dosed until this measured value is reached. If the start measured value is reached by the dosing of a start volume, then the titration starts directly.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measuring mode pH:

Input range	–20.000 to 20.000
Selection	off
Default value	off

Measuring mode U, I_{pol}:

Input range	–2000.0 to 2000.0 mV
Selection	off
Default value	off

Measuring mode U_{pol}:

Input range	–200.0 to 200.0 µA
Selection	off
Default value	off

Selection	off
Default value	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed.

Measuring mode Upol:

Input range	0.01 to 99.90 μA/min
Selection	off
Default value	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed.

Min. waiting time

The measured value is not accepted until the minimum waiting time has elapsed, even if the signal drift has already been reached. The minimum waiting time is only important for drift-controlled measurements.

Input range	0 to 999999 s
Default value	0 s

Max. waiting time

If the signal drift has been switched off or has not yet been reached, then the measured value will be accepted when the maximum waiting time has elapsed. As default value, a suitable waiting time for the signal drift is automatically calculated.

Input range	0 to 999999 s
Default value	1 s

29.2.2 Titration parameters

Under **[Titration parameters]**, the parameters influencing the run of the entire titration are defined.

Titration rate

Three predefined sets of parameters can be selected for the titration rate.

Selection	slow optimal fast user
Default value	optimal

slow

For titrations in which the finest details are also to be visible. This could however also lead to an increase in noise, which could result in unwanted equivalence points.

**NOTE**

The accuracy of the evaluation cannot be increased by using smaller increments as the measured value changes between two measuring points are then of the same order of magnitude as the noise.

Input range	0.00005 to 999.900 mL
Default value	0.10000 mL

Dosing rate

Rate at which the volume increments are dosed.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum dosing rate depends on the cylinder volume (*see chapter 32, page 443*).

Signal drift

Maximum permissible drift for the measured value acceptance, i.e. maximum change of the measured value per minute. This type of titration is often referred to as equilibrium titration. For this value, a suitable maximum waiting time is automatically calculated except you already have modified this waiting time.

**NOTE**

A constant measured value is often only reached after a certain time, as mixing and the reaction itself require a certain time. The response time of an electrode can also increase with time, i.e. reaching a constant measured value takes longer and longer. Drift-controlled measured value acceptance is particularly advisable in such cases, as the measured values are only accepted when equilibrium has almost been reached.

Measuring mode pH, U and Ipol:

Input range	0.1 to 999.0 mV/min
Default value	50.0 mV/min
Selection	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed. This can be useful when the titration reaction proceeds slowly or the electrode is slow to respond.

Measuring mode Upol:

Input range	0.01 to 99.90 μA/min
Default value	50.00 μA/min
Selection	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed. This can be useful when the titration reaction proceeds slowly or the electrode is slow to respond.

Min. waiting time

The measured value is not accepted until the minimum waiting time has elapsed, even if the signal drift has already been reached. The minimum waiting time is only important for drift-controlled measurements.

Input range	0 to 999999 s
Default value	0 s

Max. waiting time

If the signal drift has been switched off or has not yet been reached, then the measured value will be accepted when the maximum waiting time has elapsed. As default value, a suitable waiting time for the signal drift is automatically calculated.

Input range	0 to 999999 s
Default value	26 s

Table 10 Default values of the predefined titration rates for MET

	Titration rate		
	slow	optimal	fast
Volume increment	0.05000 mL	0.10000 mL	0.20000 mL
Dosing rate	maximum	maximum	maximum
Signal drift			
– pH, U and I _{pol}	20.0 mV/min	50.0 mV/min	80.0 mV/min
– U _{pol}	20.0 µA/min	50.0 µA/min	80.0 µA/min
Min. waiting time	0 s	0 s	0 s
Max. waiting time	38 s	26 s	21 s

29.2.3 Stop conditions

The conditions for canceling the titration are defined under **[Stop conditions]**.

Stop volume

The titration is stopped when the specified volume has been dosed since the start of the titration. This volume should be adjusted to the size of the titration vessel in order to prevent the contents from running over.

Input range	0.00000 to 9999.99 mL
Default value	100,000 mL
Selection	off

Stop meas. value

The titration is canceled when the specified measured value has been reached since the start of the titration.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measuring mode pH:

Input range	-20.000 to 20.000
Selection	off
Default value	off

Measuring mode U, Ipol:

Input range	-2000.0 to 2000.0 mV
Selection	off
Default value	off

Measuring mode Upol:

Input range	-200.0 to 200.0 µA
Selection	off
Default value	off

Stop EP

The titration is stopped when the specified number of equivalence points has been found.

Input range	1 to 9
Default value	9
Selection	off

Volume after EP

This volume will be added when the number of equivalence points defined under **Stop EP** has been found. The curve shape after the equivalence point can also be seen this way.

Input range	0.01000 to 9999.99 mL
Selection	off
Default value	off

Stop time

The titration is stopped when the specified time has elapsed since the start of the titration.

Input range	0 to 999999 s
Selection	off
Default value	off

Filling rate

Rate at which the dosing cylinder is filled.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum filling rate depends on the cylinder volume (*see chapter 32, page 443*).

29.2.4 Potentiometric evaluation

The parameters for the evaluation of the titration curve are defined under **[Potentiometr. evaluation]**.

EP criterion

The Equivalence point Recognition Criterion (ERC) found is compared with this value. Equivalence points whose ERC is less than the value defined here will be ignored.

Measuring mode pH:

Input range	0.10 to 9.99
Default value	0.50

Measuring mode U, Ipol:

Input range	1 to 999 mV
Default value	30 mV

Measuring mode Upol:

Input range	0.1 to 99.9 μA
Default value	2.0 μA

EP recognition

This parameter allows you to filter out only the equivalence points that are being sought.

Selection	all greatest last ascending descending off
Default value	all

all

All equivalence points will be recognized.

greatest

Only the equivalence point with the greatest ERC value, i.e. the steepest jump, will be recognized.

last

Only the last equivalence point will be recognized.

ascending

Only equivalence points with a positive slope of the titration curve will be recognized.

descending

Only equivalence points with a negative slope of the titration curve will be recognized.

off

No evaluation takes place.

Set windows

This parameter allows to recognize equivalence points only in a particular range (so-called window) of the curve. Equivalence points outside this window will not be recognized. A maximum of 9 windows can be defined. Setting windows is advisable when interference and unnecessary equivalence points are to be suppressed.

Selection	Measured value Volume off
Default value	off

Measured value

The windows are defined on the measured value axis.

Volume

The windows are defined on the volume axis.

off

The entire titration curve is being evaluated.

Setting windows

The list of defined windows is opened with **[Set window]**.

When the list is opened for the first time, a window over the entire measured value range or volume range is already defined. The windows must not overlap, they only may lie next to one another. For each window, own criteria can be defined for the equivalence point recognition.



NOTE

When you define a result variable instead of a numerical value as limit value (see chapter 32.7, page 450), then just one window can be defined

[New]

Define a new window. This is only possible when not the entire range is covered yet.

[Delete]

Delete the selected window.

[Edit]

Edit the settings of the selected window.

Lower limit

Measured value or volume for the lower limit.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measured value (measuring mode pH):

Input range	-20.000 to 20.000
Default value	-20.000

Measured value (measuring mode U, Ipol):

Input range	–2000.0 to 2000.0 mV
Default value	–2000.0 mV

Measured value (measuring mode Upol):

Input range	-200.00 to 200.00 μA
Default value	-200.00 μA

Volume:

Input range	0.00000 to 9999.99 mL
Default value	0.00000 mL

Upper limit

Measured value or volume for the upper limit.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measured value (measuring mode pH):

Input range	–20.000 to 20.000
Default value	20.000

Measured value (measuring mode U, I_{pol}):

Input range	–2000.0 to 2000.0 mV
Default value	2000.0 mV

Measured value (measuring mode Upol):

Input range	–200.00 to 200.00 µA
Default value	200.00 µA

Volume:

Input range	0.00000 to 9999.99 mL
Default value	9999.99 mL

EP criterion

The Equivalence point Recognition Criterion (ERC) found is compared with this value. Equivalence points whose ERC is less than the value defined here will be ignored.

Measuring mode pH:

Input range	0.10 to 9.99
Default value	0.50

Measuring mode U, I_{pol}:

Input range	1 to 999 mV
Default value	30 mV

Measuring mode Upol:

Input range	0.1 to 99.9 µA
Default value	2.0 µA

EP recognition

This parameter allows you to filter out only the equivalence points that are being sought.

Selection	first greatest last ascending descending
Default value	first

- V_{EP} : EP volume
 V_0 : Dosed total volume before Δ_n
 ΔV : Volume increment
 ρ : Interpolation factor according to Fortuin

For the recognition of the EPs found, the set EP criterion is compared to the ERC (Equivalence point Recognition Criterion) found. The ERC is the sum of the measured value changes before and after the jump:

$$|\Delta_{n-2}| + |\Delta_{n-1}| + |\Delta_n| + |\Delta_{n+1}| + |\Delta_{n+2}|$$

In certain cases only three or only one summand is taken into account.

EPs whose ERC is smaller than the defined EP criterion will not be recognized. The ERC is displayed in the results dialog for each discovered and recognized EP. If you adjust the EP criterion later in order to recognize more or fewer EPs, then you can reevaluate the determination (**[Recalculate]** function in the results dialog).

29.2.5 Control instrument

The control instrument the determination is carried out with is selected under **[Control device]**. Control instruments are defined under **System ► Device manager**.



NOTE

This button will only be shown when, in addition to the Ti-Touch, a Sample Processor has been configured as control instrument.

Control device

Selection of the control instrument from the list of devices. Only those devices are displayed which are able to carry out the command.

Selection	Selection of configured control instruments
-----------	---

29.2.6 Sensor

The parameters for the sensor are edited under **[Sensor]**.

Measuring input

Selection of the measuring input the sensor is connected to. The selection is not dependent on whether the control device has one or two measuring interfaces.

Selection	1 2
Default value	1

Sensor

Selection of the sensor from the sensor list. The selection depends on the measuring mode. Sensors are defined under **System ► Sensors**. You can also enter a sensor name which is not contained in the sensor list. When a determination is started there is a check whether the sensor is contained in the sensor list.

Selection	Selection of configured sensors
-----------	---------------------------------

I(pol)

The polarization current is the current that is applied to a polarizable electrode during voltametric measurement. This parameter is available only with I(pol) determinations.

The values that are actually adjustable may vary from the values specified here. The binding values can be found in the manual, in the *Technical specifications* chapter.

Input range	-125.0 to 125.0 μA (Increment: 2.5)
Default value	5.0 μA
Selection	-1.0 μA 1.0 μA

U(pol)

The polarization voltage is the voltage applied to the polarizable electrode during an amperometric measurement. This parameter is available only with U(pol) determinations.

The values that are actually adjustable may vary from the values specified here. The binding values can be found in the manual, in the *Technical specifications* chapter.

Input range	-1,250 to 1,250 mV (Increment: 25)
Default value	400 mV

Electrode check

on | off (Default value: **off**)

For the following electrodes, an electrode check can be carried out:

- pH electrodes
- Metal electrodes
- Ion-selective electrodes

A check is made that the electrode is properly connected and that no short-circuit is present. The electrode check is carried out when this command is started. Note that this electrode check has nothing to do with the electrode test (command **ELT**).

Temp. measurement

Type of temperature measurement.

Selection	continuous automatic off
Default value	automatic

continuous

A temperature sensor must be connected. The temperature is measured continuously.

automatic

If a temperature sensor is connected then the temperature will be measured continuously. Otherwise, the temperature entered manually will be used (see dialog of the titration and measuring parameters).

off

The temperature will not be measured. The temperature entered manually is used (see dialog of the titration and measuring parameters).

29.2.7 Dosing device

The parameters for the dosing device are edited under **[Dosing device]**.

Dosing device

Selection of the MSB connector to which the dosing device is connected. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4
Default value	1

Titrant

Selection of the titrant from the titrant list. We recommend always selecting the titrant. This ensures that the correct data (titer, concentration, etc.) is always used for the calculation. Titrants are defined under **Sys-tem ► Titrants**. You can also enter a name which is not contained in the titrant list. When a determination is started there is a check whether the solution is contained in the list.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct titrant has been attached and whether the type of dosing drive matches. For buret units without an integrated data chip, the cylinder volume and the type of dosing drive are checked. For the selected titrant, the validity of the titer, the working life of the titrant and the time interval for the GLP test of the buret unit are checked at the start of the determination.

Selection	Selection of configured titrants not defined
Default value	not defined

29.3 Endpoint titrations (SET)

29.3.1 Start conditions

The parameters that are carried out before the start of titration are defined under **[Start conditions]**.

Pause 1

Waiting time, e.g. for the electrode to settle down before a start volume is added.

Input range	0 to 999999 s
Default value	0 s

Start volume

Volume that is dosed prior to the start of the titration.

Input range	0.00000 to 9999.99 mL
Default value	0.00000 mL

Dosing rate

Rate at which the start volume is dosed.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum dosing rate depends on the cylinder volume (*see chapter 32, page 443*).

Pause 2

Waiting time, e.g. for the electrode to settle down after the start or a reaction time after the dosing of a start volume.

Input range	0 to 999999 s
Default value	0 s

Initial measured value

The **[Initial meas. value]** is measured before carrying out the **start conditions**. This value is used for the **titration direction determination** when the titration direction is set to **auto** and only 1 EP is defined.

Signal drift

Maximum permissible drift for the measured value acceptance, i.e. maximum change of the measured value per minute. For this value, a suitable maximum waiting time is automatically calculated except you already have modified this waiting time.

Measuring mode pH:

Input range	–20.000 to 20.000
Selection	off
Default value	off

Measuring mode U and Ipol:

Input range	–2000.0 to 2000.0 mV
Selection	off
Default value	off

Measuring mode Upol:

Input range	–200.00 to 200.00 µA
Selection	off
Default value	off

Titration rate

Three predefined sets of parameters can be selected for the titration rate.

Selection	slow optimal fast user
Default value	optimal

slow

For steep titration curves for which dosing must be carried out in small steps at the endpoint.

optimal

For all standard titrations. The parameters have been optimized for the most frequent applications.

fast

For flat titration curves for which the endpoint is reached only slowly.

user

The individual titration parameters can be modified.

The settings of the individual titration rates are listed in *table 11, page 329*.

User-defined parameters

These parameters are only accessible when **Titration rate = user**.

Dynamics

This parameter defines the control range before the specified endpoint. Individual volume steps are dosed in the control range; the dosing is finely controlled. The closer the endpoint, the slower the dosing until the dosing rate defined under **Min. rate** has been reached. The larger the control range, the slower the titration. Outside of the control range, dosing is carried out continuously; the dosing rate is defined under **Max. rate**.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measuring mode pH:

Input range	0.001 to 20.000
Default value	2.000
Selection	off

Measuring mode U and Ipol:

Input range	0.1 to 2000.0 mV
Default value	100.0 mV
Selection	off

Measuring mode Upol:

Input range	0.01 to 200.00 μA
Default value	10.00 μA
Selection	off

Max. rate

Rate at which dosing is carried out outside of the control range.

Input range	0.01 to 166.00 mL/min
Default value	10.00 mL/min
Selection	maximum

maximum

The maximum dosing rate depends on the cylinder volume (see chapter 32, page 443).

Min. rate

Rate at which dosing is carried out at the very beginning of the titration and in the control range at the end of the titration. This parameter has a decisive influence on the titration rate and thus also on the accuracy. The smaller the selected minimum rate, the slower the titration.

Input range	0.01 to 9999.00 µL/min
Default value	25.00 µL/min

Table 11 Default values of the predefined titration rates for SET

	Titration rate		
	slow	optimal	fast
Dynamics			
– pH	5.000	2.000	0.500
– U und I _{pol}	300.0 mV	100.0 mV	30.0 mV
– U _{pol}	40.00 µA	10.00 µA	5.00 µA
Max. rate	1.00 mL/min	10.00 mL/min	maximum
Min. rate	5.00 µL/min	25.00 µL/min	50.00 µL/min

Stop criterion

Stop criterion

The titration is stopped when the endpoint has been reached and this stop criterion has been fulfilled. If no stop criterion has been selected then the titration will not be stopped. The stop conditions (*see chapter 29.3.4, page 331*) always lead to a stop, even if the stop criterion has not been reached.



NOTE

In older instructions the delay time was usually defined as the stop criterion. The same delay time may however result in different stopping time points – because of different smallest possible volume increments (depending on the cylinder volume). In contrast, if the stop drift is used, the titration will always be stopped at the same curve slope dV/dt .

Selection	drift time off
Default value	drift

drift

The titration is stopped when the stop drift has been reached.

time

The titration is stopped when the endpoint has been exceeded during a certain period of time (**Delay time**).

off

The titration will not be stopped until the stop conditions have been fulfilled.

Stop drift

This parameter can only be edited with **Stop criterion = drift**.

The titration is stopped when the endpoint and the stop drift have been reached.

Input range	1 to 999 µL/min
Default value	20 µL/min

Delay time

This parameter can only be edited with **Stop criterion = time**.

When the endpoint has been reached, the specified time is allowed to elapse after the last dosing and the titration is then stopped.

Input range	0 to 999 s
Default value	10 s

Endpoint 2

The control parameters for the second endpoint are defined under **[Endpoint 2]**. The parameters and input ranges are identical with those for the first endpoint.

29.3.3 Titration parameters

Under **[Titration parameters]**, the parameters influencing the run of the entire titration are defined.

Titration direction

With the setting "auto", the titration direction is determined automatically from the initial measured value and the set endpoint. It is recommended that you specify whether the titration direction is positive or negative. If two endpoints have been set, then the titration direction will be defined automatically. In this case, the setting will be ignored.

Selection	+ - auto
Default value	auto

+

Positive measured value change by the addition of titrant, i.e. in the direction of a higher pH value, greater voltage or greater current.

—

Negative measured value change by the addition of titrant, i.e. in the direction of a lower pH value, lesser voltage or lesser current.

auto

The titration direction is determined automatically from the initial measured value and the set endpoint.

Extraction time

Minimum duration of the titration. The titration will not be stopped during the extraction time, even if the endpoint has already been reached. The titration is however stopped if a stop condition is fulfilled during this time

(see chapter 29.3.4, page 331). The entry of an extraction time may be advisable, e.g. for the titration of sparingly soluble samples.

Input range	0 to 999999 s
Default value	0 s

Temperature

Temperature entered manually. The temperature is being continuously measured when a temperature sensor is connected and when **Temp. measurement = automatic** or **continuous** is defined (see sensor dialog of the command). This value is used for temperature correction in pH measurements.

Input range	-20.0 to 150.0 °C
Default value	25.0 °C

Time interval MP

Time interval for entering a measuring point in the measuring point list. The measuring point list is limited to 1000 measuring points.

Input range	0.1 to 999999.0 s
Default value	2.0 s

29.3.4 Stop conditions

The conditions for canceling the titration are defined under **[Stop conditions]**, if this does not occur automatically. This could be the case when the endpoint set is not reached or if the stop criterion (see "Stop criterion", page 329) is not fulfilled.

Stop volume

The titration is stopped when the specified volume has been dosed since the start of the titration. This volume should be adjusted to the size of the titration vessel in order to prevent the contents from running over.

Input range	0.00000 to 9999.99 mL
Default value	100,000 mL
Selection	off

Stop time

The titration is stopped when the specified time has elapsed following the termination of the start conditions.

Input range	1 to 999999 s
Selection	off
Default value	off

Filling rate

Rate at which the dosing cylinder is filled.

Drift for manual drift correction.

Input range	0.0 to 99.9 µL/min
Default value	0.0 µL/min

Measured value display

on | off (Default value: **off**)

If this parameter is activated, the currently measured value is displayed during the conditioning.

Conditioning options [Cond. options]

Cond. stop volume

Maximum permissible volume that can be dosed during conditioning. Conditioning is stopped when the specified volume is dosed. If conditioning is continued by pressing [] once again, then the titrant volume that has already been dosed will not be taken into account; i.e. the dosing starts again at zero. The stop volume should be adjusted to the size of the titration cell in order to prevent any overflow.

Input range	0.00000 to 9999.99 mL
Default value	20.0000 mL
Selection	off

Cond. stop time

Maximum permissible time over which conditioning may take place. Conditioning is stopped when the specified time has elapsed.

Input range	1 to 999999 s
Selection	off
Default value	off

Delay reconditioning

on | off (Default value: **off**)

If this parameter is activated, the conditioning will not start before all the method commands have been processed. Otherwise, the conditioning will immediately start after the titration command.



NOTE

This parameter must be activated if commands follow which must not run simultaneously with the conditioning.

Example: Emptying the titration cell followed by adding new working medium.

29.3.6 Cell

The monitoring of the reagent used is activated or deactivated under [Cell].

Reagent monitoring

Selection of the reagent from the list of reagents. The selection depends on the titration mode. Reagents are defined under **System ► Reagents**.

Selection	Selection of configured reagents off
Default value	off

off

The reagent monitoring is not active.

29.3.7 Control instrument

The control instrument the determination is carried out with is selected under **[Control device]**. Control instruments are defined under **System ► Device manager**.



NOTE

This button will only be shown when, in addition to the Ti-Touch, a Sample Processor has been configured as control instrument.

Control device

Selection of the control instrument from the list of devices. Only those devices are displayed which are able to carry out the command.

Selection	Selection of configured control instruments
-----------	---

29.3.8 Sensor

The parameters for the sensor are edited under **[Sensor]**.

Measuring input

Selection of the measuring input the sensor is connected to. The selection is not dependent on whether the control device has one or two measuring interfaces.

Selection	1 2
Default value	1

Sensor

Selection of the sensor from the sensor list. The selection depends on the measuring mode. Sensors are defined under **System ► Sensors**. You can also enter a sensor name which is not contained in the sensor list. When a

determination is started there is a check whether the sensor is contained in the sensor list.

Selection	Selection of configured sensors
-----------	---------------------------------

I(pol)

The polarization current is the current that is applied to a polarizable electrode during voltametric measurement. This parameter is available only with I(pol) determinations.

The values that are actually adjustable may vary from the values specified here. The binding values can be found in the manual, in the *Technical specifications* chapter.

Input range	−125.0 to 125.0 μA (Increment: 2.5)
Default value	5.0 μA
Selection	−1.0 μA 1.0 μA

U(pol)

The polarization voltage is the voltage applied to the polarizable electrode during an amperometric measurement. This parameter is available only with U(pol) determinations.

The values that are actually adjustable may vary from the values specified here. The binding values can be found in the manual, in the *Technical specifications* chapter.

Input range	−1,250 to 1,250 mV (Increment: 25)
Default value	400 mV

Electrode check

on | off (Default value: **off**)

For the following electrodes, an electrode check can be carried out:

- pH electrodes
- Metal electrodes
- Ion-selective electrodes

A check is made that the electrode is properly connected and that no short-circuit is present. The electrode check is carried out when this command is started. Note that this electrode check has nothing to do with the electrode test (command **ELT**).

Temp. measurement

Type of temperature measurement.

Selection	continuous automatic off
Default value	automatic

29.3.10 Stirrer

The parameters for the stirrer are edited under **[Stirrer]**.

Stirrer

Selection of the MSB connector the stirrer is connected to. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4 off
Default value	1

off

No stirrer will be used.

Stirring rate

Setting the stirring rate. It can be set in steps of -15 to $+15$. The default setting **8** corresponds to 1000 rpm. The formula for calculating the rotational speed is specified in *chapter 32.4, page 445*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation
- "-": clockwise rotation

Input range	-15 to 15
Default value	8

Switch off automatically

on | off (Default value: **on**)

If this parameter is activated, the stirrer will be switched off automatically at the end of the titration, measurement, etc.

29.4 STAT titrations (STAT)

29.4.1 Start conditions

The parameters that are carried out before the start of titration are defined under **[Start conditions]**.

Pause 1

Waiting time, e.g. for the electrode to settle down before a start volume is added.

Input range	0 to 999999 s
Default value	0 s

Start volume

Volume that is dosed prior to the start of the titration.

Input range	0.00000 to 9999.99 mL
Default value	0.00000 mL

Dosing rate

Rate at which the start volume is dosed.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum dosing rate depends on the cylinder volume (*see chapter 32, page 443*).

Pause 2

Waiting time, e.g. for the electrode to settle down after the start or a reaction time after the dosing of a start volume.

Input range	0 to 999999 s
Default value	0 s

Initial measured value

The **[Initial meas. value]** is measured before carrying out the **start conditions**. This value is used for the **titration direction determination** when the titration direction is set to **auto**.

Signal drift

Maximum permissible drift for the measured value acceptance, i.e. maximum change of the measured value per minute. For this value, a suitable maximum waiting time is automatically calculated except you already have modified this waiting time.

Input range	0.1 to 999.0 mV/min
Selection	off
Default value	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed.

Min. waiting time

The measured value is not accepted until the minimum waiting time has elapsed, even if the signal drift has already been reached. The minimum waiting time is only important for drift-controlled measurements.

Input range	0 to 999999 s
Default value	0 s

Max. waiting time

If the signal drift has been switched off or has not yet been reached, then the measured value will be accepted when the maximum waiting time has elapsed. As default value, a suitable waiting time for the signal drift is automatically calculated.

Input range	0 to 999999 s
Default value	1 s

29.4.2 Control parameters

The parameters for the control point are defined under **[Control parameters]**.

Control point at

Measured value for the control point.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measuring mode pH:

Input range	-20.000 to 20.000
Selection	off
Default value	off

Measuring mode U:

Input range	-2000.0 to 2000.0 mV
Selection	off
Default value	off

Titration rate

Three predefined sets of parameters can be selected for the titration rate.

maximum

The maximum dosing rate depends on the cylinder volume (*see chapter 32, page 443*).

Min. rate

The rate at which dosing is carried out at the very beginning of the titration and in the control range at the end of the titration. This parameter has a decisive influence on the behavior at the control point. The minimum rate also has a decisive influence on how exactly the control point can be kept constant. The smaller the selected minimum rate, the longer it takes until the control point is reached for the first time.

The following rule of thumb should be considered in order to obtain a dosing at the control point as constant (regular) as possible:

Min. rate in $\mu\text{L/min}$ = $0.1 \cdot \text{expected reaction rate in } \mu\text{L/min}$

Input range	0.01 to 9999.00 $\mu\text{L/min}$
Default value	10.00 $\mu\text{L/min}$

Table 12 Default values of the predefined titration rates for STAT

	Reaction rate		
	50 $\mu\text{L/min}$	100 $\mu\text{L/min}$	500 $\mu\text{L/min}$
Dynamics			
– pH	1.000	1.000	1.000
– U	60.0 mV	60.0 mV	60.0 mV
Max. rate	0.25 mL/min	0.75 mL/min	2.00 mL/min
Min. rate	5.00 $\mu\text{L/min}$	10.00 $\mu\text{L/min}$	40.00 $\mu\text{L/min}$

Dosing and control range

Outside the control range the dosing rate is mainly characterized by the maximum rate. The parameters **Max. rate** and **Dynamics** should be optimized together in such a way that the titration does not overshoot too much when the control point is reached. The Dynamics should be selected in such a way that the measured value remains within the control range while the control point is maintained. A slightly larger control range should be defined for slow reactions (e.g. pH = 3, U = 180 mV). The set-point value is often set with a SET pretitration before the substrate is added. This means that the STAT titration can be started already with a small deviation from the control point.

Start meas. value

The measured values are not entered in the measuring point list until this start measured value has been reached.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measuring mode pH:

Input range	–20.000 to 20.000
Selection	off
Default value	off

Measuring mode U:

Input range	–2000.0 to 2000.0 mV
Selection	off
Default value	off

Start rate

The measured values are not entered in the measuring point list until the dosing rate has fallen below this value. Checking this parameter is not started until 10 s after the start of the command.

Input range	0.01 to 166.00 mL/min
Selection	off
Default value	off



NOTE

If several of the parameters **Start time**, **Start meas. value** or **Start rate** have been defined, then all of them must be fulfilled for measured values to be entered in the measuring point list.

29.4.4 Stop conditions

The conditions for canceling the titration are defined under **[Stop conditions]**.

Stop volume

The titration is stopped when the specified volume has been dosed since the start of the titration. This volume should be adjusted to the size of the titration vessel in order to prevent the contents from running over.

Input range	0.00000 to 9999.99 mL
Default value	100,000 mL
Selection	off

29.4.5 Monitoring

Under **[Monitoring]**, the parameters for monitoring the following quantities are defined:

- Measured value
- Dosing rate (an average rate)
- Temperature

Monitoring is not used for parameterization of the control behavior. Remote signals or RS-232 commands can be assigned to the monitored quantities (see "Communication", page 349). These remote signals or RS-232 commands are sent if a limit has been infringed. This can be used, for example, to switch a thermostat on/off if necessary.

The parameters are identical for all monitored quantities.

Monitoring

on | off (Default value: **off**)

Enabling and disabling the monitoring function. The monitoring of the dosing rate is not started until 10 s after the start of the titration.

Lower limit

When the monitored quantity falls below this value, the action defined in the following is triggered.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measuring mode pH:

Input range	–20.000 to 20.000
Default value	–20.000

Measuring mode U:

Input range	–2000.0 to 2000.0 mV
Default value	–2000.0 mV

Dosing rate:

Input range	0.00 to 166.00 mL/min
Default value	0.00 mL/min

Temperature:

Input range	–20.0 to 150.0 °C
Default value	–20.0 °C

Lower hysteresis

The lower hysteresis represents the tolerance range for the lower limit (*see figure 32, page 349*). When **Action = Wait** the dosing does not continue until the monitored quantity has once again exceeded the lower limit incl. the value entered here.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measuring mode pH:

Input range	0.000 to 20.000
Default value	0.020

Measuring mode U:

Input range	0.0 to 2000.0 mV
Default value	2.0 mV

Dosing rate:

Input range	0.00 to 166.00 mL/min
Default value	0.20 mL/min

Temperature:

Input range	0.0 to 150.0 °C
Default value	0.2 °C

Upper limit

When this value is exceeded, the action defined in the following is triggered.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measuring mode pH:

Input range	-20.000 to 20.000
Default value	20.000

Measuring mode U:

Input range	-2000.0 to 2000.0 mV
Default value	2000.0 mV

Dosing rate:

Input range	0.00 to 166.00 mL/min
Default value	166.00 mL/min

Temperature:

Input range	-20.0 to 150.0 °C
Default value	150.0 °C

Upper hysteresis

The upper hysteresis represents the tolerance range for the upper limit (see figure 32, page 349). When **Action = Wait**, the dosing does not continue until the monitored quantity has once again fallen below the upper limit incl. the value entered here.

The values that are actually measurable may vary from the values specified here. The binding measurement value ranges can be found in the manual, in the "Technical specifications" chapter.

Measuring mode pH:

Input range	0.000 to 20.000
Default value	0.020

Measuring mode U:

Input range	0.0 to 2000.0 mV
Default value	2.0 mV

Dosing rate:

Input range	0.00 to 166.00 mL/min
Default value	0.20 mL/min

Temperature:

Input range	0.0 to 150.0 °C
Default value	0.2 °C

Action

Selection of the action which is carried out when the limits of the monitored quantity are infringed.

Selection	None Exit method Cancel command Pause Wait
-----------	---

None

No action will be taken.

Exit method

The method is stopped immediately.

Cancel command

The current command is stopped und the next command in the method will be started.



Under **[Communication]**, remote signals or RS-232 commands can be assigned to the monitored quantities. These remote signals or RS-232 commands are always sent, independent of the action defined for a limit infringement (*see "Action", page 347*).

Define a new assignment.

Delete the selected assignment.

Edit the settings for the selected assignment.

Selection of the monitored quantity to which a remote signal or an RS-232 command is to be assigned.

Indicating the limit whose infringement causes the remote signal or RS-232 command defined in the following to be sent.

■■■■■■■ 349

determination is started there is a check whether the sensor is contained in the sensor list.

Selection	Selection of configured sensors
-----------	---------------------------------

Electrode check

on | off (Default value: **off**)

For the following electrodes, an electrode check can be carried out:

- pH electrodes
- Metal electrodes
- Ion-selective electrodes

A check is made that the electrode is properly connected and that no short-circuit is present. The electrode check is carried out when this command is started. Note that this electrode check has nothing to do with the electrode test (command **ELT**).

Temp. measurement

Type of temperature measurement.

Selection	continuous automatic off
Default value	automatic

continuous

A temperature sensor must be connected. The temperature is measured continuously.

automatic

If a temperature sensor is connected then the temperature will be measured continuously. Otherwise, the temperature entered manually will be used (see dialog of the titration and measuring parameters).

off

The temperature will not be measured. The temperature entered manually is used (see dialog of the titration and measuring parameters).

29.4.8 Dosing device

The parameters for the dosing device are edited under **[Dosing device]**.

Dosing device

Selection of the MSB connector to which the dosing device is connected. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4
Default value	1

Titrant

Selection of the titrant from the titrant list. We recommend always selecting the titrant. This ensures that the correct data (titer, concentration, etc.) is always used for the calculation. Titrants are defined under **System ► Titrants**. You can also enter a name which is not contained in the titrant list. When a determination is started there is a check whether the solution is contained in the list.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct titrant has been attached and whether the type of dosing drive matches. For buret units without an integrated data chip, the cylinder volume and the type of dosing drive are checked. For the selected titrant, the validity of the titer, the working life of the titrant and the time interval for the GLP test of the buret unit are checked at the start of the determination.

Selection	Selection of configured titrants not defined
Default value	not defined
	not defined
	No check takes place.

Tandem dosing

A second dosing device can be defined under **[Tandem dosing]** in order to enable uninterrupted dosing. Dosing is carried out with a combination of two dosing devices so that the second dosing device is dosing while the first one is being filled and vice versa.

Dosing device

Selection of the MSB connector to which the dosing device is connected. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4 off
Default value	off

Titrant

Selection of the titrant from the titrant list. We recommend always selecting the titrant. Titrants are defined under **System ► Titrants**. You can also enter a name which is not contained in the titrant list. When a determination is started there is a check whether the solution is contained in list.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct titrant has been attached and whether the type of dosing drive matches. For buret units without integrated data chip, the cylinder volume and the type of dosing drive are checked. For

the selected titrant, only the time interval for the GLP test of the buret unit is checked at the start of the determination.



NOTE

The concentration, validity of the titer and the working life of the titrant are ignored.

Selection	Selection of configured titrants not defined
Default value	not defined

not defined

No check takes place.

Filling rate

Rate at which the dosing cylinder is filled.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum filling rate depends on the cylinder volume (*see chapter 32, page 443*).

However, in order to ensure uninterrupted dosing, the following points must be taken into account:

- Use the highest possible filling rate in order to keep the filling times as short as possible. The filling rate should be decreased for viscous liquids.
- When you use two buret units with different dosing cylinders, the filling rate for the larger sized cylinder must be at the minimum:

$$v_{2,Fill} \geq v_{1,Fill} \cdot \frac{V_{Cyl.2}}{V_{Cyl.1}}$$

 $v_{2,Fill}$ = filling rate in mL/min for the larger sized cylinder $v_{1,Fill}$ = filling rate in mL/min for the smaller cylinder

$V_{Cyl.2}$ = cylinder volume in mL of the buret unit of the second dosing device

$V_{Cyl,1}$ = cylinder volume in mL of the buret unit of the first dosing device

Example:

Dosing device 1: volume = 20 mL, filling rate = 50 mL/min

Dosing device 2: volume = 50 mL

Filling rate $2 \geq 50 \text{ mL/min} \cdot 50 \text{ mL} / 20 \text{ mL} \geq 125 \text{ mL/min}$

- The dosing rate must not exceed 75% of the filling rate of the smaller cylinder. These values are listed in the following table, valid at maximum filling rate:

Table 13 Maximum dosing rate for different dosing cylinders

Cylinder volume	maximum dosing rate	
	Exchange unit	Dosing unit
1 mL	2.25 mL/min	–
2 mL	–	5.00 mL/min
5 mL	11.25 mL/min	12.50 mL/min
10 mL	22.50 mL/min	25.00 mL/min
20 mL	45.00 mL/min	50.00 mL/min
50 mL	112.50 mL/min	124.50 mL/min

29.4.9 Stirrer

The parameters for the stirrer are edited under **[Stirrer]**.

Stirrer

Selection of the MSB connector the stirrer is connected to. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4 off
Default value	1

off

No stirrer will be used.

Stirring rate

Setting the stirring rate. It can be set in steps of –15 to +15. The default setting **8** corresponds to 1000 rpm. The formula for calculating the rotational speed is specified in *chapter 32.4, page 445*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation
- "–": clockwise rotation

Input range	–15 to 15
Default value	8

Switch off automatically

on | off (Default value: **on**)

If this parameter is activated, the stirrer will be switched off automatically at the end of the titration, measurement, etc.

29.5 Manual titrations (MAT)

A manual titration can be carried out with a Manual Dosing Controller or with the buttons **[Dose]** and **[Cancel command]**.

In order for the manual titration to be carried out with the 6.2107.120 Manual Dosing Controller, the latter must be connected to the MSB connector of the 916 Ti-Touch via the 6.2148.010 Remote Box MSB.

Manual titration using the Manual Dosing Controller

- 1 Press the **[GO]** key to start the method.
- 2 Press the **[GO]** key during the **MAT** command for dosing.
- 3 End the MAT command with the **[FILL]** key.

Manual titration using the buttons

- 1 Press the  key to start the method.
- 2 Press the **[Dose]** button during the MAT command for dosing.
- 3 End the MAT command with the **[Cancel command]** button.



NOTE

Start of titration: Because of the technical implementation involved, a few microliters are already dosed at the start of a manual titration. The possibility of deviation from the actual endpoint volume is excluded, however. The already dosed volume is taken into account.

End of titration: A manual titration must always be ended with the **[Cancel command]** button (NOT with the **[**

29.5.1 Dosing parameters

Under **[Dosing parameters]**, the parameters influencing the run of the manual titration are defined.

Dosing ramp

The dosing ramp is a delay in the dosing rate at the time of the start of a dosing step. It is advantageous primarily for manual titrations with indicator, when small volume steps are to be dosed before the equivalence point. The selectable delay time determines after which amount of time the defined dosing rate is to be reached.

Input range	0 to 10 s
Default value	2 s

Dosing rate

The rate at which dosing takes place.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum dosing rate depends on the cylinder volume (*see chapter 32, page 443*).

Time interval MP

Time interval for entering a measuring point in the measuring point list. The measuring point list is limited to 1000 measuring points.

Input range	0.1 to 999999.0 s
Default value	5.0 s

Filling rate

Rate at which the dosing cylinder is filled.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum filling rate depends on the cylinder volume (*see chapter 32, page 443*).

Temperature

Temperature entered manually. The temperature is being continuously measured when a temperature sensor is connected and when **Temp. measurement = automatic** or **continuous** is defined (see sensor dialog

of the command). This value is used for temperature correction in pH measurements.

Input range	-20.0 to 150.0 °C
Default value	25.0 °C

29.5.2 Control instrument

The control instrument the determination is carried out with is selected under **[Control device]**. Control instruments are defined under **System ► Device manager**.



NOTE

This button will only be shown when, in addition to the Ti-Touch, a Sample Processor has been configured as control instrument.

Control device

Selection of the control instrument from the list of devices. Only those devices are displayed which are able to carry out the command.

Selection	Selection of configured control instruments
-----------	---

29.5.3 Sensor

The parameters for the sensor are edited under **[Sensor]**.

Measuring input

Selection of the measuring input the sensor is connected to. The selection is not dependent on whether the control device has one or two measuring interfaces.

Selection	1 2
Default value	1

Sensor

Selection of the sensor from the sensor list. The selection depends on the measuring mode. Sensors are defined under **System ► Sensors**. You can also enter a sensor name which is not contained in the sensor list. When a determination is started there is a check whether the sensor is contained in the sensor list.

Selection	Selection of configured sensors
-----------	---------------------------------

Temp. measurement

Type of temperature measurement.

Selection	continuous automatic off
Default value	automatic

continuous

A temperature sensor must be connected. The temperature is measured continuously.

automatic

If a temperature sensor is connected then the temperature will be measured continuously. Otherwise, the temperature entered manually will be used (see dialog of the titration and measuring parameters).

off

The temperature will not be measured. The temperature entered manually is used (see dialog of the titration and measuring parameters).

29.5.4 Dosing device

The parameters for the dosing device are edited under **[Dosing device]**.

Dosing device

Selection of the MSB connector to which the dosing device is connected. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4
Default value	1

Titrant

Selection of the titrant from the titrant list. We recommend always selecting the titrant. This ensures that the correct data (titer, concentration, etc.) is always used for the calculation. Titrants are defined under **System ► Titrants**. You can also enter a name which is not contained in the titrant list. When a determination is started there is a check whether the solution is contained in the list.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct titrant has been attached and whether the type of dosing drive matches. For buret units without an integrated data chip, the cylinder volume and the type of dosing drive are checked. For the selected titrant, the validity of the titer, the working life of the titrant and the time interval for the GLP test of the buret unit are checked at the start of the determination.

Selection	Selection of configured titrants not defined
Default value	not defined

not defined

No check takes place.

Tandem dosing

A second dosing device can be defined under **[Tandem dosing]** in order to enable uninterrupted dosing. Dosing is carried out with a combination of two dosing devices so that the second dosing device is dosing while the first one is being filled and vice versa.

Dosing device

Selection of the MSB connector to which the dosing device is connected. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4 off
Default value	off

Titrant

Selection of the titrant from the titrant list. We recommend always selecting the titrant. Titrants are defined under **System ► Titrants**. You can also enter a name which is not contained in the titrant list. When a determination is started there is a check whether the solution is contained in list.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct titrant has been attached and whether the type of dosing drive matches. For buret units without integrated data chip, the cylinder volume and the type of dosing drive are checked. For the selected titrant, only the time interval for the GLP test of the buret unit is checked at the start of the determination.



NOTE

The concentration, validity of the titer and the working life of the titrant are ignored.

Selection	Selection of configured titrants not defined
Default value	not defined

not defined

No check takes place.

Filling rate

Rate at which the dosing cylinder is filled.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum filling rate depends on the cylinder volume (*see chapter 32, page 443*).

However, in order to ensure uninterrupted dosing, the following points must be taken into account:

- Use the highest possible filling rate in order to keep the filling times as short as possible. The filling rate should be decreased for viscous liquids.
- When you use two buret units with different dosing cylinders, the filling rate for the larger sized cylinder must be at the minimum:

$$v_{2,Fill} \geq v_{1,Fill} \cdot \frac{V_{Cyl.2}}{V_{Cyl.1}}$$

$v_{2,Fill}$ = filling rate in mL/min for the larger sized cylinder

$v_{1,Fill}$ = filling rate in mL/min for the smaller cylinder

$V_{Cyl.2}$ = cylinder volume in mL of the buret unit of the second dosing device

$V_{Cyl.1}$ = cylinder volume in mL of the buret unit of the first dosing device

Example:

Dosing device 1: volume = 20 mL, filling rate = 50 mL/min

Dosing device 2: volume = 50 mL

Filling rate 2 $\geq$ 50 mL/min · 50 mL / 20 mL $\geq$ 125 mL/min

- The dosing rate must not exceed 75% of the filling rate of the smaller cylinder. These values are listed in the following table, valid at maximum filling rate:

Table 14 Maximum dosing rate for different dosing cylinders

Cylinder volume	maximum dosing rate	
	Exchange unit	Dosing unit
1 mL	2.25 mL/min	–
2 mL	–	5.00 mL/min
5 mL	11.25 mL/min	12.50 mL/min
10 mL	22.50 mL/min	25.00 mL/min
20 mL	45.00 mL/min	50.00 mL/min
50 mL	112.50 mL/min	124.50 mL/min

29.5.5 Stirrer

The parameters for the stirrer are edited under **[Stirrer]**.

Stirrer

Selection of the MSB connector the stirrer is connected to. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4 off
Default value	1

off

No stirrer will be used.

Stirring rate

Setting the stirring rate. It can be set in steps of -15 to $+15$. The default setting **8** corresponds to 1000 rpm. The formula for calculating the rotational speed is specified in *chapter 32.4, page 445*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation
- "-": clockwise rotation

Input range	-15 to 15
Default value	8

Switch off automatically

on | off (Default value: **on**)

If this parameter is activated, the stirrer will be switched off automatically at the end of the titration, measurement, etc.

29.6 Measurements (MEAS)

29.6.1 Measuring parameters

Under **[Measuring parameters]**, the parameters influencing the run of the entire measurement are defined.

Signal drift

Maximum permissible drift for the measured value acceptance, i.e. maximum change of the measured value per minute. For this value, a suitable maximum waiting time is automatically calculated except you already have modified this waiting time.



NOTE

A constant measured value is often only reached after a certain time, as mixing and the reaction itself require a certain time. The response time of an electrode can also increase with time, i.e., reaching a constant measured value takes longer and longer. Drift-controlled measured value acceptance is particularly advisable in such cases, as the measured values are only accepted when equilibrium has almost been reached.

Measuring mode pH, U and I_{pol}:

Input range	0.1 to 999.0 mV/min
Default value	10.0 mV/min
Selection	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed.

Measuring mode U_{pol}:

Input range	0.01 to 99.90 µA/min
Default value	10.00 µA/min
Selection	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed.

Measuring mode T:

Input range	0.1 to 999.0 °C/min
Default value	0.5 °C/min
Selection	off

Input range	-20.0 to 150.0 °C
Selection	off
Default value	off

This parameter is not available with the **MEAS T** command (temperature measurement).

Temperature entered manually. The temperature is being continuously measured when a temperature sensor is connected and when **Temp. measurement = automatic** or **continuous** is defined (see sensor dialog of the command). This value is used for temperature correction in pH measurements.

Input range	-20.0 to 150.0 °C
Default value	25.0 °C

Time interval for entering a measuring point in the measuring point list.
The measuring point list is limited to 1000 measuring points.

Input range	0.1 to 999999.0 s
Default value	2.0 s

The control instrument the determination is carried out with is selected under **[Control device]**. Control instruments are defined under **System ► Device manager**.



NOTE

This button will only be shown when, in addition to the Ti-Touch, a Sample Processor has been configured as control instrument.

Selection of the control instrument from the list of devices. Only those devices are displayed which are able to carry out the command.

Selection	Selection of configured control instruments
-----------	---

Electrode check

on | off (Default value: **off**)

For the following electrodes, an electrode check can be carried out:

- pH electrodes
- Metal electrodes
- Ion-selective electrodes

A check is made that the electrode is properly connected and that no short-circuit is present. The electrode check is carried out when this command is started. Note that this electrode check has nothing to do with the electrode test (command **ELT**).

Temp. measurement

Type of temperature measurement.

Selection	continuous automatic off
Default value	automatic

continuous

A temperature sensor must be connected. The temperature is measured continuously.

automatic

If a temperature sensor is connected then the temperature will be measured continuously. Otherwise, the temperature entered manually will be used (see dialog of the titration and measuring parameters).

off

The temperature will not be measured. The temperature entered manually is used (see dialog of the titration and measuring parameters).

29.6.4 Stirrer

The parameters for the stirrer are edited under **[Stirrer]**.

Stirrer

Selection of the MSB connector the stirrer is connected to. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4 off
Default value	1

off

No stirrer will be used.

Stirring rate

Setting the stirring rate. It can be set in steps of -15 to $+15$. The default setting **8** corresponds to 1000 rpm. The formula for calculating the rota-

elapsed. As default value, a suitable waiting time for the signal drift is automatically calculated.

Input range	0 to 999999 s
Default value	110 s
Selection	off

Temperature

Temperature entered manually. The temperature is being continuously measured when a temperature sensor is connected and when **Temp. measurement = automatic** or **continuous** is defined (see sensor dialog of the command).

Input range	–20.0 to 150.0 °C
Default value	25.0 °C

Sample Processor

If the calibration is carried out with a Sample Processor, the buffers or standards are automatically changed. This is why no request for the calibration temperature appears when starting the calibration. The value entered above is used.

Selection	off Remote USB
Default value	off

off

For calibrations without Sample Processor. During the calibration, a request for changing the buffer/standard is displayed.

Remote

For automatic calibrations with Sample Processors, which are connected via the Remote Box.

USB

For automatic calibrations with a Sample Processor that is connected via USB(see "*Definition of the calibration positions*", page 452).

Buffers

The buffer type and the number of buffers is defined under **[Buffers]**.

Buffer type

Selection of a predefined buffer series or definition of special buffers. In the case of predefined buffer series, the instrument automatically recognizes which buffer is involved.

Selection	Metrohm NIST DIN Fisher Fluka Basel Mettler Merck Tit. Beckman Radiometer Custom Special Merck Cer. Baker Hamilton Precisa
Default value	Metrohm

Input range	-20.000 to 20.000
Selection	off
Default value	off

Buffer 4 pH

See **Buffer 3 pH**.

Buffer 5 pH

See **Buffer 3 pH**.

Stirrer control

The parameters for the stirrer control during the calibration are edited under **[Stirrer control]**. Ensure that a stirrer has been selected in the **Edit command / Stirrer** dialog.

**NOTE**

Deactivate the **Stir solution during measurement** parameter when you carry out the calibration with a Sample Processor. The stirrer control for the tower stirrer is ineffective with these instruments.

Stir solution during measurement

on | off (Default value: **on**)

If this parameter is activated, the stirrer will be switched on and off automatically at the start or end of the measurement.

Stir before meas.

This parameter can only be edited when the parameter **Stir solution during measurement** is deactivated.

If the stirrer is switched off during the measurement, it will be stirred for the time entered before measuring.

Input range	0 to 999999 s
Default value	0 s

Pause before meas.

This parameter can only be edited when the parameter **Stir solution during measurement** is deactivated.

If it is stirred before the measurement, a waiting time can be defined here before starting the measurement. No stirring or measuring occurs during this time.

Input range	0 to 999999 s
Default value	0 s

Temp. measurement

Type of temperature measurement.

Selection	continuous automatic off
Default value	automatic

continuous

A temperature sensor must be connected. The temperature is measured continuously.

automatic

If a temperature sensor is connected then the temperature will be measured continuously. Otherwise, the calibration temperature entered manually will be used.

off

The temperature will not be measured. The calibration temperature entered manually is used. The temperature is requested after the start, except the determination is carried out with a Sample Processor.

29.7.4 Stirrer

The parameters for the stirrer are edited under **[Stirrer]**.

Stirrer

Selection of the MSB connector the stirrer is connected to. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4 off
Default value	1

off

No stirrer will be used.

Stirring rate

Setting the stirring rate. It can be set in steps of –15 to +15. The default setting **8** corresponds to 1000 rpm. The formula for calculating the rotational speed is specified in *chapter 32.4, page 445*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation
- "–": clockwise rotation

Input range	–15 to 15
Default value	8

Stirrer

A stirrer has to be connected to the titrator for the electrode test.

Temperature

Take care to ensure that the entire electrode test is carried out at a consistently uniform temperature. Work should be accomplished at room temperature wherever possible, in light of the fact that temperature exercises a great influence on response time. The preset limit values are calibrated to 25 °C.

29.8.1 Electrode test parameters

Under **[Electrode test param.]**, the parameters influencing the sequence of the electrode test are defined.

Buffer type

Selecting a predefined buffer series.

Selection	Metrohm NIST DIN Fluka Basel Mettler Merck Tit. Radiometer Merck Cer. Baker Hamilton Precisa
Default value	Metrohm

Merck Cer.

Reference temperature = 25 °C. When using Merck CertiPUR buffers (20 °C) the buffer type **Merck Tit.** must be selected.

Sample Processor

If the calibration is carried out with a Sample Processor, the buffers or standards are automatically changed. This is why no request for the calibration temperature appears when starting the calibration. The value entered above is used.

Selection	off Remote USB
Default value	off

off

For calibrations without Sample Processor. During the calibration, a request for changing the buffer/standard is displayed.

Remote

For automatic calibrations with Sample Processors, which are connected via the Remote Box.

USB

For automatic calibrations with a Sample Processor that is connected via USB(see "*Definition of the calibration positions*", page 452).

Electrode type

Selection of the electrode type.

Selection	1 2
Default value	1

Sensor

Selection of the sensor from the sensor list. The selection depends on the measuring mode. Sensors are defined under **System ► Sensors**. You can also enter a sensor name which is not contained in the sensor list. When a determination is started there is a check whether the sensor is contained in the sensor list.

Selection	Selection of configured sensors
-----------	--

Temp. measurement

Type of temperature measurement.

Selection	continuous automatic off
Default value	automatic

continuous

A temperature sensor must be connected. The temperature is measured continuously.

automatic

If a temperature sensor is connected then the temperature will be measured continuously. Otherwise, the calibration temperature entered manually will be used.

off

The temperature will not be measured. The calibration temperature entered manually is used. The temperature is requested after the start, except the determination is carried out with a Sample Processor.

29.8.4 Stirrer

The parameters for the stirrer are edited under **[Stirrer]**.

Stirrer

Selection of the MSB connector the stirrer is connected to. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4
Default value	1

Stirring rate

Setting the stirring rate. It can be set in steps of –15 to +15. The default setting **8** corresponds to 1000 rpm. The formula for calculating the rotational speed is specified in *chapter 32.4, page 445*. The optimum stirring rate can be tested in the manual control.

Table 16 Possible errors during the electrode test

Test criterion	Message	Remedy
<i>Buffer pH 9:</i> The drift in stirred solution is > 1 mV.	General problem	▪ Connect an electrode. ▪ Replace the defective electrode cable. ▪ Clean the diaphragm (see leaflet for the electrode). ▪ Replace the electrode.
<i>Buffer pH 9:</i> $-10 \text{ mV} < U(\text{not stirred}) < 10 \text{ mV}$ and The sum of the drift values after 1, 2, 3 and 4 min is < 12 mV/min.	Short circuit	Replace the electrode.
A response time does not fulfill the limit value for the rating Usable electrode .	Glass membrane / Diaphragm	Clean the diaphragm (see leaflet for the electrode).
The pH values are not defined for all the buffers at the temperature measured.	No buffer data	Repeat the electrode test at a temperature at which the pH values of all the buffers are defined.
All the slopes fulfill the limit values for the rating Usable electrode . and Uoff is outside the preset limits.	Unsuitable reference electrode	▪ Repeat the electrode test with a suitable reference system. ▪ Replace the contaminated reference electrolyte. ▪ Select Electrode type = Custom and adjust the limits for Uoff.
The streaming potential is too high.	Contaminated diaphragm	Clean the diaphragm (see leaflet for the electrode).
Two slopes do not fulfill the limit values for the rating Usable electrode .	Wrong buffer	Repeat the electrode test with buffers pH 4, 7 and 9.
All the slopes do not fulfill the limit values for the rating Usable electrode .	Partial short circuit	▪ Check the temperature sensor. ▪ Enter the correct calibration temperature.

29.9 Evaluations (EVAL)

Various **additional evaluations** (EVAL commands) can be carried out for titrations and measurements. The evaluation always refers to the last titration or measurement prior to the EVAL command. Only those evaluation commands which are available for the last titration or measurement prior to the EVAL command can be inserted in the command list. If a titration or measuring command is deleted before the EVAL command, then it will be shown in red in the list of commands, as the reference is missing.

EVAL commands can also be inserted in the method run at a later time and the determination can be reevaluated (**[Recalculate]** function in the results dialog).

Additional evaluations for DET titrations

The following additional evaluations are possible:

- **EVAL FIX-EP** (fixed endpoints)
Measured quantities:
 - Measured value
 - Time
 - Volume
- **EVAL pK/HNP** (pK value / half neutralization potential HNP)
This evaluation is only possible for the measuring modes **pH** and **U**.
- **EVAL MIN/MAX** (minimum/maximum)
- **EVAL BREAK** (break points)

Additional evaluations for MET titrations

The following additional evaluations are possible:

- **EVAL FIX-EP** (fixed endpoints)
Measured quantities:
 - Measured value
 - Time
 - Volume
- **EVAL pK/HNP** (pK value / half neutralization potential HNP)
This evaluation is only possible for the measuring modes **pH** and **U**.
- **EVAL MIN/MAX** (minimum/maximum)
- **EVAL BREAK** (break points)

Additional evaluations for SET titrations

The following additional evaluations are possible:

- **EVAL FIX-EP** (fixed endpoints)
Measured quantities:
 - Measured value
 - Time
 - Volume
- **EVAL MIN/MAX** (minimum/maximum)
- **EVAL RATE** (average dosing rate)

Additional evaluations for STAT titrations

The following additional evaluations are possible:

- **EVAL FIX-EP** (fixed endpoints)
Measured quantities:
 - Measured value
 - Time
 - Volume
- **EVAL MIN/MAX** (minimum/maximum)
- **EVAL RATE** (average dosing rate)

Additional evaluations for measurements

The following additional evaluations are possible:

- **EVAL FIX-EP** (fixed endpoints)
Measured quantities:
 - Measured value
 - Time
- **EVAL MIN/MAX** (minimum/maximum)
- **EVAL BREAK** (break points)

29.9.1 Fixed endpoint evaluation (EVAL FIX-EP)

For the fixed endpoint evaluation, the associated values are interpolated from the measuring point list for one quantity (measured value, volume, etc.). Up to nine fixed endpoints can be evaluated with one command.

Fixed quantity

Selection of the quantity to which the associated value is interpolated from the measuring point list.

Selection	Measured value Time Volume
Default value	Measured value

Fixed EP1 at

The fixed endpoint must lie between the first and the final entry in the measuring point list.

The values that are actually measurable may vary from the values specified here, depending on the hardware. The binding measurement value ranges can be found in the .

approximation, as the ionic strength is not taken into account. In order to obtain more accurate values, titrations must be carried out with decreasing ionic strength and the results extrapolated to the ionic strength zero. pK evaluation in aqueous solution is limited to the range $3.5 < \text{pK} < 10.5$ because of the leveling effect of strong acids and the lack of jumps with very weak acids. pK values of mixtures of acids and polyvalent acids can also be determined.




Figure 33 Evaluation of pK value / half neutralization potential

In non-aqueous solutions the half neutralization potential (HNP) is frequently used instead of the pK value. The HNP is evaluated in the same way as the pK value.

No parameters can be edited for the command **EVAL pK/HNP**.



NOTE

If a start volume is to be added then it must be smaller than the half endpoint volume.

29.9.3 Minimum and maximum evaluation (EVAL MIN/MAX)

For the minimum or maximum measured value the associated volume, time and temperature are interpolated from the measuring point list. The evaluation begins as soon as the slope of the curve exceeds a particular threshold value.




Figure 34 Evaluation of minimum and maximum

Evaluation

Selection whether the minimum or the maximum of the titration curve is being evaluated. With a command either the minimum or the maximum can be evaluated. If both values are needed, a second command has to be defined in the method.

Selection	Minimum Maximum
Default value	Maximum

Threshold value

The evaluation of the minimum or maximum begins as soon as the slope of the curve exceeds the set threshold value. Use a lower threshold value if the minimum or maximum is not found.

Measuring mode pH:

Input range	0.1 to 20
Default value	1.0
	With titrations, the unit is pH/mL, with measurements pH/s.

Measuring mode U, Ipol:

Input range	1.0 to 2,000.0
Default value	25.0
	With titrations the unit is mV/mL, with measurements mV/s.

Measuring mode Upol:

Input range	0.5 to 10.0
Default value	5.0
	With titrations the unit is $\mu\text{A/mL}$, with measurements $\mu\text{A/s}$.

Measuring mode T:

Input range	0.1 to 20.0 °C/s
Default value	1.0 °C/s

29.9.4 Break point evaluation (EVAL BREAK)

A break point evaluation is used to determine sharp changes of direction in the titration curve. This evaluation is primarily used for photometric and conductivity titrations. The method is based on the search for extremes in the second derivative of the curve.




Figure 35 Evaluation of a break point

EP criterion

Measure of the minimum sharpness of the break point. The smaller the EP criterion set, the more break points will be found. As this is a relative value related to the total measured value modification, even small changes in the measured value can be evaluated as a break point for a small measured value range.

Input range	0.0 to 1.0
Default value	0.3

Slope

Minimum difference between the slope before and after the break point. The smaller the difference, the more break points will be found.

Input range	0.0 to 10.0
Default value	0.9

Smoothing factor

The higher the smoothing factor, the fewer break points will be found.

Input range	2 to 20
Default value	5

Volume:

Input range	0.00000 to 9,999.99 mL
Default value	0.00000 mL

Upper limit

Measured value, time or volume for the upper limit.

The values that are actually measurable may vary from the values specified here, depending on the hardware. The binding measurement value ranges can be found in the .

Measured value (measuring mode pH):

Input range	-20.000 to 20.000
Default value	20.000

Measured value (measuring mode U, Ipol):

Input range	-2,000.0 to 2,000.0 mV
Default value	2,000.0 mV

Measured value (measuring mode Upol):

Input range	-200.00 to 200.00 µA
Default value	200.00 µA

Measured value (measuring mode T):

Input range	-20.0 to 150.0 °C
Default value	150.0 °C

Time:

Input range	0.0 to 999,999.9 s
Default value	999,999.9 s

Volume:

Input range	0.00000 to 9,999.99 mL
Default value	9,999.99 mL

29.9.5 Rate evaluation (EVAL RATE)

The rate evaluation provides the median dosing rate in a defined range (so-called window). The median dosing rate is determined by means of linear regression through a minimum of three measuring points. If the dosing is carried out with only one dosing device and if more than one cylinder volume is dosed, then the median dosing rate will be reduced as a result of the filling time.

When the command is opened for the first time, a window over the entire time range is already defined. In contrary to the potentiometric evaluation

of a DET/MET titration the individual windows may overlap. A maximum of 9 windows can be defined.

[New]

Define a new window.

[Delete]

Delete the selected window.

[Edit]

Edit the settings of the selected window.

Lower limit

Time for the lower limit.

Input range	0 to 999999 s
Default value	0 s

Upper limit

Time for the upper limit.

Input range	0 to 999999 s
Default value	999999 s

29.10 Calculations

The following calculation commands are available:

CALC	Defining calculations of a determination.
CALC LIVE	Defining the calculation whose result will be displayed in the live display (so-called live result) during a titration or measurement.

29.10.1 Calculations (CALC)

Calculations are defined with the command **CALC**. A method may contain a maximum of nine calculation commands. A maximum of nine calculations can be defined per command. A series of variables (raw data from the determination, previously calculated results, etc.) is available for the calculations.

The command contains the following options, among others:

- Save result as titer
- Save result as common variable
- Save result in result table
- Define result limits

[illegible]

The result name is specified in the list for each calculation.

[New]

Define a new calculation or load an existing template, see following chapter.

[Delete]

Delete the selected calculation from the list.

[Edit]

Edit the data of the selected calculation (see chapter 29.10.1.2, page 392).

29.10.1.1 Loading a result template

The most common calculations are already available, the so-called result templates, in order to facilitate the definition of a new calculation. You also have the option, however, of defining templates of your own (see dialog **System / Templates**).

Enter the numerical values of the wildcards

- New calculation / Load template

F1=

35.45

F6=

F2=

F7=

F3=

F8=

F4=

F9=

F5=

Cancel

Back

Next

- The editing dialog of the new calculation is displayed:

Edit command / Edit calculation

02

CALC

Calculation

Result name

Content (%)

Calc. formula R3

EP1*CONC*TITER*35.45*0.1/C00

Decimal places

2

▼

Result unit

%

▼

Note

Result variable

Result limits

Result options

The calculation can now be edited further, e.g. define result limits, see following chapter.

29.10.1.2 Editing a calculation

Result name

The result name is the text which will be shown in the result display and in the report. The default name corresponds to the result variable.

Entry	max. 24 characters
Default value	R1...R9

Calc. formula R1...R9

Shows the calculation formula. A special editor is opened for the definition (see chapter 29.10.3, page 398). The designation **R1...R9** corresponds to the result variable. Under **[Result variable]** this can later be modified.

Entry	max. 100 characters
Default value	empty

Decimal places

Number of decimal places used to display the result.

Input range	0 to 5
Default value	2

Result unit

The result unit is displayed and saved together with the result.

Entry	max. 10 characters
Selection	% mol/L mmol/L g/L mg/L mg/mL mg/100 g ppm g mg µg mL µL mg/ piece °C µg/min mL/min µL/min
Default value	%

[Note]

Entering a note on the calculation.

[Result variable]

Modifying the result variable.

[Result limits]

Defining the limits for the monitoring of the result.

[Result options]

Defining additional settings for the calculation.

Dialog "Edit calculation / Note"

In this dialog you can enter a short text, e.g. to describe the variables used.

Dialog "Edit calculation / Result variable"

When a new calculation is created, the result variable will be automatically issued. In this dialog, however, it can be modified at any time.

Result variable

To each calculation belongs an unambiguous result variable. With this result variable you can use this result in additional calculations.

Selection	R1 R2 R3 R4 R5 R6 R7 R8 R9
	The selection contains only those result variables which are not issued in this calculation command yet.

Dialog "Edit calculation / Result limits"

For each result, limit values can be defined. These result limits are monitored when the calculation is carried out. If monitoring is activated, the results are shown as follows in the result display:

- Green, if the result falls within the limit values
- Red, if the result lies outside the limit values

Monitoring result limits

on | off (Default value: **off**)

Enabling and disabling the monitoring function.

Lower limit

When the monitored quantity falls below this value, the action defined in the following is triggered.

Input range	0 to 9999999999
Default value	0

Selection	Single value Mean value off
Default value	off

Single value

The result of the current determination is saved as titer.

Mean value

If statistic calculations are carried out for the result, the current mean value of the determination series is saved.

off

The result is not saved as titer.

Save as common variable

on | off (Default value: **off**)

The calculated result can be saved as a method-independent variable, a so-called common variable. The result is then also available in other methods for calculations. In all cases, only the single value will be saved, even if the statistics function is enabled.

Variable

Selection of the common variable to which the result is being assigned.

Selection	CV01 ... CV25
-----------	----------------------

Display result

on | off (Default value: **on**)

If you deactivate this parameter, the result is neither displayed in the result dialog nor printed in the result report. This can be advisable for intermediate results.

Save result in result table

on | off (Default value: **off**)

The calculated result can be saved in the result table. This may be advisable if e.g. the results of all determinations carried out on a particular day are to be displayed clearly. A maximum of nine results from a determination can be saved in the result table.

Precision

Setting, with which accuracy the result is used in additional calculations.

Selection	Round Truncate Full precision
Default value	Round

Decimal places

Number of decimal places used to display the result.

Input range	0 to 5
Default value	2

Result unit

The result unit is displayed and saved together with the result.

Entry	max. 10 characters
Selection	% mol/L mmol/L g/L mg/L mg/mL mg/100 g ppm g mg µg mL µL mg/ piece °C µg/min mL/min µL/min
Default value	%

[Note]

Entering a note on the calculation.

Dialog "Edit calculation / Note"

In this dialog you can enter a short text, e.g. to describe the variables used.

Displaying the live result

Proceed as follows to display the live result in the live display:

1 Starting a determination

- Tap on [].

The determination starts and the live display is shown:

2 Define the display options

- Tap on the **[View]** button.
- Tap on the **[Meas. value options]** button.
A maximum of three measured values can be shown in the live display.
- Select the option **Live result** at one of the three parameters.

3 Display the live display


- Tap twice on [].

The current result is displayed,

29.10.3 Formula editor

The formulas for the calculations are entered with the formula editor. The maximum formula length is 100 characters.

The formula editor contains buttons for entering numbers, mathematical operators and variables. The variables can be selected from lists. This is advisable in order to avoid spelling mistakes. The formula editor is equipped with an automatic syntax check, which is activated when applying the formula. The standard rules of priority apply for the calculation operations.



Key	Description
C00	Sample size
CI	Sample identification CI# (# = 1 - 2)
EP	Volume of endpoint EP# (# = 1...9)
CV	Common Variable CV# (# = 1...25)
R	Result R# (# = 1 - 9)
TITER	Titer of the titrant selected in the previous command
CONC	Concentration of the titrant selected in the previous command
□	The character in front of the cursor is deleted.
Delete entry	The content of the input field is deleted completely.
⇐	The cursor within the input field is shifted to the left by one character.
⇒	The cursor within the input field is shifted to the right by one character.

If the method contains more than one data-generating command (titration, measurement etc.) before the calculation command, then you will have to indicate the command identification in front of the variable:

Variable	Description
#M.	Titration or measuring command
#C.	Calculation command
#E.	Evaluation command

"#" stands for a sequential number that you must enter manually. The number for each command identification will be counted separately and does not necessarily correspond to the number of the command line.

Example: The variable **1M.TITER** corresponds to the titer of the titrant selected in the first titration command.



CAUTION

The command identifications are not adjusted automatically when additional commands are subsequently inserted into the method. Check the formula for correctness in such cases.

[Variables]

Opening the list of additional variables (*see chapter 29.10.3.1, page 399*).

[Math. functions]

Opening the list with mathematical functions (*see chapter 29.10.3.2, page 400*).

29.10.3.1 Variables

All of the variables which can be used for calculations are divided up according to subject in the dialog **Calculation formula / Variables**:

- **System variable**
List of the system-specific variables. These variables describe the current status of the system.
- **Result and statistics variables**
List of the variables which are required for calculating results and statistics for the current determination.
- **Common variables**
List of common variables. All of the common variables are listed, no matter whether they possess a value or not.

You will find a detailed description of all of the variables in *chapter 29.10.3.3, page 400*.


Variable	Description
DD	Duration of the entire determination
Titration	
TITER	Titer of the titrant selected in the titration command
CONC	Concentration of the titrant selected in the titration command
Titration, measurements, calibrations	
EP#	Volume/quantity of endpoint EP# (# = 1 - 9)
EC#	Charge at endpoint EP# (# = 1 - 9)
EM#	Measured value of endpoint EP# (# = 1...9)
EF#	ERC of endpoint EP# (# = 1 - 9)
ET#	Temperature at endpoint EP# (# = 1...9)
ED#	Time at endpoint EP# (# = 1...9)
ESI#	Recognition of endpoint EP# (# = 1 - 9) Endpoint found = 1, no endpoint = 0
FMN	Mean gas flow during MEAS T/F ON to MEAS T/F OFF
TMN	Mean temperature during MEAS T/F ON to MEAS T/F OFF
TMI	Minimum temperature during MEAS T/F ON to MEAS T/F OFF
TMA	Maximum temperature during MEAS T/F ON to MEAS T/F OFF
MIM	Initial measured value, i.e. measured value prior to the processing of the start conditions
MIT	Initial temperature, i.e. temperature prior to the processing of the start conditions
MSA	Volume for start volume
MSP	Volume for start measured value
MSS	Volume for start slope
MSV	Volume for all the start conditions
MSD	Duration of start conditions
MSM	Start measured value, i.e. measured value after the processing of the start conditions
MST	Start temperature, i.e. temperature after the processing of the start conditions
MCD	Total duration of the titration, measurement or calibration
MTM	Type of temperature measurement (Pt1000, NTC or manually) Format = Text
MDD	Duration of effective dosing, i.e. without filling times, pauses
MDC	Drift for drift correction
DDC	Time for drift correction
MCV	End volume, i.e. total dosed volume at the end of the titration

Variable	Description
BP#	Volume of break point BP# (# = 1 - 9)
BM#	Measured value of break point BP# (# = 1 - 9)
BT#	Temperature at break point BP# (# = 1 - 9)
BD#	Time at break point BP# (# = 1 - 9)
RD#	Mean rate in time slot # (# = 1 - 9)
RDS#	Standard deviation for RD# (# = 1 - 9)
RDC#	Correlation coefficient for RD# (# = 1 - 9)
RM	Mean rate for whole range
RMS	Standard deviation for RM
RMC	Correlation coefficient for RM
Results, statistics	
R#	Result (# = 1–9)
SMN#	Mean value of result R# (# = 1–9)
SSA#	Absolute standard deviation of SMN# (# = 1 - 9)
SSR#	Relative standard deviation of SMN# in % (# = 1 - 9)
SNR#	Number of results from which the mean value SMN# was calculated (# = 1 - 9)
SSD	Number of determinations that are carried out for the statistics calculations
SNT	Statistics status (statistics activated = 1, statistics deactivated = 0)
Common variables	
CV#	Common variable (# = 1...25)
System variables	
%RN	Sample number, i.e. number that is increased by one at each start of a determination
%SC	Start counter, see dialog More determination data / Properties
%AS	Autostart status (autostart activated = 1, autostart deactivated = 0)
%AC	Autostart counter, i.e. number of autostarts already carried out
%AD	Autostart setpoint counter, i.e. number of autostarts to be carried out
%SS	Status of the sample table (sample table activated = 1, sample table deactivated = 0)
%SL	Number of the current sample table line
%SE	Last sample in the sample table processed (yes = 1, no = 0)

29.11 Reports (REPORT)

The reports to be printed out are defined with the command **REPORT**. A maximum of ten reports can be defined for each command. The following reports can be inserted, among others:

- Result report
- Calculation report
- Parameter report
- System reports (system settings, sensor list, GLP data etc.)
- Method reports (method run, method parameters etc.)



[Report options]

Defining the general settings for the report output.

[Insert]

Add a new report to the list.

[Delete]

Delete the selected report from the list.

[Edit]

Edit the settings of the selected report.

29.11.1 General report options

The settings for the report printout can be adjusted in the dialog **Edit command / Report options** to indicate, for example, whether a report header, a signature line or a frame is to be printed.

Report header

Defining the output of the report header. The report header contains general information such as device type, serial number, program version and the printing date.



Selection	off once on each page
Default value	on each page

The report header will not be printed.

The report header will only be printed on the first page.

The report header will be printed on every page.

Output of a special line for date and signature. This line will be printed at the very bottom of every page.

off

No signature line will be printed.

The signature line will only be printed on the last page.

The signature line will be printed on every page.

on | off (Default value: **on**)

If this parameter is activated, a frame is printed as lateral marking.

Fixed key [

Command REPORT: **Report** ► **Edit**

Settings can be edited for the following reports (for details, see online help):

- **Result report**
- **Curve**
Definition of the curve size, display of the individual measuring points and grid lines, etc.

- **Measuring point list**
Definition of the method command to which the measuring point list is to be printed.
- **Calculations**
Definition of the accuracy with which the variables used as well as the results calculated are printed.
- **Statistics short**
Definition of the time at which the report is to be printed (with each determination or only at the end of a sample series).
- **Statistics overview**
Definition of the time at which the report is to be printed (with each determination or only at the end of a sample series).
- **Sample table**
- **Result table**
Definition of the time at which the report is to be printed (with each determination, at the end of a sample series or at the end of the sample table).
- **Used devices**
- **Variables**
- **Monitoring**

29.11.3 List of reports

The following reports can be inserted in the command **REPORT**:



NOTE

If you are using a **CUSTOM NEOS** printer, you can print only the reports that are marked with a **#** on the display.

Report	Contents
Result report	# Report with determination properties, sample data, calculated results, etc. If there are multiple determinations, the statistics will also be printed out.
Curve	# Curve report.
Measuring point list	# Measuring point list report.
Calculations	# Details concerning the calculations carried out (parameter settings in accordance with the CALC command and the calculated results).
Used devices	The devices used for the determination, as displayed in the More determination data / Properties dialog.
Variables	All determination variables, as displayed in the results dialog.
Monitoring	Details for the monitored quantities (STAT only).
Statistics	

Report	Contents
Statistics short	# Summary of the statistics calculations. The number of determinations, the mean value, and the absolute and relative standard deviations are printed out for each result.
Statistics overview	# Detailed statistical overview. The sample data and all individual results are printed out for each determination. The number of determinations, the mean value, and the absolute and relative standard deviations are printed out for each result.
Sample data	
Sample table	List of all determinations in the sample table with the respective sample data, as entered in the sample table.
Result table	
Result table	List of all determinations in the result table with results and with the determination data, as saved in the result table.
More reports ► Method reports	
Method sequence	Method properties and the list of all method commands.
Parameters full	Method properties and options, all method commands with all parameters. All parameters which no longer have their default settings will be printed in bold . All parameters which have been modified in comparison to the stored version of the method are indicated by an *.
Titration & measuring param.	Method properties; titration, measurement and calibration commands with all parameters. All parameters which no longer have their default settings will be printed in bold . All parameters which have been modified in comparison to the stored version of the method are indicated by an *.
Modified parameters	Method properties, all method parameters which have been modified in comparison to the stored version of the method.
Non default parameters	Method properties, all method parameters which no longer have any default settings.
More reports ► System reports	
System settings	
System settings	Settings for acoustic signals and accuracy of the measured value display.
Dialog options	Settings for routine dialog and expert dialog.
Titrants	
Titrant list	List of all the titrants configured in the system.
All titrant data short	The most important titrant data of all titrants (name, concentration, titer, last titer determination).
All titrant data full	All of the titrant data of all of the titrants (working life, titer options, data concerning the exchange unit/dosing unit).

Report	Contents
Sensors	
Sensor list	List of all the sensors configured in the system.
All sensor data short	The most important sensor data for all sensors (name, calibration data).
All sensor data full	All of the sensor data for all of the sensors (working life, complete calibration data, calibration interval).
Device manager	
Device list	List of all devices configured in the system.
All device properties	Properties of all the devices configured in the system.
GLP manager	
GLP data	All data stored in the GLP manager.
Common variables	
Common variable list	List of all the common variables defined in the system, together with their most important data (name, value, status).
All common variable properties	Properties of all common variables (name, value, validity, status).
Templates	
Templates sample data	Sample identification list and sample assignment table.
Result template list	List with all of the user-generated result templates.
All result templates details	Details of the all of the user-generated result templates (calculation formula, result options, note).
Input/Output lines	List with all of the defined input and output lines at the remote interface (name, bit pattern).
Custom calibration buffers	Temperature tables for all defined custom calibration buffers.
Rack tables	
Sample rack list	List with all of the sample racks configured in the system, together with designation, number of positions and rack code.
Miscellaneous	
Form feed	If this entry is inserted between two reports, then these will each be printed out on a separate page.

29.12 Dosing and Liquid Handling

The following dosing commands are available:

Preparing the exchange unit or dosing unit (PREP)	Rinsing the cylinder and the tubings of the exchange unit / dosing unit.
Emptying the dosing unit (EMPTY)	Emptying the cylinder and the tubings of the dosing unit.
Dosing a fixed volume (ADD)	Dosing a specified volume.
Liquid Handling (LQH)	Carrying out complex dosing tasks with a Dosino.

29.12.1 Preparing an exchange or dosing unit (PREP)

The **PREP** command is used to rinse and fill air bubble-free the cylinder and tubings of the exchange or dosing unit. You should carry out this function before the first determination or once a day.

Control device

This parameter will only be shown when, in addition to the Ti-Touch, a Sample Processor has been configured as control device. Control devices are defined under **System ► Device manager**.

Selection of the control device from the list of devices.

Selection	Selection of configured control devices
-----------	---

Dosing device

Selection of the MSB connector to which the dosing device is connected. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4
Default value	1

Titrant

Selection of the titrant from the titrant list. We recommend always selecting the titrant. This ensures that the parameters defined for the titrant for the preparing as well as the tubing parameters (dosing unit only) are used. Titrants are defined under **System ► Titrants**. You can also enter a name which is not contained in the titrant list. When a determination is started there is a check whether the solution is contained in the list.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct titrant has been attached and whether

the type of dosing drive matches. For buret units without integrated data chip, the cylinder volume and the type of dosing drive are checked. For the selected titrant, the validity of the titer, the working life of the titrant and the time interval for the GLP test of the buret unit are checked at the start of the determination.

Selection	Selection of configured titrants not defined
Default value	not defined

not defined

No check takes place. The command is carried out with default parameters (see *Chapter 32.2.2, page 444* and *Chapter 32.3.2, page 445*).

Cylinder volume to

This parameter can only be edited with **Titrant = not defined**.

Selection of the port via which the titrant is ejected. This setting is only relevant for dosing units, in case of preparing an exchange unit, this setting is ignored.

Selection	Port 1 Port 2 Port 3 Port 4
Default value	Port 1

29.12.2 Emptying a dosing unit (EMPTY)

The **EMPTY** command is used to empty the cylinder and tubings of the dosing unit.

Control device

This parameter will only be shown when, in addition to the Ti-Touch, a Sample Processor has been configured as control device. Control devices are defined under **System ► Device manager**.

Selection of the control device from the list of devices.

Selection	Selection of configured control devices
-----------	---

Dosing device

Selection of the MSB connector to which the dosing device is connected. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4
Default value	1

Titrant

Selection of the titrant from the titrant list. We recommend always selecting the titrant. This ensures that the parameters defined for the titrant for the preparing as well as the tubing parameters (dosing unit only) are used.

Titration are defined under **System ► Titrants**. You can also enter a name which is not contained in the titrant list. When a determination is started there is a check whether the solution is contained in the list.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct titrant has been attached and whether the type of dosing drive matches. For buret units without integrated data chip, the cylinder volume and the type of dosing drive are checked. For the selected titrant, the validity of the titer, the working life of the titrant and the time interval for the GLP test of the buret unit are checked at the start of the determination.

Selection	Selection of configured titrants not defined
Default value	not defined

not defined

No check takes place. The command is carried out with default parameters (see Chapter 32.2.2, page 444 and Chapter 32.3.2, page 445).

Air inlet

This parameter can only be edited with **Titrant = not defined**.

Selection of the port via which air is aspirated. This setting is only relevant for dosing units, in case of preparing an exchange unit, this setting is ignored.

Selection	Port 1 Port 2 Port 3 Port 4
Default value	Port 4

29.12.3 Dosing a specified volume (ADD)

You can dose a specified volume with the command **ADD**.

29.12.3.1 Dosing parameters

The parameters for the dosing are defined under **[Dosing parameters]**.

Volume

Volume which is dosed.

Input range	0.00000 to 99999.9 mL
Default value	10.0000 mL

Dosing rate

The rate at which dosing takes place.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

titrant list. When a determination is started there is a check whether the solution is contained in the list.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct titrant has been attached and whether the type of dosing drive matches. For buret units without an integrated data chip, the cylinder volume and the type of dosing drive are checked. For the selected titrant, the validity of the titer, the working life of the titrant and the time interval for the GLP test of the buret unit are checked at the start of the determination.

Selection	Selection of configured titrants not defined
Default value	not defined

not defined

No check takes place.

Tandem dosing

A second dosing device can be defined under **[Tandem dosing]** in order to enable uninterrupted dosing. Dosing is carried out with a combination of two dosing devices so that the second dosing device is dosing while the first one is being filled and vice versa.

Dosing device

Selection of the MSB connector to which the dosing device is connected. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4 off
Default value	off

Titrant

Selection of the titrant from the titrant list. We recommend always selecting the titrant. Titrants are defined under **System ► Titrants**. You can also enter a name which is not contained in the titrant list. When a determination is started there is a check whether the solution is contained in list.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct titrant has been attached and whether the type of dosing drive matches. For buret units without integrated data chip, the cylinder volume and the type of dosing drive are checked. For the selected titrant, only the time interval for the GLP test of the buret unit is checked at the start of the determination.



NOTE

The concentration, validity of the titer and the working life of the titrant are ignored.

Selection	Selection of configured titrants not defined
Default value	not defined

not defined

No check takes place.

Filling rate

Rate at which the dosing cylinder is filled.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum filling rate depends on the cylinder volume (*see chapter 32, page 443*).

However, in order to ensure uninterrupted dosing, the following points must be taken into account:

- Use the highest possible filling rate in order to keep the filling times as short as possible. The filling rate should be decreased for viscous liquids.
- When you use two buret units with different dosing cylinders, the filling rate for the larger sized cylinder must be at the minimum:

$$v_{2,Fill} \geq v_{1,Fill} \cdot \frac{V_{Cyl.2}}{V_{Cyl.1}}$$

$$v_{2,Fill}$$
 = filling rate in mL/min for the larger sized cylinder $v_{1,Fill}$ = filling rate in mL/min for the smaller cylinder

$V_{Cyl,2}$ = cylinder volume in mL of the buret unit of the second dosing device

$V_{Cyl.1}$ = cylinder volume in mL of the buret unit of the first dosing device

Example:

Dosing device 1: volume = 20 mL, filling rate = 50 mL/min

Dosing device 2: volume = 50 mL

$$\text{Filling rate } 2 \geq 50 \text{ mL/min} \cdot 50 \text{ mL} / 20 \text{ mL} \geq 125 \text{ mL/min}$$

- The dosing rate must not exceed 75% of the filling rate of the smaller cylinder. These values are listed in the following table, valid at maximum filling rate:

Table 17 Maximum dosing rate for different dosing cylinders

Cylinder volume	maximum dosing rate	
	Exchange unit	Dosing unit
1 mL	2.25 mL/min	–
2 mL	–	5.00 mL/min
5 mL	11.25 mL/min	12.50 mL/min
10 mL	22.50 mL/min	25.00 mL/min
20 mL	45.00 mL/min	50.00 mL/min
50 mL	112.50 mL/min	124.50 mL/min

29.12.3.4 Stirrer

The parameters for the stirrer are edited under **[Stirrer]**.

Stirrer

Selection of the MSB connector the stirrer is connected to. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4 off
Default value	1

off

No stirrer will be used.

Stirring rate

Setting the stirring rate. It can be set in steps of –15 to +15. The default setting **8** corresponds to 1000 rpm. The formula for calculating the rotational speed is specified in *chapter 32.4, page 445*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation
- "–": clockwise rotation

Input range	–15 to 15
Default value	8

Switch off automatically

on | off (Default value: **on**)

If this parameter is activated, the stirrer will be switched off automatically at the end of the titration, measurement, etc.

29.12.4 Liquid Handling (LQH)

The command **LQH** can be used to carry out Liquid Handling tasks with a dosing device of the type **Dosino 800**. The manifold options of the Dosino can be exploited in full in combination with a USB Sample Processor. The four ports of the dosing units can be used at will as output or input ports. This means that not only simple dosing and filling procedures are possible. Complex Liquid Handling tasks such as pipetting or sample transfers can be carried out without difficulty. Several LQH and automation commands are required for this purpose which are best combined into subsequences.

29.12.4.1 Liquid Handling parameters

The parameters for the Liquid Handling function are defined under **[Liquid Handling param.]**.

Function

Type of Liquid Handling function.

Selection	Dose	Fill	Aspirate	Eject	Exchange position	Change port	Compensate	End volume
-----------	------	------	----------	-------	-------------------	-------------	------------	------------

Dose

The given volume is dosed. The dosing cylinder neither is automatically filled before nor afterwards.

Fill

The dosing cylinder is filled. The valve disc then remains at the selected port.

Aspirate

Liquid is aspirated. The dosing cylinder neither is automatically filled before nor afterwards. The aspirating volume has to be aspirated with a single piston stroke.

Eject

The whole cylinder content is ejected. The piston is, in contrast to the function **End volume**, lowered to the stop point, i.e. past the maximum volume mark. This function should only be carried out for preparing a dosing unit before the pipetting procedure.

Exchange position

The dosing cylinder is filled. This means that, for example, air can be aspirated via Port 4. The valve disc is then rotated to Port 2 and the dosing drive can be removed from the dosing unit.

Change port

The stopcock is moved to the given port, but no piston movement takes place.

Compensate

Because of the fact that the dosing units are interchangeable, the coupling of the Dosino piston rod (spindle) exhibits a low mechanical tolerance that can be noticed when the piston changes its direction of movement. This tolerance can be compensated with this function. A short piston movement is first made in the same direction as the previous movement, which is then followed by a piston movement in the reverse direction.

End volume

The whole cylinder content is ejected. The piston is lowered to the maximum volume mark i.e. until 10'000 pulses have been dosed. This command should be used for pipetting functions for emptying the cylinder.

In-/Outlet

Selection of the port via which the Liquid Handling command is carried out.

Selection	Port 1	Port 2	Port 3	Port 4
Port 1	Default value, when Function = Dose, Aspirate and End volume .			
Port 2	Default value, when Function = Fill, Exchange position and Change port .			
Port 4	Default value, when Function = Eject and Compensate .			

Volume

This parameter can only be edited with **Function = Dose** and **Aspirate**.
Volume which is dosed or aspirated.

Function = Dose:

Input range	0.00000 to 99999.9 mL
Default value	1.00000 mL

Function = Aspirate:

Input range	0.00000 to 50.0000 mL
Default value	1.00000 mL

Flow rate

This parameter can only be edited with **Function = Dose, Fill, Aspirate, Eject, Exchange position** and **End volume**.

Rate at which it is dosed or filled.

Input range	0.01 to 166.00 mL/min
Selection	maximum
Default value	maximum

maximum

The maximum dosing and filling rates depend on the cylinder volume (see chapter 32, page 443).

29.12.4.2 Control instrument

The control instrument the determination is carried out with is selected under **[Control device]**. Control instruments are defined under **System ► Device manager**.



NOTE

This button will only be shown when, in addition to the Ti-Touch, a Sample Processor has been configured as control instrument.

Control device

Selection of the control instrument from the list of devices. Only those devices are displayed which are able to carry out the command.

Selection	Selection of configured control instruments
-----------	---

29.12.4.3 Dosing device

The parameters for the dosing device are edited under **[Dosing device]**.

Dosing device

Selection of the MSB connector to which the dosing device is connected. The selection does not depend on the number of MSB connectors of the control instrument.

Selection	1 2 3 4
Default value	1

Titrant

Selection of the titrant from the titrant list. We recommend always selecting the titrant. This ensures that the correct data (titer, concentration, etc.) is always used for the calculation. Titrants are defined under **System ► Titrants**. You can also enter a name which is not contained in the titrant list. When a determination is started there is a check whether the solution is contained in the list.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct titrant has been attached and whether

the type of dosing drive matches. For buret units without an integrated data chip, the cylinder volume and the type of dosing drive are checked. For the selected titrant, the validity of the titer, the working life of the titrant and the time interval for the GLP test of the buret unit are checked at the start of the determination.

Selection	Selection of configured titrants not defined
Default value	not defined

not defined

No check takes place.

29.13 Communication

The following communications commands are available:

Scanning remote lines (SCAN)	Defining remote signals which are awaited before the next method command is started.
Setting remote lines (CTRL)	Defining remote signals.
Receiving RS-232 commands (SCAN RS)	Defining RS-232 commands which are awaited before the next method command is started.
Sending RS-232 commands (CTRL RS)	Defining the RS-232 commands which are sent.

29.13.1 Scanning remote lines (SCAN)

The command **SCAN** can be used to define input signals at the remote interface which are awaited before the next method command is started.

Control device

This parameter will only be shown when, in addition to the Ti-Touch, a Sample Processor has been configured as control device. Control devices are defined under **System ► Device manager**.

Selection of the control device from the list of devices.

Selection	Selection of configured control devices
-----------	--

Remote Box

Selection of the MSB connector the Remote Box is connected to. The selection does not depend on the number of MSB connectors of the control instrument. The remote signal defined is requested on this Remote Box.

A message is displayed. You can select whether you want to continue with the determination or cancel the run.

In the determination data it will be documented, that the time interval has been expired.

The determination is stopped.

The command **CTRL** can be used to define output signals at the remote interface

This parameter will only be shown when, in addition to the Ti-Touch, a Sample Processor has been configured as control device. Control devices are defined under **System ► Device manager**.

Selection of the control device from the list of devices.

Selection	Selection of configured control devices
-----------	---

Selection of the MSB connector the Remote Box is connected to. The selection does not depend on the number of MSB connectors of the control instrument. The remote signal defined is sent via this Remote Box.

Selection	1 2 3 4
Default value	1

Selection of the signal out of the templates or entering the required bit pattern. Templates are defined under **System ► Templates ► Output lines**.



A line set active is not being reset automatically, not even at the end of the determination.

Entering a bit pattern:

- 0 = line inactive
- 1 = line active
- * = retain line status
- p = set pulse

The output lines are always numbered from right to left, i.e. with the signal *******1** line 0 is set. With a pulse, the length is set to 200 ms. If it is required to set pulses with other lengths, you have to define a corresponding template.



NOTE

We recommend masking lines that are of no interest or for which no defined condition can be predicted with an asterisk (*).

Entry	Bit patterns containing exactly 14 characters or a max. of 24 characters for the name of the template
Default value	*****
Selection	Selection of the templates defined

29.13.3 Scanning the RS-232 interface (SCAN RS)

The command **SCAN RS** can be used to define RS-232 commands which are awaited before the next method command is started.

Serial port

Selection of the serial interface the peripheral device is connected to. The RS-232 command defined is requested on this interface.

Selection	COM1 COM2 COM3 COM4 COM5 COM6 COM7 COM8
Default value	COM1

Character string

Input of the RS-232 command as character string. All characters of the ASCII code page can be used. Control characters (Esc, FF etc.) have to be entered as a three-place, decimal ASCII code, starting with a slash. Each command is automatically terminated with the ASCII characters **CR** and **LF**. The * can be used for one or more characters.

Entry	ASCII string with a max. of 24 characters
Default value	!*R

Timeout

When this time interval has expired without recognizing the RS-232 command, the action defined in the following is triggered.

Input range	0 to 999 s
Default value	0 s

Action

Selection of the action which is carried out when the time interval has expired.

Selection	Display message Document message Cancel determination
Default value	Display message For all three options it is documented in the determination data (see dialog More determination data / Messages), that the time interval has been expired.

Display message

A message is displayed. You can select whether you want to continue with the determination or cancel the run.

Document message

In the determination data it will be documented, that the time interval has been expired.

Cancel determination

The determination is stopped.

29.13.4 Defining RS-232 commands (CONTROL RS)

The command **CONTROL RS** can be used to define the RS-232 commands which are sent.

Serial port

Selection of the serial interface the peripheral device is connected to. The RS-232 command defined is sent via this interface.

Selection	COM1 COM2 COM3 COM4 COM5 COM6 COM7 COM8
Default value	COM1

Character string

Input of the RS-232 command as character string. All characters of the ASCII code page can be used. Control characters (Esc, FF etc.) have to be entered as a three-place, decimal ASCII code, starting with a slash. Each command is automatically terminated with the ASCII characters **CR** and **LF**.

Entry	ASCII string with a max. of 24 characters
Default value	&M;\$G

Current sample +

Starting from the current sample (defined by the sample variable) the rack is moved forward by the number of rack positions (**1...999**) entered in the field to the right.

Current sample -

Starting from the current sample (defined by the sample variable) the rack is moved backward by the number of rack positions (**1...999**) entered in the field to the right.

Next position

Starting from the current rack position, the rack is moved forward by one position.

Prev. position

Starting from the current rack position, the rack is moved backward by one position.

Calibration pos.

For automatic calibrations with a Sample Processor (see "Definition of the calibration positions", page 452).

rotate +

Moving the rack forward by a certain increment. The rotation increment is defined in the properties of the tower.

rotate -

Moving the rack backward by a certain increment. The rotation increment is defined in the properties of the tower.

swing +

Swinging the robotic arm outward by a certain increment (towards higher angular degrees). The swing increment is defined in the properties of the Swing Head.

swing -

Swinging the robotic arm towards the center of the rack by a certain increment (towards lower angular degrees). The swing increment is defined in the properties of the Swing Head.

Beaker test action

This parameter can only be edited with **Destination = Sample, Next position** or **Prev. position**.

Selection of the action which is carried out when the beaker sensor does not detect any vessel at the sample position moved to. In addition, enable the beaker sensor in the rack table.

Selection	Rotate rack Display message
Default value	Display message

Selection	1 2
Default value	1

1

View from the front, the right tower.

2

View from the front, the left tower.

Lift position

Selecting a predefined lift position or enter any lift position. The predefined lift positions are defined in the device properties of the Sample Processor. They can, however, also be assigned in the manual control.

Input range	0 to 235 mm
Selection	Work position Shift position Rinse position Special position Home position
Default value	Work position

Home position

The lift is being located at the upper stop position, i.e. at the position 0 mm.

Lift rate

Speed at which the lift is moved.

Input range	5 to 25 mm/s
Default value	20 mm/s

29.14.3 Controlling pumps (PUMP)

The **PUMP** command is used to control the pumps mounted on or connected to the tower of the sample changer.

Tower

Selection of the tower, with which the command is carried out. Both Tower 1 and Tower 2 can always be selected, even if your Sample Processor has only one tower.

Selection	1 2
Default value	1

1

View from the front, the right tower.

2

View from the front, the left tower.

Pump

Selection of the pump.

Sample variable

Selection	= + -
Default value	+

The sample variable corresponds to the number entered in the field **Value**.

The sample variable corresponds to the current value plus the number entered in the field **Value**.

The sample variable corresponds to the current value minus the number entered in the field **Value**.

Value by which the current sample variable is to be modified.



The sample variable can only be set on rack positions which are not defined as special beakers.

Input range	1 to 999
Default value	1

It is useful, particularly with more complex tasks, e.g. rinsing sequences or Liquid Handling, to combine individual command sequences into a single subsequence. Subsequences are managed the same as method commands. A maximum of 99 commands can be combined to form a single subsequence. There is no limit to the number of subsequences within one method.

Loading a saved subsequence or saving the current subsequence.

Defining various settings which correspond to the entire subsequence.

Insert a new method command. It is inserted before the selected command.



NOTE

Not all commands are available for subsequences. Commands which cannot be inserted into subsequences are disabled.

[Delete command]

Delete the selected method command.


[Edit command]

Edit the selected method command.

29.14.6.1 Loading / saving a subsequence

[Load/ Save] is used to load stored subsequences or to save new subsequences. This means that they are available for all methods and only have to be created once.

Loading a subsequence



[Own subseq.]

Opening the list of own subsequences.

[Load]

Loading the selected subsequence.

Saving a subsequence

Custom-made subsequences are saved in the internal memory.



File name

Entry	max. 32 characters
-------	--------------------

Saving the subsequence to the internal memory.

In this dialog, the settings can be defined which correspond to the entire subsequence.

Selection	Sample sequence Start sequence End sequence Stop sequence Cond. sequence
Default value	Sample sequence

The subsequence is carried out with every determination.

This subsequence is only carried out at the start of a sample series when the autostart counter = 1.

This subsequence is only carried out at the last sample of a series when the autostart counter has reached the setpoint value (number of autostarts).

The stop sequence is only carried out when canceling a method because of one of the following cases: Manual stop with `[ ]`, stop because of an error, stop via remote signal.

This subsequence is carried out directly before the conditioning. This gives you the possibility to start e.g. a connected Polytron® high-frequency homogenizer before the conditioning as well as to set its stirring rate. Use the command **CONTROL RS** for this.

Number of directly consecutive repetitions of the subsequence.

Input range	0 to 999 s (Increment: 1)
Selection	on off
Default value	on

on

The stirrer is switched on. The stirrer runs until it is explicitly switched off. If the determination is canceled with the key [], the stirrer is also switched off.

off

The stirrer is switched off.

Stirring rate

Setting the stirring rate. It can be set in steps of -15 to $+15$. The default setting **8** corresponds to 1000 rpm. The formula for calculating the rotational speed is specified in *chapter 32.4, page 445*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation
- "-": clockwise rotation

Input range	-15 to 15
Default value	8

29.15.2 Pausing the method run (WAIT)

The method run can be paused with the command **WAIT**.

Hold sequence

on | off (Default value: **off**)

If this parameter is activated, then the method run is paused until it is continued manually.

Waiting time

This parameter can only be edited when **Hold sequence** is deactivated.

The method run is automatically continued after this waiting time.

Input range	0 to 99999 s
Default value	30 s

Message

on | off (Default value: **off**)

If this parameter is activated, the text message defined in the following is displayed during the waiting time.

Message text

This parameter can only be edited when **Message** is activated.

Text which is displayed during the waiting time.

Entry	max. 28 characters
Default value	empty

29.15.3 Scan data (REQUEST)

The following data can be scanned in the method run with the command **REQUEST**:

- Sample data
 - Sample identification 1
 - Sample identification 2
 - Sample size (value and unit)
- Common variable

You can select thereby whether the sequence should be paused or resumed in the background.

Sample ident.

Selection of the sample identification that is queried in the method run.

Selection	off Identification 1 Identification 2 ID1 & ID2
Default value	off

Sample size

on | off (Default value: **on**)

If this parameter is activated, then the value for the sample size will be requested.

Sample size unit

on | off (Default value: **off**)

If this parameter is activated, then the unit for the sample size will be requested.

Common variable

Selection of the common variable which is queried in the method run.

Selection	CV01...CV25 off
Default value	off

30.2 Maintenance

Maintenance

The electronic and mechanical functional groups of Metrohm instruments can and should be checked by specialist personnel from Metrohm as part of a regular preventive maintenance schedule. Please ask your local Metrohm representative regarding the precise terms and conditions involved in concluding a corresponding maintenance agreement.

For detailed information on this topic, please visit www.metrohm.com.

Cleaning

Cleaning the surfaces of the instrument

Prerequisites

- The instrument is disconnected from the power grid.

- 1 Clean the surfaces with a damp cloth.



NOTE

Water or ethanol can be used as a cleaning medium.



NOTE

The connectors at the rear of the instrument must only be cleaned with a dry cloth.

31.3 Results/Statistics

Problem	Cause	Remedy
No statistics are carried out for a result.	<i>Statistics is not activated.</i>	In the method options, activate the option Statistics (see chapter 16.5, page 188).
	<i>No statistics variable has been assigned to the calculation.</i>	In the result options (command CALC), define a variable for the mean value (see "Variable for mean value", page 169).
	<i>Statistics is not activated.</i>	In the control dialog, activate the Statistics option (see chapter 18, page 202).
A result has been removed from the statistics, the mean value has not been reassigned to the common variable or to the TITER variable.	<i>Assignment is not carried out automatically with retroactive modifications.</i>	Recalculate the determination manually (see chapter 24, page 230).
The result is not displayed in the result table.	<i>The column display in the result table is incorrectly configured.</i>	Modify the settings for the display accordingly (see chapter 26.1, page 253).
	<i>The parameter Save result in result table is not activated.</i>	In the result options (command CALC), activate the parameter (see "Save result in result table", page 170).

31.4 Printing

Problem	Cause	Remedy
Tapping on the fixed key [] has no effect.	<i>The fixed key is disabled.</i>	In the Dialog options / Fixed keys dialog, activate the Print option (see chapter 7.2, page 44).
	<i>A determination is running.</i>	Wait until the determination is finished.

31.5 Manual control

Problem	Cause	Remedy
Tapping on the fixed key [] has no effect.	<i>The fixed key is disabled.</i>	In the Dialog options / Fixed keys dialog, activate the Manual control option (see chapter 7.2, page 44).
A button in the manual control is disabled.	<i>The required hardware is either not connected or incorrectly connected.</i>	1. Switch off the 916 Ti-Touch. 2. Connect the hardware correctly. 3. Switch the 916 Ti-Touch back on.
	<i>The required hardware is being used in the ongoing determination.</i>	Wait until the determination is finished.

31.6 Miscellaneous

Problem	Cause	Remedy
The 916 Ti-Touch cannot be switched on.	<i>The Ti-Touch is connected to an MSB socket on the control instrument.</i>	Plug in the Ti-Touch at the Controller connector (see manual for control instrument).

31.7 SET titration

Problem	Cause	Remedy
The titration will not be finished.	<i>The minimum dosing rate is too low.</i>	Define Titration rate = user and increase the minimum rate (Min. rate) (see chapter 29.3.2, page 326).
	<i>The stop criterion is unsuitable.</i>	Adjust the control parameters (see chapter 29.3.2, page 326): ▪ Increase the stop drift. ▪ Select a short delay time.
The sample is over-titrated.	<i>The control parameters are unsuitable.</i>	Adjust the control parameters (see chapter 29.3.2, page 326): ▪ Select Titration rate = slow. ▪ Define Titration rate = user and increase the control range. ▪ Define Titration rate = user and reduce the maximum rate (Max. rate).

Problem	Cause	Remedy
		- Define Titration rate = user and reduce the minimum rate (Min. rate). - Stir faster. - Arrange the electrode and buret tip to an optimum .
	<i>The electrode responds too slowly.</i>	Replace the electrode.
The titration time is too long.	<i>The control parameters are unsuitable.</i>	Adjust the control parameters (see chapter 29.3.2, page 326): - Select Titration rate = optimal or fast. - Define Titration rate = user and reduce the control range. - Define Titration rate = user and increase the maximum rate (Max. rate). - Define Titration rate = user and increase the minimum rate (Min. rate).
The results are spread widely.	<i>The minimum dosing rate is too high.</i>	Define Titration rate = user and reduce the minimum rate (Min. rate) (see chapter 29.3.2, page 326).
	<i>The control range is too small.</i>	Define Titration rate = user and increase the control range.
	<i>The electrode responds too slowly.</i>	Replace the electrode.

31.8 STAT titration

Problem	Cause	Remedy
The control point is not kept stable. The measured value is at times too high and at times too low. The controller "oscillates".	<i>The control parameters are unsuitable.</i>	Adjust the control parameters (see chapter 29.4.2, page 339): - Reduce the Max. rate and possibly also the Min. rate. - Increase control range.
	<i>The stirring is not optimal.</i>	Check and optimize efficiency of the stirring.
	<i>The arrangement of electrode and buret tip is not optimum.</i>	Arrange the electrode and buret tip to an optimum (see manual of the instrument).

32 Appendix

32.1 Manual Dosing Controller, optional accessory for manual titration




Figure 36 Manual Dosing Controller 6.2107.120 – Overview

1 Start/Dosing key [GO]

For starting the method and for dosing during the MAT command.

2 Status LED

Indicates whether the instrument is ready.

3 Cancel command button [FILL]

For canceling the running command.

4 Connector plug

With 25-pin D-sub plug.

32.3 Exchange unit

32.3.1 Maximum dosing and filling rate

The maximum dosing rate and maximum filling rate for the exchange unit depend on the cylinder volume:

Cylinder volume	Maximum rate
1 mL	3.00 mL/min
5 mL	15.00 mL/min
10 mL	30.00 mL/min
20 mL	60.00 mL/min
50 mL	150.00 mL/min

Independent of the cylinder volume, values ranging from 0.01 to 166.00 mL/min can always be entered. When the function is carried out the rate will be, if necessary, decreased automatically to the highest possible value.

32.3.2 Default parameters for preparing (PREP)

The **PREP** command is used to rinse and fill the cylinder and tubings of the exchange unit air bubble-free. You should carry out this function before the first determination or once per day.

If no titrant is selected in the command, preparing will be carried out with the following default parameters:

- The entire cylinder volume is dosed twice at the maximum dosing rate.
- The cylinder is filled at the maximum filling rate.

32.4 Stirring rate

The stirring rate can be adjusted in steps from –15 to +15.

The approximate rotational speed for the internal magnetic stirrer (depends on the product version) can be calculated with the following formula:

$$\text{Rotational speed/min (r/min)} = 125 \cdot \text{Stirring rate}$$

Example:

Configured stirring rate: 8

Rotational speed in revolutions per minutes = $125 \cdot 8 = 1,000$




Figure 37 Rotational speed depending on the stirring rate

The information on the separately connectable 802 propeller stirrer can be found in the "802 Stirrer" manual.

32.5 Balance

The sample size and the associated unit can be sent from a connected balance. The sample size is transmitted as a number with up to ten characters (including algebraic sign and decimal point).

Sample size and unit are sent as a single character string. They are separated by a space character. The string is terminated with the ASCII characters **CR** and **LF**.

If the balance sends a negative sample size (e.g. when you are reweighing a sample), then the algebraic sign is adopted. The algebraic sign is, however, ignored for the calculations.



NOTE

With some balances, the sample identification and the method can be sent in addition to the sample size.

Make sure that the balance does not send the sample size until the end.

Mettler AX

For the Mettler AX balance, the fields that contain the sample identification or the method must be designated as follows:

- Designation for the field with the method name: **METHOD**
- Designation for the field with sample identification 1: **ID1**
- Designation for the field with sample identification 2: **ID2**

32.6 Remote interface

The 6.2148.010 Remote Box allows devices to be controlled which cannot be connected directly to the MSB interface of the Ti-Touch.




Figure 38 Connectors of the Remote Box

1 Cable

For connecting to the Ti-Touch.

2 MSB connector

Metrohm Serial Bus. For connecting external dosing devices or stirrers.

3 Remote connector

For connecting instruments with a remote interface.

32.6.1 Pin assignment of the remote interface





Figure 39 Pin assignment of remote socket and plug

The above figure of the pin assignment of a Metrohm remote interface applies not only for the Remote Box, but also for all Metrohm devices with 25-pin D-Sub remote connector.

Inputs



approx. 50 k Ω Pull-up

$t_p > 20$ ms

active = low, inactive = high

Assignment	Pin No.	Function
Output 13	20	
0 volts / GND	14	
+5 volts	15	
0 volts / GND	25	

Table 19 Explanation of the individual functions

Function	Explanation
Start	The current method is started at the time of activation. $t_{\text{pulse}} > 100 \text{ ms}$
Stop	The current method is canceled (Stop) at the time of activation. $t_{\text{pulse}} > 100 \text{ ms}$
Quit	The current command in the determination run will be canceled at the time of activation. $t_{\text{pulse}} > 100 \text{ ms}$
Ready	The instrument is ready to receive a start signal.
Conditioning OK	The line is set when Conditioning with SET titration is at OK. The line remains set until the determination is started with [START] .
Determination	The instrument performs a data-generating determination.
EOD	End of Determination. Pulse ($t_{\text{pulse}} = 200 \text{ ms}$) after a determination or after a buffer during calibration using a Sample Processor.
Error	The line is set for error message display.
Warning	The line is set for warning message display.

32.7 Result variables as parameter setting

For most of the numerical parameters, the result of a calculation can also be utilized in place of a number. This result must be calculated in a preceding calculation command (command **CALC**). A typical application is the use of a relative start volume.

The following example is used to explain how this method is used to define a sample size-dependent start volume for a titration. Proceed as follows:

1 Create a new calculation

- Insert the calculation command **CALC** in front of the titration command.
- Enter the calculation formula, e.g. $R1 = C00*3$.

Edit command / Edit calculation

01

CALC

Calculation

Result name

Relative start volume

Calc. formula R1

C00*3

Decimal places

3

▼

Result unit

mL

▼

Note

Result variable

Result limits

Result options

2 Define a result variable as parameter setting

- In the titration command, enter the variable **R1** as the start volume.

Edit command / Start conditions

02 DET pH Dynamic pH titration

Start volume mL

Dosing rate mL/min

Pause s

The selection of the result variables is displayed in the number editor.



NOTE

If the method before the command with the result variables applied contains more than one calculation commands (possibly with result variables bearing the same name), then the result variable of the preceding calculation command will always be used.

Please note that the result must be within the input range of the parameter, as otherwise the determination run will be stopped when the parameter is accessed.

32.8 Electrode calibration with Sample Processors

Basic framework for a method for automatic calibrations

If the calibration of electrodes is carried out with a Sample Processor, then the basic framework of the method will appear as follows:

1. Subsequence with the commands required for the calibration (command **SUBSEQ**)
 - a. Command for change of buffer (command **MOVE**)
 - b. Command for moving the lift (command **LIFT**)
 - c. Calibration command (command **CAL pH**)
 - d. ...
2. ...

In the interest of general comprehensibility, other commands which are required are not listed here. A method template is available for an automatic pH calibration with subsequent measurement.

Command parameterization

The following special settings are required with the individual commands:

32.9 Stored buffer series for pH calibration

The temperature-dependent pH values of several common pH buffers are stored in the Ti-Touch for automatic buffer recognition during pH calibration. In addition to the Metrohm buffer solutions, other reference buffers and technical buffers are also included in the tables.

The following tables provide an overview of the stored pH(T) series.

pH values printed in **bold** are the values for the reference temperature of the respective buffer set.

pH values printed in *italics* are interpolated or extrapolated values, all the others correspond to the manufacturer's specifications.

32.9.1 Metrohm

Table 20 Metrohm buffer solutions

Temp. (°C)	Metrohm		
	pH 4.00	pH 7.00	pH 9.00
0	3.99	7.11	9.27
5	3.99	7.08	9.18
10	3.99	7.06	9.13
15	3.99	7.04	9.08
20	3.99	7.02	9.04
25	4.00	7.00	9.00
30	4.00	6.99	8.96
35	4.01	6.98	8.93
40	4.02	6.98	8.90
45	4.03	6.97	8.87
50	4.04	6.97	8.84
55	4.06	6.97	8.81
60	4.07	6.97	8.79
65	4.09	6.98	8.76
70	4.11	6.98	8.74
75	4.13	6.99	8.73
80	4.15	7.00	8.71
85	4.18	7.00	8.70

Temp. (°C)	NIST (according to DIN standard 19266, 2015-05)				
	pH 1.679	pH 4.005	pH 6.865	pH 9.180	pH 12.454
85	1.779	4.184	6.868	8.868	-
90	1.792	4.205	6.877	8.850	-
95	1.806	4.227	6.886	8.833	-

**NOTE****Update**

The values of the individual buffers with the corresponding temperatures are kept as up to date as possible.

However, they may be changed by the respective manufacturers.

The NIST buffers are identical to the buffer solutions that are used in Chinese pharmacopoeia.

32.9.3 DIN (according to DIN standard 19267, 2012)

Table 22 *DIN buffer solutions*

Temp. (°C)	DIN (according to DIN standard 19267, 2012-08)					
	pH 1.09	pH 3.06	pH 4.65	pH 6.79	pH 9.23	pH 12.75
0	1.08	-	4.67	6.89	9.48	-
5	1.08	-	4.66	6.86	9.43	-
10	1.09	3.10	4.66	6.84	9.37	13.37
15	1.09	3.08	4.65	6.82	9.32	13.15
20	1.09	3.07	4.65	6.80	9.27	12.96
25	1.09	3.06	4.65	6.79	9.23	12.75
30	1.10	3.05	4.65	6.78	9.18	12.61
35	1.10	3.05	4.66	6.77	9.13	12.44
40	1.10	3.04	4.66	6.76	9.09	12.29
45	1.10	3.04	4.67	6.76	9.04	12.13
50	1.11	3.04	4.68	6.76	9.00	11.98
55	1.11	3.04	4.69	6.76	8.97	11.84
60	1.11	3.04	4.70	6.76	8.92	11.69
65	1.11	3.04	4.71	6.76	8.90	11.56

Temp. (°C)	Fisher			
	pH 2.00	pH 4.00	pH 7.00	pH 10.00
60	-	4.09	6.98	9.70
65	-	4.11	6.99	9.68
70	-	4.13	7.00	9.65
75	-	4.14	7.02	9.63
80	-	4.16	7.03	9.62
85	-	4.18	7.06	9.61
90	-	4.21	7.08	9.60
95	-	4.23	7.11	9.60

Table 24 Fisher buffer solutions (valid from firmware version 5.916.0044)

Temp. (°C)	Fisher			
	pH 2.00	pH 4.00	pH 7.00	pH 10.00
0	-	4.01	7.13	10.34
5	1.98	3.99	7.10	10.26
10	1.98	4.00	7.07	10.19
15	2.02	3.99	7.05	10.12
20	2.00	4.00	7.02	10.06
25	2.00	4.00	7.00	10.00
30	2.00	4.01	6.99	9.94
35	2.02	4.02	6.98	9.90
40	2.01	4.03	6.97	9.85
45	2.01	4.04	6.97	9.81
50	2.01	4.06	6.97	9.78
55	-	4.07	6.97	9.74
60	-	4.09	6.98	9.70



NOTE

Update

The values of the individual buffers with the corresponding temperatures are kept as up to date as possible.

However, they may be changed by the respective manufacturers.

32.9.5 Honeywell Fluka

Table 25 Buffer solutions Honeywell Fluka

Temp. (°C)	Fluka		
	pH	pH	pH
	4.000	7.000	9.000
0	4.030	7.130	9.240
5	4.025	7.090	9.175
10	4.020	7.050	9.110
15	4.010	7.020	9.055
20	4.000	7.000	9.000
25	4.000	6.990	8.965
30	4.000	6.980	8.930
35	4.000	6.975	8.895
40	4.000	6.970	8.860
45	4.000	6.965	8.830
50	4.000	6.960	8.800
55	4.000	6.960	8.775
60	4.000	6.960	8.750
65	4.000	6.965	8.730
70	4.000	6.970	8.710
75	4.000	6.975	8.690
80	4.000	6.980	8.670
85	4.000	6.990	8.655
90	4.000	7.000	8.640

**NOTE****Update**

The values of the individual buffers with the corresponding temperatures are kept as up to date as possible.

However, they may be changed by the respective manufacturers.

32.9.6 Mettler Toledo

Table 26 *Mettler Toledo buffer solutions*

Temp. (°C)	Mettler Toledo				
	pH 2.00	pH 4.01	pH 7.00	pH 9.21	pH 11.00
0	2.03	4.01	7.12	9.52	11.90
5	2.02	4.01	7.09	9.45	11.72
10	2.01	4.00	7.06	9.38	11.54
15	2.00	4.00	7.04	9.32	11.36
20	2.00	4.00	7.02	9.26	11.18
25	2.00	4.01	7.00	9.21	11.00
30	1.99	4.01	6.99	9.16	10.82
35	1.99	4.02	6.98	9.11	10.64
40	1.98	4.03	6.97	9.06	10.46
45	1.98	4.04	6.97	9.03	10.28
50	1.98	4.06	6.97	8.99	10.10
55	1.98	4.08	6.98	8.96	-
60	1.98	4.10	6.98	8.93	-
65	1.98	4.13	6.99	8.90	-
70	1.99	4.16	7.00	8.88	-
75	1.99	4.19	7.02	8.85	-
80	2.00	4.22	7.04	8.83	-
85	2.00	4.26	7.06	8.81	-
90	2.00	4.30	7.09	8.79	-
95	2.00	4.35	7.12	8.77	-



NOTE

Update

The values of the individual buffers with the corresponding temperatures are kept as up to date as possible.

However, they may be changed by the respective manufacturers.

32.9.7 Beckmann

Table 27 Beckmann buffer solutions (valid up to and including firm-ware version 5.916.0043)

Temp. (°C)	Beckmann		
	pH	pH	pH
	4.00	7.00	10.01
0	4.00	7.12	10.32
5	4.00	7.09	10.25
10	4.00	7.06	10.18
15	4.00	7.04	10.12
20	4.00	7.02	10.06
25	4.00	7.00	10.01
30	4.01	6.99	9.97
35	4.02	6.99	9.93
40	4.03	6.98	9.89
45	4.05	6.98	9.86
50	4.06	6.97	9.83
55	4.08	6.98	-
60	4.09	6.98	-
65	4.11	6.99	-
70	4.12	6.99	-
75	4.14	7.00	-
80	4.16	7.00	-
85	4.18	7.01	-
90	4.19	7.02	-
95	4.21	7.03	-

Table 28 Beckmann buffer solutions (valid from firmware version 5.916.0044)

Temp. (°C)	Beckmann		
	pH 4.01	pH 7.00	pH 10.01
0	4.00	7.12	10.32
5	4.00	7.09	10.25
10	4.00	7.06	10.18
15	4.00	7.04	10.12
20	4.00	7.01	10.06
25	4.01	7.00	10.01
30	4.01	6.99	9.97
35	4.02	6.99	9.93
40	4.03	6.97	9.89
45	4.05	6.97	9.86
50	4.06	6.97	9.83



NOTE

Update

The values of the individual buffers with the corresponding temperatures are kept as up to date as possible.

However, they may be changed by the respective manufacturers.

32.9.8 Radiometer Analytical

Table 29 Radiometer Analytical buffer solutions (valid up to and including firmware version 5.916.0043)

Temp. (°C)	Radiometer Analytical			
	pH 1.679	pH 4.005	pH 7.000	pH 9.180
0	1.666	4.000	7.118	9.464
5	1.668	3.998	7.087	9.395
10	1.670	3.997	7.059	9.332
15	1.672	3.998	7.036	9.276
20	1.675	4.001	7.016	9.225

Temp. (°C)	Radiometer Analytical			
	pH 1.679	pH 4.005	pH 7.000	pH 9.180
55	1.715	4.064	6.965	8.985
60	1.723	4.080	6.968	8.962
65	1.732	4.097	6.974	8.941
70	1.743	4.116	6.982	8.921
75	1.754	4.137	6.992	8.900
80	1.765	4.159	7.004	8.884
85	1.778	4.183	7.018	8.867
90	1.792	4.208	7.034	8.850

**NOTE****Update**

The values of the individual buffers with the corresponding temperatures are kept as up to date as possible.

However, they may be changed by the respective manufacturers.

32.9.9 Baker

Table 31 Baker buffer solutions

Temp. (°C)	Baker			
	pH 4.00	pH 7.00	pH 9.00	pH 10.00
0	4.00	7.13	9.23	10.30
5	4.00	7.09	9.17	10.24
10	4.00	7.05	9.10	10.17
15	4.00	7.03	9.05	10.11
20	4.00	7.00	9.00	10.05
25	4.00	6.98	8.96	10.00
30	4.01	6.98	8.91	9.96
35	4.02	6.98	8.88	9.93
40	4.03	6.97	8.84	9.89
45	4.04	6.97	8.81	9.86
50	4.05	6.96	8.78	9.82

Temp. (°C)	Baker			
	pH	pH	pH	pH
	4.00	7.00	9.00	10.00
55	4.07	6.96	8.76	9.79
60	4.08	6.96	8.73	9.76
65	4.10	6.97	8.71	9.74
70	4.12	6.97	8.69	9.72
75	4.14	6.98	8.68	9.70
80	4.16	6.98	8.66	9.68
85	4.19	6.99	8.64	9.66
90	4.21	7.00	8.62	9.64



NOTE

Update

The values of the individual buffers with the corresponding temperatures are kept as up to date as possible.

However, they may be changed by the respective manufacturers.

32.9.10 Hamilton DURACAL

Table 32 Hamilton DURACAL buffer solutions

Temp. (°C)	Hamilton DURACAL			
	pH	pH	pH	pH
	4.01	7.00	9.21	10.01
0	-	-	-	-
5	4.01	7.09	9.45	10.19
10	4.00	7.06	9.38	10.15
15	4.00	7.04	9.32	10.11
20	4.00	7.02	9.26	10.06
25	4.01	7.00	9.21	10.01
30	4.01	6.99	9.16	9.97
35	4.02	6.98	9.11	9.92
40	4.03	6.97	9.06	9.86
45	4.04	6.97	9.03	9.83
50	4.06	6.97	8.99	9.79



NOTE

Update

The values of the individual buffers with the corresponding temperatures are kept as up to date as possible.

However, they may be changed by the respective manufacturers.

32.9.11 Precisa

Table 33 Precisa buffer solutions

Temp. (°C)	Precisa		
	pH 4.00	pH 7.00	pH 9.00
0	3.99	7.11	9.27
5	3.99	7.08	9.18
10	3.99	7.06	9.13
15	3.99	7.04	9.08
20	3.99	7.02	9.04
25	4.00	7.00	9.00
30	4.00	6.99	8.96
35	4.01	6.98	8.93
40	4.02	6.98	8.90
45	4.03	6.97	8.87
50	4.04	6.97	8.84
55	4.06	6.97	8.81
60	4.07	6.97	8.79
65	4.09	6.98	8.76
70	4.11	6.98	8.74
75	4.13	6.99	8.73
80	4.15	7.00	8.71
85	4.18	7.00	8.70
90	4.20	7.01	8.68
95	4.23	7.02	8.67



NOTE

Update

The values of the individual buffers with the corresponding temperatures are kept as up to date as possible.

However, they may be changed by the respective manufacturers.

32.9.12 Merck CertiPUR 20 / Titrisol

Table 34 Merck CertiPUR 20 / Titrisol buffer solutions (valid up to and including firmware version 5.916.0043)

Temp. (°C)	Merck CertiPUR 20 / Titrisol				
	pH	pH	pH	pH	pH
	2.00	4.00	7.00	9.00	12.00
0	2.01	4.05	7.13	9.24	12.58
5	2.01	4.04	7.07	9.16	12.41
10	2.01	4.02	7.05	9.11	12.26
15	2.00	4.01	7.02	9.05	12.10
20	2.00	4.00	7.00	9.00	12.00
25	2.00	4.01	6.98	8.95	11.88
30	2.00	4.01	6.98	8.91	11.72
35	2.00	4.01	6.96	8.88	11.67
40	2.00	4.01	6.95	8.85	11.54
45	2.00	4.00	6.95	8.82	11.435
50	2.00	4.00	6.95	8.79	11.33
55	2.00	4.00	6.95	8.76	11.185
60	2.00	4.00	6.96	8.73	11.04
65	2.00	4.00	6.96	8.715	10.97
70	2.01	4.00	6.96	8.70	10.90
75	2.01	4.00	6.96	8.68	10.80
80	2.01	4.00	6.97	8.66	10.70
85	2.01	4.00	6.98	8.65	10.59
90	2.01	4.00	7.00	8.64	10.48
95	-	4.00	7.02	-	-

Table 35 Merck CertiPUR 20 / Titrisol buffer solutions (valid from firmware version 5.916.0044)

Merck CertiPUR 20 / Titrisol					
Article no.	109433	109435	109477	109476	109462
Temp. (°C)	pH 2.00	pH 4.00	pH 7.00	pH 9.00	pH 11.00
0	2.01	4.05	7.13	9.24	11.45
5	2.01	4.04	7.07	9.16	11.32
10	2.01	4.02	7.05	9.11	11.20
15	2.00	4.01	7.02	9.05	11.10
20	2.00	4.00	7.00	9.00	11.00
25	2.00	4.01	6.98	8.95	10.90
30	2.00	4.01	6.98	8.91	10.81
35	2.00	4.01	6.96	8.88	10.72
40	2.00	4.01	6.95	8.79	10.64
45	2.00	4.00	6.95	8.82	10.56
50	2.00	4.00	6.95	-	10.48



NOTE

Update

The values of the individual buffers with the corresponding temperatures are kept as up to date as possible. However, they may be changed by the respective manufacturers.

32.9.13 Merck CertiPUR 25

Table 36 Merck CertiPUR 25 buffer solutions (valid up to and including firmware version 5.916.0043)

Merck CertiPUR (25°C)				
Temp. (°C)	pH 4.00	pH 7.00	pH 9.00	pH 10.00
0	-	-	-	-
5	4.05	7.09	9.22	10.22
10	4.04	7.08	9.16	10.16
15	4.02	7.04	9.10	10.10

Temp. (°C)	Merck CertiPUR (25°C)			
	pH	pH	pH	pH
	4.00	7.00	9.00	10.00
20	4.01	7.02	9.06	10.05
25	4.00	7.00	9.00	10.00
30	3.99	6.98	8.98	9.94
35	3.98	6.98	8.93	9.90
40	3.98	6.97	8.89	9.86
45	3.98	6.97	8.86	9.80
50	3.98	6.97	8.84	9.73

Table 37 Merck CertiPUR 25 buffer solutions (valid from firmware version 5.916.0044)

Merck CertiPUR (25°C)				
Article no.	109445	109407	109408	109409
Temp. (°C)	pH 4.00	pH 7.00	pH 9.00	pH 10.00
0	-	-	-	-
5	4.05	7.09	9.22	10.22
10	4.04	7.08	9.16	10.16
15	4.02	7.04	9.10	10.10
20	4.01	7.02	9.05	10.05
25	4.00	7.00	9.00	10.00
30	3.99	6.98	8.96	9.94
35	3.98	6.98	8.93	9.90
40	3.98	6.97	8.89	9.86
45	3.98	6.97	8.86	9.80
50	3.98	6.97	8.84	9.73



NOTE

Update

The values of the individual buffers with the corresponding temperatures are kept as up to date as possible. However, they may be changed by the respective manufacturers.


32.10 Using AuditTrailViewer

When you work in accordance with the *FDA Guidance 21 CFR Part 11*, logs with detailed user activities (Audit Trail) are important. The *AuditTrailViewer* supports you in your efforts to compile, filter and export logs of this type as TXT files.

First you must install the *AuditTrailViewer* on your PC. Afterwards, you must create a backup of your Ti-Touch (see chapter 12.3, page 143). The Audit Trail is generated automatically when you create the backup. The file name of the Audit Trail is **log.madt**. In the following you will find detailed information about how you must proceed.

32.10.1 Installing AuditTrailViewer

You will find the file for installing the *AuditTrailViewer* on the USB flash drive provided.



Proceed as follows:

- 1 Double-click on the file **SetupAuditTrailViewer 1.2.3.exe**.

The dialog **Welcome to the AuditTrailViewer Setup Wizard** is displayed.

- 2 Click on **[Next >]**.

The **License Agreement** dialog is displayed.

- 3 Click on **[I Agree]** to accept the contractual conditions.

The **Choose Install Location** dialog is displayed. A memory location for the *AuditTrailViewer* is suggested in this dialog.

- 4 If you accept the suggested memory location, click on **[Next >]**.

or

Select a different memory location and click on **[Next >]**.

The **Confirm Installation** dialog is displayed.

- 5 Click on **[Install]** to install the *AuditTrailViewer*.

As soon as the installation is completed, the dialog **Installation Complete** is displayed.

- 6** Click on **[Next >]**.

The dialog **Completing the AuditTrailViewer Setup Wizard** is displayed.

- 7** If the *AuditTrailViewer* should not be started immediately, deactivate the check box **Run AuditTrailViewer**.

- 8** Click on **[Finish]**.

If the preceding instruction step has been skipped, then the *AuditTrailViewer* will start.

32.10.2 Opening the Audit Trail

Create a backup of the data and of the settings of your Ti-Touch (see *chapter 12.3, page 143*).

- 1 Open the backup on your PC.

- 2  916
-  Backup
 -  AuditTrail —  log.madt
 -  Main group
 -  RackDefinitions
 -  Files

Open the file **log.madt** with the *AuditTrailViewer*.

AuditTrailViewer dialog window

File Filter Info					
No.	Date	User	Category	Action	Details
1	2011-10-31 11:30:22	Johnson	Method	New	01 Dynamic Titration pH
2	2011-10-31 11:30:34	Johnson	Method	Delete command	New method V0 03 REPORT
3	2011-10-31 11:30:40	Johnson	Method	Edit	New method V0 01 DET pH Stirrer off
4	2011-10-31 11:30:43	Johnson	Method	Edit	New method V0 01 DET pH Titrant not defined
5	2011-10-31 11:30:43	Johnson	Method	Edit	New method V0 01 DET pH Titrant
6	2011-10-31 11:30:48	Johnson	Method	Load result template	New method V0 Content (g/L)
7	2011-10-31 11:30:52	Johnson	Method	Edit	New method V0 Calc. formula EP1*CONC*TITER*3/C00
8	2011-10-31 11:31:08	Johnson	System	Message	003-908 Number of autostarts Yes/OK
9	2011-10-31 11:31:14	Johnson	Method	Start	Start key pressed
10	2011-10-31 11:31:15	Johnson	Method	Start	New method V0 1.0 g
11	2011-10-31 11:31:28	Johnson	Method	Manual stop	New method V0
12	2011-10-31 11:31:36	Johnson	System	Message	002-907 Method modified Yes/OK
13	2011-10-31 11:31:36	Johnson	Security	Logout	
14	2011-10-31 11:31:55	Meier	Security	Login message	Wrong password
15	2011-10-31 11:31:57	Johnson	System	Message	002-102 Wrong password Yes/OK
16	2011-10-31 11:31:58	Meier	Security	Log in	
17	2011-10-31 11:32:27	Meier	Method	Load	KFT 2 V1 Internal memory
18	2011-10-31 11:32:32	Meier	Security	Logout	
19	2011-10-31 11:32:45	Meier	Security	Login message	Wrong password
20	2011-10-31 11:32:47	Meier	System	Message	002-102 Wrong password Yes/OK
21	2011-10-31 11:33:01	Meier	Security	Change password	
22	2011-10-31 11:33:03	Meier	Security	Login message	Wrong password
23	2011-10-31 11:33:05	Meier	System	Message	002-102 Wrong password Yes/OK
24	2011-10-31 11:33:08	Meier	Security	Log in	
25	2011-10-31 11:33:12	Meier	Method	Start	Start key pressed
26	2011-10-31 11:33:13	Meier	Method	Start	KFT 2 V1 1.0 g
27	2011-10-31 11:33:14	Meier	System	Message	009-111 Stirrer missing Yes/OK
28	2011-10-31 11:33:15	Meier	Security	Logout	
29	2011-10-31 11:33:27	Chang	Security	Log in	
30	2011-10-31 11:33:32	Chang	Method	New	01 Dynamic Titration pH
31	2011-10-31 11:33:38	Chang	Method	Delete command	New method V0 03 REPORT
32	2011-10-31 11:33:43	Chang	Method	Edit	New method V0 01 DET pH Start volume 5 mL
33	2011-10-31 11:33:48	Chang	Method	Load result template	New method V0 Content (g/L)

32.10.3 Contents of the Audit Trail

In the following you will find a description of the dialog window and of the menu bar of the *AuditTrailViewer*.

Table 38 Dialog window

Column	Content
No.	Each entry is numbered in sequence.

Column	Content
Symbol	Classification of the entry:  Actions that are neither relevant to safety nor alter the determination data.  Actions, e.g. changes to the login options and recalculation of determinations.  Errors that occur, e.g. entry of an incorrect password.
Date	Precise time of the event.
User	User who triggered the action.
Category	Category to which the entry belongs.
Action	Designation of the action.
Details	Details of the action.

Menu bar	
 Print	Print Audit Trail. An Audit Trail can be printed by means of the  symbol.
 Update	Update view. The list is updated automatically only at the time the dialog window is opened.
 Quick filter	Filter the Audit Trail.
 Show all	Display all entries once again.

32.10.4 Filtering the Audit Trail

The entries in the Audit Trail can be filtered. You can set the following filter criteria:

- **Date**
- **User**
- **Category**
- **Action**

Proceed as follows:

- 1 Click in a cell that contains the desired filter criterion (e.g. **Category = Method**).

- 2 Click on the  icon or the menu item **Filter ► Quick filter**.

The only entries which will be displayed are those which correspond to the desired filter criterion.

- 3 Click on the  icon or the menu item **Filter ► Show all**.

All of the entries are displayed once again.

32.10.5 Exporting the Audit Trail

You can export the Audit Trail as .txt file. This way, the Audit Trail can be archived and is available as required for inspections.

Proceed as follows:

- 1 Use the menu item **File ► Export** to open the dialog window for the export.
- 2 Enter the path and file name with the extension **.txt**.
- 3 Click on **[Open]** and then on **[Export]**.

The Audit Trail is exported.

32.11 Diagnosis

System ► Diagnosis

The electronic and mechanical functional groups of Metrohm devices can and should be checked by specialist personnel from Metrohm as part of a regular maintenance schedule. Please ask your regional Metrohm representative regarding the precise terms and conditions involved in concluding a corresponding maintenance agreement.

The following functions are available:

- **LCD test**
Check the display for faulty pixels (*see chapter 32.11.1, page 474*).
- **Format storage medium**
Format the external storage medium (*see chapter 32.11.2, page 474*).
- **Remove storage medium**
Remove the external storage medium safely (*see chapter 32.11.3, page 475*).
- **Touch adjustment**
Adjust the touch-sensitive screen (*see chapter 32.11.4, page 475*).

- **Touch screen test**
Test the function of the touch-sensitive screen (*see chapter 32.11.5, page 476*).
- **Software update**
Update the software of Ti-Touch and of the control instruments (*see chapter 32.11.6, page 477*).
- **Service**
Special functions for the service technician (*see chapter 32.11.7, page 480*).

32.11.1 LCD test

System ► Diagnosis ► LCD test

You can use the LCD test to check the display for faulty pixels. To do this, various test images are displayed one after the other.



NOTE

The fixed key is used to display the preceding test image; the test can be stopped at any time with the fixed key.

Proceed as follows:

- 1 ■ Tap on **[LCD test]**.

The entire display turns white.

2. Check the display for pixel errors and other irregularities.
 - Continue the test with [▶].

Other colors and patterns will be displayed one after the other.

- Repeat step 2 until the **System / Diagnosis** dialog is displayed again.

- 4 ■ Contact your Metrohm representative in the event of faults or irregularities.

32.11.2 Formatting an external storage medium

System ► Diagnosis ► Format storage medium

With this function, you have the option of formatting an external storage medium directly on the 916 Ti-Touch (quick formatting).

The following file systems can be selected:

- **FAT**

- FAT32
- ExFAT



CAUTION

If you format the storage medium with the **ExFAT** file system, it is possible that it will no longer be recognized by a PC. Read the pertinent information from Microsoft Support at <http://support.microsoft.com>.

Do not use this file system unless you need to be able to save more than 999 files in a single group.

32.11.3 Removing an external storage medium

System ► Diagnosis ► Remove storage medium

As long as the stored data is not being accessed (reading/saving data), you can plug in and disconnect the storage medium without difficulty at any time. The **[Remove storage medium]** function offers additional protection. This ensures that the storage medium cannot be removed while any data is being transferred.

32.11.4 Adjusting the touch screen

System ► Diagnosis ► Touch adjustment

With time the touch-sensitive screen may not register touches at their precise position. You can readjust the screen in this dialog.

Proceed as follows:

1 In the main dialog, tap on **[System]**.

2 Tap on **[Diagnosis]**.

3 Tap on **[Touch adjustment]**.

A message in English and a crosshair will appear in the center of the screen.

4 Touch a stylus (a special pen-shaped instrument for operating devices with touch screens) to the center of the crosshair and maintain this position until the cross hair moves to a new position on the screen.

5 Repeat this process as long as the crosshair continues to be repositioned on the screen.

A message in English will appear.

- 6 Tap at any location on the screen to confirm the adjustment.
- 7 Press [] to finish the adjustment.

32.11.5 Testing the touch screen

System ► Diagnosis ► Touch screen test

You can use this test to check the function of the touch-sensitive monitor screen.

A calibration is not possible with this test. If the Touch screen does not work correctly anymore, it will need to be replaced. For that please contact your Metrohm representative.



NOTE

You can use the [] fixed key to stop the test at any time.

Proceed as follows:


- 1** Tap on **[Touch screen test]**.

A red square is displayed at the upper left.

- 2** Tap on the square as close to the center as possible.

The actual detected point of contact will be displayed by a crosshair, and the deviation from the center will also be displayed on the lower edge of the image.

The next red square will also be displayed.



- 3 Repeat step 2 until all of the areas in the display have been tested.
A list with fixed keys will appear.
- 4 Tap on all of the fixed keys one after the other in any sequence.
The result is displayed in each case with a short information text, e.g.
[Print] key OK.
- 5 Finish Touch screen test with [].

32.11.6 Software update (loading program versions and language files)

System ► Diagnosis ► Software update

New program versions or language files can be loaded from an external storage medium (e.g. USB flash drive). The corresponding files must be stored on the storage medium in the **916/SwUpdates** directory (see *"Directory structure", page 142*).



NOTE

The software update to version 5.916.0046 is possible only if version 5.916.0045 is installed. If a software version older than 5.916.0045 is used, then the control software must first be updated to version 5.916.0045.



NOTE

If you wish to update a connected Sample Processor together with the Ti-Touch, then do not fail to update the control software of the Ti-Touch first and then afterwards the device program (firmware) of the Sample Processor.



NOTE

The memory will be reformatted when updating from software version 5.916.0043.

Before carrying out the software update, make a backup in order to save your data and the system settings (see *chapter 12.3, page 143*).

Use a separate USB flash drive for the software update. It must be removed from the instrument after the successful software update.

Proceed as follows:


1 Copying files to an external storage medium



NOTE

In the case of the program versions and language files for the Ti-Touch, it is imperative that both the BIN file and the associated MBIN file be copied.

- Copy the files to the "SwUpdates" directory.
- Copy also the "MetrohmTools" folder with the "MUpdateTool.exe" file to the external storage medium. Make sure that the file is stored at the top level of the external storage medium. This also applies to the "916" folder.



For the program versions 5.916.0043 to 5.916.0045, an additional folder that is on the same level as MetrohmTools with the name SRVPRF900 containing the file St900_SrPrf.cfg is required

If this structure is not maintained, then the files will not be found, because the directories are accessed directly during the update process.

- Plug in the external storage medium at the Ti-Touch.

2 Opening the update dialog

- Under **System ► Diagnosis**, tap on the button **Software update**.

3 Selecting the instrument

- Tap on the **Control device** list box and select the instrument to be updated.

32.12 Arithmetic algorithms in the Ti-Touch

Some of the algorithms and numerical formats used in the device software are explained in this chapter.

Numerical format

All floating point calculations use the IEEE 754 (1985) numerical format, either in "single precision" (32 Bit) or "double precision" (64 Bit).

Rounding off process

Measured values and results are rounded to the defined number of decimal places (commercial rounding, in accordance with the US Pharmacopeia USP). If the digit at the first dropped decimal place is **1, 2, 3 or 4**, then it will be rounded off; if this digit is **5, 6, 7, 8 or 9**, then it will be rounded up. Negative digits will be rounded in accordance with their amount, i.e. away from zero.

Examples:

2.33 yields **2.3**

2.35 yields **2.4**

2.47 yields **2.5**

-2.38 yields **-2.4**

-2.45 yields **-2.5**

Statistics

The arithmetic mean value and the absolute and relative standard deviations of results R, endpoints EP and variable C are calculated:

You can statistically evaluate a maximum of nine results ($1 \leq k \leq 9$) calculated in a determination. A statistical series can contain a maximum of 20 determinations ($1 \leq n \leq 20$).

The following convention applies to the subsequent formulas:

$1 \leq n \leq 20$ and $1 \leq k \leq 9$.

Mean value:

$$\bar{x}_k = \frac{1}{n} \cdot \sum_{i=1}^n R_{k,i}$$

Absolute standard deviation:

$$Sabs_k = + \sqrt{\frac{\sum_{i=1}^n (R_{k,i} - \bar{x}_k)^2}{n-1}}$$

Relative standard deviation (in %):

$$Srel_k = 100 \cdot \frac{Sabs_k}{\bar{x}_k}$$

33 Technical specifications

33.1 Touch screen

<i>Display</i>	VGA color display (640 pixels x 320 pixels), diagonal approx. 5.7"
<i>Touch panel</i>	resistive
<i>Resistance to chemicals</i>	resistant to the following chemicals (no visible changes after 24 h of duration of action): ▪ without protective cover: – Methanol – Toluene – Ethyl acetate – Acetone – Dichloromethane – Hydrochloric acid $c(\text{HCl}) = 2 \text{ mol/L}$ ▪ with protective cover (order number 6.2723.300): – Methanol – Propanol – Toluene – Xylene – Acetone – Chloroform – Formamide – Sulfuric acid $c(\text{H}_2\text{SO}_4) = 2 \text{ mol/L}$ – Hydrochloric acid $c(\text{HCl}) = 2 \text{ mol/L}$ – Caustic soda $c(\text{NaOH}) = 2 \text{ mol/L}$ – Composite 5 – 1-butanol – 1-hexanol – Decanol

<i>Measuring mode</i>	Determination with adjustable polarization voltage.
<i>Upol</i>	
<i>Polarization voltage</i>	–1200 to +1200 mV (increment: 10 mV) –1250 to –1210 mV / +1210 to +1250 mV: non-guaranteed values, dependent on reference voltage +2.5 V
<i>Measuring range</i>	–120 to +120 μ A
<i>Resolution</i>	0.01 μ A
<i>Measuring accuracy</i>	–

33.2.3 Temperature

A measuring input (**Temp.**) for temperature sensors of the Pt1000 or NTC type with automatic temperature compensation.

R (25 °C) and B value can be configured for NTC sensors.

<i>Measuring range</i>	
<i>Pt1000</i>	–150 to +250 °C
<i>NTC</i>	–5 to +250 °C (For an NTC sensor with R (25 °C) = 30,000 Ω and B (25/50) = 4,100 K)
<i>Resolution</i>	
<i>Pt1000</i>	0.1 °C
<i>NTC</i>	0.1 °C
<i>Measuring accuracy</i>	
<i>Pt1000</i>	± 0.2 °C (Applies for measuring range –20 to +150 °C)
<i>NTC</i>	± 0.6 °C (Applies for measuring range +10 to +40 °C)

33.3 Built-in magnetic stirrer (depending on the product version)

Rotational speed

Maximum rotational speed $\pm(1700-1900)$ rpm

Rotational speed settings ± 15 steps

Rotational speed change per step $\pm(115-125)$ rpm

Power consumption	4 W
-------------------	-----

Fuse Electronic overload protection

33.4 Interfaces

USB connector Type A, for connecting USB devices.

MSB connector For connecting dosing devices, one additional stirrer or a Remote Box.

iConnect connector For connecting an 854 iConnect with iTrode.

*Stirrer connector
(depends on prod-
uct version; only
for instruments
without built-in
magnetic stirrer)*

Ethernet connector	For connection to a data network (LAN).
--------------------	---

33.5 Energy supply

Device

Inlet

Nominal voltage 24 V DC

Power consumption (typical) 10 W

External power supply unit 6.2164.010

Inlet

Nominal voltage range 100–240 VAC ($\pm 10\%$)

Nominal frequency range 50–60 Hz ($\pm 3\%$)

Max. current 1.5 A

Outlet

Nominal voltage 24 V DC

Max. current 2.7 A

Power output 65 W

33.6 Declarations of Conformity

Up-to-date information on the *EU conformity* and the *requirements for electromagnetic compatibility* for your product can be found on the Internet under <https://www.metrohm.com/en-us/products-overview/>.

Using the article number, you can download the *Declaration of conformity for ...* under the documents of the product.

Base Steel, stainless, coated

33.10 Storage capacity

Random access memory 128 MB

Internal memory 1 GB of which 900 MB available.

Required memory 3–6 KB
per method

800 Dosino	17
801 Stirrer	18
803 Ti Stand	18
804 Ti Stand	18
805 Dosimat	17

Acoustic signals	60
ADD	411
Audit Trail	58
Archive	473
Export	473
Filter	472
Open	470
AuditTrailViewer	
Installing	469
Automatic system test	146
Automation	424
Axial distance	111

Backup	143
Interval	155
Balance	20, 129, 446
Barcode reader	
Connect	21
Beaker radius	
Sample vessel	118
Special beaker	123
Beaker sensor	
Sample vessel	119
Special beaker	123
BEEP	435
Break point	
Evaluation	385
Buffer series	453

CAL	368
CALC	388
CALC LIVE	396
Calculation	388
Add	389
Delete	389
Edit	389
Load template	389
Note	393
Result options	394
Calibration	
History	91

Calibration buffer	453
Calibration buffers	176
Calibration data	89, 234
Limit values	93
Calibration interval	94
Calibration parameters	
CAL pH	368
Calibration position 1	119
Cell	
SET	334
Command	
ADD	411
BEEP	435
CAL	368
CALC	388
CALC LIVE	396
CONTROL RS	423
CTRL	421
DET	292
ELT	374
EMPTY	410
END	435
EVAL BREAK	385
EVAL FIX-EP	381
EVAL MIN/MAX	383
EVAL pK/HNP	382
EVAL RATE	387
Insert	187
LIFT	426
LQH	416
MAT	356
MEAS	363
MET	308
MOVE	424
PREP	409
PUMP	427
RACK	428
REPORT	404
REQUEST	434
SAMPLE	428
SCAN	419
SCAN RS	422
SET	325
SIGN	435
STAT	338
STIR	432
SUBSEQ	429
WAIT	433
Command list	188

Common variable	
Delete	157
Edit	157
Properties	159
Validity	160
Communication	419
Conditioning	
SET	332
Connect	
Balance	20
Barcode reader	21
Dosing device	17
MSB devices	16
PC keyboard	21
Power supply unit	12
Printer	19
Remote Box	18
Sample Processor	22
Stirrer or titration stand	18
USB devices	19
USB hub	19
Control	202
Control instrument	104
Parameters	
305, 321, 334, 351, 358, 365,	
372, 376, 412, 418	
Control parameters	
SET	326
STAT	339
Control Remote Box	99
CONTROL RS	423
CTRL	421
Curve options	244
Curves	243
Custom logo	178

Delete	
Signature	201
Statistics data	204
DET	1, 32, 292
Determination	
Calibration data	234
Cancel	223
Carry out	221
Data	230
Define filters	241
Delete signatures	201
Endpoint details	233

List of saved determinations	239	Tubing parameters	69	F	
Load	238	E		Favorite	206
Measuring point list	233	E-mail	99	Add	207
Message	234	Electrode	85	Delete	207
Properties	235	Parameters		Edit	207
Recalculate	245	305, 321, 334, 351, 358, 366,		FDA Guidance 21 CFR Part 11	
Reevaluate	245	372, 376			47, 58
Save	243	Electrode slope		File	
Save automatically	198	Intervention limit	91	Copy	139
Send PC/LIMS report	198	Warning limit	91	Properties	140
Sign	200	Electrode test parameters	375	Rename	139
Titration data	234	Electrode type		File manager	137
Determination run	221	Custom	179	Filling	278
Device		Electrostatic charge	8	Filling rate	
Add	97	ELT	374	Maximum	444, 445
Balance	129	EMPTY	410	Filter	
Configure	97	Default parameters	444	Determination	241
Delete	96	Emptying	278	Result table	255
Edit	96	Default parameters	444	Fixed endpoint	
New	96	END	435	Evaluation	381
Printer	125	Endpoint details	233	Fixed key	10
Sample Processor	108	ERC	305, 321	Format storage medium	474
Ti-Touch titrator	104	Ethernet connector	30	Formula editor	398
USB/RS-232 adapter	131	EVAL BREAK	385	G	
Device manager	96	EVAL FIX-EP	381	GLP manager	145
Diagnosis	473	EVAL MIN/MAX	383	GLP test "Measurement" ..	147
Dialog language	42	EVAL pK/HNP	382	GLP test "Titration"	147
Dialog options	44	EVAL RATE	387	System validation	151
Dialog settings		Evaluation	380	Test tool	146
System-specific	44	Break point	385	GLP test	
Differential potentiometry	30	DET	304	Dosing unit	76
Directory structure	142	Fixed endpoint	381	Exchange unit	76
Dosing	270	HNP	382	Sensor	91
Continuous	273	Maximum	383	H	
Fixed volume	274	MET	320	Half neutralization potential	
Dosing curve	243	Minimum	383	Evaluation	382
Dosing device		pK value	382	HNP	
Connect	17	Potentiometric	300, 316	Evaluation	382
Parameters		Rate	387	Hysteresis	346
307, 323, 336, 352, 359, 412,		Set windows	301, 317, 386	I	
418		Exchange unit		ID profile	
Dosing key	443	Data	72	Create	52
Dosing parameters		GLP test	76	Identification profile	
ADD	411	Preparing parameters	74	Create	52
MAT	357	Tubing parameters	75	Initialization	436
Dosing port	67	Expiry date		Input line	
Dosing rate		Sensor	88	Template	171
Maximum	444, 445	Titration	65	Installation	
Dosing unit		External position		Set up instrument	12
Data	66	Swing angle	115	Instrument	
GLP test	76			Switch off	37
Preparing parameters	67			Switch on	36
Replace reagent	278				

- pK value
 Evaluation 382
 Potentiometric evaluation
 DET 300
 MET 316
 Power supply unit
 Connect 12
 PREP 409
 Default parameters ... 444, 445
 Parameters for the dosing unit
 67
 Parameters for the exchange
 unit 74
 Preparing 277
 Default parameters ... 444, 445
 Parameters for the dosing unit
 67
 Parameters for the exchange
 unit 74
 Printer 19, 125
 PDF settings 126
 Printing 258
 Troubleshooting 439
 Program version
 Ti-Touch 98
 Properties
 Determination 235
 File 140
 Method 196
 Swing Head 111
 Tower 110
 PUMP 427
- R**
- Rack
 Copy 116
 Delete 116
 Edit 116
 Load 116
 RACK 428
 Rack adjustment 124
 Rack name 109
 Rack offset 119
 Rate evaluation 387
 Reagent 81
 Add 81
 Delete 81
 Edit 81
 Replace 278
 Reagent monitoring 82
 Reagent replacement 83
 Reason 57
 Modification option 57
 Recalculation 245
 Reevaluation 245
 Remote
 Input 447
 Interface 447
 Output 448
 Remote Box
 Connect 18
 Pin assignment 447
 Remote lines
 Manual control 281
 Report 258
 Print 258
 REPORT 404
 Report header 178
 REQUEST 434
 Resolution
 Measured value display 60
 Result 230
 Calculation formula 392
 Limits 393
 Save as a common variable
 161
 Template 167
 Result table 252
 Define filters 255
 Delete automatically 256
 Save 257
 Result variable 393
 As parameter setting 450
 Robotic arm 113
 Max. swing angle 113
 Rotation offset 114
 Swing direction 114
 Swing increment 113
 Swing offset 113
 Swing radius 114
 Swing rate 287, 289, 426
 Rotation increment 111
 Rotation offset 114
 Routine dialog 46
 Standard configuration 46
- S**
- Safety 6
 Safety instructions 7
 SAMPLE 428
 Sample assignment table 163
 Sample data 210
 Modify live 224
 Settings 192
 Sample identification
 Modify designation 192
 Modify the designation 220
 Sample identification list 163
 Sample Processor 282
 Connect 22
 Sample rack 116
 Sample series
 Perform 222
 Sample size
 Fixed 192
 Limit value 195
 Limits 192
 Send from a balance 446
 Sample table 213
 Save as favorite 206
 Troubleshooting 438
 Save
 Determination 243
 Result table 257
 Sample table 215
 SCAN 419
 SCAN RS 422
 Sensor 85
 Add 86
 Calibration data 89
 Calibration interval 94
 Delete 86
 Edit 86
 Expiry date 88
 GLP test 91
 Parameters
 305, 321, 334, 351, 358, 366,
 372, 376
 Working life 88
 Sensor list 85
 Serial number
 Dosing cylinder 67, 73
 Dosing unit 66
 Exchange unit 73
 Service 7, 480
 SET 2, 33, 325
 Set the date 43
 Set the time 43
 Set up instrument 12
 Sign
 Determination 200
 Method 200
 SIGN 435
 Signature
 Delete 201
 Software update 477
 Special beakers
 Edit 121
 Spooler 129
 Start conditions
 DET 292
 MET 308

- SET 325
 STAT 338
 Start options 189
 STAT 2, 33, 338
 Statistics 204, 247
 Delete 204
 Troubleshooting 439
 Status LED 443
 STIR 432
 Stirrer
 Connect 18
 Parameters
 307, 324, 337, 355, 362, 367,
 373, 377, 415
 Stirring 279
 Stirring rate 445
 Stop 443
 Stop conditions
 DET 299
 MET 315
 SET 331
 STAT 343
 Stop criterion
 SET 329
 Stop options 191
 SUBSEQ 429
 Subsequence 429
 Load 430
 Note 432
 Options 431
 Save 430
 Supply voltage 8
 Swing direction 114
 Swing Head
 Properties 111
 Swing increment 113
 Swing offset 113
 Swing radius 114
 Swing rate 287, 289, 426
 Switch off 37
 Switch on 36
 System initialization 436
 System settings 42
 System test 146
 System validation 151
- ## T
- Tandem dosing 353, 360, 413
 Template
 Calibration buffers 176
 Electrode type 179
 Input line 171
 Output line 173
 Report header 178
 Result template 167
 Sample assignment table .. 163
 Sample identification list ... 163
 Templates 163
 Text input 40
 Ti-Touch 98
 Time zone 43
 Titer
 History 79
 Intervention limit 80
 Properties 78
 Warning limit 80
 Titrant 61
 Add 62
 Delete 62
 Edit 62
 Expiry date 65
 Working life 65
 Titrant data 234
 Titrant list 61
 Titration curve 243
 Titration mode
 DET 1
 MAT 2
 MET 2
 SET 2
 STAT 2
 Titration parameters
 DET 295
 MET 311
 SET 330
 STAT 342
- ## U
- Titration stand 18
 Connect 18
 Touch screen 37
 Test 476
 Tower
 Properties 110
 Tubing parameters
 Dosing unit 69
 Exchange unit 75
- ## V
- Uninterrupted dosing
 353, 360, 413
 Update 477
 USB
 Connect devices 19
 USB hub
 Connect 19
 USB/RS-232 adapter 131
 User profile
 Create 52
- ## W
- Validation interval 153
 Validity
 Calibration 94
 Common variable 160
 Titer 78
 Valve disk
 Shift direction 72