Application Note AN-NIR-088

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Fast and easy determination of cannabinoid content

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
3. DS2500CBD

2. FDS2500CBD

R ²	0.959
Standard error of calibration	0.99%
Standard error of cross-validation	1.21%

DS2500 CBD HPLC(2)11NIR

3. HPLCCBD HPLC method.

HPLC	5 min () + 40 min (HPLC)

CONTACT

143-0006 6-1-1
null 9

metrohm.jp@metrohm.jp

DS2500 Liquid Analyzer

Robust near-infrared spectroscopy for quality control, not only in laboratories but also in production environments. The DS2500 Liquid Analyzer is the tried and tested, flexible solution for routine analysis of liquids along the entire production chain. Its robust design makes the DS2500 Liquid Analyzer resistant to dust, moisture and vibrations, which means that it is eminently suited for use in harsh production environments. The DS2500 Liquid Analyzer covers the full spectral range from 400 to 2500 nm, heats samples up to 80° C and is compatible with various disposable vials and quartz cuvettes. The DS2500 Liquid Analyzer is thus adaptable to your individual sample requirements and helps you obtain accurate and reproducible results in less than one minute. The integrated sample holder detection and the self-explanatory Vision Air Software also ensure simple and safe operation by the user. In the case of larger-sized sample quantities, productivity can be considerably increased by using a flow-through cell in combination with a Metrohm sample robot.

8mmDS2500

8 mm

Vision Air 2.0 Complete

Vision Air -

Vision Air Complete

Vision Air :

-
-
- SQL

Vision Air Complete (66072208) :

-
-
-

Vision Air Complete :

- 66072207 (Vision Air Network Complete)
- 66072209 (Vision Air Pharma Complete)
- 66072210 (Vision Air Pharma Network Complete)