

1 Purpose

This document describes the fixed bugs and problems as well as improvements in the following new software version in comparison with the previous version.

- New version: 1.2.55.2
- Old version: 1.2.55.1

2 New features

None.

3 Improvements

Topic	Description
Libraries	Supports the following libraries: ▪ 66071750 MIRA Lib MCRL Organics Library Complete ▪ 66071760 MIRA Lib MCRL Organics-Halogenated ▪ 66071770 MIRA Lib MCRL Organics-Nitrogen Group ▪ 66071780 MIRA Lib MCRL Organics-Oxygen Group ▪ 66071790 MIRA Lib MCRL Organics-Sulfur Group ▪ 66071800 MIRA Lib MCRL Minerals

4 Fixed bugs

None.

5 Known issues

- No support for Arabic language.
- Some language packs have errors or missing strings in the UI.
- One cannot view XTR processed data and baselined data for the same spectrum at the same time. One cannot view baselined data of an XTR processed spectrum.