

1 Purpose

This document describes new functions, the changes, improvements and fixed bugs in the following firmware version:

- 57.1007.0016 for OS version 007

In comparison with the firmware version:

- 57.1007.0015 for OS version 007

2 New functions

The following new functions have been integrated:

- None

3 Improvements

The following improvements were made:

- Optimizations for tests in production

4 Fixed bugs

The following bugs were fixed:

- None

5 Known bugs

The following bugs are known:

- The PDF report also displays parameters that have not been activated.
- Sporadic printing errors