

Pump Unit 772 «Absaugen»

2.772.0020

Schlauchpumpe zum Absaugen von austitrierter Lösung und Suspensionen. Die Version «Absaugen» enthält umfangreiches Zubehör zum Absaugen. Die Pumpe lässt sich über eine Relay Box 731 ansteuern. Sie ist geeignet zur Förderung von Lösemitteln.

Lieferumfang 2.772.0020


Qt.	Order no.	Beschreibung	
1 STK	6.1446.160	Thermometer-Stopfen / B-14/15 / schräger Einsatz Zum Schrägstellen von Einsätzen im Titrierkopf 6.1458.010.	
1 STK	6.1543.170	Absaugspitze / Gewinde M8 Absaugspitze zu den Sample Processoren.	
1 STK	6.1621.000	Kanister / 10 L Als Spül- oder Abfallkanister in automatisierten Systemen.	

1 STK

6.1805.510

PTFE-Schlauch / M8 / 60 cm

Mit Knickschutz.



1 STK

6.1812.000

PTFE-Schlauch 4 mm / 6 mm / 4 m



2 STK

6.1820.050

Adapter, M6 / Olive

Für Pumpschlauch 6.1826.100, zu Peristaltikpumpen.



1 STK

6.1826.100

Pumpschlauch 6.4 mm / 1 m

Als Ersatz des in der 772 Pump Unit und 843 Pump Station (peristaltic) verwendeten Pumpschlauches.



1 STK

6.1828.000

PVDF-Anschlussnippel

Für Kanister 6.1621.000.



1 STK

6.2621.140

Inbusschlüssel 2.5 mm



Optionales Zubehör

Order no.	Beschreibung
6.1826.160	Pumpschlauch Viton Für Peristaltikpumpen.

A black, flexible pump hose is shown coiled into a loose 'P' shape. The hose has a smooth, slightly glossy finish. A small, faint white label is visible on the upper part of the coil, containing the text 'Heraeus' and '6.1826.160'. The hose is set against a plain white background.

