




Solid Autosampler CIC (TEI)


2.0136.0620

The Solid Autosampler CIC (TEI) is used for automated analysis of solid samples using Combustion IC. It can be used in combination with the Boat Introduction Modules (6.07311.010 or 6.07311.020). The sample rack for solids no. 1 (6.07312.010) and a cover for the sample rack are already included.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.0136.0620

Qt.	Order no.	Description
-----	-----------	-------------



Solid sample tray no. 1 for the Solid Autosampler CIC (2.0136.0620) for 20 positions.



Optional accessories

Order no.	Description
-----------	-------------

6.07312.020	Solid sample tray no. 2 Solid sample tray no. 2 for the Solid Autosampler CIC (2.0136.0620) for 20 positions.
-------------	--



6.07312.030	Solid sample tray no. 3 Solid sample tray no. 3 for the Solid Autosampler CIC (2.0136.0620) for 20 positions.
-------------	--

