



GLS Sampler CIC (TEI)


2.0136.0630

The GLS Sampler CIC (TEI) is used to analyze liquid gas or gases under increased pressure using Combustion IC. It can be used in combination with either the manual or the automated Liquid Introduction Module for the Combustion Oven (TEI).

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.0136.0630

Qt.	Order no.	Description
1 PCS	1.0136.0630	GLS Sampler CIC (TEI)



1 PCS

6.07314.010

Needle for GLS Sampler CIC (TEI), pk=3

Replacement needle for sample injection with the GLS Sampler CIC (TEI).



1 PCS

6.07314.020

Needle coupling for GLS Sampler CIC (TEI)

Coupling to connect the injection needle with the sample flow path of the GLS Sampler CIC (TEI).



Optional accessories

Order no.	Description	
6.07314.010	Needle for GLS Sampler CIC (TEI), pk=3 Replacement needle for sample injection with the GLS Sampler CIC (TEI).	
6.07314.020	Needle coupling for GLS Sampler CIC (TEI) Coupling to connect the injection needle with the sample flow path of the GLS Sampler CIC (TEI).	
6.07314.030	Quick connect, Swagelok, for GLS Sampler CIC (TEI) Quick connect for sample cartridge for connection to the GLS Sampler CIC (TEI).	
6.07314.040	Hood for GLS Sampler CIC (TEI) Acrylic glass cover for the GLS Sampler CIC (TEI) which can be connected to a fume cupboard.	