



OMNIS Coulometer with stirrer

2.1018.0020

OMNIS Coulometer with magnetic stirrer for stand-alone operation or as the heart of an OMNIS titration system. The coulometer can be enhanced with additional titration, dosing or coulometer modules.

Coulometry is the ideal method for water content determination in the trace range (10 µg to 10 mg absolute water content) in liquids, solids and gases. In addition, coulometric Karl Fischer titration is an absolute method making titer determination unnecessary.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.1018.0020

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.1018.0020

OMNIS Coulometer with stirrer

OMNIS Coulometer with magnetic stirrer



1 PCS

6.2016.050

Support rod / 300 mm



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.06001.970

Function license Coulometer

Function license for OMNIS Coulometer, includes

- Liquid Handling (LQH)
- Coulometric Karl Fischer titration (KFC)
- Coulometric determination of the bromine index according to ASTM D1492 (BRC)
- MEAS U/T/pH/CONC

MNI
OLE NEW LEVEL OF PERFORM

Optional accessories

Order no.	Description	
6.05406.000	KFC equipment for OMNIS Coulometer with generator electrode with diaphragm	
	Complete accessory set for coulometric Karl Fischer titration.	
6.05406.010	KFC equipment for OMNIS Coulometer with generator electrode without diaphragm	
	Complete accessory set for coulometric Karl Fischer titration.	
2.1009.0010	OMNIS Solvent Module	
	OMNIS Solvent Module for reliable addition of solvent and aspirating off titration vessel content without contact with the reagents. New, improved version resistant against chloroform, chlorobenzene, xylene, and glacial acetic acid (requires OMNIS 4.4.1 or higher). Manual control or actuation via OMNIS Software for maximum user-friendliness. Including accessories, solvent, and waste bottle.	
6.1446.220	Reagent Replacement Solvent Module stopper	
	Stopper for reagent replacement in a coulometric titration cell with a Solvent Module.	

6.02102.020 Cable MDL PL/SO 1 m

Connecting cable for OMNIS instruments with MDL connector (1 m)



6.02102.040 Cable MDL PL/SO 5 m

Connecting cable for OMNIS instruments with MDL connector (5 m)



6.02102.010 Cable MDL PL/SO 0.5 m

Connecting cable for OMNIS instruments with MDL connector (0.5 m)



6.02102.030 Cable MDL PL/SO 2 m

Connecting cable for OMNIS instruments with MDL connector (2 m)



6.02102.050 Cable MDL PL/SO 10 m

Connecting cable for OMNIS instruments with MDL connector (10 m)



6.06001.100 Function license DET

"Dynamic equivalence point titration" function license for one OMNIS titration system.



6.06001.110 Function license MET

"Monotonic equivalence point titration" function license for one OMNIS titration system.



6.06001.160 Function license KFT with conditioning

Function license for volumetric Karl Fischer titration with conditioning for an OMNIS titration system.



6.06001.320 Function license for 5 parallel titrations

Function license for a maximum of five parallel titrations using one OMNIS Titrator with additional Titration Modules or Dosing Modules.

