



OMNIS Coulometer Module without stirrer

2.1019.0010

Coulometer module without magnetic stirrer for connection to an OMNIS Titrator or to an OMNIS Coulometer as enhancement to include a measuring interface.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.1019.0010

Qt.	Order no.	Description
1 PCS	1.1019.0010	OMNIS Coulometer Module without stirrer
		OMNIS Coulometer Module without magnetic stirrer



Optional accessories

Order no.	Description	
6.05406.000	KFC equipment for OMNIS Coulometer with generator electrode with diaphragm	
	Complete accessory set for coulometric Karl Fischer titration.	
6.05406.010	KFC equipment for OMNIS Coulometer with generator electrode without diaphragm	
	Complete accessory set for coulometric Karl Fischer titration.	
2.1009.0010	OMNIS Solvent Module	
	OMNIS Solvent Module for reliable addition of solvent and aspirating off titration vessel content without contact with the reagents. New, improved version resistant against chloroform, chlorobenzene, xylene, and glacial acetic acid (requires OMNIS 4.4.1 or higher). Manual control or actuation via OMNIS Software for maximum user-friendliness. Including accessories, solvent, and waste bottle.	
6.1446.220	Reagent Replacement Solvent Module stopper	
	Stopper for reagent replacement in a coulometric titration cell with a Solvent Module.	

6.02102.010 Cable MDL PL/SO 0.5 m

Connecting cable for OMNIS instruments with MDL connector (0.5 m)



6.02102.020 Cable MDL PL/SO 1 m

Connecting cable for OMNIS instruments with MDL connector (1 m)



6.02102.030 Cable MDL PL/SO 2 m

Connecting cable for OMNIS instruments with MDL connector (2 m)



6.02102.040 Cable MDL PL/SO 5 m

Connecting cable for OMNIS instruments with MDL connector (5 m)



6.02102.050

Cable MDL PL/SO 10 m

Connecting cable for OMNIS instruments with MDL connector (10 m)

