



Polytron PT 1300 D

2.136.0100

Polytron PT 1300 D - Metrohm version

Homogenizer that can be controlled directly by OMNIS Software, tiamo™ or Touch Control.

The Polytron PT 1300 D is comprised of a control unit and a drive. The coupling system on the drive permits simple and rapid replacement of the aggregates without additional tools.

Solid samples can be pulverized without difficulty. However the instrument is also highly suitable for ensuring thorough mixing in viscous samples.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.136.0100

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.136.0100

Polytron PT 1300 D

Polytron PT 1300 D - Metrohm version


Homogenizer that can be controlled directly by tiamo™ or Touch Control.



1 PCS

1.136.0110

Power supply unit for Polytron PT 1300 D



Optional accessories

Order no.	Description	
6.1462.250	Swing arm for Polytron and 786 Swing Head, left swinging Swing arm for Polytron and 786 Swing Head, left swinging, with one buret tip for the addition of solvent or titrant and three spray nozzles for cleaning. The swing arm allows the usage of the Polytron for sample preparation on multi row racks.	
6.2047.110	Holder for Polytron Holder to fix the Polytron at the stand rod.	
6.9012.000	Dispersing aggregate for Polytron, 125 mm Standard dispersing aggregate for Polytron 1300 D	
6.9012.010	Dispersing aggregate for Polytron, 157 mm Special dispersing aggregate (Metrohm version) with protruding knives	

6.02008.010 Polytron holder

For fastening the Polytron in place, suitable for 601405040.



6.02008.030 Holder for homogenizer for Eco Titrations instruments

Holder for homogenizer for fastening the Polytron in place, suitable for 601405060 titration vessel lid

