



728 Magnetic Stirrer with stand

2.728.0040

For the extension of 645, 685, 665/765, 715/775, 725/776 Dosimats to Titration Stands, or supplement for Titrinos and KF Coulometers; incl. 6.2001.010 base plate with support rod and electrode holder; without power adapter, with fixed connecting cable.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.728.0040

Qt.	Order no.	Description
1 PCS	1.728.0010	728 Stirrer without stand For mounting on a stand rod 10 mm dia.; without base plate (or holder) and without stand rod; without power adapter, with fixed connecting cable.



1 PCS

6.1903.020

Stirring bar / 16 mm

Stirring Bar with magnetic core, PTFE coated, length 16 mm



1 PCS

6.1903.030

Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 PCS

6.2001.010

Support stand with support rod

For Dosimats and Titrinos



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2723.130

Cover for 655/665/765

Dust cover for 655, 665, 765 Dosimats

