

733 IC Separation Center

2.733.0120

733 IC Separation Center with two injectors

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.733.0120

Qt.	Order no.	Description	
1 PCS	6.1816.000	Silicone tubing 6 mm, 0.5 m Drainage tube for IC Separation Center chamber	
1 PCS	6.1831.010	PEEK capillary 0.25 mm i.d. / 3 m For all IC components.	
2 PCS	6.2027.030	Column holder 8.5 mm For mounting columns in column chamber of IC instruments.	

2 PCS

6.2027.040

Column holder 11.3 mm

For mounting columns in column chamber of IC instruments.



2 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



2 PCS

6.2744.020

Luer/UNF coupling

For IC instruments



4 PCS

6.2744.040

2 x UNF 10/32 coupling

For connecting 1/16 in. capillaries. For IC instruments



2 PCS


6.2816.020

Syringe 10 mL with Luer connection

For various applications in IC and VA



Optional accessories

Order no.	Description	
6.1825.210	PEEK sample loop 20 µL For injection valve, with 2 PEEK pressure screws	
6.2057.010	Holder for dialysis cell For the dialysis cell in IC column chambers	
6.5904.030	Rotor, 6/2 injector Rotor for the 6/2 injector in IC systems.	
6.5904.050	Injector with 0.25 µL internal loop Flat plate injector with 0.25 µL internal loop	

6.5904.060

Rotor (internal loop, 0.25 µL)

Rotor for 4-port flat plate injector (6.5904.050), 0.25 µL internal loop



6.5904.070

Rotor (internal loop, 1.0 µL)

Rotor for 4-port flat plate injector (6.9956.001), 1.0 µL internal loop

