

1 PCS

6.1817.000

Filling tube for MME

For filling the Multi-Mode-Electrode.




1 PCS

6.2615.030

Electrode stand for MME

Electrode stand for filling and storing the Multi-Mode Electrode.



1 PCS

6.1204.090

Stirrer tip

Stirrer tip for VA stands



1 PCS

6.1226.030

Mercury drop capillary (non-silanized) for MME

Non-silanized mercury drop capillary (10 pieces)




1 PCS

6.1244.020

Driving belt

For Rotating Disk Electrodes and stirrer in VA Stands (3 pieces).



1 PCS

6.1245.010

Electrolyte vessel made of plastic

Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.



1 PCS

6.1247.020

Needles for MME

Needles for Multi Mode Electrode (3 pieces).



1 PCS

6.1415.210

Titration vessel / 10-90 mL



7 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1801.080

Tubing for nitrogen connection

For supplying gas to voltammetry instruments



4 PCS

6.1808.000

Connection bushing M6



1 PCS

6.1824.000

Four-way micro dosing tip

With 4 x M6 thread. For VA measuring stands.



1 PCS

6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)



1 PCS

6.2406.000

Mercury drop catcher

Silver wire (to capture mercury drops by amalgamation)



1 PCS

6.2703.000

Stand ring for titration vessels



1 PCS

6.2709.080

Stopper for pipetting opening

For 694, 747 VA Stands and VA Computrace



1 PCS

6.2711.030

Tray for filling the MME

For 884 Professional VA, 797 VA Computrace, 757 VA Computrace and 647, 663, 694, 747 VA Stands



1 PCS

6.2711.040

Drip pan for VA Computrace



1 PCS

6.2730.030

Stopper

Complete with nipple and O-ring. For 6.1414.030 titration vessel lids. For KF instruments, VA measurement stands, VA Computrace and IC Dilution.



1 PCS

6.2739.000

Wrench

For tightening connectors



1 PCS

6.2816.020

Syringe 10 mL with Luer connection

For various applications in IC and VA



1 PCS

6.2816.030

Needle with Luer connector

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.



Optional accessories

Order no.	Description	
2.772.0110	772 Pump Unit Peristaltic pump for direct connection to Sample Processors.	
2.823.0010	823 Membrane Pump Unit Membrane pump for direct connection to Sample Processors.	
6.0728.000	Ag/AgCl reference electrode without electrolyte Inner reference system without electrolyte, with lock bolt, for use with 6.1245.000 electrolyte vessel made of glass.	
6.1204.130	Silver Electrode Tip Metal electrode tip, silver (Ag). Electrode disk diameter 2 mm. M3 thread. Temperature range: 0 ... 40 °C.	

6.1204.140

Gold Electrode Tip

Metal electrode tip, gold (Au). Electrode disk diameter 2 mm. M3 thread.
Temperature range: 0 ... 40 °C.



6.1204.150

Lateral Gold Electrode Tip

Metal electrode tip, gold (Au), with lateral electrode surface. Electrode disk diameter 3 mm. M3 thread. Temperature range: 0 ... 40 °C.



6.1204.170

3 mm platinum electrode tip for CVS

Metal electrode tip made of platinum (Pt), diameter of electrode disk 3 mm, polished surface, M3 thread connection. Temperature range: 0 - 40 °C.



6.1204.180

Ultra Trace electrode tip

Carbon electrode tip from Ultra Trace graphite. Electrode disk diameter 2 mm. M3 thread connection. Temperature range: 0 - 40 °C.



6.1204.210 Driving axle for rotating disk electrode (RDE)

Driving axle for rotating disk electrode (RDE) with titanium axle for VA measuring stands, M3 thread.



6.1204.220 Driving axle for rotating disk electrode (RDE) with mercury contact

Driving axle for rotating disk electrode (RDE) with titanium axle and mercury contact for VA measuring stands, M3 thread.



6.1241.020 Electrode holder

Electrode holder for VA auxiliary electrodes with lock bolt. For use with glassy carbon rod (6.1247.000) or platinum pin (6.1247.010).



6.1245.000 Electrolyte vessel made of glass

Electrolyte vessel made of glass with ground-joint diaphragm (standard ground-joint) for double junction reference electrode.



6.1247.000

Glassy carbon rod (GC)

For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.



6.1415.150

Titration vessel / 5-70 mL



6.1415.220

Titration vessel / 20-90 mL



6.1415.250

Titration vessel / 50-150 mL



6.1415.310

Titration vessel / 70-200 mL



6.1418.150

Titration vessel with thermostat jacket / 5-70 mL



6.1418.220

Titration vessel with thermostat jacket / 20-90 mL



6.1418.250

Titration vessel with thermostat jacket / 50-150 mL



6.1450.210

PFA titration vessel / 10 to 90 mL

For ultratrace analysis and for solutions containing fluoride.



6.1456.210

Measuring vessel / 10-90 mL

Measuring vessel for automated VA systems. Dilution vessel for Advanced IC Sample Processor



6.2709.040

Stopper for MME opening on VA Stand

For closing gas tubing of the Multi-Mode Electrode in VA Stands if not in use. E.301.0004 O-rings included.



6.2709.050

Stopper for stirrer opening on VA Stand

E.301.0004 O-rings included. Used for the stirrer opening of VA Stands.

