



757 VA Computrace with reduced accessories

2.757.0120

VA Computrace with USB connection and reduced accessories. Integrated PC-controlled voltammetric system for analysis and teaching. Reduced accessories without electrodes for replacement of older VA instruments (663 VA Stand or 647 VA Stand). Includes USB interface 6.5326.000 and PC program 6.6032.100 for Windows 2000 or XP Professional.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.757.0120

Qt.	Order no.	Description
-----	-----------	-------------

7 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1801.080

Tubing for nitrogen connection

For supplying gas to voltammetry instruments



4 PCS

6.1808.000

Connection bushing M6



1 PCS

6.1824.000

Four-way micro dosing tip

With 4 x M6 thread. For VA measuring stands.



1 PCS

6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)



1 PCS

6.2709.080

Stopper for pipetting opening

For 694, 747 VA Stands and VA Computrace



1 PCS

6.2711.040

Drip pan for VA Computrace



1 PCS

6.2730.030

Stopper

Complete with nipple and O-ring. For 6.1414.030 titration vessel lids. For KF instruments, VA measurement stands, VA Computrace and IC Dilution.



1 PCS


6.2739.000

Wrench

For tightening connectors



Optional accessories

Order no.	Description	
6.1817.000	Filling tube for MME For filling the Multi-Mode-Electrode.	
6.2615.030	Electrode stand for MME Electrode stand for filling and storing the Multi-Mode Electrode.	
2.772.0110	772 Pump Unit Peristaltic pump for direct connection to Sample Processors.	
2.823.0010	823 Membrane Pump Unit Membrane pump for direct connection to Sample Processors.	

6.0728.000

Ag/AgCl reference electrode without electrolyte

Inner reference system without electrolyte, with lock bolt, for use with 6.1245.000 electrolyte vessel made of glass.



6.1204.090

Stirrer tip

Stirrer tip for VA stands



6.1204.130

Silver Electrode Tip

Metal electrode tip, silver (Ag). Electrode disk diameter 2 mm. M3 thread. Temperature range: 0 ... 40 °C.



6.1204.140

Gold Electrode Tip

Metal electrode tip, gold (Au). Electrode disk diameter 2 mm. M3 thread. Temperature range: 0 ... 40 °C.



6.1204.150

Lateral Gold Electrode Tip

Metal electrode tip, gold (Au), with lateral electrode surface. Electrode disk diameter 3 mm. M3 thread. Temperature range: 0 ... 40 °C.



6.1204.170

3 mm platinum electrode tip for CVS

Metal electrode tip made of platinum (Pt), diameter of electrode disk 3 mm, polished surface, M3 thread connection. Temperature range: 0 - 40 °C.



6.1204.180

Ultra Trace electrode tip

Carbon electrode tip from Ultra Trace graphite. Electrode disk diameter 2 mm. M3 thread connection. Temperature range: 0 - 40 °C.



6.1204.210

Driving axle for rotating disk electrode (RDE)

Driving axle for rotating disk electrode (RDE) with titanium axle for VA measuring stands, M3 thread.



6.1204.220

Driving axle for rotating disk electrode (RDE) with mercury contact

Driving axle for rotating disk electrode (RDE) with titanium axle and mercury contact for VA measuring stands, M3 thread.



6.1226.030

Mercury drop capillary (non-silanized) for MME

Non-silanized mercury drop capillary (10 pieces)



6.1226.050

Mercury drop capillary (silanized) for MME

Silanized mercury drop capillary (10 pieces)



6.1241.020

Electrode holder

Electrode holder for VA auxiliary electrodes with lock bolt. For use with glassy carbon rod (6.1247.000) or platinum pin (6.1247.010).



6.1245.000

Electrolyte vessel made of glass

Electrolyte vessel made of glass with ground-joint diaphragm (standard ground-joint) for double junction reference electrode.



6.1245.010

Electrolyte vessel made of plastic

Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.



6.1247.000

Glassy carbon rod (GC)

For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.



6.1247.020

Needles for MME

Needles for Multi Mode Electrode (3 pieces).



6.1415.150 Titration vessel / 5-70 mL



6.1415.220 Titration vessel / 20-90 mL



6.1415.250 Titration vessel / 50-150 mL



6.1415.310 Titration vessel / 70-200 mL



6.1418.150

Titration vessel with thermostat jacket / 5-70 mL



6.1418.220

Titration vessel with thermostat jacket / 20-90 mL



6.1418.250

Titration vessel with thermostat jacket / 50-150 mL



6.1450.210

PFA titration vessel / 10 to 90 mL

For ultratrace analysis and for solutions containing fluoride.



6.1456.210

Measuring vessel / 10-90 mL

Measuring vessel for automated VA systems. Dilution vessel for Advanced IC Sample Processor



6.2406.000

Mercury drop catcher

Silver wire (to capture mercury drops by amalgamation)



6.2703.000

Stand ring for titration vessels



6.2709.040

Stopper for MME opening on VA Stand

For closing gas tubing of the Multi-Mode Electrode in VA Stands if not in use. E.301.0004 O-rings included.



6.2709.050

Stopper for stirrer opening on VA Stand

E.301.0004 O-rings included. Used for the stirrer opening of VA Stands.



6.2711.030

Tray for filling the MME

For 884 Professional VA, 797 VA Computrace, 757 VA Computrace and 647, 663, 694, 747 VA Stands



6.2816.020

Syringe 10 mL with Luer connection

For various applications in IC and VA



6.2816.030

Needle with Luer connector

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.

