



772 Pump Unit – aspirate

2.772.0120

Peristaltic pump for direct connection to Sample Processors. The "Aspiration" version contains extensive accessories for the aspiration of titrated-out solutions.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.772.0120

Qt.	Order no.	Description
1 PCS	1.772.0110	772 Pump Unit
		Peristaltic pump for direct connection to Sample Processors.



1 PCS

6.1446.160

Thermometer stopper, B-14/15, oblique insertion

To have the inserts of the 6.1458.010 titration head in a sloping position.



1 PCS

6.1543.170

Aspiration tip / M8 thread

Aspiration tip for the Sample Processors.



1 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



2 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



1 PCS

6.1805.510

PTFE tubing / M8 / 60 cm

With protection against kinks.



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



2 PCS

6.1820.050

M6 / Olive adapter

For 6.1826.100 pump tubing, to peristaltic pumps.



1 PCS

6.1820.070

Coupling 4 mm / 6 mm / M8 inner

For extending aspiration tubing in automated systems. Connects the M8 aspiration tubing directly to the 6.1812.000 PTFE tubing.



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



1 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



Optional accessories

Order no.	Description	
6.1826.160	Viton pump tubing For peristaltic pumps.	
6.2160.010	Adapter cable / 2 m / M8 socket / 2 x plug B Adapter cable / 2 m / M8 socket / 2 x plug B. For 772 Pump Unit.	