



786 Swing Head with transfer head left

2.786.0010

Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.786.0010

Qt.	Order no.	Description
1 PCS	1.786.0010	786 Swing Head with transfer head left
		Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.



1 PCS

6.1462.030

Robotic arm with transfer head for 786 Swing Head,
left swinging

Transfer robotic arm for liquid handling purposes



1 PCS

6.1562.100

Pipetting tubing / 10 mL

Pipetting tubing for fully automatic transfer of sample solutions.



1 PCS

6.1823.010

Guide shaft to pipetting tubing

For assembly in the transfer robotic arm on the Robotic Sample Processors



1 PCS

6.2621.120

Hexagon key 1.5 mm

