



833 Advanced IC Liquid Handling Sample Preparation Unit

2.833.0030

Metrohm inline sample preparation for the analysis of anions. Optimal for inline neutralization of strong bases. Cation exchanger for universal use. Unlimited capacity thanks to offline regeneration and rinsing.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.833.0030

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1803.020

PTFE capillary 0.97 mm i.d. / 5 m

For all IC instruments



2 PCS

6.2128.180

Cable 762/830 - 752/753/754/833

Connecting IC Interface - IC Pump Unit, IC Suppressor Module, IC Dialysis Unit or IC Liquid Handling Unit



2 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS

6.2744.030

Coupling olive/UNF 10/32, 4x

Connection between pressure screw and pump tubing. 4 pieces. For IC instruments with peristaltic pump



1 PCS

6.2744.040

2 x UNF 10/32 coupling

For connecting 1/16 in. capillaries. For IC instruments



2 PCS



6.2821.120

Inline filter 2 µm

Inline filter 2 µm



Optional accessories

Order no.	Description	
6.1803.030	PTFE capillary 0.5 mm i.d. / 3 m Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	
6.1826.060	Pump tubing (orange/yellow) For Compact IC , Personal IC , Basic IC	
6.2755.000	Tubing cartridge For the peristaltic pumps of IC instruments. Contact pressure adjustable.	
6.2821.130	Spare filter for inline filter Spare filter for inline filter.	

6.2835.000

SPM Rotor A

Rotor for inline sample preparation with cation exchange



6.2835.010

Connecting piece to MSM and SPM module

Universal connecting piece for MSM in Professional IC Vario and Compact Flex instruments and for MSM-HC and SPM module. With inlets and outlets.

