



843 VA Peristaltic Pump Station for Compact VA systems

2.843.0140

The 843 Peristaltic Pump Station has two built-in peristaltic pumps that can be controlled directly by the 797 VA Computrace via the remote signal. It is used for the automatic rinsing and emptying of the measuring vessel in partly automated and fully automated VA and CVS systems with the 797 VA Computrace.

The version for voltammetry is equipped with the complete accessories for the automatic rinsing and emptying of the measuring vessel.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.843.0140

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.843.0120

843 Pump Station (peristaltic)

The 843 Pump Station (peristaltic) has two built-in peristaltic pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button.



4 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1456.210

Measuring vessel / 10-90 mL

Measuring vessel for automated VA systems. Dilution vessel for Advanced IC Sample Processor



1 PCS

6.1602.115

Bottle attachment / S 40 / green

For direct dosing from reagent bottles with S 40 thread (Merck...).



2 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



1 PCS

6.1805.020

FEP tubing / M6 / 52 cm

With light and kink protection



1 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection



1 PCS

6.1805.180

FEP tubing / M6 / 16 cm

With light and kink protection



4 PCS 6.1805.530 FEP tubing / M6 / 2 m
With light and kink protection



1 PCS 6.1808.000 Connection bushing M6



2 PCS 6.1819.010 PTFE tube
As gas inlet in washing bottle of VA Stand



4 PCS 6.1820.020 M10/M6 adapter
Connector with M10 thread for 6.1812.XXX tubing and M6 thread for 6.1805.XXX tubing.



4 PCS

6.1820.050

M6 / Olive adapter

For 6.1826.100 pump tubing, to peristaltic pumps.



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



1 PCS

6.1828.020

Five-way M6 tubing connector for canister

For connecting up to 5 tubings



1 PCS

6.1829.020

FEP aspiration tube M6, 0.5 m

For Dosing Unit




1 PCS

6.2141.280

Connection cable 843 Pump Station to 797 VA
Computrace

Cable for the remote connection of the 797 VA Computrace with the 843
Pump Station.



Optional accessories

Order no.	Description	
6.1621.100	Cap for 10-L container Cap for 10-L container. For easy and safe transportation of filled containers.	
6.1801.120	PVC tubing / 4 mm / 6 mm / 2 m	
6.1808.010	Coupling piece 2 X M6 2 M6 outer threads. Part of the screw connections for Exchange Units and 676 Sample Changer	
6.1812.000	PTFE tubing 4 mm / 6 mm / 4 m	

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.

