



NIRS XDS RapidLiquid Module

2.921.0410

Rapid, precise analyses of liquids and suspensions of all types

The NIRS XDS RapidLiquid Analyzer enables rapid, precise analyses of liquid formulations and substances. Precise measurement results at the push of a button make the NIRS XDS RapidLiquid Analyzer an equally reliable and simple solution for quality monitoring in laboratories and processes. The samples are transferred to quartz cuvettes designed for multiple use or disposable glass vials; a tempered sample compartment ensures reproducible analysis conditions and thus accurate measurement results.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.921.0410

Qt.	Order no.	Description
-----	-----------	-------------

Rapid, precise analyses of liquids and suspensions of all types

The NIRS XDS RapidLiquid Analyzer enables rapid, precise analyses of liquid formulations and substances. Precise measurement results at the push of a button make the NIRS XDS RapidLiquid Analyzer an equally reliable and simple solution for quality monitoring in laboratories and processes. The samples are transferred to quartz cuvettes designed for multiple use or disposable glass vials; a tempered sample compartment ensures reproducible analysis conditions and thus accurate measurement results.



Optional accessories

Order no.	Description	
6.7403.030	XDS spacer for disposable 10 mm vials (6.7402.100)	
	Spacer for the disposable glass vials with 10 mm diameter.	
6.7403.020	XDS spacer for disposable 2 mm vials (6.7402.070)	
	Spacer for the disposable glass vials with 2 mm diameter.	
6.7403.000	NIRS XDS spacer for disposable 8 mm vials (6.7402.000)	
	Spacer for disposable glass vials with 8 mm diameter.	
6.7402.070	Disposable vials, 2 mm diameter, transmission	
	200 lockable disposable glass vials (borosilicate) with a diameter of 2 mm for analyses of liquids in transmission mode.	
	Compatible with:	
	• Holder OMNIS NIR, vial, 2 mm (6.07401.050) • DS2500 holder for flow cells (6.7492.000)	

6.7402.100 Disposable vials, 10 mm diameter, transmission

195 lockable disposable glass vials (borosilicate) with a diameter of 10 mm for analyses of liquids in transmission.

Compatible with:

- Holder OMNIS NIR, vial, 10 mm (6.07401.080)
- DS2500 holder for 4 mm disposable vials (6.7492.030)



6.7402.240 Disposable vial, 8 mm, transmission, qty. 100

100 disposable glass vials (borosilicate) with an optical path length of 8 mm for analyses of liquids in transmission. The disposable vials are supplied with the associated stoppers (number of pieces = 100).

Compatible with:

- Holder OMNIS NIR, vial, 8 mm (6.07401.070)
- DS2500 holder for 8 mm disposable vials (6.7492.020)



6.7402.260 Stopper for disposable vial, 4 mm, qty. 100

100 stoppers for 4 mm transmission disposable vials (6.7402.250).

