



Generator electrode with diaphragm

6.00348.100

Generator electrode with diaphragm for coulometric Karl Fischer titrations, SGJ 29/22. This electrode has a cable connection angled upwards, which simplifies cable guiding and operation when used with the OMNIS Sample Robot Oven.

This generator electrode is well-suited for:

- Samples with very low absolute water content 50 ppm
- Coulometric water content determination in oil samples
- Applications where ketone reagents are used
- Applications where low-conducting reagents are used (addition of chloroform, xylene >10%)

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.2104.120

Electrode cable / 1 m / H

With plug H. For connecting generator electrode - KF Coulometer.



6.2104.620

Electrode cable / 0.65 m / H

With plug H, for the generator electrode - coulometer connection

