



Double Au ring electrode

6.00353.100

Glass shaft electrode with two gold rings which are polarized for redox titrations (bivoltammetric titration).

This electrode is well suited for vitamin C determination and can be easily cleaned even from very sticky samples.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.00353.100

Qt.	Order no.	Description
-----	-----------	-------------

10 Stk

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



10 Stk

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Optional accessories

Order no.	Description	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.02104.040	Electrode cable plug-in head G (pol.) / plug P, 0.55 m Cable for connecting an analog polarizable electrode (plug-in head G) to an analog measuring module from OMNIS.	
6.02104.050	Electrode cable plug-in head G (pol.) / plug P, 1.5 m Cable for connecting an analog polarizable electrode (plug-in head G) to an analog measuring module from OMNIS.	