



Ag Titrode with AgBr coating

6.00430.100BR

Combined silver ring electrode with a pH glass membrane as reference electrode.

The silver ring is coated with bromide (AgBr) for a higher sensitivity and a better limit of detection.

This maintenance-free electrode is suitable for precipitation titrations when the pH value remains constant (titrant silver nitrate), for example of:

- chloride, bromide, iodide
- sulfides
- hydrogen sulfide
- mercaptans
- cyanides

This electrode is stored in distilled water.

Below, the accessories are grouped into Scope of delivery and Optional accessories.


Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.00430.100BR

Qt.	Order no.	Description
100 Stk	6.1236.050	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm, polyethylene.
		
100 Stk	6.2008.040	Storage vessel
		Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.
		

Optional accessories

Order no.	Description	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.02104.000	Electrode cable plug-in head G / plug P, 0.55 m Cable for connecting an analog electrode (plug-in head G) to an analog measuring module from OMNIS.	
6.02104.010	Electrode cable plug-in head G / plug P, 1.5 m Cable for connecting an analog electrode (plug-in head G) to an analog measuring module from OMNIS.	