



Combined F-ISE with Pt1000

6.00500.600

Combined fluoride-selective electrode with crystal membrane with integrated Pt1000 temperature sensor.

This ISE is suitable for:

- Ion measurements of F^- (10^{-6} mol/L to sat.)
- Automated ion measurements
- Titrations

Reference electrolyte: $c(KCl) = 3$ mol/L

The electrode is stored in the reference electrolyte.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.00500.600

Qt.	Order no.	Description
-----	-----------	-------------

100 Stk

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



100 Stk


6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Optional accessories

Order no.	Description	
6.2104.600	Electrode cable for plug in head U/plug F, 2x2 mm B, 1m For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).	
6.2104.610	Electrode cable for plug in head U/plug F, 2x2 mm B, 2m For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).	
6.02104.600	Electrode cable plug-in head U / plug P, 0.55 m Cable for connecting an electrode with plug-in head U to an analog measuring module from OMNIS.	
6.02104.610	Electrode cable plug-in head U / plug P, 1.5 m Cable for connecting an electrode with plug-in head U to an analog measuring module from OMNIS.	

6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)

