



Eco titration vessel lid for volumetric Karl Fischer titration

6.01405.050

Titration vessel lid for volumetric Karl Fischer titration with an Eco KF titrator

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.01405.050

Qt.	Order no.	Description
1 PCS	6.1454.130	Form seal made of FPM for 61414xxx for KF and VA instruments



Optional accessories

Order no.	Description	
6.01408.010	Set of seals for Karl Fischer titration cell Complete set of seals, suitable for 6.05609.000 and 6.05609.010.	
6.02709.010	Stopper Plastic stoppers, suitable for 2.1001.4220 OMNIS Karl Fischer Titrator.	
6.02709.020	Septum stopper Plastic septum stopper, suitable for 6014050X0 titration vessel lids.	
6.02709.030	Screw nipple Plastic screw nipple, suitable for 6014050X0 titration vessel lids.	

6.02709.040

Stopper

Stopper for titration vessel lids 6014050X0.



6.02711.000

Spoon for paste

Suitable for use in Karl Fischer water content determinations of pasty substances such as honey, cosmetics, greases or various dairy products.

The spoon for paste facilitates the addition of viscous samples. It is handled like any other spatula and enables convenient, rapid, and loss-free addition of pastes into the Karl Fischer titration vessel.

We furthermore recommend the use of a homogenizer for dissolving viscous samples quickly and completely.

Compatible with 6014050X0 titration vessel lids.

