



Spearhead electrode

6.0226.100

Maintenance-free pH electrode (gel electrolyte) for measurements in semisolid samples such as cheese, meat, or dough.

This electrode is stored in $c(\text{KCl}) = \text{sat.}$ (6.2308.000) and not suitable for ion-deficient solutions.

Thanks to the lifetime indicator, you always know early on when the electrode needs to be replaced.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.0226.100

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Optional accessories

Order no.	Description	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.2308.000	Electrolyte KCl (sat.), 250 mL Electrolyte KCl (sat.), 250 mL	