



Screw cap with septum

6.02743.090

The screw cap with septum is perfectly suitable for maintaining samples and standards over an extended time period. Screw cap with septum for bottle, PP, 50 mL (6.02743.030); PE bottle, 125 ml (6.1627.000) and PE bottle, 50 ml (6.1608.100).

In the case of the 858 Professional Sample Processor, the coated sample needle made of steel (6.2624.200) or the Pick-up needle (6.2624.100) must be mounted in order to penetrate the septum.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.1608.100

PE bottle 50 mL

Bottle for Extension Module Liquid Handling.



6.1627.000

PE bottle, 125 mL for IC

This bottle is used with fully automatic dilution.



6.02743.030

Bottle, PP, 50 mL, 12 pieces

The 50 mL bottles are made of polypropylene (PP) plastic and can be used for samples, standards, check standards and rinsing solutions. Production and packaging both take place under sterile and clean room conditions. This guarantees minimum contamination. The sample vessels are subjected to random checks for leachable anions and cations.



6.2624.200

Sample needle made of steel, coated

Alternative to the sample needle made of zirconium oxide. Particularly suitable for septum vessels. Metal-free thanks to special coating. Requires the 1/16 in. needle holder (6.2833.010) for installation on the Sample Processor.



6.2624.100

Pick-up needle MF

Coated steel needle for metalfree pick-up sampling as well as for Eppendorf sample vials.

