



Screw cap without septum, ID 14.2 mm

6.02743.100

The screw cap without septum can be used in combination with the bottle PP, 50 mL (6.02743.030). The stopper with perforation (6.2743.070) can be slid into the hole and seals the bottle. This provides a perfect sealing option if the standard needle (6.2846.010) is used on the 858 Professional Sample Processor.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.1608.100

PE bottle 50 mL

Bottle for Extension Module Liquid Handling.



6.1627.000

PE bottle, 125 mL for IC

This bottle is used with fully automatic dilution.



6.02743.030

Bottle, PP, 50 mL, 12 pieces

The 50 mL bottles are made of polypropylene (PP) plastic and can be used for samples, standards, check standards and rinsing solutions. Production and packaging both take place under sterile and clean room conditions. This guarantees minimum contamination. The sample vessels are subjected to random checks for leachable anions and cations.



6.2743.070

Stopper with perforation

2000 pieces. For IC Sample Processors



6.2743.077

Stopper with perforation

200 pieces. For IC Sample Processors

