



Adapter Luer inner / M6 outer

6.02744.050

Adapter for syringe for pre-rinsing the aspiration filter (6.2821.090).

The adapter can be screwed onto the aspiration filter. A syringe can be used to rinse the aspiration filter with ultrapure water or eluent.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.2821.090

Aspiration filter

Pore size: 27 µm. Set of 5 pieces. For aspiration tubing (6.1834.000) and filter tubes (6.1821.040) and (6.1821.050).



6.2816.020

Syringe 10 mL with Luer connection

For various applications in IC and VA

