

1 PCS

6.1446.030

Ball stopper / B-14/(15)

For the movable holding of buret tips in titration vessels



2 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



2 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection



1 PCS

6.2008.030

Buret tip storage vessel for 807 Dosing Unit

Together with the 6.2008.050 storage vessel holder, it provides a holder for the buret tip on the 807 Dosing Unit.



1 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders



1 PCS

6.2739.000

Wrench

For tightening connectors



Optional accessories

Order no.	Description	
6.05700.140	Consumable Kit, cylinder unit OMNIS	
	FEP tubing connections M6 (40 cm) and antidiffusion tip for 6.03001.xxx and 6.03004.xxx OMNIS cylinder units.	
6.05700.441	Consumable Kit, cylinder unit OMNIS 50 mL	
	Cylinder unit OMNIS 50 mL complete with cylinder, piston, cylinder base, valve disk, centering tube, and cylinder housing, including programmed data chip.	