

1 PCS

6.01508.210

Cylinder unit OMNIS special 10 mL, without accessories

Intelligent cylinder unit 10 mL for one OMNIS Titrator, Titration Module or Dosing Module. This cylinder unit is especially recommended for the following solutions:

- Aqueous alkali solutions
- Titrant 5
- Silver nitrate solutions
- Nonaqueous alkali solutions
- Permanganate solutions
- EDTA solutions

Without dosing tubing and antidiffusion tip.



1 PCS

6.1446.030

Ball stopper / B-14/(15)

For the movable holding of buret tips in titration vessels



2 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



2 PCS 6.1805.100 FEP tubing / M6 / 40 cm
With light and kink protection



1 PCS 6.2008.030 Buret tip storage vessel for 807 Dosing Unit

Together with the 6.2008.050 storage vessel holder, it provides a holder for the buret tip on the 807 Dosing Unit.



1 PCS 6.2709.070 Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders



1 PCS 6.2739.000 Wrench

For tightening connectors



Optional accessories

Order no.	Description	
6.05700.140	Consumable Kit, cylinder unit OMNIS	
	FEP tubing connections M6 (40 cm) and antidiffusion tip for 6.03001.xxx and 6.03004.xxx OMNIS cylinder units.	
6.05700.460	Consumable Kit, cylinder unit OMNIS special 10 mL	
	Cylinder unit special 10 mL, complete with cylinder, piston, cylinder base, valve disk, centering tube, housing, and data chip. For 6.03004.210 OMNIS cylinder unit.	