




## Separate Pt-wire electrode

6.0301.100

Electrode tip made of Pt wire. This electrode is used as an auxiliary electrode for differential potentiometry.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Scope of delivery 6.0301.100

| Qt.   | Order no.  | Description                             |
|-------|------------|-----------------------------------------|
| 1 PCS | 6.1236.050 | Sleeve with SGJ 14/12 mm                |
|       |            | Sleeve with SGJ 14/12 mm, polyethylene. |



Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



# Optional accessories

| Order no.  | Description                                                                                                      |                                                                                     |
|------------|------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 6.2106.020 | Strand / 1 m / 2 x B<br><br>Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B |  |
| 6.2106.060 | Strand / 2 m / 2 x B<br><br>Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B |  |