



Micro Ag Titrode

6.0433.110

Combined silver ring electrode with a pH glass membrane as reference electrode, narrowed to a diameter of 6.4 mm.

This maintenance-free electrode is well suited for precipitation titrations when the pH value remains constant (titrant silver nitrate) of, e.g.:

- chloride, bromide, iodide
- sulfides
- hydrogen sulfide
- mercaptans
- cyanides

This electrode is stored in deion. water.

Depending on the application, we recommend using an Ag Titrode with an Ag_2S coating, which can be ordered accordingly.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
