



## Micro Au Titrode (installation length 17.8 cm)

6.0435.110

Combined gold ring electrode with a pH glass membrane as reference electrode, installation length = 17.8 cm, narrowed to a diameter of 6.4 mm.

This maintenance-free electrode is well suited for redox titrations when the pH value remains constant, e.g.:

- determination of the chemical oxygen demand (COD, ferrometry)
- determination of penicillin and ampicillin
- for titrations using  $\text{Hg}(\text{NO}_3)_2$
- for titrations in presence of chromium or iron

This electrode is stored in deion. water.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Optional accessories

---

| Order no.  | Description                                                                                                               |
|------------|---------------------------------------------------------------------------------------------------------------------------|
| 6.2104.020 | Electrode cable / 1 m / F<br><br>For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F). |

