



Comb. iAg ring electrode

6.0450.300

Intelligent, combined silver ring electrode with integrated memory chip for storing of sensor data. This electrode is suitable for precipitation titrations when the pH value varies (titrant silver nitrate) of, e.g.:

- Chloride, Bromide, Iodide
- Sulfides
- Hydrogen sulfide
- Mercaptans
- Cyanides

$c(\text{KNO}_3) = 1 \text{ mol/L}$ is used as reference electrolyte and for storage.

iTodes can be used on Titrando, Ti-Touch or 913/914 meters.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.0450.300

Qt.	Order no.	Description
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm Sleeve with SGJ 14/12 mm, polyethylene.
		
1 PCS	6.2008.040	Storage vessel Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.
		

Optional accessories

Order no.	Description	
6.2310.010	Electrolyte KNO3-1mol/L 250 mL Electrolyte solution KNO ₃ 1 M (reference electrolyte for combined silver electrodes)	
2.854.0010	854 iConnect The 854 iConnect – Electrode cable and measuring amplifier for intelligent "iTrodes" electrodes (cable length 150 cm).	