



## NH<sub>3</sub>-selective gas membrane electrode (high conc.)

Combined ammonia-selective gas membrane electrode for concentrations  $>10^{-4}$  mol/L, e.g.:

- in waste water
- in samples containing oil
- for long-term measurements

An individual contaminated membrane can easily be cost-effectively replaced.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Scope of delivery

| Qt. | Order no. | Description |
|-----|-----------|-------------|
|-----|-----------|-------------|

100 Stk

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



---

100 Stk

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



# Optional accessories

| Order no.  | Description                                                                                                                  |                                                                                     |
|------------|------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 6.2316.030 | Inner filling solution to 6.0506.100 electrode<br><br>Measuring electrolyte for NH <sub>3</sub> -selective electrode (50 mL) |  |
| 6.2104.020 | Electrode cable / 1 m / F<br><br>For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).    |  |