



Combined polymer membrane electrode, K

6.0510.110

Combined potassium-selective electrode with polymer membrane.
This ISE is suitable for:

- ion measurements of K^+ (10^{-7} to 1 mol/L) in aqueous solutions
- titrations in aqueous solutions

Thanks to a robust/break-proof plastic shaft made of propylene and an impact protection for the polymer membrane, this sensor is mechanically very resistant.

The reference electrolyte used is $c(CH_3COOLi) = 1$ mol/L.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.0510.110

Qt.	Order no.	Description
-----	-----------	-------------

10 Stk

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



10 Stk

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Optional accessories

Order no.	Description	
6.2328.000	Lithium acetate electrolyte, 1 mol/L (50 mL) Lithium acetate electrolyte, 1 mol/L (reference electrolyte for combined K ⁺ -ISE 6.0510.110)	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	