




IC equipment: Inline Ultrafiltration 2 - push mode

6.05330.110

Accessory set for Inline Ultrafiltration 2 - push mode. For use with the 858 Professional Sample Processor/ 919 IC Autosampler plus.

Recommended set for the combination with Inline Dilution (MIDT).

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.05330.110

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.2744.070

Pressure screw short

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



1 PCS

6.2621.030

Hexagon key 4 mm



1 PCS

6.1831.160

PEEK capillary, 0.5 mm ID / 70 cm



1 PCS

6.2744.034

Coupling olive/UNF 10/32, 2x

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.




2 PCS

6.2744.160

Nozzle/UNF 10/32 coupling with security device

Security device for connecting pump tubing to a nozzle.



2 PCS

6.1803.040

PTFE capillary 0.5 mm i.d. / 1 m

Capillary for sample handling in IC.



1 PCS

6.02057.030

Ultrafiltration cell 2 holder

Holder for mounting the 6.02729.110 ultrafiltration cell 2 on the 858 Professional Sample Processor / 919 IC Autosampler plus.



1 PCS

6.02729.110

Ultrafiltration cell 2

For the Inline ultrafiltration of samples together with the 858 Professional Sample Processor / 919 IC Autosampler plus. As an alternative to the time-consuming manual filtration with a syringe filter.



2 PCS

6.1803.050

PTFE capillary 0.5 mm i.d. / 20 cm

Capillary for inline dialysis and filtration, for Dialysis Unit, IC Filtration Sample Processor



1 PCS

6.1826.390

Pump tubing LFL (yellow/yellow), 3 stoppers

For the sample solution in Inline Ultrafiltration.



1 PCS

6.2714.020

Filtration membrane reg. cellulose, 0.2 µm, 50 pieces

Nominal pore size: 0.2 µm. Set of 50 pieces. For inline filtration.



1 PCS

6.2744.000

Pressure screw PVDF

With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)



1 PCS

6.1826.320

Pump tubing LFL (orange/yellow) , 3-stopper

For suppressor solutions and for acceptor solutions for inline dialysis and for inline ultrafiltration.



1 PCS

6.2831.010

Plastic tweezers

For dialysis and filtration membranes



Optional accessories

Order no.	Description
6.5330.120	IC equipment: Inline dilution
	Accessory set for Inline Dilution. For use with the 858 Professional Sample Processor, 800 Dosino and 741 Magnetic Stirrer.

