




OMNIS IC-Ausr.: MIDT

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery

Order Number	Description	Quantity
6.02744.010	flanschlose Verschraubung M6, 1/8", 2 Stk.	5 PCS
		
6.01815.000	Spiralband, 3-15 mm, 0.5 m	1 PCS
		
6.05000.010	Accessory kit: Bottle OMNIS IC, PP, 2 L OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.	1 PCS
		

6.2744.070

Pressure screw short

2 PCS

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



6.2744.320

PEEK coupling M6 / UNF 10/32

2 PCS

Connecting transfer tubing (M6) with capillaries (UNF 10/32).



6.1831.050

PEEK capillary 0.5 mm i.d., 40 cm

1 PCS

For inline dialysis and for sample changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit, Compact VA Autosampler



6.2744.200

Adapter UNF 10/32 outer / M6 inner

2 PCS


Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)



6.01532.140

LQH-Transferschlauch, FEP, AD 1/8


1 PCS



6.01532.130

Transferschlauch 5 mL, FEP, AD 1/8


1 PCS



6.01805.400

FEP-Schlauch M6, 1/8" x 1.59 mm, 3 m


2 PCS



6.01831.230

PEEK-Kapillare, ID 0.75 mm, 100 cm

1 PCS



6.1831.180

PEEK Capillary 0.5 mm i.d., 3 m

1 PCS



6.1903.000

Stirring bar / 8 x 3 mm

1 PCS

Stirring bar with magnetic core, PTFE covering, 8 mm long with 3 mm diameter.

