



OMNIS IC-Ausr.: AnCat Switch

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery

Order Number	Description	Quantity
6.02744.010	flanschlose Verschraubung M6, 1/8", 2 Stk.	2 PCS
		
6.01806.000	PEEK capillary, 0.18 mm ID / 2 m For use in the high-pressure flow path, e.g. as preheating capillary for Metrosep Carb 2 columns in the 2 mm version (6.01090.2x0).	1 PCS
		
6.05000.010	Accessory kit: Bottle OMNIS IC, PP, 2 L OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.	1 PCS
		

6.2744.070

Pressure screw short

1 PCS


Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



6.01831.170

PEEK-Kapillare INJ - COL, 540/0.18 mm

1 PCS



6.2744.040

2 x UNF 10/32 coupling

1 PCS

For connecting 1/16 in. capillaries. For IC instruments



6.2744.014

Pressure screw 2x

1 PCS


With UNF 10/32 connection. For the connection of PEEK capillaries



6.01805.400

FEP-Schlauch M6, 1/8" x 1.59 mm, 3 m


2 PCS



6.02845.000

Vorwärmkartusche, ID 0.18 mm

1 PCS



6.1831.180

PEEK Capillary 0.5 mm i.d., 3 m

1 PCS

