



Consumable Kit OMNIS Pump Module (2-channel)

6.05700.010

Two pump tubings each and associated press clamps for replacement in the OMNIS pump module. One Consumable Kit can be used for re-equipping two peristaltic pumps of the pump module.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 6.05700.010

Qt.	Order no.	Description
-----	-----------	-------------

2 PCS

6.01801.020

Pharm-A-line pump tubing to the peristaltic pump module

The pump tubing comes pre-equipped with adapters that permit direct installation in the peristaltic pump of the pump module. The Pharm-A-line pump tubing is used for the transport of both aqueous and nonaqueous solutions and must be operated in connection with the white press clamp in the pump module.



2 PCS

6.02703.000

Press clamp, white, for "Peristaltic" pump module

The press clamp is used to fix and secure the pump tubing in the "Peristaltic" pump module.

