



Conductivity measuring cell $c = 0.9 \text{ cm}^{-1}$

6.0910.120

Conductivity measuring cell with cell constant $c = 0.9 \text{ cm}^{-1}$, without integrated temperature sensor. This sensor is suitable for measurements of medium conductivities (1 $\mu\text{S}/\text{cm}$ to 100 mS/cm), e.g., in:

- drinking water
- surface water
- wastewater

This sensor is suitable for automated applications (includes an SGJ 14/15).

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.2301.060 KCl ion standard 250 mL

KCl ion standard ($c(\text{ion}) = 0.1000 \pm 0.0005 \text{ mol/L}$), can also be used as standard solution for conductometry 12.87 mS/cm (25°C).



6.2324.010 Conductivity standard $100 \mu\text{S/cm}$ 250 mL

Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$.



6.2324.110 Conductivity standard $100 \mu\text{S/cm}$, 5 x 30 mL

Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$.



6.2104.020 Electrode cable / 1 m / F

For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).



6.2104.030

Electrode Cable 2 m / F

For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).

