




Multi-Mode Electrode pro

6.1246.120

Mercury electrode for voltammetry. Can be operated as DME, SMDE or HMDE.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1226.030	Mercury drop capillary (non-silanized) for MME
	Non-silanized mercury drop capillary (10 pieces)



6.1226.050 Mercury drop capillary (silanized) for MME

Silanized mercury drop capillary (10 pieces)



6.1247.020 Needles for MME

Needles for Multi Mode Electrode (3 pieces).



6.1247.110 Screw cap for the mercury refill opening for MME pro



6.1247.120 Retaining screw for capillary for MME pro



6.1247.140

Needle holder for MME pro



6.1247.150

Silicone seal, green 13.9 mm for MME pro



6.1247.210

Assembly ring for needle holder for MME



6.1247.220


Wrench for MME



6.2615.130

Electrode stand for MME pro

Electrode stand for filling and storing the Multi-Mode Electrode pro.



6.2112.130

Adapter plug contact 5.5 mm / screw contact, gold-plated

For VA electrodes

