



KF titration vessel lid with 5 openings

6.1414.030

For volumetric KF titrations.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1414.030

Qt.	Order no.	Description
1 PCS	6.1454.130	Form seal made of FPM for 61414xxx
		for KF and VA instruments



Optional accessories

Order no.	Description	
6.1244.040	Set of seals	
	Set of seals for 6.1414.030. The set is comprised of one seal and various O-rings.	
6.2730.010	Screw nipple	
	Made of PP. For 6.1448.010 Septum or for mounting 6.1403.040 Drying tube or 6.0338.100 Electrode in 6.1414.030 and 6.1414.040 KF Titration vessel lids. Used for KF titrations.	
6.2730.020	Septum stopper	
	With E.301.0041 O-ring. For 6.1448.010 Septum and 6.1414.030 Titration vessel lid. Used for KF applications.	