



Sample beaker, plastic (PP), 120 mL, 1 piece

6.1459.307

Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 1 piece

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1459.307

Qt.	Order no.	Description
0 PCS	6.1459.300	Sample beaker, plastic (PP), 120 mL, 250 pieces
		Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 250 pieces



Optional accessories

Order no.	Description	
6.1459.320	Screw cap for 61459300 / 250 pieces For safe closure of 120 mL sample beakers	
6.1459.300	Sample beaker, plastic (PP), 120 mL, 250 pieces Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 250 pieces	
6.2041.470	Sample rack, 22 x 120 mL, (PP) Sample rack for 810 Sample Changer, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 22 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307. Plastic: Polypropylene (PP)	
6.2041.840	Sample rack, 59 x 120 mL, (PVC) Sample rack for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic", suitable for 59 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307. Plastic: Polyvinyl chloride (PVC)	

6.2041.850 Sample rack, 59 x 120 mL, (PP)

Sample rack for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic", suitable for 59 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.


Plastic: Polyvinyl chloride (PVC)



6.2068.010 Sample rack, 20 x 120 mL, (PP)

Sample rack for 864 Balance Sample Processor, suitable for 20 sample beakers. The following sample beakers can be used: 6.1459.300.

Plastic: Polypropylene (PP)



6.2068.020 Sample rack, 24 x 120 mL, (PVC)

Sample rack for 815 Robotic Soliprep, suitable for 24 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.

Plastic: Polyvinyl chloride (PVC)



6.2068.030 Sample rack, 10 x 250 mL, (PVC)

Sample rack for 815 Robotic Soliprep, suitable for 10 sample beakers. The following sample beakers can be used: 6.1432.320, 6.1432.323, 6.1453.220, 6.1453.250.



6.9920.141

Sample rack, 53 x 120 mL, (PP)

Sample rack for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic", suitable for 53 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.

Plastic: Polypropylene (PP)



6.1458.040

Titration Head, 3 x SGJ 14

Titration head for 120 mL beakers. With 3 openings SGJ 14 and 4 openings 6.4 mm, 2 dosing tips included. Suitable for rinsing and aspiration equipment.



6.2037.080

Foil holder

Foil holder for 120 mL PP beaker (6.1459.300)



6.1909.050

Stirring propeller / 102 mm

Propeller for 802 Stirrer when this is used with sample racks for 120 mL beakers (6.1459.300)



6.1909.060

Stirring propeller / 95 mm

Stirring propeller with 3 wings for intensive stirring

