



Transfer tubing 15 mL, 2 x M6

6.1562.160

For liquid handling in ion chromatography. To avoid contamination of the Dosino.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1562.160

Qt.	Order no.	Description
2 PCS	6.2744.320	PEEK coupling M6 / UNF 10/32

Connecting transfer tubing (M6) with capillaries (UNF 10/32).

