



Pipette adapter 1-10 mL for transfer head

6.1808.250

Adapter with outside thread M6 and Eppendorf conus to use with transfer head 6.1462.240 for usage of 10 mL Eppendorf tips

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1808.250

Qt.	Order no.	Description
1 PCS	6.1808.000	Connection bushing M6



1 PCS

6.1831.040

PEEK capillary 0.5 mm i.d., 15 cm

For inline dialysis. For IC Dialysis Unit, IC Dialysis Sample Processor, IC
Liquid Handling Dialysis Unit



Optional accessories

Order no.	Description
6.1562.240	Pipetting tips 1-10 mL (200 pieces) Eppendorf tips for use with 864 Robotic Balance Sample Processor

